

Cassette Cross Reference

Which Hu-Friedy and Zirc cassettes fit each Dentsonic cleaner. All dimensions L × W × H, in inches.

MODEL	UC-250	UC-300	UC-400	UC-600	UC-750	UC-1000
Internal basket size	9¾ × 4 × 3	9¾ × 6¾ × 3	7¼ × 4¼ × 2⅝	10⅝ × 7½ × 3⅞	10¾ × 8⅞ × 4⅞	9½ × 8¼ × 4¼
Stainless steel basket	—	—	—	10⅞ × 6½ × 3.9	11½ × 8 × 2⅝	—

01 HU-FRIEDY

Large Cassette, 16 inst.	8 × 11 × 1.25	—	✓	—	✓	✓	✓
Large Cassette, 20 inst.	4.5 × 8 × 1	—	✓	—	✓	✓	✓
Small Cassette, 8 inst.	5.5 × 8 × 1.25	✓	✓	✓	✓	✓	✓
Small Cassette, 10 inst.	5.5 × 8 × 1.25	✓	✓	✓	✓	✓	✓
Exam Cassette, 5 inst.	3 × 8 × 1.25	✓	✓	✓	✓	✓	✓
Space Saver, 7 inst.	4.5 × 8 × 1	✓	✓	✓	✓	✓	✓
Double Decker, 8 inst.	3 × 8 × 1.5	✓	✓	✓	✓	✓	✓
Double Decker, 14 inst.	4.5 × 8 × 1.5	✓	✓	✓	✓	✓	✓
Double Decker, 20 inst.	6 × 9.7 × 1.5	—	✓	✓	✓	✓	✓
Signal Stat Cassette, 12 inst.	6.5 × 10.5 × 1.25	—	—	—	✓	✓	✓
Oral Surgery Cassette, 12 inst.	8 × 14.5 × 1.25	—	—	—	—	—	—
Radiography Cassette	4.5 × 8 × 1.5	—	✓	✓	✓	✓	✓
Handpiece Cassette	3 × 8 × 1.25	✓	✓	✓	✓	✓	✓
Ultrasonic Insert Cassette	1.5 × 8 × 1.25	✓	✓	✓	✓	✓	✓

02 ZIRC

E-Z Jett Cassette, 10-place	7.93 × 7 × 1.125	—	—	—	✓	✓	✓
E-Z Jett Cassette, 8-place	7.93 × 4.44 × 1.125	✓	✓	✓	✓	✓	✓
E-Z Jett Cassette, 5-place	8 × 3.125 × 1.19	✓	✓	✓	✓	✓	✓
Compact Cassette	7.125 × 3.875 × 0.625	✓	✓	✓	✓	✓	✓
Mini Mat	3.875 × 3.5 × 0.5	✓	✓	✓	✓	✓	✓
Instrument Mat	6.75 × 2 × 0.375	✓	✓	✓	✓	✓	✓
Standard Steri-Container	8.125 × 1.875 × 1.875	✓	✓	✓	✓	✓	✓
Compact Steri-Container	7.125 × 1.5 × 1.5	✓	✓	✓	✓	✓	✓

WARNING · NO METAL-TO-METAL CONTACT

There must be **no** metal-to-metal contact inside the tank during the ultrasonic cycle: no loose instrument-to-instrument, cassette-to-tank, cassette-to-cassette or cassette-to-basket contact. If no basket is used, make sure cassettes have standard silicone insulation and feet.

Basket sizes vary by unit. Questions? Contact your cassette manufacturer or TPC at (800) 560-8222.