

TPC 3D



Disinfection and Sterilization Guide of Scanner Tips

Intraoral Scanner User Guide

Everyone with a healthy and beautiful smile

3.1.1 Disinfection

Step 1: Prepare sealed bags, gauze, and CaviWipes, and prepare to clean the scanner tip



Step 2: Rinse the scanner tip



Step 3: Dry the scanner tip with CaviWipes



Step 4: Blow dry the scanner tip



Step 5: Fill the scanner tip with gauze



3.1.2 High-temperature and high-pressure sterilization

Step 1: Place the scanner tip filled with gauze into a sterilized bag



Step 2: Seal the aseptic bag



Step 3: Put it into the sterilization cabinet



Step 4: Choose sterilization at 134°C for 4 minutes

