

Model LED 50N

#### Features

- Low battery signal – visual and audible alarm
- Easy to read – LCD display
- High light intensity – 5 to 10 second cure
- Flexibility – Cordless light cure
- Quiet – No fan
- Sleep mode – conserves power, saves battery life

#### Specifications and Accessories

##### Switching Power Adapter

Power supply: AC 100-240V / 50-60 Hz  
Dimensions: 10"L x 2"W x 8"H




##### Handpiece

Light source: 5W light emitting diode  
Wavelength: 430-490 nm  
Light Intensity: 800 - 1000 mW/cm<sup>2</sup>  
Battery: Lithium ion  
Warranty: 2 years



#### Standard Accessories

- LED 50N Handpiece
- Charging base
- 8mm light guide
- Eye protector
- Lithium ion battery
- Power supply
- User's Manual



Shipping Dimensions: 10"L x 7"W x 6"H / 3 lbs.  
(8/cs.) 21"L x 14"W x 10"H / 16 lbs.