



The PC2635 is a portable, self contained delivery system suitable for use in-office, nursing homes, hospitals, or wherever any on-site dentistry is required. The PC2635 is a lightweight, portable system that is easily transported. The air, water and waste management are all included with the unit. The PC2635 is compatible to common brands of handpieces.

Package includes the following:

- 4 hole or 2 hole instrument tubings-specify ✓
- 3 way syringe ✓
- Self contained water bottle system ✓
- Self contained 3/4HP oil-less compressor (580W / 7L) ✓
- Foot control ✓
- Automatic handpiece flush system ✓
- Built-in carrying handle ✓

Shipping Dimensions:

30" L x 22" W x 26"H , 67 lbs