

Model LED 70N



Features

- Low battery signal – visual and audible alarm
- Easy to read – LCD display
- High light intensity – 5 to 10 second cure
- No down time
- Flexibility – corded and cordless light cure
- Quiet fan – extends life of LED
- Sleep mode – conserves power, saves battery life



Specifications and Accessories

Switching Power Adapter

Power supply: AC 100-240V / 50-60 Hz
Dimensions: 10"L x 2"W x 8"H
Weight: 24 ounces

Handpiece

Light source: 10W light emitting diode
Wavelength: 430-490 nm
Light Intensity: 2,200 mW/cm²
Battery: Lithium ion
Warranty: 2 years



Standard Accessories

- LED 70 Handpiece
- Charging base
- 8mm light guide
- Eye protector
- Lithium ion battery
- Power supply
- User's Manual

Shipping Dimensions: 10"L x 7"W x 6"H / 3 lbs.
(8/cs.) 21"L x 14"W x 10"H / 16 lbs.

