

Model 3000 Orthodontic Chair

- Cast aluminum with steel reinforcement
- Hydraulic base & back
- Synchronized back & seat movement
- Adjustable backrest speed
- 2 programmable pre-positions
- Trendelenburg position
- 60° horizontal swivel base
- Double articulating headrest
- Trip pan safety shut-off
- Automatic return/Last position presets
- Foot Control

Dental Light

Specifications:

Halogen lamp:	12 Volts, 55 Watts
Input voltage:	110 Volts, 60 Hz
Output voltage:	High 12 Volts; Low 10.5 Volts
Color Temperature:	4500° K
Bulb type:	H3
Light Intensity (Lux):	High 24,000; Low 12,000
Focal length:	27" (70 cm)
Illumination area:	3" x 7" at focal length



Model #MOP3000

Model OMC2375 Orthodontic Cabinet

Standard Cabinet Features:

- Five easy to remove and clean full extension drawers, and one pull-out writing surface
- Lockable drawers
- High quality hardware throughout
- Oversized caster wheels
- Four stationary handles for easy cart positioning
- A removable countertop plastic tray
- Four-outlet power strip
- Color: white
- 1 Year warranty

Optional Items:

- Piezo scaler
- LED curing light
- Fiber optic handpiece control

Standard Control Features:

- Automatic selection for one high speed and one low speed handpiece
- 3-way syringe
- Vacuum package with saliva ejector and HVE
- Self contained water bottle system
- Variable speed disc-type foot control
- Water coolant on/off switch with flow control
- Master control switch
- Air pressure regulator

Packing Dimensions (L x W x H)	Weight
Chair: 54" x 34" x 31"	310 lbs.
Cushion: 51" x 9" x 23"	24 lbs.
Light: 37" x 17" x 15"	28 lbs.
Cabinet: 26" x 22" x 34"	115 lbs.

