



Model #P555

Model #P557

EZ Nozzle

- Carefully designed length, angle and curve of nozzle for applications to both anterior and posterior teeth.
- Allows for accurate direction of the power jet
- Includes two prophy nozzles
- Handpiece nozzle is fully autoclavable up to 134 °C / 273 °F

Ergo Comfort Design

- Fits comfortably in large and small hands. Notched grip allows firm and accurate tip control

No Installation Cost

- Simply connects to handpiece tubing
- No extra foot control needed
- No external air or water connections needed

Cleaning Power

- Jet Polishing powder (130g) has a refreshing lemon scent

Powder Holder Covers

- Three translucent color coded holder covers included (blue, green, purple)

Quick Connection

- Simply pull back the quick release ring to release the nozzle from the handpiece

Packing dimensions: 12"L x 9"W x 4"H / 2 lbs.
(10/cs) 27"L x 12"W x 8"H / 13 lbs.



Model #P560



Package includes the following

- 1 – 130g of prophy powder
- 3 – Color coded (red, blue, green) powder holder covers
- 2 – Nozzles (360° swivel)
- 1 – Cleaning tool
- 1 – O-ring set

