




TPC



INNOVATION • RELIABILITY • QUALITY
2026-2027 Catalog U.S.A.

Large Equipment

tpcdental.com

OUR VISION

A Preferred Company

To create a quality environment where people can work together to provide simple, cost effective creative solutions for the betterment of dentistry.

OUR MISSION

Continuous Improvement

We will manufacture innovative, high quality and cost effective dental products that enable distribution partners, dentists, hygienists, teaching institutions and government organizations throughout the world to perform efficient quality dental care. TPC provides innovative solutions and quality equipment that enables dental professionals to work together for the benefit of dentistry.

TABLE OF CONTENTS

Our Vision & Our Mission	Inside Front Cover
<u>Terms & Conditions</u>	3
<u>Chair, Unit, Light Packages</u>	4-11
<u>Chairs & Delivery Units</u>	12-34
<u>Air Compressors / Wet Ring Vacuum Systems</u>	35-37
<u>Mobile & Portable Dental Systems</u>	38-43
<u>Extra Oral Suction System</u>	44
<u>Stools & X-Ray Chairs</u>	45-48
<u>Operatory Lights</u>	49-53
<u>Mobile Cabinets</u>	54
<u>Imaging Products & Accessories</u>	55-63
<u>Dental Handpiece & Parts</u>	64-77
<u>Handpiece Cleaning & Lubrication Systems</u>	78
<u>Curing & Bleaching</u>	79-87
<u>Scaling Products & Accessories</u>	88-94
<u>Ultrasonic Bath</u>	95-98
<u>Digital Amalgamator</u>	99
<u>Micro Motor & Handpiece</u>	100
<u>Digital Thermometer</u>	101
<u>Deluxe Autoclavable Valves</u>	102
<u>Parts & Components</u>	103-109
<u>Custom Length Umbilical / Utility Centers</u>	110
<u>Order Form</u>	111
<u>Corporate Staff Contacts</u>	112
<u>Upholstery Color Chart</u>	113



Ordering Information:

When placing orders, please include the following information. All orders will be processed upon receipt of a faxed, mailed, or emailed hardcopy of the purchase order.

- Model and/or Part Number
- Requested Ship Date
- Description / Voltage
- Purchase Order Number
- Quantity
- Shipping Instructions
- Cancellations must be confirmed in writing
- All credit cards are accepted unless otherwise noted

Returned Goods:

All returns require prior authorization. Returns must be authorized by TPC within 7 days of date received by the customer.

- Contact Customer Service for Return Authorization Number
- Returned products are subject to a 25% restocking charge
- All returned goods must have Return Authorization Number
- All returned goods must be shipped prepaid
- All authorized returns must be received by TPC in original packaging, accessories, documentation and in original condition for credit considerations
- All outbound repair / warranty freight charges are the responsibility of the customer.
- **No returns accepted after 14 days.**

Special Orders:

- Special order items and custom made products are not returnable for credit

Minimum Orders:

- Orders less than \$50.00 will incur an additional \$5.00 handling charge.
- All drop ship orders will incur an additional \$5.00 fee.

Terms of Sale:

- Claims: All claims for merchandise received that are not in accordance with the original purchase order must be made within five (5) working days of the receipt of the merchandise. All merchandise is quality inspected and packaged in perfect condition. The common carrier assumes full responsibility for the safe delivery of all merchandise after it is removed from the F.O.B. shipping point. All claims for shortages, loss or damage must be documented and filed promptly with the freight carrier. Concealed damage must be reported to the freight carrier immediately upon discovery, with a request for inspection. If a request is made more than 15 days after the receipt of the merchandise, reasonable evidence must be submitted indicating that damage could not have occurred after delivery.
- All pricing in catalog is subject to change. Verify all pricing prior to placing order.

Warranty Information (All stated warranties are intended for USA market ONLY):

All of our products sold are guaranteed to be free from defects in workmanship and materials for one year from date of purchase, unless otherwise stated. TPC will repair or replace any defective part at no charge. TPC will not be responsible for labor charges or shipping charges to / from the TPC facility. This guarantee does not cover normal wear or stains on surface finish. The guarantee does not cover damage resulting from improper installation, misuse or accidents incurred in shipping and handling. All claims against the freight carrier must be initiated at the time the damaged items are received. The claim is the responsibility of the customer. We are improving our products on a continuous basis. We reserve the right to make modifications without the need for prior notification.

Products:

Patient Chairs

Curing & Bleaching Systems

Imaging Products

Scaling Products

Air Compressor

Warranty Period:

5 year limited

2 years

1 year

1 year

5 year limited

Products:

Delivery Units

Seating

Mobile Cabinets

Warranty Period:

5 year limited

1 year

1 year

Extended Manufacturer Warranty (selective models):

2 Year Warranty

H6000, H6005, H6025, H7005 Lubrication System

LED Curing Lights (All Models)

Advance CAM Intraoral Camera System(s)

LED Operatory Lights

Limited Manufacturer Warranty (selective models):

6 Months Warranty

Fiber optic systems, S333, M333, QS222, QM222,

ELCA760, ELS736.

Superb Vac

Wet-Ring Vacuum



Model	WV3	WV3R	WV5	WV5R
Water consumption without water recycling:	4L/min	/	4L/min	/
Water consumption with TPC water recycling:	/	0.8L/min	/	0.8L/min
Suction Rate	900L/min	900L/min	1500L/min	1500L/min
Vacuum Level (1N Hg) Factory set to 10HG	8-15HG	8-15HG	8-15HG	8-15HG
Noise Level dB (A)	58	58	58	58
Horse Power	1.2HP	1.2HP	2.0HP	2.0HP
Voltage	110V or 220V	110V or 220V	220V	220V
Frequency	50/60Hz	50/60Hz	50/60Hz	50/60Hz
Max.users	3	3	5	5
Current (A)	8A	8A	13A	13A
Circuit Breaker (A)	20A	20A	20A	20A
Inlet Connection Size	3/4"	3/4"	3/4"	3/4"
Drain Connection Size	1.25"	1.25"	1.25"	1.25"
City Water Connection Size	3/8"	3/8"	3/8"	3/8"
Air Exhaust Connection Size	2 3/8"	2 3/8"	2 3/8"	2 3/8"
Inlet Water Pressure (PSI)	10Psi/0.7Bar	10Psi/0.7Bar	10Psi/0.7Bar	10Psi/0.7Bar
Ambient Temperature	2 C ~40 C	2 C ~40 C	2 C ~40 C	2 C ~40 C
Product Size (LxWxH) Inches	17x14x18	17x14x18	17x14x18	17x14x18
Product sit Weight	62	62	66	66
Dimension	17x18x20	21x23x23	17x18x20	21x23x23



WV3R/WV5R

Features:

- Durable aluminum-bronze construction
- Stainless steel impeller
- Computer balanced, sealed bearings in motor
- Large capacity, easy to clean filter with clear bowl for easy inspection
- Available in single and twin models to accommodate 1-5 users
- 5-year full product limited warranty
- Unique water recycling to save up to 80% water
- UL/cUL classified

Description	Model	Suggested Retail Price
Superb Vac 1.2 HP	WV3	\$3,190.00
Superb Vac 1.2 HP (with recycler)	WV3R	\$3,590.00
Superb Vac 2 HP	WV5	\$3,590.00
Superb Vac 2 HP (with recycler)	WV5R	\$3,990.00

Product includes the following:

3 Simultaneous users	✓	✓	✓	✓
5 Simultaneous users			✓	✓
Water recycler & Air/Water Separator		✓		✓



WV3/WV5

Description	Model	Suggested Retail Price
#WV-AWS	Single Pump Air/Water Separator	\$250.00
#CP-3	Control Panel – 3 Button- Vac / Air / Water	\$480.00
#CP-4	Control Panel – 4 Button- Vac1 / Vac2 / Air / Water	\$560.00
#CV22	Air Compressor and Vacuum Pump Stacking Rack	\$900.00
#CP-2 AV	Control panel - 2 Button-Air/Vac	\$400.00
#CP-3 AV	Control panel - 3 Button-Air/Vac1/Vac2	\$480.00



Superb Vac

Wet-Ring Vacuum



WV6(R)/WV10(R)

Model	WV6	WV6R	WV10	WV10R
Water consumption without water recycling:	6 L/min	/	6 L/min	/
Water consumption with TPC water recycling:	/	1.01L/min	/	1.01L/min
Suction Rate	1800L/min	1800L/min	3000L/min	3000L/min
Vacuum Level (1N Hg) Factory set to 10HG	8-15HG	8-15HG	8-15HG	8-15HG
Noise Level dB (A)	65	65	65	65
Horse Power	2.4HP	2.4HP	4.0HP	4.0HP
Voltage	220V	220V	220V	220V
Frequency	50/60Hz	50/60Hz	50/60Hz	50/60Hz
Max.users	6	6	10	10
Current (A)	8A x 2	8A x 2	13A x 2	13A x 2
Circuit Breaker (A)	20A x 2	20A x 2	20A x 2	20A x 2
Inlet Connection Size	1"	1"	1"	1"
Drain Connection Size	1.25"	1.25"	1.25"	1.25"
City Water Connection Size	3/8"	3/8"	3/8"	3/8"
Air Exhaust Connection Size	2 3/8"	2 3/8"	2 3/8"	2 3/8"
Inlet Water Pressure (PSI)	10Psi/0.7Bar	10Psi/0.7Bar	10Psi/0.7Bar	10Psi/0.7Bar
Ambient Temperature	2 C ~40 C	2 C ~40 C	2 C ~40 C	2 C ~40 C
Product Dimension (LxWxH) Inches	28x17x18	28x17x18	28x17x18	28x17x18
Product Weight (LBS)	176	176	198	198
Shipping Dimension	32x25x24	32x25x24	32x25x24	32x25x24
Shipping Weight	188Lbs	188Lbs	208Lbs	208Lbs

The TPC Superb Vac wet ring vacuum pump offers a quiet operation for small to medium size dental practices. Our twin TPC Superb Vac wet ring vacuum pump provides vacuum suction when operating in peak work periods. The twin pump offers backup protection to minimize downtime and costly interruption to work flow.

Features:

- Durable aluminum-bronze construction
- Stainless steel impeller
- Computer balanced, sealed bearings in motor
- Large capacity, easy to clean filter with clear bowl for easy inspection
- Available in single and twin models to accommodate 1-10 users
- 5-year full product limited warranty
- Unique water recycling to save up to 80% water
- UL/cUL classified
- One or Two motors can be activated to fit 1-10 chairs use

Description	Model	Suggested Retail Price
Superb Vac 2.4 HP (2x1.2HP)	WV6(Twin)	\$6,590.00
Superb Vac 2.4 HP (2x1.2HP) (with recycler)	WV6R(Twin)	\$6,990.00
Superb Vac 4 HP (2x2HP)	WV10(Twin)	\$7,390.00
Superb Vac 4 HP (2x2HP) (with recycler)	WV10R(Twin)	\$7,790.00

Product includes the following:

6 Simultaneous users	✓	✓	✓	✓
10 Simultaneous users			✓	✓
Water recycler		✓		✓
220V 50/60Hz	✓	✓	✓	✓
Adjustable vacuum level from 8-15Hg	✓	✓	✓	✓
Low voltage control panel connections	✓	✓	✓	✓
Air Water Separator	✓	✓	✓	✓

Model	Description	Suggested Retail Price
#CP-3AV	Control Panel – 3 Button- Vac1 / Vac2 / Air	\$480.00
#CP-4	Control Panel – 4 Button- Vac1 / Vac2 / Air / Water	\$560.00
#CV22	Air Compressor and Vacuum Pump Stacking Rack	\$900.00

Superb Air

Oil-less Desiccant Dryer Air Compressor



DC6224 D/A



The TPC Superb Air Oil-less Air Compressors meet the high standards of the dental professionals, providing the reliability you need from a brand you can trust. The new motor design delivers a reduced noise with a 3-year warranty. The **Desiccant air drying** system efficiently filters out compressed air contaminants and moisture from the compressed air. Purge tank removes residual humidity in the compressor following each cycle.

Features:

- New motor design to reduce noise and increase reliability
- Quiet operation
- **Desiccant air drying** system to remove moisture from compressed air
- Energy efficient thermally protected motors
- Superb Air Oil-less compressors are available in single or twin heads
- 5-year warranty

Model User Capacity	DC6112D/A 1-3users	DC6122D/A 1-3users	DC6214 D/A 3-4users	DC6224D/A 5users	DC6326D/A 6-7users	DC6228D/A 8-9users	DC6428D/A 8-9users
CDFM@80psi	4.9	4.9	9.8	9.8	14.7	19.6	19.6
# of Heads	1	1	2	2	3	2	4
Total HP	1.15	1.15	2.30	2.30	3.45	4.80	4.60
Voltage	110V	230V	110V	230V	230V	230V	230V
Amps	8	4	16	8	12	16	16
Breaker Amps	20	20	30	20	30	30	30
dBA	62	62	62	62	72	72	72
Tank Capacity	10 Gallon	10 Gallon	16 Gallon	16 Gallon	23 Gallon	30 Gallon	30 Gallon
Size (LxWxH)	16x23x27	16x23x27	26x21x29	26x21x29	34x22x29	41x22x36	41x22x36
Shipping Dimension (WxDxH)	29x23x35	29x23x35	34x26x35	34x26x35	39x25x38	46x27x41	46x27x41

Description

Model

Suggested Retail Price

Superb Air Compressor, 1.15HP, 110V	DC6112	\$3,390.00
Superb Air Compressor, 1.15HP, 230V	DC6122	\$3,390.00
Superb Air Compressor, 2 x 1.15HP(2.3HP), 110V	DC6214	\$4,990.00
Superb Air Compressor, 2 x 1.15HP(2.3HP), 230V	DC6224	\$4,990.00
Superb Air Compressor, 3 x 1.15HP(3.45HP), 230V	DC6326	\$6,390.00
Superb Air Compressor, 2 x 2.4HP(4.8HP), 230V	DC6228	\$7,790.00
Superb Air Compressor, 4 x 1.15HP(4.6HP), 230V	DC6428	\$7,940.00

Description

Model

Suggested Retail Price

Control Panel 3 button-Vac / Air / Water	CP-3	\$480.00
Control Panel 4 button-Vac1 / Vac2 / Air / Water	CP-4	\$560.00
Control Panel 2 button-Vac/Air	CP-2AV	\$400.00
Control Panel 3 button-Vac1/Vac2/Air	CP-3-AV	\$480.00
Air Compressor and Vacuum Pump Stacking Rack	CV22	\$900.00

Shipping Dimensions: L X W X H

DC6112: 29" x 23" x 35"	DC6214: 34" x 26" x 35"	DC6326: 39" x 25" x 38"	DC6428: 46" x 27" x 41"
DC6122: 29" x 23" x 35"	DC6224: 34" x 26" x 35"	DC6228: 46" x 27" x 41"	

SuperbAIR air compressors not to be used for CAD / 3D printers

Control Panels & Stacking Rack



CP-2AV



CP-3AV



CP-3



CP-4



CV-22
33 1/2" x 22" x 26"
L x W x H



Mirage2.0

Hydraulic Patient Chair with Cast Aluminum Base



4000-2.0-CB

Foot Control

Hands-free positioning of the patient.
2-programmable positions plus an entry/exit position and last position.



Cast Aluminum
Chair Base



Armrests swing-out for patient
entry and exit from the chair



Reinforced seat upholstery provides
extra comfort and high durability



Locking lever allows user to lock
chair into desired position



Articulating dual axis headrest rotation allows
for optimum positioning of the oral cavity
without the need to reposition the chair



Seat rotates 30° to left and right of center to
accommodate different working styles and to
maximize access to intraoral x-ray system



Safety shut-off plate disables seat down/backrest
recline functions to prevent injury to the operator
or damage to items beneath the plate.

Specifications:

Lifting/lowering mechanism:	Hydraulic, cantilever type
Initial height/stroke:	15.5" / 15"
Backrest reclining angle:	0° (Horizontal) – 70°
Seat tilting:	10° – 20°, synchronized with backrest
Preset position:	Microcomputer control. Three preset position, auto return, last position
Chair rotation angle:	30° left / 30° right of centerline
Headrest:	Dual axis type
Safety shut-off plate:	Under the main link
Emergency stop:	Any switch can cancel chair movement
Warranty	5 Year Limited

Standard Features:

- Synchronous backrest
- Articulating dual axis headrest
- Cast aluminum base
- Swing-away armrests
- Programmable foot control
- Long pump cover


Mirage 2.0

Hydraulic Patient Chair with Cast Aluminum Base

Description	Model	Suggested Retail Price
Mirage 2.0 Hydraulic Patient Chair with Cast Aluminum Base	4000-2.0-CB	\$8,590.00

Product includes the following:

Cast Aluminum Base
 Two programmable positions
 One preset exit position
 Electronic foot control with auto return function
 Built-in fuse / surge protection
 Toe guard protection to prevent scuffing
 All-metal construction
 Safety shut-off base plate
 Double articulating headrest, adjustable for height and wheelchair positioning
 Chair top swivel with lock, 60 degree
 Pivoting armrests for easy entry to and exit from chair
 Seamless Upholstery
 Wide chair back
 Coordinated leg rest tilt
 Input voltage 110V or 220V



Standard: PC-1000-DU
Wide backrest



Optional: TP-2005 Touch pad control installed on dental chair

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 34	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
Chair Bracket 11 x 5 x 25	35 lbs.

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

Retail

TP-2005	Touch pad control – chair mounted, 3 preset, auto return, last position	\$358.00
TPC-BRKT-2.0	Mirage M4000-2.0 post mount chair bracket	\$336.00
PC-1000-DU	Replacement Dental Chair Upholstery LK Vinyl -Wide Back	\$900.00
PC-1000-Vinyl-UC	Upcharge Non-Standard Dental Chair Upholstery Vinyl -Wide Back	\$490.00
PC-1000-ULTR-UC	Upcharge Ultra leather Dental Chair Upholstery leather -Wide Back	\$760.00
PC-1000-PR-UC	Upcharge Nauga Soft Dental Chair Upholstery Vinyl -Wide Back (PR50, 49, 47, 40, 37)	\$490.00

* Specify color and product code when ordering *

Mirage2.0

Hydraulic Patient Chair with Steel Base



4000-2.0

Foot Control

Hands-free positioning of the patient. 2-programmable positions plus an entry/exit position and last position.



Steel base with high durability plastic cover



Armrests swing-out for patient entry and exit from the chair



Reinforced seat upholstery provides extra comfort and high durability



Locking lever allows user to lock chair top into desired position



Articulating dual axis headrest rotation allows for optimum positioning of the oral cavity without the need to reposition the chair



Seat rotates 30° to left and right of center to accommodate different working styles and to maximize access to intraoral x-ray system



Safety shut-off plate disables seat down/backrest recline functions to prevent injury to the operator or damage to items beneath the plate.

Specifications:

Lifting/lowering mechanism:	Hydraulic, cantilever type
Initial height/stroke:	15.5" / 15"
Backrest reclining angle:	0° (Horizontal) – 70°
Seat tilting:	10° – 20°, synchronized with backrest
Preset position:	Microcomputer control. Three preset position, auto return, last position
Chair rotation angle:	30° left / 30° right of centerline
Headrest:	Dual axis type
Safety shut-off plate:	Under the main link
Emergency stop:	Any switch can cancel chair movement
Warranty	5 Year Limited

Standard Features:

- Synchronous backrest
- Articulating dual axis headrest
- Steel plate base
- Swing-away armrests
- Programmable foot control
- Long pump cover


Mirage 2.0

Hydraulic Patient Chair with Steel Base

Description	Model	Suggested Retail Price
Mirage 2.0 Hydraulic Patient Chair with Steel Base	4000-2.0	\$7,390.00

Product includes the following:

Two programmable positions
 One preset exit position
 Electronic foot control with auto return function
 Built-in fuse / surge protection
 Toe guard protection to prevent scuffing
 All-metal construction
 Safety shut-off base plate
 Double articulating headrest, adjustable for height and wheelchair positioning
 Chair top swivel with lock, 60 degree
 Pivoting armrests for easy entry to and exit from chair
 Seamless Upholstery
 Wide chair back
 Coordinated leg rest tilt
 Input voltage 110V or 220V



Standard: PC-1000-DU
Wide backrest



Optional: TP-2005 Touch pad control installed on dental chair

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 34	40 lbs.
Cushion: 51 x 24 x 10	35 lbs.
Chair Bracket 11 x 5 x 25	35 lbs.

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

Retail

TP-2005	Touch pad control – chair mounted, 3 preset, auto return, last position	\$358.00
TPC-BRKT-2.0	Mirage M4000-2.0 post mount chair bracket	\$336.00
PC-1000-DU	Replacement Dental Chair Upholstery LK Vinyl -Wide Back	\$900.00
PC-1000-Vinyl-UC	Upcharge Non-Standard Dental Chair Upholstery Vinyl -Wide Back	\$490.00
PC-1000-ULTR-UC	Upcharge Ultra leather Dental Chair Upholstery leather -Wide Back	\$760.00
PC-1000-PR-UC	Upcharge Nauga Soft Dental Chair Upholstery Vinyl -Wide Back (PR50, 49, 47, 40, 37)	\$490.00

* Specify color and product code when ordering *



INNOVATION•RELIABILITY•QUALITY

Mirage 2.0 Standard Upholstery Color Options



LK-03
Sand



LK-08
Taupe



LK-04
Stormy



LK-11
Grey Stone



LK-27
Orange



LK-31
Coblat



LK-34
Desert Rose



LK-49
Grey Teal



LK-59
Atlantis



LK-72
Black

TPC recommends selecting colors from an actual physical color swatch. Colors displayed in print or on digital screens may vary due to monitor, device, printer, and display settings and may not accurately represent the final product color.

Laguna 2.0

Electromechanical Patient Chair



L2000-2.0

Foot Control

Hands-free positioning of the patient.
2-programmable positions plus an entry/exit position.



Wide backrest provides extra comfort and durability



Articulating dual axis headrest rotation allows for optimum positioning of the oral cavity without the need to reposition the chair



Armrests pivot in and out for patient entry and exit from the chair



Solid steel base

Specifications:

Lifting/lowering mechanism:	Electromechanical, cantilever type
Initial height/stroke:	15.5" / 15"
Backrest reclining angle:	0° (Horizontal) – 70°
Seat tilting:	10° – 20°, synchronized with backrest
Preset position:	Microcomputer control. Three preset position, auto return, last position
Chair rotation angle:	Fixed
Headrest:	Dual axis type
Safety shut-off plate:	Under the main link
Emergency stop:	Any switch can cancel chair movement
Warranty	5 Year Limited

Standard Features:

- Independent backrest
- Articulating dual axis headrest
- Solid steel base
- Programmable foot control

Laguna 2.0

Electromechanical Patient Chair

Description	Model	Suggested Retail Price
Laguna 2.0 Electromechanical Patient Chair	L2000-2.0	\$4,590.00

Product includes the following:

Two programmable positions
 One preset exit position
 Electronic foot control with auto return function
 Built-in fuse / surge protection
 Toe guard protection to prevent scuffing
 All-metal construction
 Double articulating headrest, adjustable for height and wheelchair positioning
 Seamless Upholstery
 Input voltage 110V or 220V



Optional: TP-2005 Touch pad control installed on dental chair

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.

Model#L2000-2.0



Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."		Retail
TP-2005	Touch pad control – Chair mounted, 3 preset, auto return, last position	\$358.00
PC-5000-DU	Upcharge Dental Chair Upholstery- Wide Back LK Vinyl- Specify Color and Product Code	\$900.00
PC-5000-Vinyl-UC	Upcharge Non-Standard Dental Chair Upholstery- Wide Back Vinyl- Specify Color and Product Code	\$490.00
PC-5000-ULTR-UC	Upcharge Ultra leather Dental Chair Upholstery- Wide Back Ultra leather- Specify Color and Product Code	\$760.00
LAG.2.0-BRKT	Laguna L2000-2.0 Post mount chair bracket	\$336.00

* Specify color and product code when ordering *



INNOVATION•RELIABILITY•QUALITY

Laguna 2.0 Standard Upholstery Color Options



LK-03
Sand



LK-08
Taupe



LK-04
Stormy



LK-11
Grey Stone



LK-27
Orange



LK-31
Cobalt



LK-34
Desert Rose



LK-49
Grey Teal



LK-59
Atlantis



LK-72
Black

TPC recommends selecting colors from an actual physical color swatch. Colors displayed in print or on digital screens may vary due to monitor, device, printer, and display settings and may not accurately represent the final product color.

Ortho 1200

Fixed Base Orthodontic Patient Chair



OR1200



Optional: Pony Tail headrest

Standard Features:

- Independent backrest
- Articulating dual axis headrest
- Steel plate base
- Fixed base chair
- Foot control and last position memory

Ortho 1200

Fixed Base Orthodontic Patient Chair


Description	Model	Suggested Retail Price
Fixed Base Orthodontic Patient Chair	OR1200	\$2,700.00

Product includes the following:

Fixed base height
 One preset exit position
 Electronic foot control
 Independent backrest
 Toe guard protection to prevent scuffing
 Double articulating headrest, Steel plate base
 Seamless Upholstery



Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	440 lbs.
Cushion: 51 x 24 x 10	35 lbs.



Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."		Retail
PC-5000-DU	Upcharge Dental Chair Upholstery- Wide Back LK Vinyl- Specify Color and Product Code	\$900.00
PC-5000-Vinyl-UC	Upcharge Non-Standard Dental Chair Upholstery- Wide Back Vinyl- Specify Color and Product Code	\$490.00
PC-5000-ULTR-UC	Upcharge Ultra leather Dental Chair Upholstery- Wide Back Ultra leather- Specify Color and Product Code	\$760.00
L2001	OR-1200 Orthodontic chair bracket	\$336.00
PTH 101	Pony Tail Headrest	\$218.00
M42408C	Right or Left side Armrest with bracket	\$300.00



INNOVATION•RELIABILITY•QUALITY

QR 1200 Standard Upholstery Color Options



LK-03
Sand



LK-08
Taupe



LK-04
Stormy



LK-11
Grey Stone



LK-27
Orange



LK-31
Coblat



LK-34
Desert Rose



LK-49
Grey Teal



LK-59
Atlantis



LK-72
Black

TPC recommends selecting colors from an actual physical color swatch. Colors displayed in print or on digital screens may vary due to monitor, device, printer, and display settings and may not accurately represent the final product color.

Bellagio

*Bellagio Chair Mounted Operatory System
with Cast Aluminum Base*

BP5000-L700LED-2.0-CB



L700-LED

5000-2.0

Cast Aluminum base

4000-2.0-CB

Optional: Highspeed
Handpieces

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
Unit: 39 x 29 x 23	60 lbs.
Assist PKG: 43 x 25 x 12	40 lbs.
Dental Light Arm L600: 34 x 11 x 5	30 lbs.
Dental Light Head L910: 17 x 15 x 7	
Cuspidor 26 x 18 x 10	20 lbs.



BP5015-L700LED-2.0 CB



Model DR-1102
Optional

Model AS-1101
Optional

Bellagio

Bellagio Chair Mounted Operatory System with Cast Aluminum Base

Description	Model	Suggested Retail Price
Bellagio Operatory Package with Cuspidor	BP5000-L700LED-2.0-CB	\$15,870.00
Bellagio Operatory Package with Assistant's Instrumentation	BP5015-L700LED-2.0-CB	\$15,170.00

Package includes the following:

(with cuspidor) (without cuspidor)

4000-2.0-CB	Hydraulic Cast Aluminum Base patient chair with wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 3 programmable positions with automatic return and foot control, cast aluminum base	✓	✓
5015-2.0	Asepsis 3 Handpiece automatic control with coolant controls, flex arm with air brake, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air / water filters and regulators, master shut off valves. Includes 3 way syringe and light post with cutout for internal wires and monitor bracket Asepsis telescoping 4 position assistant's package, 3-way syringe, Hve and saliva ejector, water / air quick disconnects, solids collector clean water system, Built-in Piezo Scaler, Built-in LED Curing Light, internal timer, integrated touch pad controls	✓	✓
L700-LED	Vivid LED post mounted light with adjustable intensity, 10,000-35,000 Lux composite cure mode, motion sensor, manual toggle on/off	✓	✓
2000-C-2.0	Porcelain cuspidor with timed bowl rinse and cup filler	✓	


Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

Retail

P2335+TB	Deluxe autoclavable Hve Long valve with tubing	\$67.00
WIC800	Standard flat screen monitor bracket	\$128.00
WIC840	Mirage flat screen monitor bracket with 16" extension Arm	\$210.00
MO17	17" LCD multi media flat screen monitor	\$540.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System- ISO C (6 pin) – tubing, control box, power supply	\$270.00
53016-L-100	Disposable Solids Collector Trap, 2 1/8' (pack of 100)	\$45.00
2000-C-2.0	Mirage Cuspidor	\$700.00
EMB-901	Built in Electric Motor System	\$1,790.00
AS-1101	Mirage Assistants stool	\$818.00
DR-1102	Mirage Doctor's stool	\$650.00

Bellagio

Bellagio Chair Mounted Operatory System with Steel Base



Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
Unit: 39 x 29 x 23	60 lbs.
Assist PKG: 43 x 25 x 12	40 lbs.
Dental Light Arm L600: 34 x 11 x 5	30 lbs.
Dental Light Head L570: 17 x 15 x 7	



BP5000-L550LED-2.0



Bellagio

Bellagio Chair Mounted Operatory System with Steel Base

Description	Model	Suggested Retail Price
Bellagio Operatory Package with Cuspidor	BP5000-L550LED-2.0	\$14,380.00
Bellagio Operatory Package with Assistant's Instrumentation	BP5015-L550LED-2.0	\$13,680.00

Package includes the following:


4000-2.0	Hydraulic patient chair with steel base wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 3 programmable positions with automatic return and foot control	✓	✓
5015-2.0	Asepsis 3 Handpiece automatic control with coolant controls, flex arm with air brake, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves. Includes 3 way syringe and light post with cutout for internal wires and monitor bracket Asepsis telescoping 4 position assistant's package, 3-way syringe, Hve and saliva ejector clean water / air quick disconnects, solids collector with touch pad control • Built-in Scaler • Built-in Curing light	✓	✓
L550-LED	Lustrous LED post mounted light with adjustable intensity, 8,000-30,000 Lux composite cure mode, motion sensor, manual toggle on/off	✓	✓
2000-C-2.0	Porcelain cuspidor with timed bowl rinse and cup filler	✓	

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

		Retail
P2335+TB	Deluxe autoclavable Hve Long valve with tubing	\$67.00
WIC800	Standard flat screen monitor bracket	\$128.00
WIC840	Mirage flat screen monitor bracket with 16" extension Arm	\$210.00
MO17	17" LCD multi media flat screen monitor	\$540.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System- ISO C (6 pin) – tubing, control box, power supply	\$270.00
53016-L-100	Disposable Solids Collector Trap, 2 1/8' (pack of 100)	\$45.00
2000-C-2.0	Mirage Cuspidor	\$700.00
EMB-901	Built in Electric Motor System	\$1,790.00
AS-1101	Mirage Assistants stool	\$818.00
DR-1102	Mirage Doctor's stool	\$650.00

Bellagio

Bellagio Chair Mounted Operatory System with Steel Base



Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
Unit: 39 x 29 x 23	60 lbs.
Assist PKG: 43 x 25 x 12	40 lbs.
Dental Light Arm L-7601: 34 x 11 x 5	30 lbs.
Dental Light Head L7601: 17 x 15 x 7	



Bellagio

Bellagio Chair Mounted Operatory System with Steel Base

Description	Model	Suggested Retail Price
Bellagio Operatory Package with Cuspidor	BP5000-L-7601LED-2.0	\$15,980.00
Bellagio Operatory Package with Assistant's Instrumentation	BP5015-L-7601LED-2.0	\$15,280.00

Package includes the following:

4000-2.0	Hydraulic patient chair with steel base wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 3 programmable positions with automatic return and foot control	✓	✓
5015-2.0	Asepsis 3 Handpiece automatic control with coolant controls, flex arm with air brake, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves. Includes 3 way syringe and light post with cutout for internal wires and monitor bracket Asepsis telescoping 4 position assistant's package, 3-way syringe, Hve and saliva ejector clean water / air quick disconnects, solids collector with touch pad control • Built-in Scaler • Built-in Curing light	✓	✓
L-7601-LED	Luminous LED post mounted light with adjustable intensity, 3,000-35,000 Lux, manual toggle on/off, optional poly-block filter	✓	✓
2000-C-2.0	Porcelain cuspidor with timed bowl rinse and cup filler	✓	

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."	Retail
Deluxe autoclavable Hve Long valve with tubing	\$67.00
Standard flat screen monitor bracket	\$128.00
Mirage flat screen monitor bracket with 16" extension Arm	\$210.00
17" LCD multi media flat screen monitor	\$540.00
Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
Illumination System- ISO C (6 pin) – tubing, control box, power supply	\$270.00
Disposable Solids Collector Trap, 2 1/8' (pack of 100)	\$45.00
Mirage Cuspidor	\$700.00
Built in Electric Motor System	\$1,790.00
Mirage Assistants stool	\$818.00
Mirage Doctor's stool	\$650.00
LU-PB Antipolymerization Filter "Poli-Block"	\$130.00

Bellagio

Swing Mount Operatory System

Model #3550-2.0



MSP5500-2.0-CB

4000-2.0-CB

Model #5500-2.0



Optional: TP-2005 touch pad control

Optional: Post LP100 & L700 LED Light



Optional: Model #L795-LED Ceiling mounted LED operatory light



4000-2.0

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
5500: 36 x 28 x 39	100 lbs.
3550: 29 x 12 x 24	25 lbs.
LP100: 45 x 26 x 44	30 lbs.

Bellagio

Swing Mount Operatory System

Description	Model	Suggested Retail Price
Bellagio Swing Mount with Cast Aluminum Base	MSP5500-2.0-CB	\$14,570.00
Bellagio Swing Mount with Steel Base	MSP5500-2.0	\$13,370.00

Package includes the following:

4000-2.0-CB	Hydraulic cast aluminum base patient chair with wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 4 programmable positions with automatic return and foot control	✓	
5500-2.0	Rotating asepsis 3 Handpiece automatic control with coolant controls, flex arm with air brake, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves. Includes 3 way syringe and clean water system, integrated Piezo Scaler, Built-in LED Curing Light internal timer, integrated touch pad controls	✓	✓
3550-2.0	Rear radius asepsis assistant's arm system, Hve and saliva, 3 Way Syringe	✓	✓
4000-2.0	Hydraulic patient chair with steel base, wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 3 programmable positions with automatic return and foot control		✓

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

Retail

5500-2.0	Bellagio swing mount unit with chair bracket (for TPC chair only)	\$4,790.00
3550-2.0	Mirage rear swing mount suction package (for TPC chair only)	\$1,190.00
LP-100	Mirage swing mount light post (ordered w/unit)	\$770.00
LP-105	Mirage swing mount monitor post with monitor mounting bracket (ordered w/unit)	\$960.00
L700-LED	Vivid LED Operatory Light	\$1,590.00
L795-LED	Vivid LED Ceiling Mounted Operatory Light	\$2,360.00
L550-LED	Lustrous LED Operatory Light	\$1,300.00
L7601-LED	Luminous LED Operatory Light	\$2,900.00
P2335+TB	Deluxe autoclavable Hve Long valve with tubing	\$67.00
WIC800	Standard flat screen monitor bracket	\$128.00
WIC840	Mirage flat screen monitor bracket with 16" extension	\$210.00
MO17	17" LCD multi media flat screen monitor	\$540.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply	\$270.00
53016-L	Disposable Solids Collector Trap 2 1/8" (pack of 100)	\$45.00
EMB-901	Built in Electric Motor System	\$1,790.00
AS-1101	Mirage Assistants stool	\$818.00
DR1102	Mirage Doctor's Stool	\$650.00

Mirage2.0

*Mirage Chair Mounted Operatory System
with Cast Aluminum Base*

MP2000-L550LED-2.0-CB



Cast Aluminum base

4000-2.0-CB

Optional:
Highspeed
Handpieces



Optional:
Articulating Tray

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
Unit: 39 x 29 x 23	60 lbs.
Assist PKG: 43 x 25 x 12	40 lbs.
Dental Light Arm L600: 34 x 11 x 5	30 lbs.
Dental Light Head L910: 17 x 15 x 7	
Cuspidor 26 x 18 x 10	20 lbs.



MP2015-L550LED-2.0 CB



Model DR-1102
Optional

Model AS-1101
Optional

Mirage2.0

Mirage Chair Mounted Operatory System with Cast Aluminum Base

Description	Model	Suggested Retail Price
Mirage 2.0 Operatory Package with Cuspidor	MP2000-L550LED-2.0-CB	\$14,780.00
Mirage 2.0 Operatory Package with Assistant's Instrumentation	MP2015-L550LED-2.0-CB	\$14,080.00

Package includes the following:

		(with cuspidor)	(without cuspidor)
4000-2.0-CB	Hydraulic Cast Aluminum Base patient chair with wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 3 programmable positions with automatic return and foot control, cast aluminium base	✓	✓
2015-2.0	Asepsis 3 Handpiece automatic control with coolant controls, flex arm with air brake, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves. Includes 3 way syringe and light post with cutout for internal wires and monitor bracket Asepsis telescoping 4 position assistant's package, Hve and saliva ejector, water / air quick disconnects, solids collector clean water system Touch pad Control on Both Dr. & Assistant instrument	✓	✓
L550-LED	Lustrous LED post mounted light with adjustable intensity, 8,000-30,000 Lux composite cure mode, motion sensor, manual toggle on/off	✓	✓
2000-C	Porcelain cuspidor with timed bowl rinse and cup filler	✓	

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

	Retail
TPC3WAY-PRT	Additional autoclavable 3 way syringe with tubing \$78.00
P2335+TB	Deluxe autoclavable Hve Long valve with tubing \$67.00
WIC800	Standard flat screen monitor bracket \$128.00
WIC840	Mirage flat screen monitor bracket with 16" extension Arm \$210.00
MO17	17" LCD multi media flat screen monitor \$540.00
A800N-1	Advance 800 piezo scaler – built in \$570.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply \$310.00
H7660	Illumination System- ISO C (6 pin) – tubing, control box, power supply \$270.00
ALED55BN-A	LED 5W curing light – built in to dental unit \$438.00
53016-L-100	Disposable Solids Collector Trap, 2 1/8' (pack of 100) \$45.00
TPE-2011	Touch pad Control - Chair & Electric Motor System touch pad \$460.00
EMB-901	Built in Electric Motor System \$1,790.00
XV100	Super thin X-ray film viewer \$140.00
XV150-R	X-ray film viewer bracket (curved) for TPC Mirage unit only \$80.00
AS-1101	Mirage Assistant's stool \$818.00
DR-1102	Mirage Doctor's stool \$650.00

Mirage2.0

Mirage Chair Mounted Operatory System with Steel Base

MP2015-L550LED-2.0



L550-LED

2015-2.0

4000-2.0

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
Unit: 39 x 29 x 23	60 lbs.
Assist PKG: 43 x 25 x 12	40 lbs.
Dental Light Arm L600: 34 x 11 x 5	30 lbs.
Dental Light Head L570: 17 x 15 x 7	



Model DR-1102
Optional



Model AS-1101
Optional



MP2000-L550LED-2.0

Mirage2.0

Mirage Chair Mounted Operatory System with Steel Base

Description	Model	Suggested Retail Price
Mirage 2.0 Operatory Package with Cuspidor	MP2000-L550LED-2.0	\$13,580.00
Mirage 2.0 Operatory Package with Assistant's Instrumentation	MP2015-L550LED-2.0	\$12,880.00

Package includes the following:

4000-2.0	Hydraulic patient chair with steel base, wide back and full lumbar support, rotating armrests left or right, double articulating head-rest with sliding height adjustment, 30 degree swivel from left and right of center, 3 programmable positions with automatic return and foot control	✓	✓
2015-2.0	Asepsis 3 Handpiece automatic control with coolant controls, flex arm with air brake, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves. Includes 3 way syringe and light post with cutout for internal wires and monitor bracket Asepsis telescoping 4 position assistant's package, Hve and saliva ejector, water / air quick disconnects, solids collector clean water system, Touch pad Control on Both Dr. & Assistant instrument	✓	✓
L550-LED	Lustrous LED post mounted light with adjustable intensity, 8,000-30,000 Lux composite cure mode, motion sensor, manual toggle on/off	✓	✓
2000-C-2.0	Porcelain cuspidor with timed bowl rinse and cup filler	✓	

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

	Retail
TPC3WAY-PRT Additional autoclavable 3 way syringe with tubing	\$78.00
P2335+TB Deluxe autoclavable Hve Long valve with tubing	\$67.00
WIC800 Standard flat screen monitor bracket	\$128.00
WIC840 Mirage flat screen monitor bracket with 16" extension Arm	\$210.00
MO17 17" LCD multi media flat screen monitor	\$540.00
A800N-1 Advance 800 piezo scaler – built in	\$570.00
H7000 Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660 Illumination System- ISO C (6 pin) – tubing, control box, power supply	\$270.00
ALED55BN-A LED 5W curing light – built in to dental unit	\$438.00
53016-L-100 Disposable Solids Collector Trap, 2 1/8' (pack of 100)	\$45.00
TPE-2011 Touch pad Control - Chair & Electric Motor System touch pad	\$460.00
EMB-901 Built in Electric Motor System	\$1,790.00
XV100 Super thin X-ray film viewer	\$140.00
XV150-R X-ray film viewer bracket (curved) for TPC Mirage unit only	\$80.00
AS-1101 Mirage Assistant's stool	\$818.00
DR-1102 Mirage Doctor's stool	\$650.00

Mirage 2.0

Swing Mount Operatory System

Model #3550-2.0

#MSP3500-2.0



Model #3500-2.0

4000-2.0

Optional: LP100 swing mount light post with optional L550-LED Operatory Light



4000-2.0 CB



Optional: Model #L570-LED Ceiling mounted LED operatory light

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
5500: 36 x 28 x 39	100 lbs.
3550: 29 x 12 x 24	25 lbs.
LP100: 45 x 26 x 44	30 lbs.

Mirage 2.0

Swing Mount Operatory System

Description	Model	Suggested Retail Price	
Mirage 2.0 Swing Mount Operatory System with Steel Base	MSP3500-2.0		\$12,570.00
Mirage 2.0 Swing Mount Operatory System with Cast Aluminum Base		MSP3550-2.0-CB	\$13,770.00
Package includes the following:			
4000-2.0	Hydraulic patient chair with steel base, wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 4 programmable positions with automatic return and foot control	✓	
3500-2.0	Rotating asepsis 3 Handpiece automatic control with coolant controls, flex arm with air brake, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves. Includes 3 way syringe and clean water system, with touch pad control	✓	✓
3550-2.0	Rear radius asepsis assistant's arm system, Hve and saliva ejector, water quick disconnect and solids collector, 3 Way Syringe	✓	✓
4000-2.0-CB	Hydraulic Cast Aluminum Base patient chair with wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 3 programmable positions with automatic return and foot control, cast aluminium base		✓

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

Retail

3500-2.0	Mirage swing mount unit with chair bracket (for TPC chair only)	\$3,990.00
3550-2.0	Mirage rear swing mount suction package (for TPC chair only)	\$1,190.00
LP-100	Mirage swing mount light post (ordered w/unit)	\$770.00
LP-105	Mirage swing mount monitor post with monitor mounting bracket (ordered w/unit)	\$960.00
L700-LED	Vivid LED Operatory Light	\$1,590.00
L795-LED	Vivid LED Ceiling Mounted Operatory Light	\$2,360.00
L550-LED	Lustrous LED Operatory Light	\$1,300.00
L7601-LED	Luminous LED Operatory Light	\$2,900.00
P2335+TB	Deluxe autoclavable Hve Long valve with tubing	\$67.00
WIC800	Standard flat screen monitor bracket	\$128.00
WIC840	Mirage flat screen monitor bracket with 16" extension	\$210.00
MO17	17" LCD multi media flat screen monitor	\$540.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply	\$270.00
53016-L	Disposable Solids Collector Trap 2 1/8" (pack of 100)	\$45.00
EMB-901	Built in Electric Motor System	\$1,790.00
AS-1101	Mirage Assistants stool	\$818.00
DR1102	Mirage Doctor's Stool	\$650.00

Mirage 2.0

Hygiene Operatory System with Steel Base



#H1001

MHP1001-2.0

4000-2.0



W/Chair Mount light



Optional: Model #L570-LED Ceiling mounted LED operatory light



W/Cuspidor & Post Mount light

Mirage 2.0

Hygiene Operatory System with Steel Base

Description	Model	Suggested Retail Price
Mirage 2.0 Hygiene Operatory System (without light)	MHP1001	\$11,090.00

Package includes the following:

4000-2.0	Hydraulic patient chair with steel base, wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 4 programmable positions with automatic return and foot control	✓
H1001	Rear radius asepsis 3 Handpiece automatic control with coolant controls, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves, Unit Head Height Adjustment from 27" to 39" Includes 3 way syringe and clean water system, with touch pad control	✓

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."		Retail
LP-100	Mirage swing mount light post (ordered w/unit)	\$770.00
LP-105	Mirage swing mount monitor post with monitor mounting bracket (ordered w/unit)	\$960.00
L570-LED	Lustrous LED Ceiling Mounted Operatory Light	\$1,960.00
L550-LED	Lustrous LED Operatory Light	\$1,300.00
LP48	Lustrous LED Operatory Light	\$210.00
WIC800	Standard flat screen monitor bracket	\$128.00
WIC840	Mirage flat screen monitor bracket with 16" extension	\$210.00
MO17	17" LCD multi media flat screen monitor	\$540.00
A800N-1	Advance 800 piezo scaler – built in	\$570.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply	\$270.00
ALED55BN-A	LED 5W curing light – built in to dental unit	\$438.00
TPCRDSBRKT	TPC Radius Bracket	\$460.00
TPE-2011	Touch pad Control - Chair & Electric Motor System Touch pad	\$460.00
EMB-901	Built in Electric Motor System	\$1,790.00
AS-1101	Mirage Assistants stool	\$818.00
DR1102	Mirage Doctor's Stool	\$650.00
TPC-BRKT-2.0	Mirage chair Bracket	\$336.00
2000-PMU	PMU with cuspidor & cup filler	\$1,590.00

Mirage 2.0

Orthodontic Package with Steel Base

Standard Features:

- Synchronous backrest
- Articulating dual axis headrest
- Laser cut steel base
- Long pump cover
- Hydraulic pump for quiet operation
- Programmable foot control

MOP3000-L550LED-2375



Optional: Pony Tail headrest

#LP48

L550-LED

Standard Features:

- Composite non cure mode
- 24V DC power Supply
- Bulb type: 4 White LEDs 2 Yellow
- Light intensity (Lux): High 30,000 Low 8,000
- Focal length: 27" (70 cm)
- Illumination area: 100mm x 180mm
- Sensor mode manual toggle

Standard Control Features:

- 3 HP Automatic selection
- 3-way Syringe
- Vacuum package with saliva ejector and hve
- Self contained water bottle system
- Variable speed disc-type foot control
- Water coolant on/off switch with flow control
- 1 Year warranty

Requires central air and vacuum for plumbing

Optional


Self contained Compressor & Vacuum

(No plumbing required)



OMC2375

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
Light Arm L600 26 x 14 x 8	30 lbs.
Light Head L570 17 x 15 x 7	
Cabinet 32 x 28 x 41	140 lbs.



Counter top dimensions



Rear cabinet door view



Pull out writing surface



Optional: OMC-2200 LCD Bracket

Mirage 2.0

Orthodontic Package with Steel Base

Description	Model	Suggested Retail Price
Mirage 2.0 Orthodontic Package	MOP3000-550LED-2375	\$11,840.00
Mirage 2.0 Orthodontic Package	MOP3000-550LED-2375CV	\$13,400.00
Mirage 2.0 Orthodontic Package	MOP3000-550LED-2375CV-LS	\$14,300.00

Package includes the following:

		2375	2375CV	2375CV-LS
3000-2.0	Hydraulic patient chair with steel base, wide back and full lumbar support, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 4 programmable positions with automatic return and foot control	✓	✓	✓
OMC2375	This orthodontic mobile delivery cabinet offers a pullout writing surface, two instrument drawers, one medium and two large storage drawers for materials and supplies that fully extend and quickly remove for easy cleaning. In addition, these cabinets provide convenient access from the front. Oversized locking dual-caster wheel, central plumbing with junction box	✓		
OMC2375CV	OMC2375+Compressor+ Vacuum (Self contained)		✓	
OMC2375CV-LS	OMC2375+Compressor+ Vacuum+LED Curing light+Scaler (Self contained)			✓
L550-LED	Lustrous LED operatory light with Poliblock mode 8,000 - 30,000 LUX , dual axis light head rotation	✓	✓	✓

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."


		Retail
WIC840	Mirage flat screen monitor bracket with 16" extension	\$210.00
MO17	17" LCD multi media flat screen monitor	\$540.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply	\$270.00
ALED55BN-A	LED 5W curing light – built in to dental unit	\$438.00
TP-2005	Touch pad control – dental chair mounted 3 preset, auto return, last position	\$358.00
OMC2200	LCD monitor bracket	\$330.00
OMC-ST81	Optional Stone Surface on OMC 2375, Color Artic	\$418.00
OMC-ST82	Optional Stone Surface on OMC 2375, Color White	\$418.00
OMC-ST85	Optional Stone Surface on OMC 2375, Color Cream	\$418.00
OMC-ST86	Optional Stone Surface on OMC 2375, Color Bamboo	\$418.00
M4000-09	Optional Magnetic Headrest	\$420.00
AS-1101	Mirage Assistants stool	\$818.00
DR-1102	Mirage Doctor's stool	\$650.00
PTH 101	Pony Tail Headrest	\$218.00



Optional: Stone Counter top
OMC-ST81



Optional: Stone Counter top
OMC-ST82



Optional: Stone Counter top
OMC-ST85



Optional: Stone Counter top
OMC-ST86

Laguna 2.0

Chair Mounted Operatory System

LP2100-550LED-2.0

Optional
MO17 LCD Monitor

L550-LED

2000-2.0

L2000-2.0

Optional



DR-5045



AS-5047



LP2115-550LED-2.0

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
Unit: 39 x 29 x 23	60 lbs.
Unit: 43 x 25 x 12	40 lbs.
Dental Light: 34 x 11 x 5 Light	30 lbs.
Head: 12 x 15 x 7	
Cuspidor 20 x 13 x 14	20 lbs.

Laguna 2.0

Chair Mounted Operatory System

Description	Model	Suggested Retail Price
Laguna 2.0 Operatory Package with Cuspidor	LP2100-550LED-2.0	\$10,780.00
Laguna 2.0 Operatory Package with Assistant's Instrumentation	LP2115-550LED-2.0	\$10,080.00
Package includes the following:	(with cuspidor) (without cuspidor)	
L2000-2.0 Electromechanical patient chair with wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 3 programmable positions with automatic return and foot control	✓	✓
2015-2.0 Asepsis 3 Handpiece automatic control with coolant controls, flex arm with air brake, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves. Includes 3 way syringe and light post with cutout for internal wires and monitor bracket Asepsis rotating 4 position assistant's package, Hve and saliva ejector, water disconnect solids collector and clean water system. With Touch pad control. on both Dr. & Assistant instrument	✓	✓
L550-LED Lustrous LED post mounted light with adjustable intensity, 8,000-30,000 Lux	✓	✓
2000-C-2.0 Porcelain cuspidor with timed bowl rinse and cup filler	✓	

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

	Retail
TPC3WAY-PRT Additional autoclavable 3 way syringe with tubing	\$78.00
P2335+TB Deluxe autoclavable Hve Long valve with tubing	\$67.00
WIC800 Standard flat screen monitor bracket	\$128.00
WIC840 Mirage flat screen monitor bracket with 16" extension Arm	\$210.00
MO17 17" LCD multi media flat screen monitor	\$540.00
A800N-1 Advance 800 piezo scaler – built in	\$570.00
H7000 Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660 Illumination System- ISO C (6 pin) – tubing, control box, power supply	\$270.00
ALED55BN-A LED 5W curing light – built in to dental unit	\$438.00
53016-L-100 Disposable Solids Collector Trap, 2 1/8' (pack of 100)	\$45.00
TPE-2011 Touch pad Control - Chair & Electric Motor System touch pad	\$460.00
EMB-901 Built in Electric Motor System	\$1,790.00
XV100 Super thin X-ray film viewer	\$140.00
XV150-R X-ray film viewer bracket (curved) for TPC Mirage unit only	\$80.00
AS-5047 Laguna Assistant's stool	\$590.00
DR-5045 Laguna Doctor's stool	\$480.00

Laguna 2.0

Swing Mount Operatory System

Model #3550-2.0

#LSP3500-2.0



Model #3500-2.0

L2000-2.0

Optional: LP100 swing mount light post with optional L550-LED Operatory Light



LP100+L550 LED

L2000-2.0



Optional: Model #L570-LED Ceiling mounted LED operatory light

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
5500: 36 x 28 x 39	100 lbs.
3550: 29 x 12 x 24	25 lbs.
LP100: 45 x 26 x 44	30 lbs.

Laguna 2.0

Swing Mount Operatory System

Description	Model	Suggested Retail Price
Laguna 2.0 Swing Mount Operatory System with Steel Base	LSP3500-2.0	\$9,770.00

Package includes the following:

L2000-2.0	Electromechanical patient chair with wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 4 programmable positions with automatic return and foot control	✓
3500-2.0	Rotating asepsis 3 Handpiece automatic control with coolant controls, flex arm with air brake, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves. Includes 3 way syringe and clean water system, with touch pad control	✓
3550-2.0	Rear radius asepsis assistant's arm system, Hve and saliva ejector, water quick disconnect and solids collector, 3 Way Syringe	✓

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."		Retail
3500-2.0	Mirage swing mount unit with chair bracket (for TPC chair only)	\$3,390.00
3550-2.0	Mirage rear swing mount suction package (for TPC chair only)	\$1,190.00
LP-100	Mirage swing mount light post (ordered w/unit)	\$770.00
LP-105	Mirage swing mount monitor post with monitor mounting bracket (ordered w/unit)	\$960.00
L795-LED	Vivid LED Ceiling Mounted Operatory Light	\$2,360.00
L570-LED	Lustrous LED Ceiling Mounted Operatory Light	\$1,960.00
L550-LED	Lustrous LED Operatory Lights	\$1,300.00
P2335+TB	Deluxe autoclavable Hve Long valve with tubing	\$67.00
WIC800	Standard flat screen monitor bracket	\$128.00
WIC840	Mirage flat screen monitor bracket with 16" extension	\$210.00
MO17	17" LCD multi media flat screen monitor	\$540.00
A800N-1	Advance 800 piezo scaler – built in	\$570.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply	\$270.00
ALED55BN-A	LED 5W curing light – built in to dental unit	\$438.00
53016-L	Disposable Solids Collector Trap 2 1/8" (pack of 100)	\$45.00
TPE-2011	Touch pad Control - Chair & Electric Motor System Touch pad	\$460.00
EMB-901	Built in Electric Motor System	\$1,790.00
AS-5047	Laguna Assistants stool	\$590.00
DR-5045	Laguna doctor's stool	\$480.00

Laguna 2.0

Hygiene operatory package

#LHP1001-2.0

H1001



L2000-2.0



W/Chair Mount light



W/Ceiling
Mount light

Dimension
Chair
H1001
Cushion



W/Cuspidor & Post Mount light

Laguna 2.0

Hygiene Operatory System

Description	Model	Suggested Retail Price
Laguna 2.0 Hygiene Operatory System (without light)	LHP1001-2.0	\$8,290.00

Package includes the following:

L2000-2.0	Electromechanical patient chair with wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 4 programmable positions with automatic return and foot control	✓
H1001	Rear radius asepsis 3 Handpiece automatic control with coolant controls, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves, Unit Head Height Adjustment from 27" to 39" Includes 3 way syringe and clean water system, with touch pad control.	✓
	Hve and saliva ejector, and solids collector, piezo scaler with light	✓

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."		Retail
LP-100	Mirage swing mount light post (ordered w/unit)	\$770.00
LP-105	Mirage swing mount monitor post with monitor mounting bracket (ordered w/unit)	\$960.00
L570-LED	Lustrous LED Ceiling Mounted Operatory Light	\$1,960.00
L550-LED	Lustrous LED Operatory Light	\$1,300.00
LP48	Lustrous LED Operatory Light	\$210.00
WIC800	Standard flat screen monitor bracket	\$128.00
WIC840	Mirage flat screen monitor bracket with 16" extension	\$210.00
MO17	17" LCD multi media flat screen monitor	\$540.00
A800N-1	Advance 800 piezo scaler – built in	\$570.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply	\$270.00
ALED55BN-A	LED 5W curing light – built in to dental unit	\$438.00
TPCRDSBRKT	TPC Radius Bracket	\$460.00
TPE-2011	Touch pad Control - Chair & Electric Motor System Touch pad	\$460.00
EMB-901	Built in Electric Motor System	\$1,790.00
AS-5047	Laguna Assistants stool	\$590.00
DR5045	Laguna Doctor's Stool	\$480.00
M4001	Mirage chair Bracket	\$336.00
2000-PMU	PMU with cuspidor & cup filler	\$1,590.00

Laguna 2.0

Orthodontic Package

Standard Features:

- Independent backrest
- Articulating dual axis headrest
- Solid steel base
- Programmable foot control
- Built in fuse/surge protection
- Electric motor operated on chair base and backrest

Standard Cabinet Features:

- Five easy to remove and clean full extension drawers, and one pull-out writing surface.
- High quality hardware throughout
- Oversized caster wheels with locks
- A removable counter - top plastic tray
- Standard 12" X 12" Junction Box
- 7' length of umbilical



Standard Features:

- Composite cure mode
- 24V DC power Supply
- Bulb type: 4 White LEDs 2 Yellow
- Light intensity (Lux): High 30,000 Low 8,000
- Focal length: 27" (70 cm)
- Illumination area: 100mm x 180mm
- Sensor mode, manual toggle on/off

Standard Control Features:

- 3 HP Automatic selection
- 3-way Syringe
- Vacuum package with saliva ejector and hve
- Self contained water bottle system
- wet / dry foot control
- 1 Year warranty
- **Requires central Air / vacuum systems**
- **Includes utility center and junction box**

Optional:

Self contained compressor & Vacuum
(No plumbing required)



OMC2375

Counter top dimensions



Optional: OMC-
2200 LCD Bracket

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
Light Arm L600: 26 x 14 x 10	
Light head L550 17 x 15 x 7	30 lbs.
Cabinet 32 x 28 x 41	140 lbs.

Laguna 2.0

Orthodontic Package

Description	Model	Suggested Retail Price		
Laguna 2.0 Orthodontic Package	LOP2000-550LED-2375	\$9,090.00		
Laguna 2.0 Orthodontic Package	LOP2000-550LED-2375CV	\$10,600.00		
Laguna 2.0 Orthodontic Package	LOP2000-550LED-2375CV-SL	\$11,500.00		
Package includes the following:		2375	2375CV	2375CV-SL

L2000-2.0	Electromechanical patient chair with wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 4 programmable positions with automatic return and foot control	✓	✓	✓
OMC-2375	This orthodontic mobile delivery cabinet offers a pullout writing surface, two instrument drawers, one medium and two large storage drawers for materials and supplies that fully extend and quickly remove for easy cleaning. In addition, these cabinets provide a large open space for accessories with convenient access from the front. Oversized locking dual-casters, central plumbing with junction Box	✓		
OMC2375CV	OMC2375+Compressor+Vacuum(Self contained)		✓	
OMC2375CV-SL	OMC2375+Compressor+Vacuum+LED Curing light+Scaler(Self contained)			✓
L550-LED	Lustrous LED post mounted light with dual intensity, 8,000 - 30,000 LUX, dual axis light head rotation	✓	✓	✓
L2001-LP48	The Laguna chair bracket and LP48" Light post are included at no charge (value: \$470) when purchasing the Orthodontic Package	✓	✓	✓


Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."		Retail
WIC840	Mirage flat screen monitor bracket with 16" extension	\$210.00
MO17	17" LCD multi media flat screen monitor	\$540.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply	\$270.00
OMC2200	LCD monitor bracket	\$330.00
OMC-ST81	Optional Stone Surface on OMC 2375, Color Artic	\$418.00
OMC-ST82	Optional Stone Surface on OMC 2375, Color White	\$418.00
OMC-ST85	Optional Stone Surface on OMC 2375, Color Cream	\$418.00
OMC-ST86	Optional Stone Surface on OMC 2375, Color Bamboo	\$418.00
L2001	Laguna L2000-2.0 Chair Bracket	\$336.00
AS-5047	Laguna Assistants stool	\$590.00
DR-5045	Laguna Doctor's stool	\$480.00
PTH 101	Pony Tail Headrest	\$218.00



Optional: Stone Counter top
OMC-ST81



Optional: Stone Counter top
OMC-ST82



Optional: Stone Counter top
OMC-ST85



Optional: Stone Counter top
OMC-ST86

OR-1200

Orthodontic Package

Standard Features:

- Independent backrest
- Articulating dual axis headrest
- Fixed base
- Built in fuse/surge protection
- Last position memory foot control
- Short pump cover
- Electric motor on backrest

Standard Cabinet Features:

- Five easy to remove and clean full extension drawers, and one pull-out writing surface.
- High quality hardware throughout
- Oversized caster wheels with locks
- A removable counter - top plastic tray
- Standard 12" X 12" Junction Box
- 7' length of umbilical



Standard Features:

- Composite cure mode
- 24V DC power Supply
- Bulb type: 4 White LEDs 2 Yellow
- Light intensity (Lux): High 30,000 Low 8,000
- Focal length: 27" (70 cm)
- Illumination area: 100mm x 180mm
- Sensor mode, manual toggle on/off

Standard Control Features:

- 3 HP Automatic selection
- 3-way Syringe
- Vacuum package with saliva ejector and hve
- Self contained water bottle system
- wet / dry foot control
- 1 Year warranty
- **Requires central Air / vacuum systems**
- **Includes utility center and junction box**

Optional:

Self contained compressor & Vacuum
(No plumbing required)



L550-LED


OR 1200

ORP1200-550LED



OMC2375

Optional OMC 2200
Monitor Bracket



Counter top dimensions

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
Light Arm: L600 26 x 14 x 8 Light Head: L550 17 x 15 x 7	30 lbs.
Cabinet 32 x 28 x 41	140 lbs.



Optional: Pony Tail
headrest



Pull out writing surface



Optional: OMC-
2200 LCD Bracket

OR-1200

Orthodontic Package

Description	Model	Suggested Retail Price
ORTHO 1200 Orthodontic Package	ORP1200-550LED-2375	\$7,736.00
	ORP1200-550LED-2375CV	\$9,264.00
	ORP1200-550LED-2375CV-SL	\$10,146.00

Package includes the following:

		2375	2375-CV	2375CV-SL
OR 1200	Fixed base electromechanical patient chair with wide back and full lumbar support, double articulating headrest with sliding height adjustment, One Armrest included foot control & last position Memory	✓	✓	✓
OMC-2375	This orthodontic mobile delivery cabinet offers a pullout writing surface, two instrument drawers, one medium and two large storage drawers for materials and supplies that fully extend and quickly remove for easy cleaning. In addition, these cabinets provide convenient access from the front. Oversized locking dual-casters.	✓		
OMC2375CV	OMC2375+Compressor+Vacuum(Self contained)		✓	
OMC2375CV-SL	OMC2375+Compressor+Vacuum+LED Curing light+Scaler(Self contained)			✓
L550-LED	Lustrous LED post mounted light with dual intensity, 8,000 - 30,000 LUX, dual axis light head rotation	✓	✓	✓

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

Retail


WIC840	Mirage flat screen monitor bracket with 16" extension	\$210.00
MO17	17" LCD multi media flat screen monitor	\$540.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply	\$270.00
OMC2200	LCD monitor bracket	\$330.00
OMC-ST81	Optional Stone Surface on OMC 2375, Color Artic	\$418.00
OMC-ST82	Optional Stone Surface on OMC 2375, Color White	\$418.00
OMC-ST85	Optional Stone Surface on OMC 2375, Color Cream	\$418.00
OMC-ST86	Optional Stone Surface on OMC 2375, Color Bamboo	\$418.00
L2001	Laguna L2000-2.0 Chair Bracket	\$336.00
AS-5047	Laguna Assistants stool	\$590.00
DR-5045	Laguna Doctor's stool	\$480.00
PTH 101	Pony Tail Headrest	\$218.00



Optional: Stone Counter top
OMC-ST81



Optional: Stone Counter top
OMC-ST82



Optional: Stone Counter top
OMC-ST85



Optional: Stone Counter top
OMC-ST86

Bellagio

Chair Mounted Delivery System



Standard Features

- Asepsis automatic control for three handpieces
- Individual handpiece pressure adjustments
- Non-retracting water coolant with on/off
- Individual water coolant flow controls
- Air coolant flow control
- Handpiece tubing flush system
- Ultra soft handpiece tubing
- 5-Hole straight gray asepsis tubing with Midwest connector.
- Disc foot control wet/dry toggle with 7' of smooth tubing
- Equipment utility center with asepsis housing
- Vacuum canister with solids collector screen
- Water quick disconnect

- Self contained water system with city/bottle selector switch
- Straight light post for unit mount dental light
- Pre-wiring for fiber optic light source
- Air-actuated counterbalanced arm lock
- 3 way syringe on unit head
- Integrated touch pad control
- Integrated piezo scaler
- Built in LED Curing Light
- Integrated timer

Model #5000-2.0



Electric Drain Cuspidor

- Contoured bowl is made of fine porcelain for easy cleaning and infection control.
- The cuspidor operates with a timed bowl rinse and cup filler
- Swing away cuspidor bowl
- Electric solenoid water flow controlled

Assistant's Instrumentation on Rotating arm

- 1 autoclavable HVE with asepsis tubing
- 1 autoclavable saliva ejector with asepsis tubing
- Assistant's side 3-way syringe
- Touch pad control

Options

- Built In Electric Handpiece
- H.P. Illumination light system
- Post-mounted utility tray
- 17" LCD Monitor



Control Head: Easy access controls for simple adjustment. Ultra soft handpiece tubing minimizing hand and wrist stress.



Cuspidor Bowl: The deep contoured bowl is made of fine porcelain for easy cleaning and infection control.



Optional: Multi media flat screen monitor for chair side intraoral imaging, patient education and DVD viewing.



Built in self contained 1500ml water bottle system including city/bottle selector switch.



Assistant's instrument



More space available for use

Bellagio

Chair Mounted Delivery System

Model #5015-2.0



Description	Model	Suggested Retail Price
Bellagio Chair Mounted Delivery System with Cuspidor	5000-2.0	\$5,690.00
Bellagio Chair Mounted Delivery System	5015-2.0	\$4,990.00

Package includes the following:

(with cuspidor)

(without cuspidor)

Asepsis 3 handpiece automatic control with coolant controls, flex arm with dual air brake, tri block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves. Includes 3 way syringe and light post with cutout for internal wires and monitor bracket, 7' umbilical with touch pad control



Porcelain cuspidor with sensor controlled bowl rinse and cup filler, rotating 4 position assistant's package, Hve and saliva ejector, water quick disconnect, solids collector and clean water system with touch pad control



Rotating 4 position assistant's package, Hve and saliva ejector,



water quick disconnect, solids collector and clean water system, light Post



Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

Retail

WIC800	Standard flat screen monitor bracket	\$128.00
WIC840	Mirage flat screen monitor bracket with 16" extension	\$210.00
#MO17	17" LCD multi media flat screen monitor	\$540.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply	\$270.00
L700-LED	Vivid LED Operatory Light	\$1,590.00
L550-LED	Lustrous LED Operatory Light "with motion sensor, with manual toggle"	\$1,300.00
BELBRKT	BLMT chair mount bracket (Bel 20/Bel 50)	\$238.00
53016 -L-100	Disposable Solids Collector Trap, 2 1/8" (pack of 100)	\$45.00
EMB-901	Built In Electric Motor system	\$1,790.00
XV100	Super Thin X-ray Viewer	\$140.00
XV150-R	Xray viewer bracket (round cut-out) – mount to Mirage unit head (Round Arm)	\$80.00

Packing Dimensions (L x W x H)	Weight
Unit: 39 x 29 x 23	60 lbs.
Unit PMU: 43 x 25 x 12	40 lbs.
Cuspidor: 26 x 18 x 10	20 lbs.

Mirage 2.0

Chair Mounted Delivery System



Standard Features

- Asepsis automatic control for three handpieces
- Individual handpiece pressure adjustments
- Non-retracting water coolant with on/off
- Individual water coolant flow controls
- Air coolant flow control
- Handpiece tubing flush system
- Ultra soft handpiece tubing
- 6-pin straight gray asepsis tubing with Midwest connector. (5-Hole tubing optional)
- Disc foot control wet/dry toggle with 7' of smooth tubing
- Equipment utility center with asepsis housing
- Vacuum canister with solids collector screen
- Water quick disconnect
- Self contained water system with city/bottle selector switch
- Straight light post for unit mount dental light
- Pre-wiring for fiber optic light source
- Air-actuated counterbalanced arm lock
- 3 way syringe on unit head
- Touch pad controls



Model #2000-2.0

Electric Drain Cuspidor

- Contoured bowl is made of fine porcelain for easy cleaning and infection control.
- The cuspidor operates with a timed bowl rinse and cup filler
- Swing away cuspidor bowl
- Electric solenoid water flow controlled
- Touch pad control

Assistant's Instrumentation on Rotating Arm

- 1 autoclavable HVE with asepsis tubing
- 1 autoclavable saliva ejector with asepsis tubing
- Optional Assistant's side 3-way syringe
- Touch pad control

Options

- Built In Electric Handpiece
- Built-in Led curing light
- Built-in scaler
- Illumination light source
- Post-mounted utility tray
- 17" LCD Monitor
- LCD monitor bracket



Control Head: Easy access controls for simple adjustments. Ultra soft handpiece tubing minimizing hand and wrist stress.



Cuspidor Bowl: The deep contoured bowl is made of fine porcelain for easy cleaning and infection control.



Optional: Multi media flat screen monitor for chair side intraoral imaging, patient education and DVD viewing.



Built in self contained 1500ml water bottle system including city/bottle selector switch.



Optional: Assistant's side 3-way syringe



Built in water quick disconnect for ultrasonic scaling.

Mirage 2.0

Chair Mounted Delivery System

Model #2015-2.0

Optional: Assistant's side
3-way syringe shown



Optional: TP-2005
Touchpad control



Optional: TPE-2011
Touchpad control

Packing Dimensions (L x W x H)	Weight
Unit: 39 x 29 x 23	60 lbs.
Unit PMU: 43 x 25 x 12	40 lbs.
Cuspidor: 26 x 18 x 10	20 lbs.

Description

Mirage 2.0 Chair Mounted Delivery System with Cuspidor
Mirage 2.0 Chair Mounted Delivery System

Model

Suggested Retail Price

2000-2.0

2015-2.0

\$4,890.00

\$4,190.00

Package includes the following:

(with cuspidor)

(without cuspidor)

Asepsis 3 handpiece automatic control with coolant controls, flex arm with dual air brake, tri block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves. Includes 3 way syringe and light post with cutout for internal wires and monitor bracket, 7' umbilical with touch pad control



Porcelain cuspidor with sensor controlled bowl rinse and cup filler, rotating 4 position assistant's package, Hve and saliva ejector, water quick disconnect, solids collector and clean water system with touch pad control



Rotating 4 position assistant's package, Hve and saliva ejector,



water quick disconnect, solids collector and clean water system, light Post



Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

Retail

TPC3WAY-PRT	Additional autoclavable 3 way syringe with tubing	\$78.00
P2335+TB	Deluxe autoclavable Hve valve with tubing	\$67.00
WIC800	Standard flat screen monitor bracket	\$128.00
WIC840	Mirage flat screen monitor bracket with 16" extension	\$210.00
#MO17	17" LCD multi media flat screen monitor	\$540.00
A800N-1	Advance 800 piezo scaler – built in	\$570.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply	\$270.00
ALED55BN-A	LED 5W curing light – built in	\$438.00
L700-LED	Vivid LED Operatory Light	\$1,590.00
L600-LED	Lucent LED Operatory Light	\$1,500.00
L550-LED	Lustrous LED Operatory Light "with motion sensor, with manual toggle"	\$1,300.00
BELBRKT	BLMT chair mount bracket BLMT chair mount bracket (20 / 50)	\$238.00
TPCBRKT	TPC chair mount bracket	\$290.00
TP-2005	Touch pad control – unit head or dental chair mounted, 3 preset, auto return, last position	\$358.00
TPE-2011	Touch pad Control - Chari & Electric Motor System touch pad	\$460.00
53016 -L-100	Disposable Solids Collector Trap, 2 1/8" (pack of 100)	\$45.00
2000-C-2.0	Mirage Cuspidor	\$700.00
EMB-901	Built In Electric Motor system	\$1,790.00
XV100	Super Thin X-ray Viewer	\$140.00
XV150-R	Xray viewer bracket (round cut-out) – mount to Mirage unit head (Round Arm)	\$80.00

Mirage 2.0 – BLMT

Swing Mount Delivery System



Model #2450-2.0



Model #2400-2.0

Compatible with BLMT 7, 10, 20, 20N, 50 patient chairs



Optional: LP100 swing mount light post with optional L700-LED Operatory Light

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
3500: 36 x 28 x 39	100 lbs.
3550: 29 x 12 x 24	25 lbs.
LP100: 45 x 26 x 44	30 lbs.

Mirage 2.0 – BLMT

Swing Mount Delivery System

Description	Model	Suggested Retail Price
Mirage 2.0 BLMT Swing Mount Delivery System (Unit + Asst pkg)	MSP2400-2.0	\$6,200.00

Packages include the following:

2400-2.0	Rotating asepsis 3 Handpiece automatic control with coolant controls, flex arm with air brake, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves. Includes 3 way syringe and clean water system,	* Fits 7, 10, 20, 20N, 50
2450-2.0	Rear radius asepsis assistant's arm system, Hve and saliva ejector, water quick disconnect and solids collector	* Specify * Fits 7, 10, 20, 20N, 50 chair at time of order

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."		Retail
2400-2.0	BLMT– TPC Swing Mount Delivery System- 7, 10, 20, 20N, 50	\$4,900.00
2450-2.0-O	BLMT- TPC Rear Radius Vacuum Pkg- 20N, 50	\$1,300.00
2450-2.0-U	BBLMT-TPC Rear Radius Vacuum Pkg- 7, 10, 20	\$1,300.00
2450-BRKT-BEL-N	BLMT-TPC Rear Radius Bracket ONLY- 20N, 50	\$252.00
2450-BRKT-BEL	BLMT-TPC Rear Radius Bracket ONLY- 7,10,20	\$252.00
BELBRKT	BLMT Chair Mount Bracket	\$238.00
2400-BRKT-BEL	BLMT- TPC Radius Chair Bracket ONLY- 7, 10, 20, 20N, 50	\$550.00
BEL20N-FC	BLMT Foot Control- 20N, 50	\$390.00
BEL20-FC	BLMT Foot Control- 7,10,20	\$390.00
LP-100	Mirage swing mount light post (ordered w/unit)	\$770.00
LP-105	Swing mount monitor post with monitor mounting bracket (ordered w/unit)	\$960.00
L600-LED	Lucent LED Operatory Light	\$1,500.00
L690-LED	Lucent LED Ceiling Mounted Operatory Light	\$1,590.00
L550-LED	Lustrous LED Operatory Light	\$1,300.00
P2335+TB	Deluxe autoclavable Hve Long valve with tubing	\$67.00
WIC800	Standard flat screen monitor bracket	\$128.00
WIC840	Mirage flat screen monitor bracket with 16" extension	\$210.00
MO17	17" LCD multi media flat screen monitor	\$540.00
A800N-1	Advance 800 piezo scaler – built in	\$570.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply	\$270.00
ALED55BN-A	LED 5W curing light – built in to dental unit	\$438.00
XV100	Super thin X-ray film viewer	\$140.00
XV150-R	X-ray film viewer bracket (curved) for TPC Mirage Unit only	\$80.00
AS-1101	Mirage Assistants stool	\$818.00
DR-1102	Mirage doctor's stool	\$650.00

BLMT Compatible Components



#2450-BRKT-BEL
 • Fits BLMT Chairs 7, 10, 20
 • Includes bolt / washer set



#2450-BRKT-BEL-N
 • Fits BLMT Chairs 20N, 50
 • Includes bolts / washer kit



#BELBRKT
 • Fits BLMT Chairs
 • Includes bolts / washer set



#2400-BRKT-BEL
 • Fits BLMT Chair 7, 10, 20, 20N, 50
 • Includes bearing / washer kit
 • Used to mount TPC swing unit to belmont chair



Square Connector



#BEL20-FC
 • Fits BLMT Chairs 7, 10, 20
 • Includes wire harness



Flat Connector



#BEL20N-FC
 • Fits BLMT Chairs 20N, 50
 • Includes wire harness

Description

BLMT-TPC Radius Chair Bracket. Fits 7, 10, 20, 20N, Bel 50
 BLMT-Chair Mount Bracket
 BLMT-Chair Foot Control. Fits 7, 10, 20
 BLMT-Chair Foot Control. Fits 20N, 50
 BLMT-TPC Rear Radius Bracket only 7, 10, 20
 BLMT-TPC Rear Radius Bracket only 20N, 50
 BLMT-TPC Rear Radius Vacuum Pkg- 7, 10, 20
 BLMT-TPC Rear Radius Vacuum Pkg- 20N, 50

Model

2400-BRKT-BEL
 BELBRKT
 BEL20-FC
 BEL20N-FC
 2450-BRKT-BEL
 2450-BRKT-BEL-N
 2450-2.0-U
 2450-2.0-O

Suggested Retail Price


\$550.00
 \$238.00
 \$390.00
 \$390.00
 \$252.00
 \$252.00
 \$1,300.00
 \$1,300.00

Mirage 2.0

Mirage 2.0 Mobile Doctors Cart



Model #MC-501-2.0



Model #MC-500-2.0



Description	Model		Suggested Retail Price
Mirage 2.0 Mobile Doctors Cart	MC-500-2.0		\$2,500.00
Mirage 2.0 Mobile Doctors Cart with Vacuum Package		MC-501-2.0	\$2,800.00
Asepsis 3 handpiece automatic selection with bottle system	✓	✓	
3 way syringe with quick change autoclavable tip	✓	✓	
Master on/off toggle	✓	✓	
Air / Water pressure adjustment for each handpiece	✓	✓	
Pressure gauge	✓	✓	
Non-retracting water / air coolant adjustment and purge	✓	✓	
Asepsis handpiece tubing	✓	✓	
Painted H-frame with height adjustment, 27" to 39"	✓	✓	
Vacuum canister with removable screen 2 1/8"		✓	
1 autoclavable Hve and saliva ejector and tubing		✓	
Variable speed foot control	✓	✓	
10 feet of umbilical	✓	✓	
Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."			Retail
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply		\$310.00
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply		\$270.00
A800N-1	Advance 800 N piezo scaler – built in		\$570.00
53016-L-100	Disposable Solids Collector Trap, 2 1/8" Pk of 100		\$45.00
MC500-UMB-FT	Additional Umbilical Per Foot		\$30.00
MC501-UMB-FT	Additional Umbilical Per Foot		\$30.00

Model #MC-501-2.0

Model #MC-500-2.0

Packing Dimensions: 29" x 23" x 21" / 50 lbs


Packing Dimensions: 29" x 23" x 21" / 50 lbs

Mirage 2.0


Wall / Cabinet Mounted System




Model #MS3000-2.0



Model #MR4000-2.0



Model #MVP5000-2.0



Description

Suggested Retail Price

Mirage Cabinet Mounted Side Delivery Unit	MS-3000	\$3,300.00
Mirage Cabinet Mounted Rear Delivery Unit	MR-4000.2.0	\$2,200.00
Mirage Cabinet Mounted Vacuum Package	MVP-5000.2.0	\$950.00
Asepsis 3 handpiece automatic selection with bottle system	✓	✓
3 way syringe with quick change autoclavable tip	✓	✓
Master on/off toggle	✓	✓
Air / Water pressure adjustment for each handpiece	✓	✓
Pressure gauge	✓	✓
Non-retracting water / air coolant adjustment and purge	✓	✓
Asepsis handpiece tubing	✓	✓
Asepsis flex arm system, w/ 9" rigid arm, pneumatic brake	✓	
2" diameter wall bracket	✓	
Pivoting arm mounted vacuum canister w/ screen		✓
3 way syringe with quick change autoclavable tip		✓
Water quick disconnect		✓
Vacuum canister with removable screen		✓
4 position holder, 1 autoclavable Hve and saliva ejector & tubing		✓
Variable speed foot control	✓	✓

Options:

Retail

H7660	Illumination system – ISO C (6 pin) – tubing, control box	\$270.00
H7000	Illumination system – ISO B (5 hole) – tubing, control box, bulb	\$310.00
A800N	Advance 800 piezo scaler – built in	\$570.00
P2335+TB	Deluxe autoclavable Hve valve with tubing	\$67.00
53016-L-100	Disposable Solids Collector Trap, 2 1/8" (Pack of 100)	\$45.00
Packing Dimensions: 17" x 21" x 51" / 60 lbs	Packing Dimensions: 19" x 15" x 20" / 25 lbs	Packing Dimensions: 13" x 24" x 11" / 15 lbs
MS-3000	MR-4000	MVP-5000

Mirage


Panel & Cabinet Mount Systems



Model #MP-2600
Mirage Panel Mount
Rear Delivery System

Model #MR-4500 L/R
Mirage Rear Cabinet Mount
Delivery System with Assistants Package

Model #MP-2601
Mirage Wall / Cabinet Mount
Rear Assistants package



Standard Features

- 2 Handpiece Tubing Positions
- 3 Position HP Holder
- TPC Press ring style syringe
- Individual drive air & water coolant controls
- Foot control
- Control Dimensions 8 1/2" x 1 3/4" x 5"
- Cut-out dimensions 7 1/2" x 11 1/2"
- Syringe has adjustable air & water controls
- Air / Water master control with Regulator.

Standard Features

- Asepsis 3 handpiece automatic selection with bottle system
- 3 way syringe with quick change autoclavable tip
- Master on/off toggle
- Air / Water pressure adjustment for each handpiece
- Pressure gauge
- Non-retracting water / air coolant adjustment and purge
- Asepsis handpiece tubing
- Self contained bottle system
- Stainless Steel tray
- Vacuum canister with removable screen
- 4 position holder, 1 autoclavable Hve and saliva ejector & tubing
- Variable speed foot control

Standard Features

- 2 HVE Suciton Valves
- 1 SE Suction Valve
- 3 Position Suction Holder
- Vacuum canister with removable screen

Description	Model	Suggested Retail Price
Mirage Panel Mount Rear Delivery System (2 HP manual select, 3-Way Syringe, 3 position HP holder, Utility Valves)	MP-2600	\$990.00
Mirage Wall & Cabinet Mount Assistance Package (2 HVE, 1 SE, 3 position holder)	MP-2601	\$350.00
Mirage Rear Cabinet Mount Delivery System w/ Assistance Package (2 HVE, 1 SE, 3 position holder)	MR-4500 L/R	\$2,990.00

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

		Retail
H7660	Illumination system – ISO C (6 pin) – tubing, control box	\$270.00
H7000	Illumination system – ISO B (5 hole) – tubing, control box, bulb	\$310.00
A800N	Advance 800 piezo scaler – built in	\$570.00

Packing Dimensions: 9" x 17" x 15" / 9 lbs
MP-2600

Packing Dimensions: 8" x 17" x 15" / 5 lbs
MP-2601

Packing Dimensions: 19" x 15" x 20" / 30 lbs
MP-4500 L/R

OMC 2375

Orthodontic Mobile Delivery Cabinet



Optional OMC 2200
Monitor Bracket



Rear compartment
4" W X 12"D X 20" H



OMC2375

Writing surface with
optional solid surface



Artic

OMC-ST81



White

OMC-ST82



Cream

OMC-ST85



Bamboo

OMC-ST86

Optional hard surface stone top

Function:


This orthodontic mobile delivery cabinet offers a pullout breadboard, two instrument drawers, one medium and two large storage drawers for materials and supplies that fully extend and quickly remove for easy cleaning. In addition, these cabinets provide a large open space for accessories with convenient access from the front. The oversized dual-caster wheels, makes positioning of the cabinet simple and easy on all floor surfaces, including carpet.

Standard Cabinet Features:

- Five easy to remove and clean full extension drawers, and one pull-out writing surface.
- High quality hardware throughout
- Oversized caster wheels with LOCKS
- Four stationary handles for easy cart positioning
- A removable counter - top plastic tray
- Standard 12" X 12" Junction Box
- 7' Length of umbilical

Standard Control Features:

- 3 HP Automatic selection
- 3-way syringe
- Vacuum Package with Saliva Ejector and HVE
- Self contained water bottle system
- Variable speed disc-type foot control
- Water coolant on/off switch with flow control
- 1 Year warranty
- Includes Utility center + junction box
- Requires central air and vacuum connections



Description

Model

Suggested Retail Price

Orthodontic Mobile Delivery Cabinet	OMC2375	\$3,190.00
-------------------------------------	---------	------------

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

	Retail
A800N-1 Advance 800N Piezo Scaler - Built in	\$530.00
ALED55BN-A LED 55BN 5W LED Curing Light -Built in	\$438.00
2375-UMB-FT Umbilical Per Foot	\$30.00
OMC2200 LCD Monitor Bracket W/ Post	\$418.00
OMC-ST81 Optional Stone Surface on OMC 2375, Color Artic	\$418.00
OMC-ST82 Optional Stone Surface on OMC 2375, Color White	\$418.00

Packing Dimensions: 41" x 28" x 32" / 140 lbs

OMC 2375CV-SL

Orthodontic Mobile Delivery Cabinet



Optional OMC 2200
Monitor Bracket



OMC2375CVSL



Rear compartment
4" W X 12" D X 20" H

Function:


This self-contained orthodontic mobile delivery cabinet offers a pullout breadboard, two instrument drawers, one medium and two large storage drawers for materials and supplies that fully extend and quickly remove for easy cleaning. In addition, this cabinets provides a large open space for accessories with convenient access from the front. The oversized dual-caster wheels, makes positioning of the cabinet simple and easy on all floor surfaces, including carpet.

Standard Cabinet Features:

- Five easy to remove and clean full extension drawers, and one pull-out writing surface.
- High quality hardware throughout
- Oversized caster wheels with LOCKS
- Four stationary handles for easy cart positioning
- A removable counter - top plastic tray

Standard Control Features:


- 3 HP Automatic selection
- 3-way syringe
- Vacuum Package with Saliva Ejector and HVE
- Self contained water bottle system
- Variable speed disc-type foot control
- Water coolant on/off switch with flow control
- 1 Year warranty
- **Self contained dual compressor**



OMC-ST81



OMC-ST82



OMC-ST85



OMC-ST86

Optional hard surface stone top

Description

Orthodontic Mobile Delivery Cabinet (Self contained)
Orthodontic Mobile Delivery Cabinet (Self contained)

Model

OMC2375CV-SL
OMC2375CV(No Scaler or LED curing light)

Suggested Retail Price

\$5,600.00
\$4,700.00

Options:

OMC2200	LCD Monitor Bracket W/ Post	\$330.00
OMC-ST81	Optional Stone Surface on OMC 2375, Color Artic	\$418.00
OMC-ST82	Optional Stone Surface on OMC 2375, Color White	\$418.00
OMC-ST85	Optional Stone Surface on OMC 2375, Color Cream	\$418.00
OMC-ST86	Optional Stone Surface on OMC 2375, Color Bamboo	\$418.00

Packing Dimensions: 32" x 28" x 41" / 180 lbs CV-SL model



MC3600CV(SL)

Mirage Mobile Delivery Cart

- Moves with ease
- Efficient
- Simple maintenance
- No plumbing required
- Cost effective
- Compact



Clean Water Bottle



Waste Bottle



MC3600CVSL



Optional:

- #OMC-2260 Monitor Bracket
- #MO17 LCD Monitor
- #UC100 Utility Tray
- #AIC5888 Intraoral camera

Specifications

Electrical: (2) X 550W head (115)VAC, 50/60 Hz 6.0 Amps

Sound: 68 db

Motor: Oil less air compressor

Construction: Metal Case

Vacuum: 80" of water (6" of Hg)

Vacuum: 80" of water (6" of Hg)

Max psi: 115 PSI / 40L per min

Waste Canister: 750ml

Shipping: Freight only

Physical Dimensions

31"H x 20"W x 16.5"D

Weight: 96 lbs

The MC3600CV (SL) is a self contained mobile delivery cart suitable for use by the dentist in a location when it is not feasible to connect to central plumbing

Description	Model	Suggested Retail Price
Mirage Mobile Delivery Cart	MC3600CV	\$4,270.00
Mirage Mobile Delivery Cart (with built in scaler and led curing light)	MC3600CV-SL	\$5,170.00

Package Includes the following

- 3 way syringe, saliva ejector, hve
- 2 handpiece positions
- Variable speed disc-type foot control
- Self contained (dual compressor)
- Clean water Bottle
- Waste water Bottle
- Advance 800 Built in Piezo Scaler
- LED 39 Built in Curing Light
- Dual 110v outlet (rear of cart)
- 1 year warranty

✓
✓
✓
✓
✓
✓
✓
✓
✓
✓

✓
✓
✓
✓
✓
✓
✓
✓
✓
✓

Shipping Dimensions:

Packing Dimensions: 39" x 27" x 23" / 160 lbs

Options:

OMC-2260	LCD Monitor Bracket W/ Post
AIC3N1N1W	17" LCD Monitor W/AIC5888 Camera
US100	Utility Tray

Suggested Retail Price

\$330.00
\$1,500.00
\$70.00

PC 2635

Portable Delivery System



PC-2635

The **PC2635** is a portable, self contained delivery system suitable for use in-office, nursing homes, hospitals, or wherever any off-site dentistry is required. The PC2635 is a lightweight, portable system that is easily transported. Air and water supply are included on the PC2635 portable dental system. Optional Tornado High and Slow speed handpieces are available.

Description	Model	Suggested Retail Price
Portable Dental System	PC2635	\$1,790.00

Package Includes the following

• 4 hole or 2 hole instrument tubings (specify)	✓
• 3 way syringe	✓
• Self contained oil less compressor (580W / 7L)	✓
• Foot control	✓
• Automatic handpiece flush System	✓
• Built in carrying handle	✓
• Self contained water bottle system	✓
• 1 year warranty	✓

Shipping Dimensions:

Packing Dimensions: 23" x 17" x 22" / 75 lbs

PC 2630

Portable Delivery System

- PC2930: Includes
- LED Curing light
 - Piezo Scaler



PC-2630



PC-2930



The **PC2630** is a fully self contained air dental system. Quick and easy to setup and operate. Ready for travel. Tornado HP available upon request. Ask for details.

Description	Model	Suggested Retail Price
Portable Dental System	PC2630	\$2,390.00
Portable Dental System W/ Scaler & LED Curing Light	PC2930	\$3,300.00

Package Includes the following

• 4 hole or 2 hole instrument tubings (specify)	✓	✓
• 3 way syringe	✓	✓
• Saliva ejector	✓	✓
• Self contained water bottle system	✓	✓
• Self contained oil less compressor (580W / 7L)	✓	✓
• Foot control	✓	✓
• Air venturi vacuum system	✓	✓
• Transportation case with pull out handle.	✓	✓
• 1 year warranty	✓	✓
• Built In Scaler		✓
• Built In Curing Light		✓

Options:

	Retail
A800N-1 Advance 800 piezo scaler – built in	\$570.00
INSTALL-1 LABOR (To Install Build-In Piezo Scaler with Additional Holder on PC2630)	\$278.00
53012-INSTL Female Quick disconnect (factory installed)	\$160.00

Shipping Dimensions: pc2630

Packing Dimensions: 16" x 20" x 28" / 60 lbs

Shipping Dimensions: pc2930

Packing Dimensions: 16" x 20" x 28" / 65 lbs

PC 2700

Portable Dental Chair Package



Description	Model	Suggested Retail Price
Portable Dental Chair Package	PC 2700	\$1,996.00

(Includes Patient Chair, Light, Stool, Tray, Spittoon & heavy duty Carry bags)

1 Year Warranty

Options	Model	Retail
PC 2720	Portable Patient Chair with Carrying Bag	\$1,090.00
PC 2740	Portable Dr. Stool with Carrying Bag	\$480.00
PC 2760	Portable Operator Light chair mounted	\$260.00
PC 2780	Stainless Steel tray and spit cup	\$290.00
PC 2790	Portable chair heavy duty carry bag	\$190.00
PC 2795	Stool heavy duty carry bag	\$100.00
PC2750	Portable LED Operator Light with Tripod Base with Carrying Bag	\$378.00
PC2760-PS	Replacement Power supply for portable LED operator light	\$38.00
PC2798	Carry bag for portable LED operator light	\$50.00

Shipping Dimensions:	PC 2700	PC 2740	PC 2750
Packing Dimensions:	49" x 12" x 24" / 73 lbs	49" x 12" x 24" / 17 lbs	49" x 12" x 24" / 10 lbs.

PC 2530

Mobile Vacuum System



PC2530

- Moves with Ease
- Efficient
- Simple Maintenance
- Perfect Backup
- Cost Effective
- Compact

Specifications

Electrical:	7/8 Horsepower 115 VAC, 50/60 Hz 230 VAC, 50/60 Hz 4.0 Amps
Sound:	70 db
Motor:	Variable speed control
Construction:	Metal Case
Vacuum:	80" of water (6" of Hg)
Airflow:	75 cfm
Canister:	2 Liter
Dimensions:	14" L x 12" W x 31" H
Shipping:	Freight

The PC 2530 Mobile Vac System is a portable, self contained vacuum system suitable for use by the dentist in a temporary location or when it is not feasible to install a central vacuum system.

Description	Model	Suggested Retail Price
Mobile Vacuum System	PC2530	\$2,100.00

Package Includes the following

- 2 HVE and 1 SE Valves
- Drain Line
- 2 SE suction tips
- Power cord
- Solids trap
- 1 year warranty



Shipping Dimensions:

Packing Dimensions: 17" x 25" x 37" / 95 lbs

ProClean

Extra Oral Suction System



Benefits:

- Captures and destroys viruses, germs and microorganism aerosols and droplets generated from high speed drilling and ultrasonic scaling at chairside
- Captures mercury vapors from amalgam removal
- Captures BPA particulates from composite material
- Neutralizes noxious chemical disinfectants odors and VOC gases
- Neutralizes smoke plum during dental laser and electrosurgery procedures
- Captures air abrasion and micro-etching dust

Filtration Process:

- Moisture Absorption Filter: Absorption of medium sized water droplets, moisture, aerosols and oils
- 2 Layer HEPA Filter: H13 Fine particle filter. Filtration of fine and ultra-fine dust, bacteria, viruses, fungi and prions. Efficiency: >99.97% @ >0.3 microns
- UV-C Germicidal Lamp: Disinfection and sanitation of air beyond filtration. Power: 10 watts, Wavelength: 254.7 nm
- Cold Plasma Ionizer: Eliminates particles as small as 0.01 microns including viruses, odor molecules and cigarette smoke

Technical Specifications:

Power Requirements:	100-120V AC, 50/60 Hz (International: 220-240V AC, 50/50 Hz)
Power	1200 Watts
Suction Pressure:	29 Kpa
Suction Flow:	4,000L/min (5Kpa)
Noise Level:	<65dB
New Weight:	64 lbs
Air Intake:	Segmented rotary arm with removable hood- 71.5" length x 2.5" diameter
UC-C Lamp Sterilization:	254.7 nm, 10 Watt
Plasma Disinfection:	5,000V
HEPA Grade H13:	0.3 micron, > 99.97%

Features:

- Kills 99.97% of viruses, germs and microorganisms
- Adjustable suction power
- Segmented rotary arm with replaceable wide-mouth suction hood
- Medical grade filtration system
- 2-year warranty

Description	Model	Suggested Retail Price
ProClean Extra Oral Suction System	PC2800	\$1,590.00

Package Includes the following

• Power Cord	✓
• Operation Manual	✓
• 2 year warranty	✓

Options:

	Retail
#PC2800-SYNFLTR	Replacement washable synthetic filter
#PC2800-HEPA13SML	Replacement HEPA-13 Small Filter
#PC2800-UVLAMP	Replacement UV Lamp
#PC2800-PLSMAGEN	Replacement Plasma Generator
#PC2800-HEPA13LG	Replacement HEPA-13 Large filter
#PC2800-Motor	Replacement Suction Motor
	\$9.00
	\$69.00
	\$69.00
	\$69.00
	\$95.00
	\$390.00


Shipping Dimensions:

ProClean : 29" x 19" x 38" / 85 lbs



Vivid LED

Operatory Light



Model #L700-LED

Model L795-LED
Ceiling mounted

Packing Dimensions (L x W x H)	Weight
L900-08: 15 x 3 x 3	3 lbs.
L900-09: 30 x 3 x 3	5 lbs.
L900-10: 44 x 4 x 4	7 lbs.

#LP48 Light Post must be ordered if dental light is ordered separately

Specifications:

Dimmer mode: From 3,000-35,000 Lux

Input voltage: 110V/220V, 60/50 Hz

Output voltage: 24 Volts AC / DC

Light intensity (Lux): 3,000 - 35,000 Lux @ 700mm

Cooling System: Static Heat sink (no fan)

Illumination area: 100mm x 180mm @ 700mm

Description

Vivid LED Post Mount Operatory Light With Motion Sensor + Toggle

Vivid LED Ceiling Mounted Light With Motion Sensor + Toggle
(specify 8', 9' or 10' ceiling)

Model

L700-LED

Suggested Retail Price

\$1,590.00

L795-LED

\$2,360.00

Package includes the following:

Input voltage 24 Vac 50/60Hz

2 year warranty

Manual toggle / Motion sensor standard

PoliBlock Filter mode

✓

✓

✓

✓

✓

✓

✓

✓

Options:

RLL-PS 110 - 240 V Power supply in 24 V out

LP-48 48" light post

L690-08 Ceiling post for 8' ceiling height

L690-09 Ceiling post for 9' ceiling height

L690-010 Ceiling post for 10' ceiling height

L690-011 Ceiling post for 11' ceiling height (special order)

Retail

\$70.00

\$210.00

\$220.00

\$220.00

\$220.00

\$276.00

Packing Dimensions (L x W x H)	Weight	Packing Dimensions (L x W x H)	Weight
L700-LED: 38 x 17 x 15	10 lbs.	Light Head: 14 x 9 x 14	7 lbs.
L795-LED: 38 x 17 x 15	18 lbs.		



Lustrous LED


Operatory Light



Model #L550-LED



Model L570-LED
Ceiling mounted



Packing Dimensions (L x W x H)	Weight
L900-08: 15 x 3 x 3	3 lbs.
L900-09: 30 x 3 x 3	5 lbs.
L900-10: 44 x 4 x 4	7 lbs.

#LP48 Light Post must be ordered if dental light is ordered separately

Specifications:

Dimmer mode:	From 8,000-30,000 Lux	Light intensity (Lux):	8,000 - 30,000 Lux @ 700mm 5,000 Kelvin
Input voltage:	110V/220V, 60/50 Hz	Cooling System:	Static Heat sink (no fan)
Output voltage:	24 Volts AC / DC	Illumination area:	100mm x 180mm @ 700mm

Description	Model	Suggested Retail Price
Lustrous LED Post Mount Operatory Light	L550-LED	\$1,300.00
Lustrous LED Ceiling Mounted Light (specify 8', 9' or 10' ceiling)	L570-LED	\$1,960.00

Package includes the following:

Input voltage 24 Vac 50/60Hz	✓	✓
2 year warranty	✓	✓
Manual toggle / Motion sensor standard	✓	✓
PoliBlock Filter mode	✓	✓

Options:

		Retail
RLL-PS	110 - 240 V Power supply in 24 V out	\$70.00
LP-48	48" light post	\$210.00
L690-08	Ceiling post for 8' ceiling height	\$220.00
L690-09	Ceiling post for 9' ceiling height	\$220.00
L690-010	Ceiling post for 10' ceiling height	\$220.00
L690-011	Ceiling post for 11' ceiling height (special order)	\$276.00

Packing Dimensions (L x W x H)	Weight
L550: 17 x 15 x 7	7 lbs.
L570: 17 x 15 x 7	7 lbs.
L600 Arm: 34 x 11 x 5	10 lbs.
L690 Arm: 32 x 18 x 6	18 lbs.




Lucent LED

Operatory Light



Model #L600-LED



Model L690-LED
Ceiling mounted
Lucent LED
operatory light.

Packing Dimensions (L x W x H)	Weight
L900-08: 15 x 3 x 3	3 lbs.
L900-09: 30 x 3 x 3	5 lbs.
L900-10: 44 x 4 x 4	7 lbs.

#LP48 Light Post must be ordered if dental light is ordered separately
Specifications:

Absorbed Power:	Max 8vA	Light intensity (Lux):	3,000 - 35,000 Lux @ 700mm
Input voltage:	110V/220V, 60/50 Hz	Cooling System:	Static Heat sink (no fan)
Output voltage:	24 Volts AC / DC	Illumination area:	100mm x 180mm @ 700mm

Description	Model	Suggested Retail Price
Lucent LED Post Mount Operatory Light	L600-LED	\$1,500.00
Lucent LED Post Mount Operatory Light with motion sensor	L625-LED	\$1,700.00
Lucent LED Ceiling Mounted Light (specify 8', 9' or 10' ceiling)	L690-LED	\$2,160.00
Lucent LED Ceiling Mounted Light (specify 8', 9' or 10' ceiling) with motion sensor	L695-LED	\$2,360.00

Package includes the following:

Input voltage 110V	✓	✓	✓	✓
2 year warranty	✓	✓	✓	✓
Dimmer Mode	✓	✓	✓	✓

Options:

RLL-H04	Replacement LED light optical bulb group	\$410.00
RLL-H06	Replacement light handle (Left and Right)	\$127.00
RLL-H01	Replacement shield	\$130.00
RLL-H03	Replacement reflector / set of 2	\$150.00
RLL-PS	110 - 240V Power supply in 24V out	\$70.00
LP-48	48" light post	\$210.00
L690-08	Ceiling post for 8' ceiling height	\$220.00
L690-09	Ceiling post for 9' ceiling height	\$220.00
L690-010	Ceiling post for 10' ceiling height	\$220.00
L690-011	Ceiling post for 11' ceiling height (special order)	\$276.00

Packing Dimensions (L x W x H)	Weight
L600: 34 x 11 x 5	10 lbs.
L690: 32 x 18 x 6	18 lbs.
Packing Dimensions (L x W x H)	Weight
Light Head: 14 x 9 x 14	7 lbs.



Radiant LED

Operatory Light



Radiant light
handles are removable
and sterilizable

R6101-LED

#LP48-RI Light Post must be ordered if dental light is ordered separately

Specifications:

Absorbed Power: Max 8vA

Input voltage: 110V/220V, 60/50 Hz

Output voltage: 24 Volts AC / DC

Bulb type: LED

Light intensity (Lux): 3,000 - 35,000 Lux @ 700mm

Cooling System: Static Heat sink (no fan)

Illumination area: 100mm x 180mm @ 700mm

Description

Model

Suggested Retail Price

Radiant LED Operatory Light

R6101-LED

\$1,990.00

Radiant LED Operatory Light with motion sensor

R6105-LED

\$2,190.00

Package includes the following:

Input voltage 110V - 240 Vac

✓

✓

2 year warranty

✓

✓

Dimmer Mode

✓

✓

Options:

Retail

RLL-H04 Replacement LED light optical bulb group

\$410.00

RLL-H06 Replacement light handle (Left and Right)

\$127.00

RLL-H01 Replacement shield

\$130.00

RLL-H03 Replacement reflector / set of 2

\$150.00

RLL-PS 110 - 240 V Power supply in 24 V out

\$70.00

LP-48-RI 48" light post

\$210.00

408500 Wall Mount Bracket

\$700.00

Packing Dimensions (L x W x H)

Weight

Arm: 17 x 6 x 36

14 lbs.

Head: 14 x 14 x 9

8 lbs.

Nylon Bushing

#LP-RI-S

• OD 1 7/8"

• ID 1 3/8"

– Page 60 –





Luminous LED

LED Operatory Light

with PoliBlock Filter



Luminous light handles are removable and sterilizable

#L7601-LED

#LP48-RI Light Post must be ordered if dental light is ordered separately
Specifications:

Absorbed Power:	Max 8vA
Input voltage:	110V/220V, 60/50 Hz
Output voltage:	24 Volts AC / DC
Bulb type:	LED

Light intensity (Lux):	3,000 - 35,000 Lux @ 700mm
Cooling System:	Static Heat sink (no fan)
Illumination area:	100mm x 180mm @ 700mm

Description	Model	Suggested Retail Price
Luminous LED Operatory Light	L-7601-LED	\$2,900.00
Luminous LED Operatory Light with motion sensor	L-7605-LED	\$3,190.00
Luminous LED Operatory Light Ceiling Mount	L-7901-LED	\$4,390.00

Package includes the following:

Input voltage 110V - 240 Vac



2 year warranty




Options:

LU-PB	Antipolymerization Filter "Poli-Block"
408500	Wall Mount Bracket

Retail

\$130.00
\$700.00



* Ceiling light extension post must be cut at the time of installation.

Nylon Bushing

#LP-RI-S	Included with models
• OD 1 7/8"	R6101-LED, R6105-LED
• ID 1 3/8"	



Packing Dimensions (L x W x H)	Weight
Light Arm: 17" x 6" x 36"	18 lbs.
Light Head: 14 x 14 x 9	

Mobile Cabinet Systems



Utility Mobile Cabinet
#TMC-180

- 14.5" W x 16" D x 32" H
- designed for ultrasonic scalers, curing lights, electrosurge, intraoral cameras, etc.
- 1 instrument drawer
- 2 adjustable shelves
- 1 year warranty
- Shipping dimensions: 20" x 22" x 34" / 65 lbs.



Doctor's Mobile Cabinet
#TMC-140

- 21.5" W x 19" D x 32" H
- 4 instrument drawers
- 3 storage drawers
- 1 year warranty
- Shipping dimensions: 33" x 22" x 24" / 110 lbs.



Assistant's "Alabama" Mobile Cabinet
#TMC-160

- 21.5" W x 19" D x 32" H
- Sliding counter top with amalgamator well
- 1 instrument drawer
- 2 storage drawers
- Side to side sliding counter
- 1 year warranty
- Shipping dimensions: 33" x 22" x 24" / 110 lbs.



Assistant's "North Carolina" Mobile Cabinet
#TMC-120

- 21.5" W x 19" D x 32" H
- Solid top with pull handle
- Amalgamator well drop leaf door
- 3 instrument drawers
- 1 storage drawer
- 1 year warranty
- Shipping dimensions: 33" x 22" x 24" / 110 lbs.



- White will be shipped if color choice not indicated on purchase order
- Mobile cabinets ship on pallet via freight carrier

Description	Model (light grey)	Model (white)	Suggested Retail Price
Utility Mobile Cabinet	TMC-180-G	TMC-180-W	\$478.00
Doctor's Mobile Cabinet	TMC-140-G	TMC-140-W	\$990.00
Assistant's "Alabama" Mobile Cabinet	TMC-160-G	TMC-160-W	\$990.00
Assistant's "North Carolina" Mobile Cabinet	TMC-120-G	TMC-120-W	\$990.00
Assistant's "North Carolina" Mobile Cabinet RT pull out		TMC-120R-W	\$990.00

Package includes the following:

External surfaces of high pressure laminate	✓	✓
with complimentary white interior surface	✓	✓
for easy cleaning. Smooth flowing ball bearing glides. Durable casters.	✓	✓
Color choice: white or light grey.	✓	✓
1 year warranty	✓	✓

Mirage

Operatory Stools

Model SS-5105



Model DR-1102



Model DR-7000



Dimensions:

Vertical adjustment range: 0" – 3.5"
 Seat width: 16"
 Seat height ranges: 24.5" – 31.8"
 5-Leg base diameter: 24"
 Shipping dimensions: 22" x 21" x 15" / 23 lbs.

Dimensions:

Backrest vertical adjustment range: 0" – 3.25"
 Seat width: 18"
 Seat height ranges: 19.5" – 24.8"
 5-Leg base diameter: 24"
 Shipping dimensions: 23" x 23" x 17" / 33 lbs.

Dimensions:

Backrest vertical adjustment range: 0" – 3.25"
 Seat width: 18.3"
 Seat height ranges: 19.5" – 24.8"
 5-Leg base diameter: 24"
 Shipping dimensions: 22" x 23" x 17" / 33 lbs.

Description

Model

Suggested Retail Price

Mirage saddle stool	SS-5105		\$690.00
Mirage doctor's stool		DR-1102	\$650.00
Solace doctor's stool		DR-7000	\$818.00

Stool Package includes the following:

Five-leg, ultra stable base, 24" diameter	✓	✓	✓
Gas cylinder lift for seat height adjustment	✓	✓	✓
Adjustable for right or left hand use	✓	✓	✓
Large, firm cushion seat	✓	✓	✓
Adjustable full ring foot rest			
Dual wheel casters	✓	✓	✓
Adjustable anatomical body support bar			
Height adjustment	23" - 31"	19" - 23"	19" - 23"
16 Upholstery color choices	✓	✓	✓
Seat tilt	✓	✓	✓
Optional adjustable armrest		✓	✓
1 year warranty	✓	✓	✓

Options *Specify Color And Product Code*

	Retail
DR1102-9 Replacement DR1102 Seat Cushion ONLY Duration Vinyl	\$178.00
DR1102-11 Replacement DR1102 Chair Back Cushion ONLY Duration Vinyl	\$140.00
DR1102-Vinyl-UC Upcharge Non Standard Vinyl DR1102 Upholstery upcharge	\$148.00
DR-1102-ULTR-UC Upcharge Ultra leather DR1102 Upholstery upcharge	\$280.00
DR-7000-Vinyl-UC Upcharge Non Standard Vinyl DR7000 Upholstery upcharge	\$148.00
DR-7000-ULTR-UC Upcharge Ultra leather DR-7000 Upholstery upcharge	\$290.00
SS-5105-Vinyl-UC Upcharge Non Standard Vinyl SS5105 Upholstery upcharge	\$160.00
SS-5105-ULTR-UC Upcharge Ultra leather SS5105 Upholstery upcharge	\$300.00

Mirage

Operatory Stool

Model AS-1101



Dimensions:

Backrest vertical adjustment range: 0" – 3.25"
 Seat width: 18.5"
 Seat height ranges: 23" – 30"
 5-Leg base diameter: 24"
 Packing dimensions: 24" x 23" x 21" / 42 lbs.

Model HY-5107



Dimensions:

Backrest vertical adjustment range: 0" – 3.25"
 Seat width: 17"
 Seat height ranges: 20" – 25"
 5-Leg base diameter: 24"
 Packing dimensions: 22" x 21" x 15" / 30 lbs.

Model DRSA-7100-DU100



Dimensions:

Backrest vertical adjustment range: 0" – 3.25"
 Seat width: 18" x 15.7"
 Seat height ranges: 20" – 25"
 5-Leg base diameter: 24"
 Packing dimensions: 24" x 23" x 22" / 45 lbs.
 Only Available in black DU-100

Description

Model

Suggested Retail Price

Mirage Assistant's stool	AS-1101		\$818.00
Mirage Hygiene stool		HY-5107	\$690.00
Mirage Mica Surgeon Stool		DRSA-7100-DU100	\$1,150.00

Stool Package includes the following:

Five-leg, ultra stable base, 24" diameter	✓	✓	✓
Gas cylinder lift for seat height adjustment	✓	✓	✓
Large, firm cushion seat	✓	✓	✓
Durable tile / carpet Casters	✓	✓	✓
Adjustable full ring foot rest	✓		
Adjustable for right or left handed use	✓		
Adjustable anatomical body support bar	✓		
Seat tilt		✓	✓
Height adjustment	22" - 29"	19" - 23"	19"-25"
Adjustable armrest			✓
Available in 4 colors (DU-18, DU-100, DU-116, DU-220)		✓	DU-100 Black only
1 year warranty	✓	✓	✓

Options *Specify Color And Product Code*

AS1101-16	Replacement AS1101 Seat Cushion ONLY Duration Vinyl	\$178.00
AS1101-21	Replacement AS1101 Chair Back Cushion ONLY Duration Vinyl	\$140.00
AS1101-18	Replacement AS1101 Abdominal Bar Cushion ONLY Duration Vinyl	\$156.00
AS1101-Vinyl-UC	Upcharge Non Standard Vinyl AS1101 Upholstery Upcharge	\$170.00
AS1101-ULTR-UC	Upcharge Ultra leather AS1101 Upholstery Upcharge	\$320.00

Laguna

Operatory Stool

Model DR-5045



Dimensions:
 Backrest vertical adjustment range: 0" – 3.25"
 Seat width: 18" x 15.7"
 Seat height ranges: 20" – 25"
 5-Leg base diameter: 24"
 Packing dimensions: 24" x 23" x 20" / 33 lbs.

Model AS-5047



Dimensions:
 Seat width: 18.5"
 Seat height ranges: 23" – 30"
 5-Leg base diameter: 24"
 Packing dimensions: 24" x 23" x 20" / 42 lbs.

Description

Model

Suggested Retail Price

Laguna Assistants stool

AS-5047

\$590.00

Laguna Hygiene stool

DR-5045

\$480.00

Upholstery Selections:



Sand
LK-03



Stormy
LK-04



Taupe
LK-08



Greystone
LK-11



Orange
LK-27



Cobalt
LK-31



Desert Rose
LK-34



Grey Teal
LK-49



Atlantis
LK-59



Black
LK-72

Stool Package includes the following:

Five-leg, ultra stable base, 24" diameter
 Gas cylinder lift for seat height adjustment
 Large, firm cushion seat
 Durable tile / carpet Casters
 Adjustable full ring foot rest
 Adjustable for right or left handed use
 Adjustable anatomical body support bar
 Seat tilt
 Height adjustment
 Adjustable armrest
 1 year warranty



22" - 29"



19" - 23"





Model XR-6101

Mirage

X-ray Chair

Color Chart:



Sand
LK-03



Stormy
LK-04



Taupe
LK-08



Greystone
LK-11



Orange
LK-27



Cobalt
LK-31



Desert Rose
LK-34



Grey Teal
LK-49



Atlantis
LK-59



Black
LK-72

Dimensions:

Seat width: 21"

Seat height ranges: 21.5" – 27.5"

Base diameter: 24"

Packing dimensions: 31" x 20" x 11" / 45 lbs.

8" x 8" x 18" / 15 lbs.

24" x 24" x 4" / 28 lbs.

Description

Mirage Xray chair

Model

XR-6101

Retail

\$1,550

X-ray stool Package includes the following:

Adjustable headrest



Gas cylinder lift for seat height adjustment



Pneumatic pump on base



Large, firm cushion seat



Large diameter base



Height adjustment



16 Upholstery color choices



Arm rest



Leg rest



1 year warranty



Options

XR-6101-8-C

Replacement XR6101 Duration Vinyl Seat Cushion ONLY

\$370.00

XR-6101-9-C

Replacement XR6101 Duration Vinyl Chair Back Cushion ONLY

\$370.00

XR-6101-10-C

Replacement XR6101 Duration Vinyl Headrest Cushion ONLY

\$120.00

XR-6101-Vinyl-UC

Upcharge Non Standard Vinyl XR6101 Upholstery Upcharge

\$340.00

XR-6101-ULTR-UC

Upcharge Ultra leather XR6101 Upholstery Upcharge

\$540.00

Suggested Retail Price

