

TPC



INNOVATION • RELIABILITY • QUALITY
2026-2027 Catalog U.S.A.

tpcdental.com

OUR VISION

A Preferred Company

To create a quality environment where people can work together to provide simple, cost effective creative solutions for the betterment of dentistry.

OUR MISSION

Continuous Improvement

We will manufacture innovative, high quality and cost effective dental products that enable distribution partners, dentists, hygienists, teaching institutions and government organizations throughout the world to perform efficient quality dental care. TPC provides innovative solutions and quality equipment that enables dental professionals to work together for the benefit of dentistry.

TABLE OF CONTENTS

Our Vision & Our Mission	Inside Front Cover
<u>Terms & Conditions</u>	3
<u>Chair, Unit, Light Packages</u>	4-11
<u>Chairs & Delivery Units</u>	12-34
<u>Air Compressors / Wet Ring Vacuum Systems</u>	35-37
<u>Mobile & Portable Dental Systems</u>	38-43
<u>Extra Oral Suction System</u>	44
<u>Stools & X-Ray Chairs</u>	45-48
<u>Operatory Lights</u>	49-53
<u>Mobile Cabinets</u>	54
<u>Imaging Products & Accessories</u>	55-63
<u>Dental Handpiece & Parts</u>	64-77
<u>Handpiece Cleaning & Lubrication Systems</u>	78
<u>Curing & Bleaching</u>	79-87
<u>Scaling Products & Accessories</u>	88-94
<u>Ultrasonic Bath</u>	95-98
<u>Digital Amalgamator</u>	99
<u>Micro Motor & Handpiece</u>	100
<u>Digital Thermometer</u>	101
<u>Deluxe Autoclavable Valves</u>	102
<u>Parts & Components</u>	103-109
<u>Custom Length Umbilical / Utility Centers</u>	110
<u>Order Form</u>	111
<u>Corporate Staff Contacts</u>	112
<u>Upholstery Color Chart</u>	113



Ordering Information:

When placing orders, please include the following information. All orders will be processed upon receipt of a faxed, mailed, or emailed hardcopy of the purchase order.

- Model and/or Part Number
- Requested Ship Date
- Description / Voltage
- Purchase Order Number
- Quantity
- Shipping Instructions
- Cancellations must be confirmed in writing
- All credit cards are accepted unless otherwise noted

Returned Goods:

All returns require prior authorization. Returns must be authorized by TPC within 7 days of date received by the customer.

- Contact Customer Service for Return Authorization Number
- Returned products are subject to a 25% restocking charge
- All returned goods must have Return Authorization Number
- All returned goods must be shipped prepaid
- All authorized returns must be received by TPC in original packaging, accessories, documentation and in original condition for credit considerations
- All outbound repair / warranty freight charges are the responsibility of the customer.
- **No returns accepted after 14 days.**

Special Orders:

- Special order items and custom made products are not returnable for credit

Minimum Orders:

- Orders less than \$50.00 will incur an additional \$5.00 handling charge.
- All drop ship orders will incur an additional \$5.00 fee.

Terms of Sale:

- Claims: All claims for merchandise received that are not in accordance with the original purchase order must be made within five (5) working days of the receipt of the merchandise. All merchandise is quality inspected and packaged in perfect condition. The common carrier assumes full responsibility for the safe delivery of all merchandise after it is removed from the F.O.B. shipping point. All claims for shortages, loss or damage must be documented and filed promptly with the freight carrier. Concealed damage must be reported to the freight carrier immediately upon discovery, with a request for inspection. If a request is made more than 15 days after the receipt of the merchandise, reasonable evidence must be submitted indicating that damage could not have occurred after delivery.
- All pricing in catalog is subject to change. Verify all pricing prior to placing order.

Warranty Information (All stated warranties are intended for USA market ONLY):

All of our products sold are guaranteed to be free from defects in workmanship and materials for one year from date of purchase, unless otherwise stated. TPC will repair or replace any defective part at no charge. TPC will not be responsible for labor charges or shipping charges to / from the TPC facility. This guarantee does not cover normal wear or stains on surface finish. The guarantee does not cover damage resulting from improper installation, misuse or accidents incurred in shipping and handling. All claims against the freight carrier must be initiated at the time the damaged items are received. The claim is the responsibility of the customer. We are improving our products on a continuous basis. We reserve the right to make modifications without the need for prior notification.

Products:

Patient Chairs

Curing & Bleaching Systems

Imaging Products

Scaling Products

Air Compressor

Warranty Period:

5 year limited

2 years

1 year

1 year

5 year limited

Products:

Delivery Units

Seating

Mobile Cabinets

Warranty Period:

5 year limited

1 year

1 year

Extended Manufacturer Warranty (selective models):

2 Year Warranty

H6000, H6005, H6025, H7005 Lubrication System

LED Curing Lights (All Models)

Advance CAM Intraoral Camera System(s)

LED Operatory Lights

Limited Manufacturer Warranty (selective models):

6 Months Warranty

Fiber optic systems, S333, M333, QS222, QM222,

ELCA760, ELS736.

Superb Vac

Wet-Ring Vacuum



Model	WV3	WV3R	WV5	WV5R
Water consumption without water recycling:	4L/min	/	4L/min	/
Water consumption with TPC water recycling:	/	0.8L/min	/	0.8L/min
Suction Rate	900L/min	900L/min	1500L/min	1500L/min
Vacuum Level (1N Hg) Factory set to 10HG	8-15HG	8-15HG	8-15HG	8-15HG
Noise Level dB (A)	58	58	58	58
Horse Power	1.2HP	1.2HP	2.0HP	2.0HP
Voltage	110V or 220V	110V or 220V	220V	220V
Frequency	50/60Hz	50/60Hz	50/60Hz	50/60Hz
Max.users	3	3	5	5
Current (A)	8A	8A	13A	13A
Circuit Breaker (A)	20A	20A	20A	20A
Inlet Connection Size	3/4"	3/4"	3/4"	3/4"
Drain Connection Size	1.25"	1.25"	1.25"	1.25"
City Water Connection Size	3/8"	3/8"	3/8"	3/8"
Air Exhaust Connection Size	2 3/8"	2 3/8"	2 3/8"	2 3/8"
Inlet Water Pressure (PSI)	10Psi/0.7Bar	10Psi/0.7Bar	10Psi/0.7Bar	10Psi/0.7Bar
Ambient Temperature	2 C ~40 C	2 C ~40 C	2 C ~40 C	2 C ~40 C
Product Size (LxWxH) Inches	17x14x18	17x14x18	17x14x18	17x14x18
Product sit Weight	62	62	66	66
Dimension	17x18x20	21x23x23	17x18x20	21x23x23

Features:

- Durable aluminum-bronze construction
- Stainless steel impeller
- Computer balanced, sealed bearings in motor
- Large capacity, easy to clean filter with clear bowl for easy inspection
- Available in single and twin models to accommodate 1-5 users
- 5-year full product limited warranty
- Unique water recycling to save up to 80% water
- UL/cUL classified



WV3R/WV5R

Description	Model	Suggested Retail Price
Superb Vac 1.2 HP	WV3	\$3,190.00
Superb Vac 1.2 HP (with recycler)	WV3R	\$3,590.00
Superb Vac 2 HP	WV5	\$3,590.00
Superb Vac 2 HP (with recycler)	WV5R	\$3,990.00

Product includes the following:

3 Simultaneous users	✓	✓	✓	✓
5 Simultaneous users			✓	✓
Water recycler & Air/Water Separator		✓		✓



WV3/WV5

Description	Model	Suggested Retail Price
#WV-AWS	Single Pump Air/Water Separator	\$250.00
#CP-3	Control Panel – 3 Button- Vac / Air / Water	\$480.00
#CP-4	Control Panel – 4 Button- Vac1 / Vac2 / Air / Water	\$560.00
#CV22	Air Compressor and Vacuum Pump Stacking Rack	\$900.00
#CP-2 AV	Control panel - 2 Button-Air/Vac	\$400.00
#CP-3 AV	Control panel - 3 Button-Air/Vac1/Vac2	\$480.00



Superb Vac

Wet-Ring Vacuum



WV6(R)/WV10(R)

Model	WV6	WV6R	WV10	WV10R
Water consumption without water recycling:	6 L/min	/	6 L/min	/
Water consumption with TPC water recycling:	/	1.01L/min	/	1.01L/min
Suction Rate	1800L/min	1800L/min	3000L/min	3000L/min
Vacuum Level (1N Hg) Factory set to 10HG	8-15HG	8-15HG	8-15HG	8-15HG
Noise Level dB (A)	65	65	65	65
Horse Power	2.4HP	2.4HP	4.0HP	4.0HP
Voltage	220V	220V	220V	220V
Frequency	50/60Hz	50/60Hz	50/60Hz	50/60Hz
Max.users	6	6	10	10
Current (A)	8A x 2	8A x 2	13A x 2	13A x 2
Circuit Breaker (A)	20A x 2	20A x 2	20A x 2	20A x 2
Inlet Connection Size	1"	1"	1"	1"
Drain Connection Size	1.25"	1.25"	1.25"	1.25"
City Water Connection Size	3/8"	3/8"	3/8"	3/8"
Air Exhaust Connection Size	2 3/8"	2 3/8"	2 3/8"	2 3/8"
Inlet Water Pressure (PSI)	10Psi/0.7Bar	10Psi/0.7Bar	10Psi/0.7Bar	10Psi/0.7Bar
Ambient Temperature	2 C ~40 C	2 C ~40 C	2 C ~40 C	2 C ~40 C
Product Dimension (LxWxH) Inches	28x17x18	28x17x18	28x17x18	28x17x18
Product Weight (LBS)	176	176	198	198
Shipping Dimension	32x25x24	32x25x24	32x25x24	32x25x24
Shipping Weight	188Lbs	188Lbs	208Lbs	208Lbs

The TPC Superb Vac wet ring vacuum pump offers a quiet operation for small to medium size dental practices. Our twin TPC Superb Vac wet ring vacuum pump provides vacuum suction when operating in peak work periods. The twin pump offers backup protection to minimize downtime and costly interruption to work flow.

Features:

- Durable aluminum-bronze construction
- Stainless steel impeller
- Computer balanced, sealed bearings in motor
- Large capacity, easy to clean filter with clear bowl for easy inspection
- Available in single and twin models to accommodate 1-10 users
- 5-year full product limited warranty
- Unique water recycling to save up to 80% water
- UL/cUL classified
- One or Two motors can be activated to fit 1-10 chairs use

Description	Model	Suggested Retail Price
Superb Vac 2.4 HP (2x1.2HP)	WV6(Twin)	\$6,590.00
Superb Vac 2.4 HP (2x1.2HP) (with recycler)	WV6R(Twin)	\$6,990.00
Superb Vac 4 HP (2x2HP)	WV10(Twin)	\$7,390.00
Superb Vac 4 HP (2x2HP) (with recycler)	WV10R(Twin)	\$7,790.00

Product includes the following:

6 Simultaneous users	✓	✓	✓	✓
10 Simultaneous users			✓	✓
Water recycler		✓		✓
220V 50/60Hz	✓	✓	✓	✓
Adjustable vacuum level from 8-15Hg	✓	✓	✓	✓
Low voltage control panel connections	✓	✓	✓	✓
Air Water Separator	✓	✓	✓	✓

Model	Description	Suggested Retail Price
#CP-3AV	Control Panel – 3 Button- Vac1 / Vac2 / Air	\$480.00
#CP-4	Control Panel – 4 Button- Vac1 / Vac2 / Air / Water	\$560.00
#CV22	Air Compressor and Vacuum Pump Stacking Rack	\$900.00

Superb Air

Oil-less Desiccant Dryer Air Compressor



DC6224 D/A



The TPC Superb Air Oil-less Air Compressors meet the high standards of the dental professionals, providing the reliability you need from a brand you can trust. The new motor design delivers a reduced noise with a 3-year warranty. The **Desiccant air drying** system efficiently filters out compressed air contaminants and moisture from the compressed air. Purge tank removes residual humidity in the compressor following each cycle.

Features:

- New motor design to reduce noise and increase reliability
- Quiet operation
- **Desiccant air drying** system to remove moisture from compressed air
- Energy efficient thermally protected motors
- Superb Air Oil-less compressors are available in single or twin heads
- 5-year warranty

Model User Capacity	DC6112D/A 1-3users	DC6122D/A 1-3users	DC6214 D/A 3-4users	DC6224D/A 5users	DC6326D/A 6-7users	DC6228D/A 8-9users	DC6428D/A 8-9users
CDFM@80psi	4.9	4.9	9.8	9.8	14.7	19.6	19.6
# of Heads	1	1	2	2	3	2	4
Total HP	1.15	1.15	2.30	2.30	3.45	4.80	4.60
Voltage	110V	230V	110V	230V	230V	230V	230V
Amps	8	4	16	8	12	16	16
Breaker Amps	20	20	30	20	30	30	30
dBA	62	62	62	62	72	72	72
Tank Capacity	10 Gallon	10 Gallon	16 Gallon	16 Gallon	23 Gallon	30 Gallon	30 Gallon
Size (LxWxH)	16x23x27	16x23x27	26x21x29	26x21x29	34x22x29	41x22x36	41x22x36
Shipping Dimension (WxDxH)	29x23x35	29x23x35	34x26x35	34x26x35	39x25x38	46x27x41	46x27x41

Description

Model

Suggested Retail Price

Superb Air Compressor, 1.15HP, 110V	DC6112	\$3,390.00
Superb Air Compressor, 1.15HP, 230V	DC6122	\$3,390.00
Superb Air Compressor, 2 x 1.15HP(2.3HP), 110V	DC6214	\$4,990.00
Superb Air Compressor, 2 x 1.15HP(2.3HP), 230V	DC6224	\$4,990.00
Superb Air Compressor, 3 x 1.15HP(3.45HP), 230V	DC6326	\$6,390.00
Superb Air Compressor, 2 x 2.4HP(4.8HP), 230V	DC6228	\$7,790.00
Superb Air Compressor, 4 x 1.15HP(4.6HP), 230V	DC6428	\$7,940.00

Description

Model

Suggested Retail Price

Control Panel 3 button-Vac / Air / Water	CP-3	\$480.00
Control Panel 4 button-Vac1 / Vac2 / Air / Water	CP-4	\$560.00
Control Panel 2 button-Vac/Air	CP-2AV	\$400.00
Control Panel 3 button-Vac1/Vac2/Air	CP-3-AV	\$480.00
Air Compressor and Vacuum Pump Stacking Rack	CV22	\$900.00

Shipping Dimensions: L X W X H

DC6112: 29" x 23" x 35"	DC6214: 34" x 26" x 35"	DC6326: 39" x 25" x 38"	DC6428: 46" x 27" x 41"
DC6122: 29" x 23" x 35"	DC6224: 34" x 26" x 35"	DC6228: 46" x 27" x 41"	

SuperbAIR air compressors not to be used for CAD / 3D printers

Control Panels & Stacking Rack



CP-2AV



CP-3AV



CP-3



CP-4



CV-22
33 1/2" x 22" x 26'
L x W x H



Mirage2.0

Hydraulic Patient Chair with Cast Aluminum Base



4000-2.0-CB

Foot Control

Hands-free positioning of the patient. 2-programmable positions plus an entry/exit position and last position.



Cast Aluminum Chair Base



Armrests swing-out for patient entry and exit from the chair



Reinforced seat upholstery provides extra comfort and high durability



Locking lever allows user to lock chair into desired position



Articulating dual axis headrest rotation allows for optimum positioning of the oral cavity without the need to reposition the chair



Seat rotates 30° to left and right of center to accommodate different working styles and to maximize access to intraoral x-ray system



Safety shut-off plate disables seat down/backrest recline functions to prevent injury to the operator or damage to items beneath the plate.

Specifications:

Lifting/lowering mechanism:	Hydraulic, cantilever type
Initial height/stroke:	15.5" / 15"
Backrest reclining angle:	0° (Horizontal) – 70°
Seat tilting:	10° – 20°, synchronized with backrest
Preset position:	Microcomputer control. Three preset position, auto return, last position
Chair rotation angle:	30° left / 30° right of centerline
Headrest:	Dual axis type
Safety shut-off plate:	Under the main link
Emergency stop:	Any switch can cancel chair movement
Warranty	5 Year Limited

Standard Features:

- Synchronous backrest
- Articulating dual axis headrest
- Cast aluminum base
- Swing-away armrests
- Programmable foot control
- Long pump cover

Mirage 2.0


Hydraulic Patient Chair with Cast Aluminum Base

Description	Model	Suggested Retail Price
-------------	-------	------------------------

Mirage 2.0 Hydraulic Patient Chair with Cast Aluminum Base	4000-2.0-CB	\$8,590.00
--	-------------	------------

Product includes the following:

Cast Aluminum Base
 Two programmable positions
 One preset exit position
 Electronic foot control with auto return function
 Built-in fuse / surge protection
 Toe guard protection to prevent scuffing
 All-metal construction
 Safety shut-off base plate
 Double articulating headrest, adjustable for height and wheelchair positioning
 Chair top swivel with lock, 60 degree
 Pivoting armrests for easy entry to and exit from chair
 Seamless Upholstery
 Wide chair back
 Coordinated leg rest tilt
 Input voltage 110V or 220V



Standard: PC-1000-DU
 Wide backrest



Optional: TP-2005 Touch pad
 control installed on dental chair

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 34	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
Chair Bracket 11 x 5 x 25	35 lbs.

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

Retail

TP-2005	Touch pad control – chair mounted, 3 preset, auto return, last position	\$358.00
TPC-BRKT-2.0	Mirage M4000-2.0 post mount chair bracket	\$336.00
PC-1000-DU	Replacement Dental Chair Upholstery LK Vinyl -Wide Back	\$900.00
PC-1000-Vinyl-UC	Upcharge Non-Standard Dental Chair Upholstery Vinyl -Wide Back	\$490.00
PC-1000-ULTR-UC	Upcharge Ultra leather Dental Chair Upholstery leather -Wide Back	\$760.00
PC-1000-PR-UC	Upcharge Nauga Soft Dental Chair Upholstery Vinyl -Wide Back (PR50, 49, 47, 40, 37)	\$490.00

* Specify color and product code when ordering *

Mirage2.0

Hydraulic Patient Chair with Steel Base



4000-2.0

Foot Control

Hands-free positioning of the patient. 2-programmable positions plus an entry/exit position and last position.



Steel base with high durability plastic cover



Armrests swing-out for patient entry and exit from the chair



Reinforced seat upholstery provides extra comfort and high durability



Locking lever allows user to lock chair top into desired position



Articulating dual axis headrest rotation allows for optimum positioning of the oral cavity without the need to reposition the chair



Seat rotates 30° to left and right of center to accommodate different working styles and to maximize access to intraoral x-ray system



Safety shut-off plate disables seat down/backrest recline functions to prevent injury to the operator or damage to items beneath the plate.

Specifications:

Lifting/lowering mechanism:	Hydraulic, cantilever type
Initial height/stroke:	15.5" / 15"
Backrest reclining angle:	0° (Horizontal) – 70°
Seat tilting:	10° – 20°, synchronized with backrest
Preset position:	Microcomputer control. Three preset position, auto return, last position
Chair rotation angle:	30° left / 30° right of centerline
Headrest:	Dual axis type
Safety shut-off plate:	Under the main link
Emergency stop:	Any switch can cancel chair movement
Warranty	5 Year Limited

Standard Features:

- Synchronous backrest
- Articulating dual axis headrest
- Steel plate base
- Swing-away armrests
- Programmable foot control
- Long pump cover


Mirage 2.0

Hydraulic Patient Chair with Steel Base

Description	Model	Suggested Retail Price
Mirage 2.0 Hydraulic Patient Chair with Steel Base	4000-2.0	\$7,390.00

Product includes the following:

Two programmable positions
 One preset exit position
 Electronic foot control with auto return function
 Built-in fuse / surge protection
 Toe guard protection to prevent scuffing
 All-metal construction
 Safety shut-off base plate
 Double articulating headrest, adjustable for height and wheelchair positioning
 Chair top swivel with lock, 60 degree
 Pivoting armrests for easy entry to and exit from chair
 Seamless Upholstery
 Wide chair back
 Coordinated leg rest tilt
 Input voltage 110V or 220V



Standard: PC-1000-DU
Wide backrest



Optional: TP-2005 Touch pad control installed on dental chair

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 34	40 lbs.
Cushion: 51 x 24 x 10	35 lbs.
Chair Bracket 11 x 5 x 25	35 lbs.

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

Retail

TP-2005	Touch pad control – chair mounted, 3 preset, auto return, last position	\$358.00
TPC-BRKT-2.0	Mirage M4000-2.0 post mount chair bracket	\$336.00
PC-1000-DU	Replacement Dental Chair Upholstery LK Vinyl -Wide Back	\$900.00
PC-1000-Vinyl-UC	Upcharge Non-Standard Dental Chair Upholstery Vinyl -Wide Back	\$490.00
PC-1000-ULTR-UC	Upcharge Ultra leather Dental Chair Upholstery leather -Wide Back	\$760.00
PC-1000-PR-UC	Upcharge Nauga Soft Dental Chair Upholstery Vinyl -Wide Back (PR50, 49, 47, 40, 37)	\$490.00

* Specify color and product code when ordering *

Laguna 2.0

Electromechanical Patient Chair



L2000-2.0

Foot Control

Hands-free positioning of the patient.
2-programmable positions plus an entry/exit position.



Wide backrest provides extra comfort and durability



Articulating dual axis headrest rotation allows for optimum positioning of the oral cavity without the need to reposition the chair



Armrests pivot in and out for patient entry and exit from the chair

Specifications:

Lifting/lowering mechanism:	Electromechanical, cantilever type
Initial height/stroke:	15.5" / 15"
Backrest reclining angle:	0° (Horizontal) – 70°
Seat tilting:	10° – 20°, synchronized with backrest
Preset position:	Microcomputer control. Three preset position, auto return, last position
Chair rotation angle:	Fixed
Headrest:	Dual axis type
Safety shut-off plate:	Under the main link
Emergency stop:	Any switch can cancel chair movement
Warranty	5 Year Limited



Solid steel base

Standard Features:

- Independent backrest
- Articulating dual axis headrest
- Solid steel base
- Programmable foot control

Laguna 2.0

Electromechanical Patient Chair

Description	Model	Suggested Retail Price
Laguna 2.0 Electromechanical Patient Chair	L2000-2.0	\$4,590.00

Product includes the following:

Two programmable positions
 One preset exit position
 Electronic foot control with auto return function
 Built-in fuse / surge protection
 Toe guard protection to prevent scuffing
 All-metal construction
 Double articulating headrest, adjustable for height and wheelchair positioning
 Seamless Upholstery
 Input voltage 110V or 220V



Optional: TP-2005 Touch pad control installed on dental chair

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.

Model#L2000-2.0



Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."		Retail
TP-2005	Touch pad control – Chair mounted, 3 preset, auto return, last position	\$358.00
PC-5000-DU	Upcharge Dental Chair Upholstery- Wide Back LK Vinyl- Specify Color and Product Code	\$900.00
PC-5000-Vinyl-UC	Upcharge Non-Standard Dental Chair Upholstery- Wide Back Vinyl- Specify Color and Product Code	\$490.00
PC-5000-ULTR-UC	Upcharge Ultra leather Dental Chair Upholstery- Wide Back Ultra leather- Specify Color and Product Code	\$760.00
LAG.2.0-BRKT	Laguna L2000-2.0 Post mount chair bracket	\$336.00

* Specify color and product code when ordering *

Ortho 1200

Fixed Base Orthodontic Patient Chair



OR1200



Optional: Pony Tail headrest



Standard Features:

- Independent backrest
- Articulating dual axis headrest
- Steel plate base
- Fixed base chair
- Foot control and last position memory

Ortho 1200

Fixed Base Orthodontic Patient Chair


Description	Model	Suggested Retail Price
Fixed Base Orthodontic Patient Chair	OR1200	\$2,700.00

Product includes the following:

Fixed base height
 One preset exit position
 Electronic foot control
 Independent backrest
 Toe guard protection to prevent scuffing
 Double articulating headrest, Steel plate base
 Seamless Upholstery



Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	440 lbs.
Cushion: 51 x 24 x 10	35 lbs.



Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."		Retail
PC-5000-DU	Upcharge Dental Chair Upholstery- Wide Back LK Vinyl- Specify Color and Product Code	\$900.00
PC-5000-Vinyl-UC	Upcharge Non-Standard Dental Chair Upholstery- Wide Back Vinyl- Specify Color and Product Code	\$490.00
PC-5000-ULTR-UC	Upcharge Ultra leather Dental Chair Upholstery- Wide Back Ultra leather- Specify Color and Product Code	\$760.00
L2001	OR-1200 Orthodontic chair bracket	\$336.00
PTH 101	Pony Tail Headrest	\$218.00
M42408C	Right or Left side Armrest with bracket	\$300.00

Bellagio

*Bellagio Chair Mounted Operatory System
with Cast Aluminum Base*

BP5000-L700LED-2.0-CB



L700-LED

5000-2.0

Cast Aluminum base

4000-2.0-CB

Optional: Highspeed
Handpieces

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
Unit: 39 x 29 x 23	60 lbs.
Assist PKG: 43 x 25 x 12	40 lbs.
Dental Light Arm L600: 34 x 11 x 5	30 lbs.
Dental Light Head L910: 17 x 15 x 7	
Cuspidor 26 x 18 x 10	20 lbs.



BP5015-L700LED-2.0 CB



Model DR-1102
Optional

Model AS-1101
Optional

Bellagio

Bellagio Chair Mounted Operatory System with Cast Aluminum Base

Description	Model	Suggested Retail Price	
Bellagio Operatory Package with Cuspidor	BP5000-L700LED-2.0-CB	\$15,870.00	
Bellagio Operatory Package with Assistant's Instrumentation	BP5015-L700LED-2.0-CB	\$15,170.00	
Package includes the following:		(with cuspidor)	(without cuspidor)
4000-2.0-CB	Hydraulic Cast Aluminum Base patient chair with wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 3 programmable positions with automatic return and foot control, cast aluminum base	✓	✓
5015-2.0	Asepsis 3 Handpiece automatic control with coolant controls, flex arm with air brake, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air / water filters and regulators, master shut off valves. Includes 3 way syringe and light post with cutout for internal wires and monitor bracket Asepsis telescoping 4 position assistant's package, 3-way syringe, Hve and saliva ejector, water / air quick disconnects, solids collector clean water system, Built-in Piezo Scaler, Built-in LED Curing Light, internal timer, integrated touch pad controls	✓	✓
L700-LED	Vivid LED post mounted light with adjustable intensity, 10,000-35,000 Lux composite cure mode, motion sensor, manual toggle on/off	✓	✓
2000-C-2.0	Porcelain cuspidor with timed bowl rinse and cup filler	✓	

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

Retail

P2335+TB	Deluxe autoclavable Hve Long valve with tubing	\$67.00
WIC800	Standard flat screen monitor bracket	\$128.00
WIC840	Mirage flat screen monitor bracket with 16" extension Arm	\$210.00
MO17	17" LCD multi media flat screen monitor	\$540.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System- ISO C (6 pin) – tubing, control box, power supply	\$270.00
53016-L-100	Disposable Solids Collector Trap, 2 1/8' (pack of 100)	\$45.00
2000-C-2.0	Mirage Cuspidor	\$700.00
EMB-901	Built in Electric Motor System	\$1,790.00
AS-1101	Mirage Assistants stool	\$818.00
DR-1102	Mirage Doctor's stool	\$650.00

Bellagio

Bellagio Chair Mounted Operatory System with Steel Base



Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
Unit: 39 x 29 x 23	60 lbs.
Assist PKG: 43 x 25 x 12	40 lbs.
Dental Light Arm L600: 34 x 11 x 5	30 lbs.
Dental Light Head L570: 17 x 15 x 7	



Model DR-1102
Optional

Model AS-1101
Optional

BP5000-L550LED-2.0

Bellagio

Bellagio Chair Mounted Operatory System with Steel Base

Description	Model	Suggested Retail Price
Bellagio Operatory Package with Cuspidor	BP5000-L550LED-2.0	\$14,380.00
Bellagio Operatory Package with Assistant's Instrumentation	BP5015-L550LED-2.0	\$13,680.00

Package includes the following:

4000-2.0	Hydraulic patient chair with steel base wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 3 programmable positions with automatic return and foot control	✓	✓
5015-2.0	Asepsis 3 Handpiece automatic control with coolant controls, flex arm with air brake, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves. Includes 3 way syringe and light post with cutout for internal wires and monitor bracket Asepsis telescoping 4 position assistant's package, 3-way syringe, Hve and saliva ejector clean water / air quick disconnects, solids collector with touch pad control • Built-in Scaler • Built-in Curing light	✓	✓
L550-LED	Lustrous LED post mounted light with adjustable intensity, 8,000-30,000 Lux composite cure mode, motion sensor, manual toggle on/off	✓	✓
2000-C-2.0	Porcelain cuspidor with timed bowl rinse and cup filler	✓	

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

		Retail
P2335+TB	Deluxe autoclavable Hve Long valve with tubing	\$67.00
WIC800	Standard flat screen monitor bracket	\$128.00
WIC840	Mirage flat screen monitor bracket with 16" extension Arm	\$210.00
MO17	17" LCD multi media flat screen monitor	\$540.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System- ISO C (6 pin) – tubing, control box, power supply	\$270.00
53016-L-100	Disposable Solids Collector Trap, 2 1/8' (pack of 100)	\$45.00
2000-C-2.0	Mirage Cuspidor	\$700.00
EMB-901	Built in Electric Motor System	\$1,790.00
AS-1101	Mirage Assistants stool	\$818.00
DR-1102	Mirage Doctor's stool	\$650.00

Bellagio

Swing Mount Operatory System

Model #3550-2.0



MSP5500-2.0-CB

4000-2.0-CB

Model #5500-2.0



Optional: TP-2005 touch pad control



Optional: Post LP100 & L700 LED Light



Optional: Model #L795-LED Ceiling mounted LED operatory light

4000-2.0

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
5500: 36 x 28 x 39	100 lbs.
3550: 29 x 12 x 24	25 lbs.
LP100: 45 x 26 x 44	30 lbs.

Bellagio

Swing Mount Operatory System

Description	Model	Suggested Retail Price
Bellagio Swing Mount with Cast Aluminum Base	MSP5500-2.0-CB	\$14,570.00
Bellagio Swing Mount with Steel Base	MSP5500-2.0	\$13,370.00

Package includes the following:

4000-2.0-CB	Hydraulic cast aluminum base patient chair with wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 4 programmable positions with automatic return and foot control	✓	
5500-2.0	Rotating asepsis 3 Handpiece automatic control with coolant controls, flex arm with air brake, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves. Includes 3 way syringe and clean water system, integrated Piezo Scaler, Built-in LED Curing Light internal timer, integrated touch pad controls	✓	✓
3550-2.0	Rear radius asepsis assistant's arm system, Hve and saliva, 3 Way Syringe	✓	✓
4000-2.0	Hydraulic patient chair with steel base, wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 3 programmable positions with automatic return and foot control		✓

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

Retail

5500-2.0	Bellagio swing mount unit with chair bracket (for TPC chair only)	\$4,790.00
3550-2.0	Mirage rear swing mount suction package (for TPC chair only)	\$1,190.00
LP-100	Mirage swing mount light post (ordered w/unit)	\$770.00
LP-105	Mirage swing mount monitor post with monitor mounting bracket (ordered w/unit)	\$960.00
L700-LED	Vivid LED Operatory Light	\$1,590.00
L795-LED	Vivid LED Ceiling Mounted Operatory Light	\$2,360.00
L550-LED	Lustrous LED Operatory Light	\$1,300.00
L7601-LED	Luminous LED Operatory Light	\$2,900.00
P2335+TB	Deluxe autoclavable Hve Long valve with tubing	\$67.00
WIC800	Standard flat screen monitor bracket	\$128.00
WIC840	Mirage flat screen monitor bracket with 16" extension	\$210.00
MO17	17" LCD multi media flat screen monitor	\$540.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply	\$270.00
53016-L	Disposable Solids Collector Trap 2 1/8" (pack of 100)	\$45.00
EMB-901	Built in Electric Motor System	\$1,790.00
AS-1101	Mirage Assistants stool	\$818.00
DR1102	Mirage Doctor's Stool	\$650.00

Mirage2.0

**Mirage Chair Mounted Operatory System
with Cast Aluminum Base**



MP2000-L550LED-2.0-CB

L550-LED

2000-2.0

Optional:
Highspeed
Handpieces

Cast Aluminum base

4000-2.0-CB



Optional:
Articulating Tray

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
Unit: 39 x 29 x 23	60 lbs.
Assist PKG: 43 x 25 x 12	40 lbs.
Dental Light Arm L600: 34 x 11 x 5	30 lbs.
Dental Light Head L910: 17 x 15 x 7	
Cuspidor 26 x 18 x 10	20 lbs.



MP2015-L550LED-2.0 CB



Model DR-1102
Optional

Model AS-1101
Optional

Mirage2.0

Mirage Chair Mounted Operatory System with Cast Aluminum Base

Description	Model	Suggested Retail Price
Mirage 2.0 Operatory Package with Cuspidor	MP2000-L550LED-2.0-CB	\$14,780.00
Mirage 2.0 Operatory Package with Assistant's Instrumentation	MP2015-L550LED-2.0-CB	\$14,080.00

Package includes the following:

		(with cuspidor)	(without cuspidor)
4000-2.0-CB	Hydraulic Cast Aluminum Base patient chair with wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 3 programmable positions with automatic return and foot control, cast aluminium base	✓	✓
2015-2.0	Asepsis 3 Handpiece automatic control with coolant controls, flex arm with air brake, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves. Includes 3 way syringe and light post with cutout for internal wires and monitor bracket Asepsis telescoping 4 position assistant's package, Hve and saliva ejector, water / air quick disconnects, solids collector clean water system Touch pad Control on Both Dr. & Assistant instrument	✓	✓
L550-LED	Lustrous LED post mounted light with adjustable intensity, 8,000-30,000 Lux composite cure mode, motion sensor, manual toggle on/off	✓	✓
2000-C	Porcelain cuspidor with timed bowl rinse and cup filler	✓	

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

	Retail
TPC3WAY-PRT	Additional autoclavable 3 way syringe with tubing \$78.00
P2335+TB	Deluxe autoclavable Hve Long valve with tubing \$67.00
WIC800	Standard flat screen monitor bracket \$128.00
WIC840	Mirage flat screen monitor bracket with 16" extension Arm \$210.00
MO17	17" LCD multi media flat screen monitor \$540.00
A800N-1	Advance 800 piezo scaler – built in \$570.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply \$310.00
H7660	Illumination System- ISO C (6 pin) – tubing, control box, power supply \$270.00
ALED55BN-A	LED 5W curing light – built in to dental unit \$438.00
53016-L-100	Disposable Solids Collector Trap, 2 1/8' (pack of 100) \$45.00
TPE-2011	Touch pad Control - Chair & Electric Motor System touch pad \$460.00
EMB-901	Built in Electric Motor System \$1,790.00
XV100	Super thin X-ray film viewer \$140.00
XV150-R	X-ray film viewer bracket (curved) for TPC Mirage unit only \$80.00
AS-1101	Mirage Assistant's stool \$818.00
DR-1102	Mirage Doctor's stool \$650.00

Mirage2.0

Mirage Chair Mounted Operatory System with Steel Base

MP2015-L550LED-2.0



4000-2.0

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
Unit: 39 x 29 x 23	60 lbs.
Assist PKG: 43 x 25 x 12	40 lbs.
Dental Light Arm L600: 34 x 11 x 5	30 lbs.
Dental Light Head L570: 17 x 15 x 7	



Model DR-1102
Optional



Model AS-1101
Optional



MP2000-L550LED-2.0

Mirage2.0

Mirage Chair Mounted Operatory System with Steel Base

Description	Model	Suggested Retail Price
Mirage 2.0 Operatory Package with Cuspidor	MP2000-L550LED-2.0	\$13,580.00
Mirage 2.0 Operatory Package with Assistant's Instrumentation	MP2015-L550LED-2.0	\$12,880.00

Package includes the following:

4000-2.0	Hydraulic patient chair with steel base, wide back and full lumbar support, rotating armrests left or right, double articulating head-rest with sliding height adjustment, 30 degree swivel from left and right of center, 3 programmable positions with automatic return and foot control	✓	✓
2015-2.0	Asepsis 3 Handpiece automatic control with coolant controls, flex arm with air brake, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves. Includes 3 way syringe and light post with cutout for internal wires and monitor bracket Asepsis telescoping 4 position assistant's package, Hve and saliva ejector, water / air quick disconnects, solids collector clean water system, Touch pad Control on Both Dr. & Assistant instrument	✓	✓
L550-LED	Lustrous LED post mounted light with adjustable intensity, 8,000-30,000 Lux composite cure mode, motion sensor, manual toggle on/off	✓	✓
2000-C-2.0	Porcelain cuspidor with timed bowl rinse and cup filler	✓	

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

	Retail
TPC3WAY-PRT Additional autoclavable 3 way syringe with tubing	\$78.00
P2335+TB Deluxe autoclavable Hve Long valve with tubing	\$67.00
WIC800 Standard flat screen monitor bracket	\$128.00
WIC840 Mirage flat screen monitor bracket with 16" extension Arm	\$210.00
MO17 17" LCD multi media flat screen monitor	\$540.00
A800N-1 Advance 800 piezo scaler – built in	\$570.00
H7000 Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660 Illumination System- ISO C (6 pin) – tubing, control box, power supply	\$270.00
ALED55BN-A LED 5W curing light – built in to dental unit	\$438.00
53016-L-100 Disposable Solids Collector Trap, 2 1/8' (pack of 100)	\$45.00
TPE-2011 Touch pad Control - Chair & Electric Motor System touch pad	\$460.00
EMB-901 Built in Electric Motor System	\$1,790.00
XV100 Super thin X-ray film viewer	\$140.00
XV150-R X-ray film viewer bracket (curved) for TPC Mirage unit only	\$80.00
AS-1101 Mirage Assistant's stool	\$818.00
DR-1102 Mirage Doctor's stool	\$650.00

Mirage 2.0

Swing Mount Operatory System

Model #3550-2.0

#MSP3500-2.0



Model #3500-2.0

4000-2.0

Optional: LP100 swing mount light post with optional L550-LED Operatory Light



4000-2.0 CB



Optional: Model #L570-LED Ceiling mounted LED operatory light

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
5500: 36 x 28 x 39	100 lbs.
3550: 29 x 12 x 24	25 lbs.
LP100: 45 x 26 x 44	30 lbs.

Mirage 2.0

Swing Mount Operatory System

Description	Model	Suggested Retail Price	
Mirage 2.0 Swing Mount Operatory System with Steel Base	MSP3500-2.0		\$12,570.00
Mirage 2.0 Swing Mount Operatory System with Cast Aluminum Base		MSP3550-2.0-CB	\$13,770.00
Package includes the following:			
4000-2.0	Hydraulic patient chair with steel base, wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 4 programmable positions with automatic return and foot control	✓	
3500-2.0	Rotating asepsis 3 Handpiece automatic control with coolant controls, flex arm with air brake, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves. Includes 3 way syringe and clean water system, with touch pad control	✓	✓
3550-2.0	Rear radius asepsis assistant's arm system, Hve and saliva ejector, water quick disconnect and solids collector, 3 Way Syringe	✓	✓
4000-2.0-CB	Hydraulic Cast Aluminum Base patient chair with wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 3 programmable positions with automatic return and foot control, cast aluminium base		✓

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

Retail

3500-2.0	Mirage swing mount unit with chair bracket (for TPC chair only)	\$3,990.00
3550-2.0	Mirage rear swing mount suction package (for TPC chair only)	\$1,190.00
LP-100	Mirage swing mount light post (ordered w/unit)	\$770.00
LP-105	Mirage swing mount monitor post with monitor mounting bracket (ordered w/unit)	\$960.00
L700-LED	Vivid LED Operatory Light	\$1,590.00
L795-LED	Vivid LED Ceiling Mounted Operatory Light	\$2,360.00
L550-LED	Lustrous LED Operatory Light	\$1,300.00
L7601-LED	Luminous LED Operatory Light	\$2,900.00
P2335+TB	Deluxe autoclavable Hve Long valve with tubing	\$67.00
WIC800	Standard flat screen monitor bracket	\$128.00
WIC840	Mirage flat screen monitor bracket with 16" extension	\$210.00
MO17	17" LCD multi media flat screen monitor	\$540.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply	\$270.00
53016-L	Disposable Solids Collector Trap 2 1/8" (pack of 100)	\$45.00
EMB-901	Built in Electric Motor System	\$1,790.00
AS-1101	Mirage Assistants stool	\$818.00
DR1102	Mirage Doctor's Stool	\$650.00

Mirage 2.0

Hygiene Operatory System with Steel Base



Mirage 2.0

Hygiene Operatory System with Steel Base

Description	Model	Suggested Retail Price
Mirage 2.0 Hygiene Operatory System (without light)	MHP1001	\$11,090.00

Package includes the following:

4000-2.0	Hydraulic patient chair with steel base, wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 4 programmable positions with automatic return and foot control	✓
H1001	Rear radius asepsis 3 Handpiece automatic control with coolant controls, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves, Unit Head Height Adjustment from 27" to 39" Includes 3 way syringe and clean water system, with touch pad control	✓

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."		Retail
LP-100	Mirage swing mount light post (ordered w/unit)	\$770.00
LP-105	Mirage swing mount monitor post with monitor mounting bracket (ordered w/unit)	\$960.00
L570-LED	Lustrous LED Ceiling Mounted Operatory Light	\$1,960.00
L550-LED	Lustrous LED Operatory Light	\$1,300.00
LP48	Lustrous LED Operatory Light	\$210.00
WIC800	Standard flat screen monitor bracket	\$128.00
WIC840	Mirage flat screen monitor bracket with 16" extension	\$210.00
MO17	17" LCD multi media flat screen monitor	\$540.00
A800N-1	Advance 800 piezo scaler – built in	\$570.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply	\$270.00
ALED55BN-A	LED 5W curing light – built in to dental unit	\$438.00
TPCRDSBRKT	TPC Radius Bracket	\$460.00
TPE-2011	Touch pad Control - Chair & Electric Motor System Touch pad	\$460.00
EMB-901	Built in Electric Motor System	\$1,790.00
AS-1101	Mirage Assistants stool	\$818.00
DR1102	Mirage Doctor's Stool	\$650.00
TPC-BRKT-2.0	Mirage chair Bracket	\$336.00
2000-PMU	PMU with cuspidor & cup filler	\$1,590.00

Mirage 2.0

Orthodontic Package with Steel Base

Standard Features:

- Synchronous backrest
- Articulating dual axis headrest
- Laser cut steel base
- Long pump cover
- Hydraulic pump for quiet operation
- Programmable foot control

MOP3000-L550LED-2375



Optional: Pony Tail headrest

#LP48



L550-LED



3000-2.0

Standard Features:

- Composite non cure mode
- 24V DC power Supply
- Bulb type: 4 White LEDs 2 Yellow
- Light intensity (Lux): High 30,000 Low 8,000
- Focal length: 27" (70 cm)
- Illumination area: 100mm x 180mm
- Sensor mode manual toggle

Standard Control Features:

- 3 HP Automatic selection
- 3-way Syringe
- Vacuum package with saliva ejector and hve
- Self contained water bottle system
- Variable speed disc-type foot control
- Water coolant on/off switch with flow control
- 1 Year warranty

Requires central air and vacuum for plumbing

Optional


Self contained Compressor & Vacuum

(No plumbing required)



OMC2375

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
Light Arm L600 26 x 14 x 8	30 lbs.
Light Head L570 17 x 15 x 7	
Cabinet 32 x 28 x 41	140 lbs.



Counter top dimensions



Rear cabinet door view



Pull out writing surface



Optional: OMC-2200 LCD Bracket

Mirage 2.0

Orthodontic Package with Steel Base

Description	Model	Suggested Retail Price
Mirage 2.0 Orthodontic Package	MOP3000-550LED-2375	\$11,840.00
Mirage 2.0 Orthodontic Package	MOP3000-550LED-2375CV	\$13,400.00
Mirage 2.0 Orthodontic Package	MOP3000-550LED-2375CV-LS	\$14,300.00

Package includes the following:

		2375	2375CV	2375CV-LS
3000-2.0	Hydraulic patient chair with steel base, wide back and full lumbar support, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 4 programmable positions with automatic return and foot control	✓	✓	✓
OMC2375	This orthodontic mobile delivery cabinet offers a pullout writing surface, two instrument drawers, one medium and two large storage drawers for materials and supplies that fully extend and quickly remove for easy cleaning. In addition, these cabinets provide convenient access from the front. Oversized locking dual-caster wheel, central plumbing with junction box	✓		
OMC2375CV	OMC2375+Compressor+ Vacuum (Self contained)		✓	
OMC2375CV-LS	OMC2375+Compressor+ Vacuum+LED Curing light+Scaler (Self contained)			✓
L550-LED	Lustrous LED operatory light with Poliblock mode 8,000 - 30,000 LUX , dual axis light head rotation	✓	✓	✓

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."


		Retail
WIC840	Mirage flat screen monitor bracket with 16" extension	\$210.00
MO17	17" LCD multi media flat screen monitor	\$540.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply	\$270.00
ALED55BN-A	LED 5W curing light – built in to dental unit	\$438.00
TP-2005	Touch pad control – dental chair mounted 3 preset, auto return, last position	\$358.00
OMC2200	LCD monitor bracket	\$330.00
OMC-ST81	Optional Stone Surface on OMC 2375, Color Artic	\$418.00
OMC-ST82	Optional Stone Surface on OMC 2375, Color White	\$418.00
OMC-ST85	Optional Stone Surface on OMC 2375, Color Cream	\$418.00
OMC-ST86	Optional Stone Surface on OMC 2375, Color Bamboo	\$418.00
M4000-09	Optional Magnetic Headrest	\$420.00
AS-1101	Mirage Assistants stool	\$818.00
DR-1102	Mirage Doctor's stool	\$650.00
PTH 101	Pony Tail Headrest	\$218.00




Optional: Stone Counter top
OMC-ST81



Optional: Stone Counter top
OMC-ST82



Optional: Stone Counter top
OMC-ST85



Optional: Stone Counter top
OMC-ST86

Laguna 2.0

Chair Mounted Operatory System

LP2100-550LED-2.0

Optional
MO17 LCD Monitor

L550-LED

2000-2.0

L2000-2.0

Optional



DR-5045



AS-5047



LP2115-550LED-2.0

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
Unit: 39 x 29 x 23	60 lbs.
Unit: 43 x 25 x 12	40 lbs.
Dental Light: 34 x 11 x 5 Light	30 lbs.
Head: 12 x 15 x 7	
Cuspidor 20 x 13 x 14	20 lbs.

Laguna 2.0

Chair Mounted Operatory System

Description	Model	Suggested Retail Price
Laguna 2.0 Operatory Package with Cuspidor	LP2100-550LED-2.0	\$10,780.00
Laguna 2.0 Operatory Package with Assistant's Instrumentation	LP2115-550LED-2.0	\$10,080.00
Package includes the following:	(with cuspidor) (without cuspidor)	
L2000-2.0 Electromechanical patient chair with wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 3 programmable positions with automatic return and foot control	✓	✓
2015-2.0 Asepsis 3 Handpiece automatic control with coolant controls, flex arm with air brake, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves. Includes 3 way syringe and light post with cutout for internal wires and monitor bracket Asepsis rotating 4 position assistant's package, Hve and saliva ejector, water disconnect solids collector and clean water system. With Touch pad control. on both Dr. & Assistant instrument	✓	✓
L550-LED Lustrous LED post mounted light with adjustable intensity, 8,000-30,000 Lux	✓	✓
2000-C-2.0 Porcelain cuspidor with timed bowl rinse and cup filler	✓	

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

	Retail
TPC3WAY-PRT Additional autoclavable 3 way syringe with tubing	\$78.00
P2335+TB Deluxe autoclavable Hve Long valve with tubing	\$67.00
WIC800 Standard flat screen monitor bracket	\$128.00
WIC840 Mirage flat screen monitor bracket with 16" extension Arm	\$210.00
MO17 17" LCD multi media flat screen monitor	\$540.00
A800N-1 Advance 800 piezo scaler – built in	\$570.00
H7000 Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660 Illumination System- ISO C (6 pin) – tubing, control box, power supply	\$270.00
ALED55BN-A LED 5W curing light – built in to dental unit	\$438.00
53016-L-100 Disposable Solids Collector Trap, 2 1/8' (pack of 100)	\$45.00
TPE-2011 Touch pad Control - Chair & Electric Motor System touch pad	\$460.00
EMB-901 Built in Electric Motor System	\$1,790.00
XV100 Super thin X-ray film viewer	\$140.00
XV150-R X-ray film viewer bracket (curved) for TPC Mirage unit only	\$80.00
AS-5047 Laguna Assistant's stool	\$590.00
DR-5045 Laguna Doctor's stool	\$480.00

Laguna 2.0

Swing Mount Operatory System

Model #3550-2.0

#LSP3500-2.0



Model #3500-2.0

L2000-2.0

Optional: LP100 swing mount light post with optional L550-LED Operatory Light



L2000-2.0

LP100+L550 LED



Optional: Model #L570-LED Ceiling mounted LED operatory light

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
5500: 36 x 28 x 39	100 lbs.
3550: 29 x 12 x 24	25 lbs.
LP100: 45 x 26 x 44	30 lbs.

Laguna 2.0

Swing Mount Operatory System

Description	Model	Suggested Retail Price
Laguna 2.0 Swing Mount Operatory System with Steel Base	LSP3500-2.0	\$9,770.00

Package includes the following:

L2000-2.0	Electromechanical patient chair with wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 4 programmable positions with automatic return and foot control	✓
3500-2.0	Rotating asepsis 3 Handpiece automatic control with coolant controls, flex arm with air brake, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves. Includes 3 way syringe and clean water system, with touch pad control	✓
3550-2.0	Rear radius asepsis assistant's arm system, Hve and saliva ejector, water quick disconnect and solids collector, 3 Way Syringe	✓

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."		Retail
3500-2.0	Mirage swing mount unit with chair bracket (for TPC chair only)	\$3,390.00
3550-2.0	Mirage rear swing mount suction package (for TPC chair only)	\$1,190.00
LP-100	Mirage swing mount light post (ordered w/unit)	\$770.00
LP-105	Mirage swing mount monitor post with monitor mounting bracket (ordered w/unit)	\$960.00
L795-LED	Vivid LED Ceiling Mounted Operatory Light	\$2,360.00
L570-LED	Lustrous LED Ceiling Mounted Operatory Light	\$1,960.00
L550-LED	Lustrous LED Operatory Lights	\$1,300.00
P2335+TB	Deluxe autoclavable Hve Long valve with tubing	\$67.00
WIC800	Standard flat screen monitor bracket	\$128.00
WIC840	Mirage flat screen monitor bracket with 16" extension	\$210.00
MO17	17" LCD multi media flat screen monitor	\$540.00
A800N-1	Advance 800 piezo scaler – built in	\$570.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply	\$270.00
ALED55BN-A	LED 5W curing light – built in to dental unit	\$438.00
53016-L	Disposable Solids Collector Trap 2 1/8" (pack of 100)	\$45.00
TPE-2011	Touch pad Control - Chair & Electric Motor System Touch pad	\$460.00
EMB-901	Built in Electric Motor System	\$1,790.00
AS-5047	Laguna Assistants stool	\$590.00
DR-5045	Laguna doctor's stool	\$480.00

Laguna 2.0

Hygiene operatory package

#LHP1001-2.0

H1001



L2000-2.0



W/Chair Mount light



W/Ceiling Mount light



W/Cuspidor & Post Mount light

Dimension
Chair
H1001
Cushion

Laguna 2.0

Hygiene Operatory System

Description	Model	Suggested Retail Price
Laguna 2.0 Hygiene Operatory System (without light)	LHP1001-2.0	\$8,290.00

Package includes the following:

L2000-2.0	Electromechanical patient chair with wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 4 programmable positions with automatic return and foot control	✓
H1001	Rear radius asepsis 3 Handpiece automatic control with coolant controls, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves, Unit Head Height Adjustment from 27" to 39" Includes 3 way syringe and clean water system, with touch pad control.	✓
	Hve and saliva ejector, and solids collector, piezo scaler with light	✓

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."		Retail
LP-100	Mirage swing mount light post (ordered w/unit)	\$770.00
LP-105	Mirage swing mount monitor post with monitor mounting bracket (ordered w/unit)	\$960.00
L570-LED	Lustrous LED Ceiling Mounted Operatory Light	\$1,960.00
L550-LED	Lustrous LED Operatory Light	\$1,300.00
LP48	Lustrous LED Operatory Light	\$210.00
WIC800	Standard flat screen monitor bracket	\$128.00
WIC840	Mirage flat screen monitor bracket with 16" extension	\$210.00
MO17	17" LCD multi media flat screen monitor	\$540.00
A800N-1	Advance 800 piezo scaler – built in	\$570.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply	\$270.00
ALED55BN-A	LED 5W curing light – built in to dental unit	\$438.00
TPCRDSBRKT	TPC Radius Bracket	\$460.00
TPE-2011	Touch pad Control - Chair & Electric Motor System Touch pad	\$460.00
EMB-901	Built in Electric Motor System	\$1,790.00
AS-5047	Laguna Assistants stool	\$590.00
DR5045	Laguna Doctor's Stool	\$480.00
M4001	Mirage chair Bracket	\$336.00
2000-PMU	PMU with cuspidor & cup filler	\$1,590.00

Laguna 2.0

Orthodontic Package

Standard Features:

- Independent backrest
- Articulating dual axis headrest
- Solid steel base
- Programmable foot control
- Built in fuse/surge protection
- Electric motor operated on chair base and backrest

Standard Cabinet Features:

- Five easy to remove and clean full extension drawers, and one pull-out writing surface.
- High quality hardware throughout
- Oversized caster wheels with locks
- A removable counter - top plastic tray
- Standard 12" X 12" Junction Box
- 7' length of umbilical



Standard Features:

- Composite cure mode
- 24V DC power Supply
- Bulb type: 4 White LEDs 2 Yellow
- Light intensity (Lux): High 30,000 Low 8,000
- Focal length: 27" (70 cm)
- Illumination area: 100mm x 180mm
- Sensor mode, manual toggle on/off

Standard Control Features:

- 3 HP Automatic selection
- 3-way Syringe
- Vacuum package with saliva ejector and hve
- Self contained water bottle system
- wet / dry foot control
- 1 Year warranty
- **Requires central Air / vacuum systems**
- **Includes utility center and junction box**

Optional:

Self contained compressor & Vacuum
(No plumbing required)



OMC2375

Counter top dimensions



Optional: OMC-2200 LCD Bracket

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
Light Arm L600: 26 x 14 x 10	30 lbs.
Light head L550 17 x 15 x 7	
Cabinet 32 x 28 x 41	140 lbs.

Laguna 2.0

Orthodontic Package

Description	Model	Suggested Retail Price		
Laguna 2.0 Orthodontic Package	LOP2000-550LED-2375	\$9,090.00		
Laguna 2.0 Orthodontic Package	LOP2000-550LED-2375CV	\$10,600.00		
Laguna 2.0 Orthodontic Package	LOP2000-550LED-2375CV-SL	\$11,500.00		
Package includes the following:		2375	2375CV	2375CV-SL

L2000-2.0	Electromechanical patient chair with wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 4 programmable positions with automatic return and foot control	✓	✓	✓
OMC-2375	This orthodontic mobile delivery cabinet offers a pullout writing surface, two instrument drawers, one medium and two large storage drawers for materials and supplies that fully extend and quickly remove for easy cleaning. In addition, these cabinets provide a large open space for accessories with convenient access from the front. Oversized locking dual-casters, central plumbing with junction Box	✓		
OMC2375CV	OMC2375+Compressor+Vacuum(Self contained)		✓	
OMC2375CV-SL	OMC2375+Compressor+Vacuum+LED Curing light+Scaler(Self contained)			✓
L550-LED	Lustrous LED post mounted light with dual intensity, 8,000 - 30,000 LUX, dual axis light head rotation	✓	✓	✓
L2001-LP48	The Laguna chair bracket and LP48" Light post are included at no charge (value: \$470) when purchasing the Orthodontic Package	✓	✓	✓


Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."		Retail
WIC840	Mirage flat screen monitor bracket with 16" extension	\$210.00
MO17	17" LCD multi media flat screen monitor	\$540.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply	\$270.00
OMC2200	LCD monitor bracket	\$330.00
OMC-ST81	Optional Stone Surface on OMC 2375, Color Artic	\$418.00
OMC-ST82	Optional Stone Surface on OMC 2375, Color White	\$418.00
OMC-ST85	Optional Stone Surface on OMC 2375, Color Cream	\$418.00
OMC-ST86	Optional Stone Surface on OMC 2375, Color Bamboo	\$418.00
L2001	Laguna L2000-2.0 Chair Bracket	\$336.00
AS-5047	Laguna Assistants stool	\$590.00
DR-5045	Laguna Doctor's stool	\$480.00
PTH 101	Pony Tail Headrest	\$218.00



Optional: Stone Counter top
OMC-ST81



Optional: Stone Counter top
OMC-ST82



Optional: Stone Counter top
OMC-ST85



Optional: Stone Counter top
OMC-ST86

OR-1200

Orthodontic Package

Standard Features:

- Independent backrest
- Articulating dual axis headrest
- Fixed base
- Built in fuse/surge protection
- Last position memory foot control
- Short pump cover
- Electric motor on backrest

Standard Cabinet Features:

- Five easy to remove and clean full extension drawers, and one pull-out writing surface.
- High quality hardware throughout
- Oversized caster wheels with locks
- A removable counter - top plastic tray
- Standard 12" X 12" Junction Box
- 7' length of umbilical



Standard Features:

- Composite cure mode
- 24V DC power Supply
- Bulb type: 4 White LEDs 2 Yellow
- Light intensity (Lux): High 30,000 Low 8,000
- Focal length: 27" (70 cm)
- Illumination area: 100mm x 180mm
- Sensor mode, manual toggle on/off

Standard Control Features:

- 3 HP Automatic selection
- 3-way Syringe
- Vacuum package with saliva ejector and hve
- Self contained water bottle system
- wet / dry foot control
- 1 Year warranty
- **Requires central Air / vacuum systems**
- **Includes utility center and junction box**

Optional:

Self contained compressor & Vacuum
(No plumbing required)



L550-LED


OR 1200

ORP1200-550LED



OMC2375

Optional OMC 2200
Monitor Bracket



Counter top dimensions

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
Light Arm: L600 26 x 14 x 8 Light Head: L550 17 x 15 x 7	30 lbs.
Cabinet 32 x 28 x 41	140 lbs.



Optional: Pony Tail
headrest



Pull out writing surface



Optional: OMC-
2200 LCD Bracket

OR-1200

Orthodontic Package

Description	Model	Suggested Retail Price
ORTHO 1200 Orthodontic Package	ORP1200-550LED-2375	\$7,736.00
	ORP1200-550LED-2375CV	\$9,264.00
	ORP1200-550LED-2375CV-SL	\$10,146.00

Package includes the following:

		2375	2375-CV	2375CV-SL
OR 1200	Fixed base electromechanical patient chair with wide back and full lumbar support, double articulating headrest with sliding height adjustment, One Armrest included foot control & last position Memory	✓	✓	✓
OMC-2375	This orthodontic mobile delivery cabinet offers a pullout writing surface, two instrument drawers, one medium and two large storage drawers for materials and supplies that fully extend and quickly remove for easy cleaning. In addition, these cabinets provide convenient access from the front. Oversized locking dual-casters.	✓		
OMC2375CV	OMC2375+Compressor+Vacuum(Self contained)		✓	
OMC2375CV-SL	OMC2375+Compressor+Vacuum+LED Curing light+Scaler(Self contained)			✓
L550-LED	Lustrous LED post mounted light with dual intensity, 8,000 - 30,000 LUX, dual axis light head rotation	✓	✓	✓

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."


		Retail
WIC840	Mirage flat screen monitor bracket with 16" extension	\$210.00
MO17	17" LCD multi media flat screen monitor	\$540.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply	\$270.00
OMC2200	LCD monitor bracket	\$330.00
OMC-ST81	Optional Stone Surface on OMC 2375, Color Artic	\$418.00
OMC-ST82	Optional Stone Surface on OMC 2375, Color White	\$418.00
OMC-ST85	Optional Stone Surface on OMC 2375, Color Cream	\$418.00
OMC-ST86	Optional Stone Surface on OMC 2375, Color Bamboo	\$418.00
L2001	Laguna L2000-2.0 Chair Bracket	\$336.00
AS-5047	Laguna Assistants stool	\$590.00
DR-5045	Laguna Doctor's stool	\$480.00
PTH 101	Pony Tail Headrest	\$218.00




Optional: Stone Counter top
OMC-ST81



Optional: Stone Counter top
OMC-ST82



Optional: Stone Counter top
OMC-ST85



Optional: Stone Counter top
OMC-ST86

Bellagio

Chair Mounted Delivery System



Standard Features

- Asepsis automatic control for three handpieces
- Individual handpiece pressure adjustments
- Non-retracting water coolant with on/off
- Individual water coolant flow controls
- Air coolant flow control
- Handpiece tubing flush system
- Ultra soft handpiece tubing
- 5-Hole straight gray asepsis tubing with Midwest connector.
- Disc foot control wet/dry toggle with 7' of smooth tubing
- Equipment utility center with asepsis housing
- Vacuum canister with solids collector screen
- Water quick disconnect

- Self contained water system with city/bottle selector switch
- Straight light post for unit mount dental light
- Pre-wiring for fiber optic light source
- Air-actuated counterbalanced arm lock
- 3 way syringe on unit head
- Integrated touch pad control
- Integrated piezo scaler
- Built in LED Curing Light
- Integrated timer

Model #5000-2.0



Electric Drain Cuspidor

- Contoured bowl is made of fine porcelain for easy cleaning and infection control.
- The cuspidor operates with a timed bowl rinse and cup filler
- Swing away cuspidor bowl
- Electric solenoid water flow controlled

Assistant's Instrumentation on Rotating arm

- 1 autoclavable HVE with asepsis tubing
- 1 autoclavable saliva ejector with asepsis tubing
- Assistant's side 3-way syringe
- Touch pad control

Options

- Built In Electric Handpiece
- H.P. Illumination light system
- Post-mounted utility tray
- 17" LCD Monitor



Control Head: Easy access controls for simple adjustment. Ultra soft handpiece tubing minimizing hand and wrist stress.



Cuspidor Bowl: The deep contoured bowl is made of fine porcelain for easy cleaning and infection control.



Optional: Multi media flat screen monitor for chair side intraoral imaging, patient education and DVD viewing.



Built in self contained 1500ml water bottle system including city/bottle selector switch.



Assistant's instrument



More space available for use

Bellagio

Chair Mounted Delivery System

Model #5015-2.0



Description	Model	Suggested Retail Price
Bellagio Chair Mounted Delivery System with Cuspidor	5000-2.0	\$5,690.00
Bellagio Chair Mounted Delivery System	5015-2.0	\$4,990.00

Package includes the following:

(with cuspidor)

(without cuspidor)

Asepsis 3 handpiece automatic control with coolant controls, flex arm with dual air brake, tri block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves. Includes 3 way syringe and light post with cutout for internal wires and monitor bracket, 7' umbilical with touch pad control



Porcelain cuspidor with sensor controlled bowl rinse and cup filler, rotating 4 position assistant's package, Hve and saliva ejector, water quick disconnect, solids collector and clean water system with touch pad control



Rotating 4 position assistant's package, Hve and saliva ejector,



water quick disconnect, solids collector and clean water system, light Post



Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

Retail

WIC800	Standard flat screen monitor bracket	\$128.00
WIC840	Mirage flat screen monitor bracket with 16" extension	\$210.00
#MO17	17" LCD multi media flat screen monitor	\$540.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply	\$270.00
L700-LED	Vivid LED Operatory Light	\$1,590.00
L550-LED	Lustrous LED Operatory Light "with motion sensor, with manual toggle"	\$1,300.00
BELBRKT	BLMT chair mount bracket (Bel 20/Bel 50)	\$238.00
53016 -L-100	Disposable Solids Collector Trap, 2 1/8" (pack of 100)	\$45.00
EMB-901	Built In Electric Motor system	\$1,790.00
XV100	Super Thin X-ray Viewer	\$140.00
XV150-R	Xray viewer bracket (round cut-out) – mount to Mirage unit head (Round Arm)	\$80.00

Packing Dimensions (L x W x H)	Weight
Unit: 39 x 29 x 23	60 lbs.
Unit PMU: 43 x 25 x 12	40 lbs.
Cuspidor: 26 x 18 x 10	20 lbs.

Mirage 2.0

Chair Mounted Delivery System



Standard Features

- Asepsis automatic control for three handpieces
- Individual handpiece pressure adjustments
- Non-retracting water coolant with on/off
- Individual water coolant flow controls
- Air coolant flow control
- Handpiece tubing flush system
- Ultra soft handpiece tubing
- 6-pin straight gray asepsis tubing with Midwest connector. (5-Hole tubing optional)
- Disc foot control wet/dry toggle with 7' of smooth tubing
- Equipment utility center with asepsis housing
- Vacuum canister with solids collector screen
- Water quick disconnect
- Self contained water system with city/bottle selector switch
- Straight light post for unit mount dental light
- Pre-wiring for fiber optic light source
- Air-actuated counterbalanced arm lock
- 3 way syringe on unit head
- Touch pad controls



Model #2000-2.0

Electric Drain Cuspidor

- Contoured bowl is made of fine porcelain for easy cleaning and infection control.
- The cuspidor operates with a timed bowl rinse and cup filler
- Swing away cuspidor bowl
- Electric solenoid water flow controlled
- Touch pad control

Assistant's Instrumentation on Rotating Arm

- 1 autoclavable HVE with asepsis tubing
- 1 autoclavable saliva ejector with asepsis tubing
- Optional Assistant's side 3-way syringe
- Touch pad control

Options

- Built In Electric Handpiece
- Built-in Led curing light
- Built-in scaler
- Illumination light source
- Post-mounted utility tray
- 17" LCD Monitor
- LCD monitor bracket



Control Head: Easy access controls for simple adjustments. Ultra soft handpiece tubing minimizing hand and wrist stress.



Cuspidor Bowl: The deep contoured bowl is made of fine porcelain for easy cleaning and infection control.



Optional: Multi media flat screen monitor for chair side intraoral imaging, patient education and DVD viewing.



Built in self contained 1500ml water bottle system including city/bottle selector switch.



Optional: Assistant's side 3-way syringe



Built in water quick disconnect for ultrasonic scaling.

Mirage 2.0

Chair Mounted Delivery System

Model #2015-2.0

Optional: Assistant's side
3-way syringe shown



Optional: TP-2005
Touchpad control



Optional: TPE-2011
Touchpad control

Packing Dimensions (L x W x H)	Weight
Unit: 39 x 29 x 23	60 lbs.
Unit PMU: 43 x 25 x 12	40 lbs.
Cuspidor: 26 x 18 x 10	20 lbs.

Description

Mirage 2.0 Chair Mounted Delivery System with Cuspidor
Mirage 2.0 Chair Mounted Delivery System

Model

Suggested Retail Price

2000-2.0
2015-2.0

\$4,890.00
\$4,190.00

Package includes the following:

(with cuspidor)

(without cuspidor)

Asepsis 3 handpiece automatic control with coolant controls, flex arm with dual air brake, tri block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves. Includes 3 way syringe and light post with cutout for internal wires and monitor bracket, 7' umbilical with touch pad control



Porcelain cuspidor with sensor controlled bowl rinse and cup filler, rotating 4 position assistant's package, Hve and saliva ejector, water quick disconnect, solids collector and clean water system with touch pad control



Rotating 4 position assistant's package, Hve and saliva ejector,



water quick disconnect, solids collector and clean water system, light Post



Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

Retail

TPC3WAY-PRT	Additional autoclavable 3 way syringe with tubing	\$78.00
P2335+TB	Deluxe autoclavable Hve valve with tubing	\$67.00
WIC800	Standard flat screen monitor bracket	\$128.00
WIC840	Mirage flat screen monitor bracket with 16" extension	\$210.00
#MO17	17" LCD multi media flat screen monitor	\$540.00
A800N-1	Advance 800 piezo scaler – built in	\$570.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply	\$270.00
ALED55BN-A	LED 5W curing light – built in	\$438.00
L700-LED	Vivid LED Operatory Light	\$1,590.00
L600-LED	Lucent LED Operatory Light	\$1,500.00
L550-LED	Lustrous LED Operatory Light "with motion sensor, with manual toggle"	\$1,300.00
BELBRKT	BLMT chair mount bracket BLMT chair mount bracket (20 / 50)	\$238.00
TPCBRKT	TPC chair mount bracket	\$290.00
TP-2005	Touch pad control – unit head or dental chair mounted, 3 preset, auto return, last position	\$358.00
TPE-2011	Touch pad Control - Chari & Electric Motor System touch pad	\$460.00
53016 -L-100	Disposable Solids Collector Trap, 2 1/8" (pack of 100)	\$45.00
2000-C-2.0	Mirage Cuspidor	\$700.00
EMB-901	Built In Electric Motor system	\$1,790.00
XV100	Super Thin X-ray Viewer	\$140.00
XV150-R	Xray viewer bracket (round cut-out) – mount to Mirage unit head (Round Arm)	\$80.00

Mirage 2.0 – BLMT

Swing Mount Delivery System



Model #2450-2.0



Model #2400-2.0

Compatible with BLMT 7, 10, 20, 20N, 50 patient chairs



Optional: LP100 swing mount light post with optional L700-LED Operatory Light

Packing Dimensions (L x W x H)	Weight
Chair: 44 x 29 x 33	340 lbs.
Cushion: 51 x 24 x 10	35 lbs.
3500: 36 x 28 x 39	100 lbs.
3550: 29 x 12 x 24	25 lbs.
LP100: 45 x 26 x 44	30 lbs.

Mirage 2.0 – BLMT

Swing Mount Delivery System

Description	Model	Suggested Retail Price
Mirage 2.0 BLMT Swing Mount Delivery System (Unit + Asst pkg)	MSP2400-2.0	\$6,200.00

Packages include the following:

2400-2.0	Rotating asepsis 3 Handpiece automatic control with coolant controls, flex arm with air brake, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves. Includes 3 way syringe and clean water system,	* Fits 7, 10, 20, 20N, 50
2450-2.0	Rear radius asepsis assistant's arm system, Hve and saliva ejector, water quick disconnect and solids collector	* Specify * Fits 7, 10, 20, 20N, 50 chair at time of order

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."		Retail
2400-2.0	BLMT- TPC Swing Mount Delivery System- 7, 10, 20, 20N, 50	\$4,900.00
2450-2.0-O	BLMT- TPC Rear Radius Vacuum Pkg- 20N, 50	\$1,300.00
2450-2.0-U	BBLMT-TPC Rear Radius Vacuum Pkg- 7, 10, 20	\$1,300.00
2450-BRKT-BEL-N	BLMT-TPC Rear Radius Bracket ONLY- 20N, 50	\$252.00
2450-BRKT-BEL	BLMT-TPC Rear Radius Bracket ONLY- 7,10,20	\$252.00
BELBRKT	BLMT Chair Mount Bracket	\$238.00
2400-BRKT-BEL	BLMT- TPC Radius Chair Bracket ONLY- 7, 10, 20, 20N, 50	\$550.00
BEL20N-FC	BLMT Foot Control- 20N, 50	\$390.00
BEL20-FC	BLMT Foot Control- 7,10,20	\$390.00
LP-100	Mirage swing mount light post (ordered w/unit)	\$770.00
LP-105	Swing mount monitor post with monitor mounting bracket (ordered w/unit)	\$960.00
L600-LED	Lucent LED Operatory Light	\$1,500.00
L690-LED	Lucent LED Ceiling Mounted Operatory Light	\$1,590.00
L550-LED	Lustrous LED Operatory Light	\$1,300.00
P2335+TB	Deluxe autoclavable Hve Long valve with tubing	\$67.00
WIC800	Standard flat screen monitor bracket	\$128.00
WIC840	Mirage flat screen monitor bracket with 16" extension	\$210.00
MO17	17" LCD multi media flat screen monitor	\$540.00
A800N-1	Advance 800 piezo scaler – built in	\$570.00
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply	\$310.00
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply	\$270.00
ALED55BN-A	LED 5W curing light – built in to dental unit	\$438.00
XV100	Super thin X-ray film viewer	\$140.00
XV150-R	X-ray film viewer bracket (curved) for TPC Mirage Unit only	\$80.00
AS-1101	Mirage Assistants stool	\$818.00
DR-1102	Mirage doctor's stool	\$650.00

BLMT Compatible Components



#2450-BRKT-BEL
• Fits BLMT Chairs 7, 10, 20
• Includes bolt / washer set



#2450-BRKT-BEL-N
• Fits BLMT Chairs 20N, 50
• Includes bolts / washer kit



#BELBRKT
• Fits BLMT Chairs
• Includes bolts / washer set



#2400-BRKT-BEL
• Fits BLMT Chair 7, 10, 20, 20N, 50
• Includes bearing / washer kit
• Used to mount TPC swing unit to belmont chair



Square Connector

#BEL20-FC
• Fits BLMT Chairs 7, 10, 20
• Includes wire harness



Flat Connector

#BEL20N-FC
• Fits BLMT Chairs 20N, 50
• Includes wire harness

Description

BLMT-TPC Radius Chair Bracket. Fits 7, 10, 20, 20N, Bel 50
BLMT-Chair Mount Bracket
BLMT-Chair Foot Control. Fits 7,10, 20
BLMT-Chair Foot Control. Fits 20N, 50
BLMT-TPC Rear Radius Bracket only 7, 10, 20
BLMT-TPC Rear Radius Bracket only 20N, 50
BLMT-TPC Rear Radius Vacuum Pkg- 7, 10, 20
BLMT-TPC Rear Radius Vacuum Pkg- 20N, 50

Model

Suggested Retail Price


2400-BRKT-BEL	\$550.00
BELBRKT	\$238.00
BEL20-FC	\$390.00
BEL20N-FC	\$390.00
2450-BRKT-BEL	\$252.00
2450-BRKT-BEL-N	\$252.00
2450-2.0-U	\$1,300.00
2450-2.0-O	\$1,300.00

Mirage 2.0

Mirage 2.0 Mobile Doctors Cart



Model #MC-501-2.0



Model #MC-500-2.0



Description	Model		Suggested Retail Price
Mirage 2.0 Mobile Doctors Cart	MC-500-2.0		\$2,500.00
Mirage 2.0 Mobile Doctors Cart with Vacuum Package		MC-501-2.0	\$2,800.00
Asepsis 3 handpiece automatic selection with bottle system	✓	✓	
3 way syringe with quick change autoclavable tip	✓	✓	
Master on/off toggle	✓	✓	
Air / Water pressure adjustment for each handpiece	✓	✓	
Pressure gauge	✓	✓	
Non-retracting water / air coolant adjustment and purge	✓	✓	
Asepsis handpiece tubing	✓	✓	
Painted H-frame with height adjustment, 27" to 39"	✓	✓	
Vacuum canister with removable screen 2 1/8"		✓	
1 autoclavable Hve and saliva ejector and tubing		✓	
Variable speed foot control	✓	✓	
10 feet of umbilical	✓	✓	
Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."			Retail
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply		\$310.00
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply		\$270.00
A800N-1	Advance 800 N piezo scaler – built in		\$570.00
53016-L-100	Disposable Solids Collector Trap, 2 1/8" Pk of 100		\$45.00
MC500-UMB-FT	Additional Umbilical Per Foot		\$30.00
MC501-UMB-FT	Additional Umbilical Per Foot		\$30.00

Model #MC-501-2.0

Model #MC-500-2.0

Packing Dimensions: 29" x 23" x 21" / 50 lbs


Packing Dimensions: 29" x 23" x 21" / 50 lbs

Mirage 2.0


Wall / Cabinet Mounted System




Model #MS3000-2.0



Model #MR4000-2.0



Model #MVP5000-2.0



Description

Suggested Retail Price

Mirage Cabinet Mounted Side Delivery Unit	MS-3000	\$3,300.00
Mirage Cabinet Mounted Rear Delivery Unit	MR-4000.2.0	\$2,200.00
Mirage Cabinet Mounted Vacuum Package	MVP-5000.2.0	\$950.00
Asepsis 3 handpiece automatic selection with bottle system	✓	✓
3 way syringe with quick change autoclavable tip	✓	✓
Master on/off toggle	✓	✓
Air / Water pressure adjustment for each handpiece	✓	✓
Pressure gauge	✓	✓
Non-retracting water / air coolant adjustment and purge	✓	✓
Asepsis handpiece tubing	✓	✓
Asepsis flex arm system, w/ 9" rigid arm, pneumatic brake	✓	
2" diameter wall bracket	✓	
Pivoting arm mounted vacuum canister w/ screen		✓
3 way syringe with quick change autoclavable tip		✓
Water quick disconnect		✓
Vacuum canister with removable screen		✓
4 position holder, 1 autoclavable Hve and saliva ejector & tubing		✓
Variable speed foot control	✓	✓

Options:

Retail

H7660	Illumination system – ISO C (6 pin) – tubing, control box	\$270.00
H7000	Illumination system – ISO B (5 hole) – tubing, control box, bulb	\$310.00
A800N	Advance 800 piezo scaler – built in	\$570.00
P2335+TB	Deluxe autoclavable Hve valve with tubing	\$67.00
53016-L-100	Disposable Solids Collector Trap, 2 1/8" (Pack of 100)	\$45.00
Packing Dimensions: 17" x 21" x 51" / 60 lbs	Packing Dimensions: 19" x 15" x 20" / 25 lbs	Packing Dimensions: 13" x 24" x 11" / 15 lbs
MS-3000	MR-4000	MVP-5000

Mirage


Panel & Cabinet Mount Systems



Model #MP-2600
Mirage Panel Mount
Rear Delivery System

Model #MR-4500 L/R
Mirage Rear Cabinet Mount
Delivery System with Assistants Package

Model #MP-2601
Mirage Wall / Cabinet Mount
Rear Assistants package



Standard Features

- 2 Handpiece Tubing Positions
- 3 Position HP Holder
- TPC Press ring style syringe
- Individual drive air & water coolant controls
- Foot control
- Control Dimensions 8 1/2" x 1 3/4" x 5"
- Cut-out dimensions 7 1/2" x 11 1/2"
- Syringe has adjustable air & water controls
- Air / Water master control with Regulator.

Standard Features

- Asepsis 3 handpiece automatic selection with bottle system
- 3 way syringe with quick change autoclavable tip
- Master on/off toggle
- Air / Water pressure adjustment for each handpiece
- Pressure gauge
- Non-retracting water / air coolant adjustment and purge
- Asepsis handpiece tubing
- Self contained bottle system
- Stainless Steel tray
- Vacuum canister with removable screen
- 4 position holder, 1 autoclavable Hve and saliva ejector & tubing
- Variable speed foot control

Standard Features

- 2 HVE Suciton Valves
- 1 SE Suction Valve
- 3 Position Suction Holder
- Vacuum canister with removable screen

Description	Model	Suggested Retail Price
Mirage Panel Mount Rear Delivery System (2 HP manual select, 3-Way Syringe, 3 position HP holder, Utility Valves)	MP-2600	\$990.00
Mirage Wall & Cabinet Mount Assistance Package (2 HVE, 1 SE, 3 position holder)	MP-2601	\$350.00
Mirage Rear Cabinet Mount Delivery System w/ Assistance Package (2 HVE, 1 SE, 3 position holder)	MR-4500 L/R	\$2,990.00

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

		Retail
H7660	Illumination system – ISO C (6 pin) – tubing, control box	\$270.00
H7000	Illumination system – ISO B (5 hole) – tubing, control box, bulb	\$310.00
A800N	Advance 800 piezo scaler – built in	\$570.00

Packing Dimensions: 9" x 17" x 15" / 9 lbs
MP-2600

Packing Dimensions: 8" x 17" x 15" / 5 lbs
MP-2601

Packing Dimensions: 19" x 15" x 20" / 30 lbs
MP-4500 L/R

OMC 2375

Orthodontic Mobile Delivery Cabinet



Optional OMC 2200
Monitor Bracket




Rear compartment
4" W X 12"D X 20" H



OMC2375

Writing surface with
optional solid surface




Artic

OMC-ST81



White

OMC-ST82



Cream

OMC-ST85



Bamboo

OMC-ST86

Optional hard surface stone top

Function:


This orthodontic mobile delivery cabinet offers a pullout breadboard, two instrument drawers, one medium and two large storage drawers for materials and supplies that fully extend and quickly remove for easy cleaning. In addition, these cabinets provide a large open space for accessories with convenient access from the front. The oversized dual-caster wheels, makes positioning of the cabinet simple and easy on all floor surfaces, including carpet.

Standard Cabinet Features:

- Five easy to remove and clean full extension drawers, and one pull-out writing surface.
- High quality hardware throughout
- Oversized caster wheels with LOCKS
- Four stationary handles for easy cart positioning
- A removable counter - top plastic tray
- Standard 12" X 12" Junction Box
- 7' Length of umbilical

Standard Control Features:

- 3 HP Automatic selection
- 3-way syringe
- Vacuum Package with Saliva Ejector and HVE
- Self contained water bottle system
- Variable speed disc-type foot control
- Water coolant on/off switch with flow control
- 1 Year warranty
- Includes Utility center + junction box
- Requires central air and vacuum connections



Description

Model

Suggested Retail Price

Orthodontic Mobile Delivery Cabinet	OMC2375	\$3,190.00
-------------------------------------	---------	------------

Options: "Options are not factory-installed unless specified on the purchase order. Additional charges apply."

	Retail
A800N-1 Advance 800N Piezo Scaler - Built in	\$530.00
ALED55BN-A LED 55BN 5W LED Curing Light -Built in	\$438.00
2375-UMB-FT Umbilical Per Foot	\$30.00
OMC2200 LCD Monitor Bracket W/ Post	\$418.00
OMC-ST81 Optional Stone Surface on OMC 2375, Color Artic	\$418.00
OMC-ST82 Optional Stone Surface on OMC 2375, Color White	\$418.00

Packing Dimensions: 41" x 28" x 32" / 140 lbs

OMC 2375CV-SL

Orthodontic Mobile Delivery Cabinet



Optional OMC 2200
Monitor Bracket



OMC2375CVSL



Rear compartment
4" W X 12" D X 20" H

Function:


This self-contained orthodontic mobile delivery cabinet offers a pullout breadboard, two instrument drawers, one medium and two large storage drawers for materials and supplies that fully extend and quickly remove for easy cleaning. In addition, this cabinets provides a large open space for accessories with convenient access from the front. The oversized dual-caster wheels, makes positioning of the cabinet simple and easy on all floor surfaces, including carpet.

Standard Cabinet Features:

- Five easy to remove and clean full extension drawers, and one pull-out writing surface.
- High quality hardware throughout
- Oversized caster wheels with LOCKS
- Four stationary handles for easy cart positioning
- A removable counter - top plastic tray

Standard Control Features:


- 3 HP Automatic selection
- 3-way syringe
- Vacuum Package with Saliva Ejector and HVE
- Self contained water bottle system
- Variable speed disc-type foot control
- Water coolant on/off switch with flow control
- 1 Year warranty
- **Self contained dual compressor**



OMC-ST81



OMC-ST82



OMC-ST85



OMC-ST86

Optional hard surface stone top

Description

Orthodontic Mobile Delivery Cabinet (Self contained)
Orthodontic Mobile Delivery Cabinet (Self contained)

Model

OMC2375CV-SL
OMC2375CV(No Scaler or LED curing light)

Suggested Retail Price

\$5,600.00
\$4,700.00

Options:

OMC2200	LCD Monitor Bracket W/ Post	\$330.00
OMC-ST81	Optional Stone Surface on OMC 2375, Color Artic	\$418.00
OMC-ST82	Optional Stone Surface on OMC 2375, Color White	\$418.00
OMC-ST85	Optional Stone Surface on OMC 2375, Color Cream	\$418.00
OMC-ST86	Optional Stone Surface on OMC 2375, Color Bamboo	\$418.00

Packing Dimensions: 32" x 28" x 41" / 180 lbs CV-SL model



MC3600CV(SL)

Mirage Mobile Delivery Cart

- Moves with ease
- Efficient
- Simple maintenance
- No plumbing required
- Cost effective
- Compact



Clean Water Bottle



Waste Bottle



MC3600CVSL



Optional:

#OMC-2260 Monitor Bracket
#MO17 LCD Monitor
#UC100 Utility Tray
#AIC5888 Intraoral camera

Specifications

Electrical: (2) X 550W head (115)VAC, 50/60 Hz 6.0 Amps

Sound: 68 db

Motor: Oil less air compressor

Construction: Metal Case

Vacuum: 80" of water (6" of Hg)

Vacuum: 80" of water (6" of Hg)

Max psi: 115 PSI / 40L per min

Waste Canister: 750ml

Shipping: Freight only

Physical Dimensions

31"H x 20"W x 16.5"D

Weight: 96 lbs

The MC3600CV (SL) is a self contained mobile delivery cart suitable for use by the dentist in a location when it is not feasible to connect to central plumbing

Description	Model	Suggested Retail Price
-------------	-------	------------------------

Mirage Mobile Delivery Cart	MC3600CV	\$4,270.00
-----------------------------	----------	------------

Mirage Mobile Delivery Cart (with built in scaler and led curing light)	MC3600CV-SL	\$5,170.00
---	-------------	------------

Package Includes the following

- 3 way syringe, saliva ejector, hve
- 2 handpiece positions
- Variable speed disc-type foot control
- Self contained (dual compressor)
- Clean water Bottle
- Waste water Bottle
- Advance 800 Built in Piezo Scaler
- LED 39 Built in Curing Light
- Dual 110v outlet (rear of cart)
- 1 year warranty

✓
✓
✓
✓
✓
✓
✓
✓
✓
✓

✓
✓
✓
✓
✓
✓
✓
✓
✓
✓

Shipping Dimensions:

Packing Dimensions: 39" x 27" x 23" / 160 lbs

Options:

OMC-2260	LCD Monitor Bracket W/ Post
AIC3N1N1W	17" LCD Monitor W/AIC5888 Camera
US100	Utility Tray

Suggested Retail Price

\$330.00
\$1,500.00
\$70.00

PC 2635

Portable Delivery System



PC-2635

The **PC2635** is a portable, self contained delivery system suitable for use in-office, nursing homes, hospitals, or wherever any off-site dentistry is required. The PC2635 is a lightweight, portable system that is easily transported. Air and water supply are included on the PC2635 portable dental system. Optional Tornado High and Slow speed handpieces are available.

Description	Model	Suggested Retail Price
Portable Dental System	PC2635	\$1,790.00

Package Includes the following

• 4 hole or 2 hole instrument tubings (specify)	✓
• 3 way syringe	✓
• Self contained oil less compressor (580W / 7L)	✓
• Foot control	✓
• Automatic handpiece flush System	✓
• Built in carrying handle	✓
• Self contained water bottle system	✓
• 1 year warranty	✓

Shipping Dimensions:

Packing Dimensions: 23" x 17" x 22" / 75 lbs

PC 2630

Portable Delivery System

- PC2930: Includes
- LED Curing light
 - Piezo Scaler



PC-2630



PC-2930



The **PC2630** is a fully self contained air dental system. Quick and easy to setup and operate. Ready for travel. Tornado HP available upon request. Ask for details.

Description	Model	Suggested Retail Price
Portable Dental System	PC2630	\$2,390.00
Portable Dental System W/ Scaler & LED Curing Light	PC2930	\$3,300.00

Package Includes the following

• 4 hole or 2 hole instrument tubings (specify)	✓	✓
• 3 way syringe	✓	✓
• Saliva ejector	✓	✓
• Self contained water bottle system	✓	✓
• Self contained oil less compressor (580W / 7L)	✓	✓
• Foot control	✓	✓
• Air venturi vacuum system	✓	✓
• Transportation case with pull out handle.	✓	✓
• 1 year warranty	✓	✓
• Built In Scaler		✓
• Built In Curing Light		✓

Options:

A800N-1	Advance 800 piezo scaler – built in	\$570.00
INSTALL-1	LABOR (To Install Build-In Piezo Scaler with Additional Holder on PC2630)	\$278.00
53012-INSTL	Female Quick disconnect (factory installed)	\$160.00

Shipping Dimensions: pc2630

Packing Dimensions: 16" x 20" x 28" / 60 lbs

Shipping Dimensions: pc2930

Packing Dimensions: 16" x 20" x 28" / 65 lbs

PC 2700

Portable Dental Chair Package



Description	Model	Suggested Retail Price
Portable Dental Chair Package	PC 2700	\$1,996.00
(Includes Patient Chair, Light, Stool, Tray, Spittoon & heavy duty Carry bags)		
1 Year Warranty		
Options	Model	Retail
PC 2720	Portable Patient Chair with Carrying Bag	\$1,090.00
PC 2740	Portable Dr. Stool with Carrying Bag	\$480.00
PC 2760	Portable Operator Light chair mounted	\$260.00
PC 2780	Stainless Steel tray and spit cup	\$290.00
PC 2790	Portable chair heavy duty carry bag	\$190.00
PC 2795	Stool heavy duty carry bag	\$100.00
PC2750	Portable LED Operator Light with Tripod Base with Carrying Bag	\$378.00
PC2760-PS	Replacement Power supply for portable LED operator light	\$38.00
PC2798	Carry bag for portable LED operator light	\$50.00

Shipping Dimensions:	PC 2700	PC 2740	PC 2750
Packing Dimensions:	49" x 12" x 24" / 73 lbs	49" x 12" x 24" / 17 lbs	49" x 12" x 24" / 10 lbs.

PC 2530

Mobile Vacuum System



PC2530

- Moves with Ease
- Efficient
- Simple Maintenance
- Perfect Backup
- Cost Effective
- Compact

Specifications

Electrical:	7/8 Horsepower 115 VAC, 50/60 Hz 230 VAC, 50/60 Hz 4.0 Amps
Sound:	70 db
Motor:	Variable speed control
Construction:	Metal Case
Vacuum:	80" of water (6" of Hg)
Airflow:	75 cfm
Canister:	2 Liter
Dimensions:	14" L x 12" W x 31" H
Shipping:	Freight

The PC 2530 Mobile Vac System is a portable, self contained vacuum system suitable for use by the dentist in a temporary location or when it is not feasible to install a central vacuum system.

Description	Model	Suggested Retail Price
Mobile Vacuum System	PC2530	\$2,100.00

Package Includes the following

- 2 HVE and 1 SE Valves
- Drain Line
- 2 SE suction tips
- Power cord
- Solids trap
- 1 year warranty



Shipping Dimensions:

Packing Dimensions: 17" x 25" x 37" / 95 lbs

ProClean

Extra Oral Suction System



Benefits:

- Captures and destroys viruses, germs and microorganism aerosols and droplets generated from high speed drilling and ultrasonic scaling at chairside
- Captures mercury vapors from amalgam removal
- Captures BPA particulates from composite material
- Neutralizes noxious chemical disinfectants odors and VOC gases
- Neutralizes smoke plum during dental laser and electrosurgery procedures
- Captures air abrasion and micro-etching dust

Filtration Process:

- Moisture Absorption Filter: Absorption of medium sized water droplets, moisture, aerosols and oils
- 2 Layer HEPA Filter: H13 Fine particle filter. Filtration of fine and ultra-fine dust, bacteria, viruses, fungi and prions. Efficiency: >99.97% @ >0.3 microns
- UV-C Germicidal Lamp: Disinfection and sanitation of air beyond filtration. Power: 10 watts, Wavelength: 254.7 nm
- Cold Plasma Ionizer: Eliminates particles as small as 0.01 microns including viruses, odor molecules and cigarette smoke

Technical Specifications:

Power Requirements:	100-120V AC, 50/60 Hz (International: 220-240V AC, 50/50 Hz)
Power	1200 Watts
Suction Pressure:	29 Kpa
Suction Flow:	4,000L/min (5Kpa)
Noise Level:	<65dB
New Weight:	64 lbs
Air Intake:	Segmented rotary arm with removable hood- 71.5" length x 2.5" diameter
UC-C Lamp Sterilization:	254.7 nm, 10 Watt
Plasma Disinfection:	5,000V
HEPA Grade H13:	0.3 micron, > 99.97%

Features:

- Kills 99.97% of viruses, germs and microorganisms
- Adjustable suction power
- Segmented rotary arm with replaceable wide-mouth suction hood
- Medical grade filtration system
- 2-year warranty

Description	Model	Suggested Retail Price
ProClean Extra Oral Suction System	PC2800	\$1,590.00

Package Includes the following

• Power Cord	✓
• Operation Manual	✓
• 2 year warranty	✓

Options:

	Retail
#PC2800-SYNFLTR	Replacement washable synthetic filter
#PC2800-HEPA13SML	Replacement HEPA-13 Small Filter
#PC2800-UVLAMP	Replacement UV Lamp
#PC2800-PLSMAGEN	Replacement Plasma Generator
#PC2800-HEPA13LG	Replacement HEPA-13 Large filter
#PC2800-Motor	Replacement Suction Motor
	\$9.00
	\$69.00
	\$69.00
	\$69.00
	\$95.00
	\$390.00

Shipping Dimensions:

ProClean : 29" x 19" x 38" / 85 lbs



Vivid LED


Operatory Light



Model #L700-LED



Model L795-LED
Ceiling mounted



Packing Dimensions (L x W x H)	Weight
L900-08: 15 x 3 x 3	3 lbs.
L900-09: 30 x 3 x 3	5 lbs.
L900-10: 44 x 4 x 4	7 lbs.

#LP48 Light Post must be ordered if dental light is ordered separately
Specifications:

Dimmer mode:	From 3,000-35,000 Lux	Light intensity (Lux):	3,000 - 35,000 Lux @ 700mm
Input voltage:	110V/220V, 60/50 Hz	Cooling System:	Static Heat sink (no fan)
Output voltage:	24 Volts AC / DC	Illumination area:	100mm x 180mm @ 700mm

Description	Model	Suggested Retail Price
Vivid LED Post Mount Operatory Light With Motion Sensor + Toggle	L700-LED	\$1,590.00
Vivid LED Ceiling Mounted Light With Motion Sensor + Toggle (specify 8', 9' or 10' ceiling)	L795-LED	\$2,360.00

Package includes the following:

Input voltage 24 Vac 50/60Hz	✓	✓
2 year warranty	✓	✓
Manual toggle / Motion sensor standard	✓	✓
PoliBlock Filter mode	✓	✓

Options:	Retail
RLL-PS 110 - 240 V Power supply in 24 V out	\$70.00
LP-48 48" light post	\$210.00
L690-08 Ceiling post for 8' ceiling height	\$220.00
L690-09 Ceiling post for 9' ceiling height	\$220.00
L690-010 Ceiling post for 10' ceiling height	\$220.00
L690-011 Ceiling post for 11' ceiling height (special order)	\$276.00

Packing Dimensions (L x W x H)	Weight	Packing Dimensions (L x W x H)	Weight
L700-LED: 38 x 17 x 15	10 lbs.	Light Head: 14 x 9 x 14	7 lbs.
L795-LED: 38 x 17 x 15	18 lbs.		



Lustrous LED


Operatory Light



Model #L550-LED



Model L570-LED
Ceiling mounted



Packing Dimensions (L x W x H)	Weight
L900-08: 15 x 3 x 3	3 lbs.
L900-09: 30 x 3 x 3	5 lbs.
L900-10: 44 x 4 x 4	7 lbs.

#LP48 Light Post must be ordered if dental light is ordered separately

Specifications:

Dimmer mode:	From 8,000-30,000 Lux	Light intensity (Lux):	8,000 - 30,000 Lux @ 700mm 5,000 Kelvin
Input voltage:	110V/220V, 60/50 Hz	Cooling System:	Static Heat sink (no fan)
Output voltage:	24 Volts AC / DC	Illumination area:	100mm x 180mm @ 700mm

Description

Description	Model	Suggested Retail Price
Lustrous LED Post Mount Operatory Light	L550-LED	\$1,300.00
Lustrous LED Ceiling Mounted Light (specify 8', 9' or 10' ceiling)	L570-LED	\$1,960.00

Package includes the following:

Input voltage 24 Vac 50/60Hz	✓	✓
2 year warranty	✓	✓
Manual toggle / Motion sensor standard	✓	✓
PoliBlock Filter mode	✓	✓

Options:

Options:	Retail
RLL-PS 110 - 240 V Power supply in 24 V out	\$70.00
LP-48 48" light post	\$210.00
L690-08 Ceiling post for 8' ceiling height	\$220.00
L690-09 Ceiling post for 9' ceiling height	\$220.00
L690-010 Ceiling post for 10' ceiling height	\$220.00
L690-011 Ceiling post for 11' ceiling height (special order)	\$276.00

Packing Dimensions (L x W x H)	Weight
L550: 17 x 15 x 7	7 lbs.
L570: 17 x 15 x 7	7 lbs.
L600 Arm: 34 x 11 x 5	10 lbs.
L690 Arm: 32 x 18 x 6	18 lbs.




Lucent LED

Operatory Light



Model #L600-LED



Model L690-LED
Ceiling mounted
Lucent LED
operatory light.

Packing Dimensions (L x W x H)	Weight
L900-08: 15 x 3 x 3	3 lbs.
L900-09: 30 x 3 x 3	5 lbs.
L900-10: 44 x 4 x 4	7 lbs.

#LP48 Light Post must be ordered if dental light is ordered separately
Specifications:

Absorbed Power:	Max 8vA	Light intensity (Lux):	3,000 - 35,000 Lux @ 700mm
Input voltage:	110V/220V, 60/50 Hz	Cooling System:	Static Heat sink (no fan)
Output voltage:	24 Volts AC / DC	Illumination area:	100mm x 180mm @ 700mm

Description	Model	Suggested Retail Price
Lucent LED Post Mount Operatory Light	L600-LED	\$1,500.00
Lucent LED Post Mount Operatory Light with motion sensor	L625-LED	\$1,700.00
Lucent LED Ceiling Mounted Light (specify 8', 9' or 10' ceiling)	L690-LED	\$2,160.00
Lucent LED Ceiling Mounted Light (specify 8', 9' or 10' ceiling) with motion sensor	L695-LED	\$2,360.00

Package includes the following:

Input voltage 110V	✓	✓	✓	✓
2 year warranty	✓	✓	✓	✓
Dimmer Mode	✓	✓	✓	✓

Options:

RLL-H04	Replacement LED light optical bulb group	\$410.00
RLL-H06	Replacement light handle (Left and Right)	\$127.00
RLL-H01	Replacement shield	\$130.00
RLL-H03	Replacement reflector / set of 2	\$150.00
RLL-PS	110 - 240V Power supply in 24V out	\$70.00
LP-48	48" light post	\$210.00
L690-08	Ceiling post for 8' ceiling height	\$220.00
L690-09	Ceiling post for 9' ceiling height	\$220.00
L690-010	Ceiling post for 10' ceiling height	\$220.00
L690-011	Ceiling post for 11' ceiling height (special order)	\$276.00

Packing Dimensions (L x W x H)	Weight
L600: 34 x 11 x 5	10 lbs.
L690: 32 x 18 x 6	18 lbs.
Packing Dimensions (L x W x H)	Weight
Light Head: 14 x 9 x 14	7 lbs.



Radiant LED

Operatory Light



Radiant light handles are removable and sterilizable

R6101-LED

#LP48-RI Light Post must be ordered if dental light is ordered separately

Specifications:

Absorbed Power: Max 8vA

Input voltage: 110V/220V, 60/50 Hz

Output voltage: 24 Volts AC / DC

Bulb type: LED

Light intensity (Lux): 3,000 - 35,000 Lux @ 700mm

Cooling System: Static Heat sink (no fan)

Illumination area: 100mm x 180mm @ 700mm

Description

Model

Suggested Retail Price

Radiant LED Operatory Light

R6101-LED

\$1,990.00

Radiant LED Operatory Light with motion sensor

R6105-LED

\$2,190.00

Package includes the following:

Input voltage 110V - 240 Vac

✓

✓

2 year warranty

✓

✓

Dimmer Mode

✓

✓

Options:

Retail

RLL-H04 Replacement LED light optical bulb group

\$410.00

RLL-H06 Replacement light handle (Left and Right)

\$127.00

RLL-H01 Replacement shield

\$130.00

RLL-H03 Replacement reflector / set of 2

\$150.00

RLL-PS 110 - 240 V Power supply in 24 V out

\$70.00

LP-48-RI 48" light post

\$210.00

408500 Wall Mount Bracket

\$700.00

Packing Dimensions (L x W x H)

Weight

Arm: 17 x 6 x 36

14 lbs.

Head: 14 x 14 x 9

8 lbs.

Nylon Bushing

#LP-RI-S

• OD 1 7/8"

• ID 1 3/8"

– Page 60 –





Luminous LED

LED Operatory Light

with PoliBlock Filter



Luminous light handles are removable and sterilizable

#L7601-LED

#LP48-RI Light Post must be ordered if dental light is ordered separately
Specifications:

Absorbed Power:	Max 8vA
Input voltage:	110V/220V, 60/50 Hz
Output voltage:	24 Volts AC / DC
Bulb type:	LED

Light intensity (Lux):	3,000 - 35,000 Lux @ 700mm
Cooling System:	Static Heat sink (no fan)
Illumination area:	100mm x 180mm @ 700mm

Description	Model	Suggested Retail Price
Luminous LED Operatory Light	L-7601-LED	\$2,900.00
Luminous LED Operatory Light with motion sensor	L-7605-LED	\$3,190.00
Luminous LED Operatory Light Ceiling Mount	L-7901-LED	\$4,390.00

Package includes the following:

Input voltage 110V - 240 Vac



2 year warranty




Options:

LU-PB	Antipolymerization Filter "Poli-Block"
408500	Wall Mount Bracket

Retail

\$130.00
\$700.00



*** Ceiling light extension post must be cut at the time of installation.**

Nylon Bushing

#LP-RI-S	Included with models
• OD 1 7/8"	R6101-LED, R6105-LED
• ID 1 3/8"	



Packing Dimensions (L x W x H)	Weight
Light Arm: 17" x 6" x 36"	18 lbs.
Light Head: 14 x 14 x 9	

Mobile Cabinet Systems



Utility Mobile Cabinet
#TMC-180

- 14.5" W x 16" D x 32" H
- designed for ultrasonic scalers, curing lights, electrosurge, intraoral cameras, etc.
- 1 instrument drawer
- 2 adjustable shelves
- 1 year warranty
- Shipping dimensions: 20" x 22" x 34" / 65 lbs.



Doctor's Mobile Cabinet
#TMC-140

- 21.5" W x 19" D x 32" H
- 4 instrument drawers
- 3 storage drawers
- 1 year warranty
- Shipping dimensions: 33" x 22" x 24" / 110 lbs.



Assistant's "Alabama" Mobile Cabinet
#TMC-160

- 21.5" W x 19" D x 32" H
- Sliding counter top with amalgamator well
- 1 instrument drawer
- 2 storage drawers
- Side to side sliding counter
- 1 year warranty
- Shipping dimensions: 33" x 22" x 24" / 110 lbs.



Assistant's "North Carolina" Mobile Cabinet
#TMC-120

- 21.5" W x 19" D x 32" H
- Solid top with pull handle
- Amalgamator well drop leaf door
- 3 instrument drawers
- 1 storage drawer
- 1 year warranty
- Shipping dimensions: 33" x 22" x 24" / 110 lbs.



- White will be shipped if color choice not indicated on purchase order
- Mobile cabinets ship on pallet via freight carrier

Description	Model (light grey)	Model (white)	Suggested Retail Price
Utility Mobile Cabinet	TMC-180-G	TMC-180-W	\$478.00
Doctor's Mobile Cabinet	TMC-140-G	TMC-140-W	\$990.00
Assistant's "Alabama" Mobile Cabinet	TMC-160-G	TMC-160-W	\$990.00
Assistant's "North Carolina" Mobile Cabinet	TMC-120-G	TMC-120-W	\$990.00
Assistant's "North Carolina" Mobile Cabinet RT pull out		TMC-120R-W	\$990.00

Package includes the following:

External surfaces of high pressure laminate	✓	✓
with complimentary white interior surface	✓	✓
for easy cleaning. Smooth flowing ball bearing glides. Durable casters.	✓	✓
Color choice: white or light grey.	✓	✓
1 year warranty	✓	✓

Mirage

Operatory Stools

Model SS-5105



Model DR-1102



Model DR-7000



Dimensions:

Vertical adjustment range: 0" – 3.5"
 Seat width: 16"
 Seat height ranges: 24.5" – 31.8"
 5-Leg base diameter: 24"
 Shipping dimensions: 22" x 21" x 15" / 23 lbs.

Dimensions:

Backrest vertical adjustment range: 0" – 3.25"
 Seat width: 18"
 Seat height ranges: 19.5" – 24.8"
 5-Leg base diameter: 24"
 Shipping dimensions: 23" x 23" x 17" / 33 lbs.

Dimensions:

Backrest vertical adjustment range: 0" – 3.25"
 Seat width: 18.3"
 Seat height ranges: 19.5" – 24.8"
 5-Leg base diameter: 24"
 Shipping dimensions: 22" x 23" x 17" / 33 lbs.

Description

Mirage saddle stool
 Mirage doctor's stool
 Solace doctor's stool

Model

SS-5105

DR-1102

Suggested Retail Price

\$690.00

\$650.00

DR-7000

\$818.00

Stool Package includes the following:

Five-leg, ultra stable base, 24" diameter
 Gas cylinder lift for seat height adjustment
 Adjustable for right or left hand use
 Large, firm cushion seat
 Adjustable full ring foot rest
 Dual wheel casters
 Adjustable anatomical body support bar
 Height adjustment
 16 Upholstery color choices
 Seat tilt
 Optional adjustable armrest
 1 year warranty

✓

✓

✓

✓

✓

✓

23" - 31"

✓

✓

✓

✓

✓

✓

✓

✓

✓

✓

19" - 23"

✓

✓

✓

✓

✓

✓

✓

✓

✓

✓

19" - 23"

✓

✓

✓

✓

Options *Specify Color And Product Code*

		Retail
DR1102-9	Replacement DR1102 Seat Cushion ONLY Duration Vinyl	\$178.00
DR1102-11	Replacement DR1102 Chair Back Cushion ONLY Duration Vinyl	\$140.00
DR1102-Vinyl-UC	Upcharge Non Standard Vinyl DR1102 Upholstery upcharge	\$148.00
DR-1102-ULTR-UC	Upcharge Ultra leather DR1102 Upholstery upcharge	\$280.00
DR-7000-Vinyl-UC	Upcharge Non Standard Vinyl DR7000 Upholstery upcharge	\$148.00
DR-7000-ULTR-UC	Upcharge Ultra leather DR-7000 Upholstery upcharge	\$290.00
SS-5105-Vinyl-UC	Upcharge Non Standard Vinyl SS5105 Upholstery upcharge	\$160.00
SS-5105-ULTR-UC	Upcharge Ultra leather SS5105 Upholstery upcharge	\$300.00

Mirage

Operatory Stool

Model AS-1101



Dimensions:

Backrest vertical adjustment range: 0" – 3.25"
 Seat width: 18.5"
 Seat height ranges: 23" – 30"
 5-Leg base diameter: 24"
 Packing dimensions: 24" x 23" x 21" / 42 lbs.

Model HY-5107



Dimensions:

Backrest vertical adjustment range: 0" – 3.25"
 Seat width: 17"
 Seat height ranges: 20" – 25"
 5-Leg base diameter: 24"
 Packing dimensions: 22" x 21" x 15" / 30 lbs.

Model DRSA-7100-DU100



Dimensions:

Backrest vertical adjustment range: 0" – 3.25"
 Seat width: 18" x 15.7"
 Seat height ranges: 20" – 25"
 5-Leg base diameter: 24"
 Packing dimensions: 24" x 23" x 22" / 45 lbs.
 Only Available in black DU-100

Description

Model

Suggested Retail Price

Mirage Assistant's stool	AS-1101	\$818.00
Mirage Hygiene stool	HY-5107	\$690.00
Mirage Mica Surgeon Stool	DRSA-7100-DU100	\$1,150.00

Stool Package includes the following:

Five-leg, ultra stable base, 24" diameter	✓	✓	✓
Gas cylinder lift for seat height adjustment	✓	✓	✓
Large, firm cushion seat	✓	✓	✓
Durable tile / carpet Casters	✓	✓	✓
Adjustable full ring foot rest	✓		
Adjustable for right or left handed use	✓		
Adjustable anatomical body support bar	✓		
Seat tilt		✓	✓
Height adjustment	22" - 29"	19" - 23"	19"-25"
Adjustable armrest			✓
Available in 4 colors (DU-18, DU-100, DU-116, DU-220)		✓	DU-100 Black only
1 year warranty	✓	✓	✓

Options *Specify Color And Product Code*

AS1101-16	Replacement AS1101 Seat Cushion ONLY Duration Vinyl	\$178.00
AS1101-21	Replacement AS1101 Chair Back Cushion ONLY Duration Vinyl	\$140.00
AS1101-18	Replacement AS1101 Abdominal Bar Cushion ONLY Duration Vinyl	\$156.00
AS1101-Vinyl-UC	Upcharge Non Standard Vinyl AS1101 Upholstery Upcharge	\$170.00
AS1101-ULTR-UC	Upcharge Ultra leather AS1101 Upholstery Upcharge	\$320.00

Laguna

Operatory Stool

Model DR-5045

Model AS-5047



Dimensions:

Seat width: 18.5"
 Seat height ranges: 23" – 30"
 5-Leg base diameter: 24"
 Packing dimensions: 24" x 23" x 20" / 42 lbs.



Dimensions:

Backrest vertical adjustment range: 0" – 3.25"
 Seat width: 18" x 15.7"
 Seat height ranges: 20" – 25"
 5-Leg base diameter: 24"
 Packing dimensions: 24" x 23" x 20" / 33 lbs.

Description

Model

Suggested Retail Price

Laguna Assistant's stool

AS-5047

\$590.00

Laguna Hygiene stool

DR-5045

\$480.00

Upholstery Selections:



Sand
LK-03



Stormy
LK-04



Taupe
LK-08



Greystone
LK-11



Orange
LK-27



Cobalt
LK-31



Desert Rose
LK-34



Grey Teal
LK-49



Atlantis
LK-59



Black
LK-72

Stool Package includes the following:

Five-leg, ultra stable base, 24" diameter
 Gas cylinder lift for seat height adjustment
 Large, firm cushion seat
 Durable tile / carpet Casters
 Adjustable full ring foot rest
 Adjustable for right or left handed use
 Adjustable anatomical body support bar
 Seat tilt
 Height adjustment
 Adjustable armrest
 1 year warranty



22" - 29"



19" - 23"






Model XR-6101

Mirage

X-ray Chair

Color Chart:



Dimensions:

Seat width: 21"
 Seat height ranges: 21.5" – 27.5"
 Base diameter: 24"
 Packing dimensions: 31" x 20" x 11" / 45 lbs.
 8" x 8" x 18" / 15 lbs.
 24" x 24" x 4" / 28 lbs.

Description	Model	Retail
Mirage Xray chair	XR-6101	\$1,550

X-ray stool Package includes the following:

Adjustable headrest	✓	
Gas cylinder lift for seat height adjustment	✓	
Pneumatic pump on base	✓	
Large, firm cushion seat	✓	
Large diameter base	✓	
Height adjustment	✓	21.5" - 27.5"
16 Upholstery color choices	✓	
Arm rest	✓	
Leg rest	✓	
1 year warranty	✓	

Options

Suggested Retail Price

XR-6101-8-C	Replacement XR6101 Duration Vinyl Seat Cushion ONLY	\$370.00
XR-6101-9-C	Replacement XR6101 Duration Vinyl Chair Back Cushion ONLY	\$370.00
XR-6101-10-C	Replacement XR6101 Duration Vinyl Headrest Cushion ONLY	\$120.00
XR-6101-Vinyl-UC	Upcharge Non Standard Vinyl XR6101 Upholstery Upcharge	\$340.00
XR-6101-ULTR-UC	Upcharge Ultra leather XR6101 Upholstery Upcharge	\$540.00

17"LCD Monitor

Multimedia Monitor



Model #MO17



Standard Features:

- Multimedia Monitor
- Glass panel
- 4:3 Screen
- DVD, Intraoral Camera READY!!
- Full 17" viewable TFT Active Matrix flat panel display
- Auto adjustment function
- Built-in dual speakers
- On screen display controls (OSDC)
- Recommended resolution: 1280 x 1024 (17")
- Universal 100-240 VAC, 50/60 Hz
- Gross weight: 11.6 lbs
- Includes: AC adapter, remote control, monitor handle, camera holder
- Connections: Video-in, HDMI, RGA-in, audio-in, DC-in, Flash Drive Media Player
- 1 year warranty



Intraoral Camera Port

Multimedia Input Port

Description	Model	Suggested Retail Price
17" LCD multi media monitor	MO17	\$540.00

Shipping Dimensions 17" Monitor: 18" x 6" x 21" / 12 lbs.
(2/cs.) 19" x 12" x 22" / 27 lbs.



Advance Cam Intraoral Camera

Corded Direct USB Camera System

AIC5899



Specifications

- Imager Sony 1/4" CCD
- Capture Size 640 X 480
- Light Source 6 White LEDs
- Power USB 5VDC
- Frame Rate 30 fps MAX
- PC connection USB 2.0 Port

System Requirements

- Intel Pentium 4 .15 GHZ or Higher
- 512 MB RAM or higher
- Windows 10 or higher
- USB 2.0 port



Description

Corded Direct USB Camera System

Model

AIC5899

Suggested Retail Price

\$790

Options:

- | Options: | Retail |
|---|----------|
| AIC500 Advance Cam Barriers, 500 / Box | \$196.00 |
| AIC930 Advance Cam Camera Holder | \$30.00 |
| AIC928-USB Advance Cam USB Camera Cable | \$63.00 |

Retail

Advance Cam Intraoral Camera

Corded Camera & Docking Station System



AIC5888+AIC5925A

AIC5888

System Requirements

Intel Pentium 4 .15 GHZ or Higher

512 MB RAM or higher

Windows 10 or higher

USB 2.0 port

Standard Features:

- 1080P high definition macro lens
- Built-in 6 high-brightness white LEDs imported from the U.S.
- The docking station is equipped with AV, USB, VGA and HDMI image output
- The docking station comes with a standard 16G USB flash drive for storage
- Built-in Wi-Fi transmission function enables wireless transfer of images to a computer for storage.

Specifications:

- Camera: 1080P CMOS image sensor
- Pixels: 1920×1080
- Frame rate: 30 fps MAX
- Light source: 6 white LEDs



Description

Corded Camera & Docking Station System

Model

AIC5888+AIC5925A

Suggested Retail Price

\$1,130.00

Advance Cam Intraoral Camera

Wireless Camera & Docking Station System



Function:

- 1080P ultra-high definition macro lens.
- Built-in 6 high-brightness white LEDs imported from the US.
- Uses 5G Wi-Fi transmission for stable and lag-free image quality.
- Supports Type-C fast charging and is equipped with a large capacity lithium battery.
- The docking station is equipped with AV, USB, VGA and HDMI image output interfaces.
- The docking station comes with a standard 16G USB flash drive for instant storage to prevent loss of captured images.
- Images can be wirelessly transferred to a computer for storage via Wi-Fi.

Specifications

- Camera: 1080P CMOS image sensor
- Pixels: 1920×1080
- Frame rate: 30 fps MAX
- Battery capacity: 1200 mAh (supports 120 minutes of continuous use)
- Camera charging method: Type-C
- Docking station power: DC 9V
- Light source: 6 white LEDs

System Requirements (for USB Connection)

- Windows 10 or higher
- USB 2.0 ports
- 4 GB RAM



Description

Wireless Camera & Docking Station System

Model

AIC5855+AIC5925B

Suggested Retail Price

\$1,130.00

17"LCD Monitor

Multimedia Monitor / Camera Combo (Corded)



MO17



AIC5888

WIFI image transmission



res:

- Camera features a 1/4" High Definition image sensor
- 17" Monitor features a High Definition resolution of 1280 X 1024
- Full USB port allows for easy image transfer
- Camera buttons are easily accessible to control capture and playback functions
- Camera is illuminated with six bright LEDs; giving the best contrast, definition and resolution
- Camera features an automatic gain control to automatically adjust the light intensity to the CCD
- Input supports: HDMI and AV.
- Medial playback format includes: MP3 audio, MP4 video and JPG images
- Adjustments on monitor include: Contrast, brightness and color
- Monitor comes equipped with a medical grade tempered glass front panel to minimize scratches
- Monitor Language supported: English, Chinese, Spanish, German, Portuguese and Arabic
- 16 gb memory card (included) allows for easy image transfer and patient file management
- Built-in Wi-Fi transmission function enables wireless transfer of images to a computer for storage.
- Easy integration of LCD monitor with camera on the dental delivery unit due to no additional cabling needed

Description	Model	Suggested Retail Price
17"LCDmonitorwithintegratedIOcameracord- ed Includes: MO17DS, AIC5888	AIC3IN1W	\$1,500.00
Options:		Retail
AIC500 AdvanceCamBarriers,500/Box		\$196.00

Shipping Dimensions 17" Monitor: 18" x 6" x 21" / 12 lbs.
(2/cs.) 19" x 12" x 22" / 27 lbs.

17"LCD Monitor

Multimedia Monitor / Camera Combo (Wireless)



MO17



•Medial
playback



Features:

- Camera features a 1/4" High Definition image sensor
- 17" Monitor features a High Definition resolution of 1280 X 1024
- Full USB port allows for easy image transfer
- High Definition real time video input allows for quick playback and easy patient diagnosis
- Camera buttons are easily accessible to control capture and playback functions
- Camera is illuminated with six bright LEDs; giving the best contrast, definition and resolution
- Camera features an automatic gain control to automatically adjust the light intensity to the CCD
- Input supports: HDMI and AV.
- Built in WIFI transmitter WIFI image transmissionAIC5855

format includes: MP3 audio, MP4 video and JPG images

- Adjustments on monitor include: Contrast, brightness and color
- Monitor comes equipped with a medical grade tempered glass front panel to minimize scratches
- Monitor Language supported: English, Chinese, Spanish, German, Portuguese and Arabic
- 16gb memory card (included) allows for easy image transfer and patient file management
- Easy integration of LCD monitor with camera on the dental delivery unit due to no additional cabling needed
- Wireless camera: built in 1200mAh. Lithium battery, 120min continuous use, Type C fast charging.

Description	Model	Suggested Retail Price
17"LCDmonitorwithintegratedIOcamera(wireless) Includes: MO17DS, AIC5855	AIC3IN1WL	\$1,780.00

Options:

	Retail
AIC500 AdvanceCamBarriers,500/Box	\$196.00

Shipping Dimensions 17" Monitor: 18" x 6" x 21" / 12 lbs.
(2/cs.) 19" x 12" x 22" / 27 lbs.

LCD Monitor Brackets

Monitor Brackets



WIC800

Standard WIC800 LCD Monitor Bracket

- Easily position the angle of your monitor for maximum comfort and productivity
- 4-way pivot adjustment
- Standard 75mm spacing
- Fit 45-50mm post
- 1 Year warranty

Packing Dimensions:
11" x 7" x 6" 6 lbs.



WIC840/WIC855

Mirage Monitor Bracket

- 16" extension arm
- Horizontal travel
- 4-way pivot adjustment
- Standard 75mm spacing
- WIC840 fit 45-50mm post
- WIC855 fit 55-60mm post
- 1 Year warranty

Packing Dimensions:
15" x 10" x 5" 9 lbs.



WIC875



Monitor Adapter for AdvanceCAM Intraoral Camera & Docking Station



WIC840
WIC800
Clamp Size:
1.77" - 1.96"
45mm - 50mm

WIC855
Clamp Size:
2.16" - 2.36"
55mm - 60mm



M100mm
Expansion Plate
for 100mm Spacing
• 1 year warranty

WIC840

Description	Model	Suggested Retail Price
Standard LCD Monitor Bracket fit 45-50mm post	WIC800	\$128.00
Mirage LCD Monitor Bracket with 16" Extension Arm fit 45-50mm post	WIC840	\$210.00
Mirage LCD Monitor Bracket with 16" Extension Arm fit 55-60mm post	WIC855	\$240.00
Expansion Plate for 100 mm spacing pattern	M100MM	\$24.00
AdvanceCam Intraoral Camera & Docking Station Monitor Adapter	WIC875	\$85.00

Super Thin

LED X-Ray Viwer



- Viewing Area: 8.5" H x 11" W x 0.5" Thick
- Countertop, wall, or unit head mounted
- AC power source
- Sealed LED lights
- Instant illumination
- No grey in viewing area
- Viewing area accepts panoramic, cephalometric and regular size film
- 0.5" super thin viewer
- 1 year warranty



XV150-R
*Mount to TPC
Mirage Unit (Round Arm)



XV150-F
*Mount to Flat Arm

Description	Model	Suggested Retail Price
Super Thin X-ray Viewer	XV100	\$140
Options:		Retail
XV150-R Xray viewer bracket (round cut-out) – mount to Mirage unit head (Round Arm)		\$80.00
XV150-F Xray viewer bracket (flat cut-out) – mount to unit head (Flat Arm)		\$80.00
XV130 Replacement LED light strip		\$30.00
XV100-PSW Viewer replacement power supply		\$20.00
XV100-CLIP Viewer replacement Clip (3 Pack)		\$16.00

Package Dimensions: 17" x 11" x 2" / 3 lbs.
(10/cs.) 17" x 11" x 19" / 28 lbs.

EMC - 900

eTornado Electric Motor System



eTornado Electric Motor Features:

- Fiber optic brushless electric motor
- Two programmable settings for each mode
- Push button control for easy operation
- Light weight, compact - ergonomic design
- Counter top system is easy to install and integrate
- E-type ISO standard connection
- Speed range from 100 - 200,000 rpm

Specifications

- 2 General / Basic Endo Modes
- E-type attachment / coupling
- Manual reverse
- Unit dimensions
4.5"W X 5.25" D X 2.5"H
- Motor weight: 87 grams
- Motor Length: 50mm
- Input Voltage: 100 - 240 Vac



Easy to install and integrate with any dental unit.

Description	Model	Suggested Retail Price
e Tornado Electric Motor System - Counter top	EMC-900	\$1,900.00

Package includes the following:

e Tornado Electric Motor System - Counter top	✓
Power Pack	✓
HP Tubing	✓
E600 Electric Motor up to 40,000 rpm	✓

Options	Description	Retail Price
EMA150-FG	1:5 Friction Grip Attachment	\$1,160.00
EMA110-CA	1:1 Contra Angle Attachment	\$1,060.00
EMP30	Power Supply 30Vac, 4A for control box	\$78.00
EMC-600-T	Replacement HP tubing	\$530.00
EM-600	Replacement Motor	\$1,050.00

Shipping Dimensions: 11" x 11" x 5" / 5 lbs.

EMB - 901

eTornado Electric Motor System



Electric Motor With Tubing



eTornado
Control Box



eTornado
Touchpad

EMT-915

eTornado Electric Motor Features:

- Fiber optic brushless electric motor
- Two programmable settings for each mode
- Push button control for easy operation
- Light weight, compact - ergonomical design
- Built in system pre wired for easy installation
- Universal e-type connector compatible with most electric attachments
- Speed range from 100 - 200,000 rpm

Specifications

- 2 General / Basic Endo Modes
- e-type attachment / coupling
- Manual reverse
- Unit dimensions
4.5"W X 5.25" D X 2.5"H
- Motor weight: 87 grams
- Motor Length: 50mm
- Input Voltage: 100 - 240 Vac

Description	Model	Suggested Retail Price
eTornado Built in Electric Motor System	EMB-901	\$1,790.00

Package includes the following:

e Tornado Electric Motor System - Built in	✓
Power supply	✓
HP Tubing	✓
E600 Electric Motor up to 40,000 rpm	✓

Options

EMA150-FG	1:5 Friction Grip Attachment	\$1,160.00
EMA110-CA	1:1 Contra Angle Attachment	\$1,060.00
EMB-901-PCB	Replacement Main PCB Board	\$390.00
EMB-600-T	Replacement HP tubing	\$530.00
EM-600	Replacement Motor	\$1,050.00

Shipping Dimensions:

EMA-150FG

High Speed Contra Angle Attachment



EMA-110 CA

EMA-150FG

Description	Model	Suggested Retail Price
1:5 Friction Grip Attachment	EMA-150FG	\$1,160.00
1:1 Contra Angle Attachment	EMA-110CA	\$1,060.00

Package includes the following:

Push Button Auto Chuck	✓
Quatro Water Spray	✓
Speed Range up to 200,000 rpm	✓
Attachment to electric motor w/ E-type connection	✓
Max Lux of 32,000	✓

Shipping Dimensions:

High Speed Handpieces Overview

All High Speed Handpieces Use Ceramic Bearings

Fiber Optic Handpiece with Triple Water Spray			
	LQS(M) 666K - Kavo Type Handpiece		
			
Non-Fiber Optic Handpiece with Triple Water Spray	LQS(M) 666N - NSK Type Handpiece		
		Star Fiber Optic Swivel	 or 
	LQS(M) 666S - Star Type Handpiece		
Non-Fiber Optic Handpiece With Triple Water Spray Fixed Back End			
	QS(M) 555K-Kavo Type Handpiece		
		Star Non-Fiber Optic Swivel	
Non-Fiber Optic Handpiece With Triple Water Spray Fixed Back End	QS(M) 555N-NSK Type Handpiece		
			
	QS(M) 555S-Star Type Handpiece		
Non-Fiber Optic Handpiece With Single Water Spray		4 Hole	
	S(M)506 M4	2 Hole	
			
Non-Fiber Optic Handpiece with Single Water Spray	S(M)506 B2		
			
	QS(M) 222- TPC Type Handpiece	Fixed Back Handpiece - Does not need connector	
Non-Fiber Optic Handpiece with Single Water Spray		Fixed Back Handpiece - Does not need connector	
	S(M) 333M4- TPC Type Handpiece		
			
	S(M) 333B2- TPC Type Handpiece		

Fiber Optic Handpieces

Triple Water Spray / Push Button / Auto Chuck / Ceramic Bearing

Premius Series



Fiber Optic
w/Triple
Water Spray



LQS 666 (K,N,S)



LQM 666 (K,N,S)



Standard Head



Mini Head

FIBER OPTIC QUICK CONNECTORS

(Swivel)



FSK6 Kavo Type 6 Pin



FSK5 Kavo Type 5 Hole



FSN6 NSK Type 6 Pin

Description

Model

Retail Price

Kavo Type Fiber Optic HP Standard Head w/Triple Water Spray
Fits Kavo Connector

LQS666K

\$750.00

Kavo Type Fiber Optic HP Mini Head w/Triple Water Spray
Fits Kavo Connector

LQM666K

\$750.00

6 Hole (ISO-C) Fiber Optic Connector
Fits Handpiece with KAVO "Multiflex" Quick Connector System

FSK6 -LED

\$260.00

5 Hole (Midwest) Fiber Optic Connector
Fits Handpiece with KAVO "Multiflex" Quick Connector System

FSK5

\$260.00

NSK Type Fiber Optic HP Standard Head w/Triple Water Spray
Fits NSK Connector

LQS666N

\$750.00

NSK Type Fiber Optic HP Mini Head w/Triple Water Spray
Fits NSK Connector

LQM666N

\$750.00

6 Hole (ISO-C) Fiber Optic Connector
Fits Handpiece with NSK Quick Connector System

FSN6

\$260.00

Star Type Fiber Optic HP Standard Head w/Triple Water Spray
Fits Star Connector

LQS666S

\$750.00

Star Type Fiber Optic HP Mini Head w/Triple Water Spray
Fits Star Connector

LQM666S

\$750.00

Non-Fiber Optic Handpieces

Triple Water Spray / Push Button / Auto Chuck / Ceramic Bearing

Optimus Series



Non-Fiber Optic
w/Triple
Water Spray



QS 555 (K,N,S)



Standard Head



QM 555 (K,N,S)



Mini Head

NON-FIBER OPTIC QUICK CONNECTORS

(Swivel)



SK4 Kavo Type 4 Hole



SN4 NSK Type 4 Hole

Description	Model	Suggested Retail Price
Kavo Type Non-Fiber Optic HP <u>Standard</u> Head w/Triple Water Spray Fits Kavo Connector	QS555K	\$570.00
Kavo Type Non-Fiber Optic HP <u>Mini</u> Head w/Triple Water Spray Fits Kavo Connector	QM555K	\$570.00
4 Hole Non-Fiber Optic Connector Fits Handpiece with KAVO "Multiflex" Quick Connector System	SK4	\$134.00
NSK Type Non-Fiber Optic HP <u>Standard</u> Head w/Triple Water Spray Fits NSK Connector	QS555N	\$570.00
NSK Type Non-Fiber Optic HP <u>Mini</u> Head w/Triple Water Spray Fits NSK Connector	QM555N	\$570.00
4 Hole Non-Fiber Optic Connector Fits Handpiece with NSK Quick Connector System	SN4	\$134.00
Star Type Non-Fiber Optic HP <u>Standard</u> Head w/Triple Water Spray Fits Star Connector	QS555S	\$570.00
Star Type Non-Fiber Optic HP <u>Mini</u> Head w/Triple Water Spray Fits Star Connector	QM555S	\$570.00

Non-Fiber Optic Handpieces (Standard Back End)

Triple Water Spray / Push Button / Auto Chuck / Ceramic Bearing

Triton Series



Non-Fiber Optic
w/Triple
Water Spray



S504 M4 / S506 M4



Standard Head



M504 M4 / M506 M4



Mini Head

Description	Model	Suggested Retail Price
Non-Fiber Optic HP <u>Standard</u> Head w/Triple Water Spray Standard Back End-4 Hole	S504M4 / S506M4	\$530.00
Non-Fiber Optic HP <u>Mini</u> Head w/Triple Water Spray Standard Back End-4 Hole	M504M4 / M506M4	\$530.00
Non-Fiber Optic HP <u>Standard</u> Head w/Triple Water Spray Standard Back End-2 Hole	S502B2	\$530.00
Non-Fiber Optic HP <u>Mini</u> Head w/Triple Water Spray Standard Back End-2 Hole	M502B2	\$530.00

Non-Fiber Optic Handpieces

Single Water Spray / Push Button / Auto Chuck / Ceramic Bearing

Nexius Series



Non-Fiber Optic
w/Single
Water Spray



QS 222 Standard Head



Standard Head



QM 222 Mini Head



Mini Head

NON-FIBER OPTIC QUICK CONNECTORS (Swivel)



ST4 TPC Type 4 Hole

Description	Model	Suggested Retail Price
Non-Fiber Optic HP (Standard Head) with Single Water Spray Fits TPC Connector	QS222	\$460.00
Non-Fiber Optic HP (Mini Head) with Single Water Spray Fits TPC Connector	QM222	\$460.00
4 Hole Non-Fiber Optic Connector Fits Handpiece with TPC Quick Connector System	ST4	\$100.00

Non-Fiber Optic Handpieces (Standard Back End)

Single Water Spray / Push Button / Auto Chuck / Ceramic Bearing

Nexus Series



Non-Fiber Optic
w/Single
Water Spray



S333M4 Standard Head



Standard Head



M333M4 Mini Head



Mini Head



P333M4 Pedo Head



Pedo Head

Description	Model	Suggested Retail Price
Push Button Autochuck (Standard Head) with Single Water Spray Standard Back End -4 Hole	S333M4	\$340.00
Push Button Autochuck (Mini Head) with Single Water Spray Standard Back End -4 Hole	M333M4	\$340.00
Push Button Autochuck (Standard Head) with Single Water Spray Standard Back End -2 Hole	S333B2	\$340.00
Push Button Autochuck (Mini Head) with Single Water Spray Standard Back End -2 Hole	M333B2	\$340.00
Push Button Autochuck (Pedo Head) Standard Back End -4 Hole	P333M4	\$340.00

Low Speed Handpieces

External Water Coolant

Internal Water Coolant

Exomus Series



ELM 744: Air Motor, Speed range up to 20,000 rpm, forward / reverse speed adjustment



ELS 736: Transmission 1:1, speed range up to 20,000 rpm



ELCA 760: Transmission 1:1, push button autochuck, for latch burs



ELCA 762: Transmission 1:1, latch type, for latch burs

Inximus Series



ILM 844: Air Motor, Speed range up to 20,000 rpm, forward / reverse speed adjustment, internal air/water coolant



ILS 836: Transmission 1:1, internal water spray



ILCA 860: Transmission 1:1, internal water spray, push button auto chuck for latch burs



Description

Model

Suggested Retail Price

External Water Low Speed Complete package includes (air motor, straight nose cone HP, Push button contra angle HP)	ELK700-M4	\$810.00
External Water Low Speed Complete package includes (air motor, straight nose cone HP, Latch type contra angle HP)	ELK702-M4	\$670.00
Low Speed Motor, External Water, 20,000 rpm (4 hole)	ELM744-M4	\$350.00
Low Speed Motor, External Water, 20,000 rpm (2 hole)	ELM722-B2	\$350.00
Straight nose cone HP, 1:1	ELS736	\$192.00
Push button contra angle HP, Complete 1:1	ELCA-760	\$347.00
Latch type contra angle HP, Complete 1:1	ELCA-762	\$192.00
Internal Water Low Speed Complete package includes (air motor, straight nose cone HP, push button contra angle HP)	ILK800-4	\$1,470.00
Low Speed Motor (Internal Water) 20,000 rpm-4 hole	ILM844-4	\$540.00
Straight nose cone HP, 1:1	ILS836	\$540.00
Push Button Contra Angle HP, Complete 1:1	ILCA860	\$520.00

Fiber Optic Light Source

Light Source Package



H 7070
Fiber Optic Light Bulb
Cartridge



H 7040
Fiber Optic 5-Hole Tubing



H 7020
Control Board



H 7045
Fiber Optic 6-Pin Tubing



H 7090
Fiber Optic Transform



H 7660 - 6 Pin



H 7000 - 5 Hole

Package includes: H7020+H7045+H7090

Package includes: H7020+H7040+H7070+H7090

Description	Model	Suggested Retail Price
TPC Fiber Optic Light Source System- 6 Pin Package includes: Control board, Fiber Optic HP Tubing, Transformer	H7660	\$270.00
TPC Fiber Optic Light Source System- 5 Hole Package includes: Control board, Fiber Optic HP Tubing, Transformer, Bulb Cartridge	H7000	\$310.00
TPC Fiber Optic Control Board (accepts two Fiber Optic HP Tubing)	H7020	\$142.00
TPC Fiber Optic HP Tubing-6 Pin, Grey	H7045	\$120.00
TPC Fiber Optic HP Tubing- 5 Hole, Grey	H7040	\$120.00
TPC Fiber Optic Light Bulb Cartridge LED	H7070-LED	\$58.00
TPC Fiber Optic Light Source Transformer	H7090	\$26.00

Power Optic Quick Connects

Handpiece Quick Connect Couplers



FSK6 Kavo Type 6 Pin



FSK5 Kavo Type 5 Hole



FSN6 NSK Type 6 pin



FSN5 NSK Coupler 5 Hole



TRP213
Halogen Bulb



H7070
Halogen Bulb



TRP215
Halogen Bulb



H7070
Halogen Bulb



TRP213-LED
LED Bulb for T-Series



TRP215-LED
LED Bulb for NSK

Description	Model	Suggested Retail Price
Kavo Type 5 Hole power optic connector	FSK5	\$260.00
Kavo Type 6 Pin power optic connector	FSK6	\$260.00
NSK Type 6 Pin power optic connector	FSN6	\$260.00
NSK Type 5 Hole power optic connector	FSN5	\$260.00
Halogen light bulb for FSK6 Kavo type t-Series	TRP-213	\$27.00
LED light bulb (black casing) for FSK6 Kavo type t-series	TRP-213-LED	\$27.00
Halogen light bulb for FSN6 NSK type	TRP-215	\$27.00
LED light bulb for FSN6 NSK type	TRP-215-LED	\$27.00
LED light bulb for FSK6 Kavo type	TRP-218	\$27.90
LED light bulb black cartridge type	H7070-LED	\$58.00

Replacement Turbines & Parts

Handpiece Replacement Parts



TRP 240

REPLACEMENT TURBINES

Item #	Description	Retail
TRT 110S	Replacement turbine fits Tornado HP LQS 666, QS 555, Standard Head	\$258.00
TRT 115M	Replacement turbine fits Tornado HP LQM 666, QM 555, Mini Head	\$258.00
TRT 120S	Replacement turbine fits Tornado HP S333, Standard Head	\$225.00
TRT 122S	Replacement turbine fits Tornado HP QS222, Standard Head	\$225.00
TRT 125M	Replacement turbine fits Tornado HP QM222, M333, Mini Head	\$225.00
TRT 130	Push button turbine fits Kavo Super Torque 630 B/C and 640 B/C	\$258.00
TRT 135	Push button turbine fit for Kavo Super Torque 647B and 649 B	\$258.00
TRT 140	Push button turbine fit for Kavo Bella Torque 639B/C and 642 B/C	\$258.00
TRT 145	Push button turbine fit for Kavo Mira 635 B	\$258.00
TRT 150	Push button turbine fit for Star 430 Series Handpiece	\$258.00
TRT 155	Turbine fit for Midwest Quiet Air L F.O & Non F.O. HP	\$258.00
TRT 160	Turbine fit for Midwest Tradition L F.O & Non F.O. HP	\$258.00
TRT 165	Turbine fit for Midwest Tradition F.O & Non F.O. HP	\$138.00
TRT 170	Turbine fit for Midwest XGT/Tradition PB F.O & Non F.O. HP	\$258.00

REPLACEMENT PARTS

Item #	Description	Retail
TRP 213	Replacement Light Bulb for FSK6 (Kavo Type Quick Connector)	\$27.00
TRP 215	Replacement Light Bulb for FSN6 (NSK Type Quick Connector)	\$27.00
TRP 220	Back Cap Removal Tool for LQS 666, QS 555	\$18.00
TRP 222	Back Cap Removal Tool for LQM 666, QM 555	\$18.00
TRP 224	Back Cap Removal Tool for QS 222	\$18.00
TRP 226	Back Cap Removal Tool for QM 222, M 333	\$18.00
TRP 228	Back Cap Removal Tool for S 333	\$7.00
TRP 231	Replacement Push Button Cap for LQS 666, QS 555	\$46.00
TRP 233	Replacement Push Button Cap for LQM 666, QM 555	\$46.00
TRP 235	Replacement Push Button Cap for QS 222	\$46.00
TRP 237	Replacement Push Button Cap for QM 222, M 333	\$46.00
TRP 239	Replacement Push Button Cap for S 333	\$46.00
TRP 240	Replacement O Ring for FSK6, FSK5, Sk4 (Pack of 5 (2 white, 3 black))	\$10.00

H6000

Handpiece Lubrication System

Product Features

- Economical, efficient and time-saving
- A fail-safe maintenance solution for your office
- Filtered air used to rotate handpieces for safety
- Care for most handpieces and all handpieces brands
- Disposable mist filtration system
- Solution level indicator
- Detachable door for easy cleaning of unit
- Easy replenishment of solution tank
- Air pressure gauge
- Oil flow control valve



H6000
(3 high speed positions)



• Saves Time

• Simple to Operate

• Optimizes Performance



H6005

- 2 high speed positions
- 1 low speed E-type motor



H6025

- 1 high speed positions
- 2 low speed E-type motor



H6150-Kavo
kavo maintenance coupler



H6155-NSK
NSK maintenance coupler



H6165-W&H
W&H maintenance coupler



H6145-MW
Quick Connect for Midwest Slow Speed



H6160-BienAir
BienAir maintenance coupler



H6175-E-Type
4 Hole to E-Type Adaptor

Description

Model

Suggested Retail Price

H6000 Handpiece cleaning and lubrication system	H6000				\$1,490.00
H6005 Handpiece cleaning and lubrication system		H6005			\$1,550.00
H6025 Handpiece cleaning and lubrication system			H6025		\$1,590.00
H7005 Handpiece cleaning and librication system				H7005	\$1,598.00

Package includes the following:

3 high speed handpiece positions	✓			
2 high speed / 1 low speed E-type motor		✓		
1 high speed / 2 low speed E-type motor			✓	✓
Mist filters (pack of 12)	✓	✓	✓	✓
Lubrication fluid	✓	✓	✓	✓
Air filter	✓	✓	✓	✓
Front slide out oil tray				✓

Options:

		Retail
H6125	Lubrication fluid	\$64.00
H6130	Mist filter (pack of 12)	\$23.00
H6135	Oil absorbent door strips (pack of 20)	\$18.00
H6150-Kavo	Kavo type HP adapter	\$110.00
H6155-NSK	NSK type HP adapter	\$110.00
H6160-BienAir	Bien Air type HP apapter	\$110.00
H6165-W&H	W&H type HP adapter	\$110.00
H6145-MW	Quick connect for Midwest Slow Speed	\$110.00
H6175-E-Type	4 Hole to E-Type Adaptor	\$110.00

Packing Dimensions: 17"L x 15"W x 18"H / 24 lbs.



Advance 2505

LED Bleaching System

Features

- Quiet fan operation
- Multi-arch whitening application
- Three high intensity LED's
- Digital display with audio feedback
- High power Blue and UV light emitting diodes
- Compact mobile stand for easy storage
- Adjustable bleaching time with microprocessor-controlled timer
- Smart auto-select power: 110-240 Volts AC, 50/60 Hz
- LED bleaching system includes illuminator, mobile base, adjustable arm, power cord & user manual
- 2 year warranty



Specifications:

Broad-spectrum:	395–490 nm
Light intensity:	1,500 mW/cm ²
Weight:	14 lbs.
Dimensions:	60" (H) X 20" (W)
Voltage output:	110-240 V / 50-60 Hz
Packing Dimensions:	21" x 21" x 6" / 14lbs.

Description	Model	Suggested Retail Price
Advance 2505 LED Bleaching Light	ALED 2505	\$1,450.00
Package includes the following:		
Advance 2505 LED Bleaching Light	✓	
Mobile stand	✓	
AC power supply	✓	
2 Year warranty	✓	



LED 70N

Cordless Curing Light System

Specifications and Accessories

Switching Power Adapter

Power supply: AC 100-240V / 50-60 Hz
Dimensions: 10"L x 2"W x 8"H

Hand Piece

Light source: 10W light
Emitting diode
Wavelength: 430-490 nm
Light Intensity: 2,000 mW/cm²
Battery: Lithium ion






ALED-70

Features

- Low battery signal – visual and audible alarm
- Easy to read – LCD display
- High light intensity – 5 to 10 second cure
- No down time
- Flexibility – corded and cordless light cure
- Quiet fan – extends life of LED
- Sleep mode – conserves power, saves battery life
- Fast mode – 3 output modes 2,000 – 1,500 – 1,000
- Built-in radiometer

Standard Accessories

- LED 70 Handpiece
- Charging base
- 8mm light guide
- Eye protector
- Lithium ion battery
- Power supply
- User's manual



Description

LED 70 cordless curing light

Model

ALED-70

Suggested Retail Price

\$720.00

Package includes the following:

LED 70 cordless curing light
Charging base
Protective eye-shield
Countertop stand with power indicator
8mm Light Guide
Built-in radiometer
AC power supply



Options:

P109-LG Replacement light guide (8mm)
P208 Replacement protective eye-shield
ALED 70N-BATT Replacement rechargeable battery for LED 70
P117 11mm light guide

Retail

\$86.00
\$10.50
\$90.00
\$116.00

Shipping Dimensions: 11"L x 8"W x 4"H / 3 lbs.
(10/cs.) 18"L x 12"W x 16"H / 25 lbs.



LED 60N

Cordless Curing Light System

Specifications and Accessories

Switching Power Adapter

Power supply: AC 100-240V / 50-60 Hz

Dimensions: 10"L x 2"W x 8"H

Hand Piece

Light source: 5W light
Emitting diode

Wavelength: 430-490 nm

Light Intensity: 1,400 mW/cm²

Battery: Lithium ion






ALED-60



Features

- Low battery signal – visual and audible alarm
- Easy to read – LCD display
- High light intensity – 5 to 10 second cure
- No down time
- Flexibility – corded and cordless light cure
- Quiet fan – extends life of LED
- Sleep mode – conserves power, saves battery life
- Built-in radiometer



Standard Accessories

- LED 60 Handpiece
- Charging base
- 8mm light guide
- Eye protector
- Lithium ion battery
- Power supply
- User's Manual

Description

LED 60N cordless curing light

Model

ALED-60N

Suggested Retail Price

\$550.00

Package includes the following:

LED 60N cordless curing light

Charging base

Protective eye-shield

Countertop stand with power indicator

8mm light guide

Built-in power radiometer

AC power supply



Options:

P109-LG	Replacement light guide (8mm)
P208	Replacement protective eye-shield
ALED 60N-BATT	Replacement rechargeable battery for LED 60
P117	11mm light guide

Retail

\$86.00

\$10.50

\$70.00

\$116.00

Shipping Dimensions: 11"L x 8"W x 4"H / 3 lbs.
(10/cs.) 18"L x 12"W x 17"H / 24 lbs.

LED 55N

Cordless LED Curing Light



Specifications

- Battery: Lithium Ion, DC 3.6-4.2V
- Power adapter: 90-250V AC, 50-60 Hz
- Charging base: 9V DC
- Power consumption: 5 W LED
- Wavelength range: 430-490 nm
- Light Intensity: 1,250 mW / cm²
- Dimensions: 7"L x 1"W x 1"H
- Weight: 160 g
- Warranty: 2 year

Benefits

- Ergonomic, light weight design
- High intensity output up to 1 hour before recharging
- Cordless for ultimate convenience
- Cures most materials in 5-10 seconds
- Curing time: 20 seconds per cycle with beep every 10 seconds
- 2 Batteries included.

Charging Stand Light Indicators

LED color	Light Strength
Green	Normal
Orange	Light strength insufficient

Curing Status Light Indicators

LED color	Status Indication
Green	Normal
Red	Battery Voltage Low
Red Flashing	LED curing light handpiece is over temperature

Curing Mode Light Indicators

LED color	Working Mode
Green	Fast curing: Time setting 20 seconds (every 10 second, it sounds once)
Red	Pulse curing: interrupt radiating (every 10 seconds, it sounds once)
Orange	Step curing: soft start for 10 second from 0 to full power, then full power during the later 10 second



Charging Stand






ALED55N-1-G



ALED55N-1-R



ALED55N-1-B



LM-300

Description

LED 55N Cordless Curing Light Silver
LED 55N Cordless Curing Light Blue
LED 55N Cordless Curing Light Red

Model

ALED55N-1-G Silver
ALED55N-1-B Blue
ALED55N-1-R Red

Suggested Retail Price

\$538.00
\$538.00
\$538.00

Package includes the following:

LED 55N cordless curing light
Charging stand with power indicator
Protective eye-shield
AC power supply



Options:

P109-LG Replacement light guide for LED 55 (8mm)
ALED55N-BATT Replacement battery
P208 Replacement protective eye-shield
P117 11mm light guide
LM300 Curing light meter

Retail

\$86.00
\$90.00
\$10.50
\$116.00
\$172.00

Shipping Dimensions: 10"L x 7"W x 4"H / 3 lbs.
(10/cs.) 21"L x 11"W x 14"H / 24 lbs.



Specifications and Accessories

Switching Power Adapter

Power supply: AC 100-240V / 50-60 Hz
Dimensions: 10"L x 2"W x 8"H

Hand Piece

Light source: 5W light emitting diode
Wavelength: 430-490 nm
Light Intensity: 800 - 1000 mW/cm²
Battery: Lithium ion

Standard Accessories

- LED 50N handpiece
- charging base
- 8mm light guide
- Eye protector
- Lithium ion battery
- Power supply
- User's manual



LM-300






ALED-50N

LED 50N

Cordless Curing Light System

Features

- Low battery signal – visual and audible alarm
- Easy to read – LCD display
- High light intensity – 5 to 10 second cure
- Flexibility – Cordless light cure
- Sleep mode – conserves power, saves battery life



Description

Model

Suggested Retail Price

LED 50 Cordless Curing Light
Package includes the following:

ALED-50N

\$438.00

LED 50 cordless curing light
Charging base
Protective eye-shield
Counter top stand with power indicator
8mm light guide
Built-in power transformer
AC power supply

✓
✓
✓
✓
✓
✓
✓

Options:

Retail

P109-LG Replacement light guide (8mm)
P208 Replacement protective eye-shield
LED 50-BATT Replacement rechargeable battery for LED 50
P117 11mm light guide
LM300 Curing light meter

\$86.00
\$10.50
\$70.00
\$116.00
\$172.00

Shipping Dimensions: 10"L x 8"W x 5"H / 2 lbs.
(10/cs.) 24"L x 10"W x 16"H / 20 lbs.

LED 55BN

Built-in Curing Light System

Specifications




- Wavelength range: 430-490 nm
- Light Intensity: 1,250 mW / cm²
- Warranty: 2 year
- Power consumption: 5 W LED



Model ALED55-BN

Curing Mode Light Indicators

LED color	Working Mode
Green	Fast curing: Time setting 20 seconds (every 10 second, it sounds once)
Red	Pulse curing: interrupt radiating (every 10 seconds, it sounds once)
Orange	Step curing: soft start for 10 second from 0 to full power, then full power during the later 10 second



Curing Status Light Indicators

LED color	Status Indication
Green	Normal
Red Flashing	LED curing light handpiece is over temperature

Description

Model

Suggested Retail Price

LED 55BN Built-in Curing Light (LED 5W installed in dental unit)	ALED55BN-A (110V - 240V)	\$438.00
LED 55BN Built-in Curing Light (LED 5W installed in dental unit)	ALED55BN-E (24V)	\$438.00

Package includes the following:

Turbo Light Guide	✓
Protective Eye-Shield	✓

Options:

Retail

P109-LG	Replacement light guide for LED 55B	\$86.00
P208	Replacement protective eye-shield	\$10.50
LED 55BN-PS110	Replacement power supply 110V - 240V	\$60.00
LED 55BN-PS24	Replacement power supply 24V	\$60.00
P117	11mm light guide	\$116.00
LM300	Curing light meter	\$172.00
Shipping Dimensions: 10"L x 6"W x 2"H / 2 lbs.		
(20/cs.) 21"L x 11"W x 12"H / 27 lbs.		

LED 39N

Curing Light System



Specifications

- Class: Class II (ICE601-1)
- Wavelength range: 430-490 nm
- Light Intensity: 1,000 mW / cm²
- Warranty: 2 year

Standard Accessories

- LED 39 handpiece
- Stand
- 8mm light guide
- Eye protector
- Power supply
- User Manual






Curing Mode Light Indicators

LED color	Working Mode
Green	Fast curing: Time setting 20 seconds (every 10 second, it sounds once)
Red	Pulse curing: interrupt radiating (every 10 seconds, it sounds once)
Orange	Step curing: soft start for 10 second from 0 to full power, then full power during the later 10 second

Curing Status Light Indicators

LED color	Status Indication
Green	Normal
Red Flashing	LED curing light handpiece is over temperature



Description	Model	Suggested Retail Price
LED 39N Curing Light	ALED 39N	\$378.00
LED 39BN Curing Light	ALED 39BN (24V)	\$378.00

Package includes the following:

LED 39 curing light	✓
8mm light guide	✓
Protective eye-shield	✓
Stand	✓
Built-in power transformer	✓

Options:

	Retail
P109-LG Replacement light guide for LED 39N	\$86.00
P208 Replacement protective eye-shield	\$10.50
ALED 39N-PS Replacement power supply	\$42.00

Shipping Dimensions: 13"L x 7"W x 3"H / 2 lbs.
(20/cs.) 24"L x 14"W x 14"H / 30 lbs.

LM300

Curing Light Meter

LM-300



Description	Model	Suggested Retail Price
LED Light Meter	LM-300	\$172.00
Package includes the following:		
LED Light Meter	✓	
Battery	✓	
1 Year warranty	✓	
Shipping Dimensions: 8"L x 4"W x 2"H / 1 lbs.		

Advance 850LED

LED Piezo Scaler w/Bottle System



A850-LED

Standard Scaler Features

- Compact unit housing for easy placement anywhere in the operatory
- Extremely gentle on patient, causes little or no sensitivity
- No heat build up in handpiece tip with extended use
- Autoclavable handpiece
- Extensive power range to allow widest choice of clinical applications includes General, Perio and Endo
- Fine water spray eliminates excessive overspray onto patient
- Bottle water or External water source
- 5 tips included with unit G1, G2, G4, P4, P5
- Endodontic file holder & files included
- Easy and rapid operating for root canal treatment
- Achieve enlargement of canal quickly
- 1 year warranty



5 Included tips
G1, G2, G4, P4, P5

P3 Optional

Technical Data

- Absorbed Power: Max 15W
- Scaler: Max 10W, 24 -28kHz, automatic tuning
- Water Pressure: 30 - 145 PSI
- Water Consumption: 10 - 50 cc/min
- Approvals: ETL, according to UL - 2601 - 1
- Classification according to Medical Device
- Directive 93/42 EEC: Class 11a



Description

Model

Suggested Retail Price

Advance 850-LED Piezo Ultrasonic Scaler & Endodontic Unit

A850-LED

\$938.00

Package includes the following:

Indicate EMS or Sateltec style tips at time of order

Advance 850-LED piezo ultrasonic scaler with handpiece
G1, G2, G4, P4, P5 tips included
Foot control
Tip wrench
Water Bottle
Endo file holder + files

✓
✓
✓
✓
✓
✓

Options		Retail	Replacement Parts		Retail
A760	Endo files and holder (21mm #25)	\$85.00	A746	Piezo scaler handpiece w/ LED	\$298.00
A761	Piezo scaler tip #G1	\$32.00	A776	Tip wrench	\$25.00
A762	Piezo scaler tip #G2	\$32.00	A749	Piezo scaler O-rings (5 pack)	\$9.00
A763	Piezo scaler tip #P3 Slim	\$32.00	A778	Replacement water line 7.5' w/ Q.C.	\$28.00
A767	Piezo scaler tip #G4	\$32.00	A775	Replacement scaler foot control	\$60.00
A764	Piezo scaler tip #P4	\$32.00	A850-RP-BTL	Replacement water bottle	\$25.00
A766	Piezo scaler tip #P5	\$32.00	A781	Replacement LED HP tubing	\$136.00
			A746-RP101	Replacement LED light	\$50.00

Shipping Dimensions: 18"L x 14"W x 5"H / 6 lbs.
(4/cs.) 18"L x 15"W x 21"H / 26 lbs.

Advance 750LED

LED Piezo Scaler

Standard Scaler Features

- Compact unit housing for easy placement anywhere in the operatory
- intense LED Light Source
- Extremely gentle on patient, causes little or no sensitivity
- No heat build up in handpiece tip with extended use
- Autoclavable handpiece
- Extensive power range to allow widest choice of clinical applications
- Fine water spray eliminates excessive overspray onto patient
- Thin subgingival tips available
- 5 tips included with unit G1, G2, G4, P4, P5

Optional Endodontics Unit Features

- Endodontic file holder
- Easy and rapid operating for root canal treatment
- Achieve enlargement of canal quickly

Specifications

- Input rating: 115/230V, 50/60 Hz
- Dimensions: 160mm(W) x 155mm(D) x 70mm(H)
- Weight: 2.1 kg
- Frequency: 32 KHz
- Water pressure: 1-4 kg/cm²
- Fuses: 230V, 2 per unit ; F0.5A 250V-CQ-S500
- Fuses: 110V, 1 per unit ; F1A 5A 250V-CQ-S500
- 5 tips included with unit G1, G2, G4, P4, P5
- 1 year warranty (tips excluded)



5 Included tips
G1, G2, G4, P4, P5

P3 Optional

Description

Model

Suggested Retail Price

Advance 750N Piezo Ultrasonic Scaler & Endodontic Unit	A750N	\$650.00
Advance 750N Piezo Ultrasonic Scaler & Endodontic Unit W/LED	A750LED	\$750.00

Package includes the following:

Indicate EMS or Sateltec style tips at time of order

Advance 850-LED piezo ultrasonic scaler with handpiece
G1, G2, G4, P4, P5 tips included
Foot control
Tip wrench
Water bottle
Endo file holder + files

✓
✓
✓
✓
✓
✓

Options	Retail	Replacement Parts	Retail
A760 Endo files and holder (21mm #25)	\$85.00	A746 Piezo scaler handpiece w/ LED	\$298.00
A761 Piezo scaler tip #G1	\$32.00	A776 Tip wrench	\$25.00
A762 Piezo scaler tip #G2	\$32.00	A749 Piezo scaler O-rings (5 pack)	\$9.00
A763 Piezo scaler tip #P3 Slim	\$32.00	A778 Replacement water line 7.5' w/ Q.C.	\$28.00
A767 Piezo scaler tip #G4	\$32.00	A745 Replacement Scaler HP non LED	\$216.00
A764 Piezo scaler tip #P4	\$32.00	A782 Replacement HP tubing W LED source	\$136.00
A766 Piezo scaler tip #P5	\$32.00	A780 Replacement HP tubing non LED	\$100.00
		A746-RP101 Replacement LED light	\$50.00

Shipping Dimensions: 13"L x 10"W x 5"H / 5 lbs.
(8/cs.) 21"L x 14"W x 18"H / 42 lbs.

Advance 800N

Piezo Ultrasonic Built-In Scaler

A800N



A800LED



Features

- Extremely gentle on patient, causes little or no sensitivity
- No heat build up in handpiece tip with extended use
- Autoclavable handpiece
- Extensive power range to allow widest choice of clinical application
- Fine water spray eliminates excessive overspray onto patient
- Thin subgingival tips available
- 5 tips included with unit G1, G2, G4, P4, P5
- 1 year warranty (tips excluded)

Additional Benefits

- The Advance 800N eliminates the need for external water line and power cord connections
- The Advance 800N operates from unit foot control - no additional foot control needed
- The Advance 800N allows for a clean appearing operatory eliminating all visible tubes and power lines.

Specifications

- Transformer: 24V, 1 amp
- Activated by unit foot control – no additional foot control needed
- Small electronic box, 2-3/8"L x 1.5"W x 1.25"H
- Easy tip system allows tips to be changed easily
- Absorbed power: max 15 watt
- Scaler: max 10 watt, 24-28 kHz, automatic tuning
- Water pressure: 30-145 PSI
- Water consumption: 10-50 cc/min
- No external water line / power cord
- 1 year warranty (tips excluded)
- Indicate EMS or Sateltec Tips on all orders

5 Included tips
G1, G2, G4, P4, P5



P3 Optional

Description	Model	Suggested Retail Price
Advance 800N Piezo Ultrasonic Built In Scaler	A800N	\$570.00
Advance 800N Piezo Ultrasonic Built In Scaler W/LED	A800LED	\$658.00

Package includes the following:

Advance 800N piezo ultrasonic scaler with handpiece
#G1, #G2, #G4, #P4, #P5 tip
Tip wrench
Transformer
Control box
Advance 800 LED handpiece

Indicate EMS or Sateltec style tips at time of order

✓	
✓	✓
✓	✓
✓	✓
✓	✓

Options	Retail	Replacement Parts	Retail
A760 Endo files and holder (21mm #25)	\$85.00	A746 Piezo scaler handpiece w/ LED	\$298.00
A761 Piezo scaler tip #G1	\$32.00	A776 Tip wrench	\$25.00
A762 Piezo scaler tip #G2	\$32.00	A749 Piezo scaler O-rings (5 pack)	\$9.00
A763 Piezo scaler tip #P3 Slim	\$32.00		
A767 Piezo scaler tip #G4	\$32.00	A745 Replacement Scaler HP non LED	\$216.00
A764 Piezo scaler tip #P4	\$32.00	A782 Replacement HP tubing W LED source	\$124.00
A766 Piezo scaler tip #P5	\$32.00	A780 Replacement HP tubing non LED	\$100.00
		A746-RP101 Replacement LED light	\$50.00
		A800N-PCB Main PCB	\$157.50
		A800N-WH-HARNESS	\$32.00

Shipping Dimensions: 9"L x 7"W x 4"H / 4 lbs.
(10/cs.) 14"L x 10"W x 18"H / 40 lbs.

PowerMAX25

Ultrasonic Scaling System



Features

- Available in 25Khz or 30Khz
- Precision automatic tuning circuitry
- Linear power control for exact power setting
- Accepts Cavitron®, P, SLI, TFI, and FSI scaler inserts
- Compact and lightweight
- Auto-ranging 110-240 VAC input
- Includes one ultrasonic scaler insert and foot control
- 1 year warranty

Specifications

Weight: 6 lbs

Dimensions: 8"(W) x 8"(D) x 3"(H)

Output power to tip: Operatory adjust approximately 10-40 watts

Handpiece cable: 6'

Power supply to unit: 110-240 VAC Input, 50/60 Hz

Water flow and supply pressure: 18 cc to 30 cc per minute flow from insert tip

Supply pressure to scaler is 25-35 PSI city water or bottle system



Description	Model	Suggested Retail Price
PowerMAX25 Ultrasonic Scaler with insert	PM25	\$790.00
PowerMAX30 Ultrasonic Scaler with insert	PM30	\$790.00

Package includes the following:

PowerMAX ultrasonic scaler	✓	✓
25K universal scaler insert	✓	
30K universal scaler insert		✓
Foot control	✓	✓
1 year warranty	✓	✓

Options:

		Retail
S1025	25K universal insert-external water flow	\$78.00
S1030	30K universal insert- external water flow	\$78.00
S2025	25K universal slim insert- external water flow	\$78.00
S2030	30K universal slim insert- external water flow	\$78.00
S3025	25K universal insert- plastic handle- internal water flow	\$78.00
S3030	30K universal insert- plastic handle- internal water flow	\$78.00
S4025	25K universal slim insert- plastic handle- internal water flow	\$78.00
S4030	30K universal slim insert- plastic handle- internal water flow	\$78.00
S8008	30K Replacement handpiece w/ Tubing	\$120.00
S8001	25K Replacement handpiece w/ Tubing	\$120.00
PM775	Power Max replacement foot control	\$87.00
A778	Replacement water line 6' long with male quick disconnect (with shut off)	\$28.00

Shipping Dimensions: 12"L x 12"W x 4"H / 8 lbs.
(5/cs.) 18"L x 13"W x 14"H / 38 lbs.

Ultrasonic Scaler Inserts



#S1025
25K Universal Scaler Insert
(metal handle)

#S1030
30K Universal Scaler Insert
(metal handle)

Features

- External water
- Interchangeable for use with other brands of scalers
- 90 day warranty against defects arising from faulty materials or workmanship



#S2025
25K Universal Slim Scaler Insert
(metal handle)

#S2030
30K Universal Slim Scaler Insert
(metal handle)

Features

- Slim tip design
- External water
- Interchangeable for use with other brands of scalers
- 90 day warranty against defects arising from faulty materials or workmanship



#S3025
25K Universal Scaler Insert
(black plastic handle)

#S3030
30K Universal Scaler Insert
(black plastic handle)

Features

- Internal water delivery system delivers water directly to the insert tip and focuses it where it is needed most
- Large, round, non-slip plastic grip for greater tactile sensitivity, rotational control and reduced hand fatigue
- Better water management and less water usage means greater patient comfort
- Interchangeable for use with other brands of scalers
- 90 day warranty against defects arising from faulty materials or workmanship
- Ergonomic grip design



#S4025
25K Universal Slim Scaler Insert
(blue plastic handle)

#S4030
30K Universal Slim Scaler Insert
(blue plastic handle)

Features

- Internal water delivery system delivers water directly to the insert tip and focuses it where it is needed most
- Large, round, non-slip plastic grip for greater tactile sensitivity, rotational control and reduced hand fatigue
- Better water management and less water usage means greater patient comfort
- Interchangeable for use with other brands of scalers
- 90 day warranty against defects arising from faulty materials or workmanship
- Ergonomic grip design

Model	Description	Retail
S1025	25K universal insert—external water flow	\$78.00
S1030	30K universal insert- external water flow	\$78.00
S2025	25K universal slim insert- external water flow	\$78.00
S2030	30K universal slim insert- external water flow	\$78.00
S3025	25K universal insert- black plastic handle- internal water flow	\$78.00
S3030	30K universal insert- black plastic handle- internal water flow	\$78.00
S4025	25K universal slim insert- blue plastic handle- internal water flow	\$78.00
S4030	30K universal slim insert- blue plastic handle- internal water flow	\$78.00
S8003	Ultrasonic scaler insert O-rings (5 pack)	\$10.00

Shipping Dimensions: 9"L x 2"W x 1"H / 1 lbs.

AS 5000

Sonic Air Scaler

Features

- Includes 3 tips AS1, AS2, AS3.
- Tip wrench included.
- Light weight
- Portable
- Saves time
- Autoclavable
- Available in 2 or 4 hole.

Specifications

- Air pressure 0.25 - 0.35 Mpa
- Water flow 50ml/min
- Noise ≤ 70bd



AS5000 M4



AS5000 B2



AS5000-K



Wrench



Sk4 Adaptor



AS1

AS2

AS3

Description

Model

Suggested Retail Price

Air Sonic Scaler -4 Hole Kit	AS5000-M4	\$310.00
Air Sonic Scaler -2 Hole Kit	AS5000-B2	\$310.00
Air Sonic Scaler -4 Hole Kit (Quick Connect(KAVO Type) SK4 adaptor not included)	AS5000-K	\$310.00

Package includes the following:

AS5000 Air sonic scaler handpiece	✓	✓	✓
AS1 Universal Tip	✓	✓	✓
AS2 Sickle Tip	✓	✓	✓
AS3 Perio Tip	✓	✓	✓
A776 Tip Wrench	✓	✓	✓
6 Month warranty	✓	✓	✓
Quick Connect(KAVO Type) SK4 adaptor not included			✓

Options:

Retail

AS1 Universal Tip	\$27.50
AS2 Sickle Tip	\$27.50
AS3 Perio Tip	\$27.50
A776 Scaler Tip Wrench Only	\$25.00
SK4 Kavo Type Quick Connector	\$67.00

Shipping Dimensions: 9"L x 4"W x 2"H / 1 lbs.

Air Polishing System

Polishing System

EZ Nozzle

- Carefully designed length, angle and curve of nozzle for applications to both anterior and posterior teeth. Allows for accurate direction of the power jet
- Includes two prophyl nozzles
- Handpiece nozzle is fully autoclavable up to 134 °C / 273 °F

Ergo Comfort Design

- Fits comfortably in large and small hands. Notched grip allows firm and accurate tip control

Lightweight & Flexible

No Installation Cost

- Simply connects to handpiece tubing
- No extra foot control needed
- No external air or water connections needed

Cleaning Power

- Jet Polishing powder (130g) has a refreshing mint scent

Powder Holder Covers

- Three translucent color coded holder covers are included (blue, green, purple)

Quick Connection

- Simply pull back the quick release ring to release the nozzle from the handpiece



P560

Description	Model		Suggested Retail Price
Air Polishing System	P555-4H		\$270.00
Air Polishing System		P557-4H	\$310.00
Package includes the following:			
1 – 130g of prophyl powder	✓	✓	
3 – color coded (red, blue, green) powder holder covers	✓		
2 – nozzles (360 degree swivel)	✓	✓	
1 – cleaning tool	✓	✓	
1 – O-ring set	✓	✓	
1 year warranty	✓	✓	
Options:			Retail
P560	130g of prophyl powder		\$14.90
P660	Nozzle for P555		\$56.00
P665	Nozzle for P557		\$68.00

Shipping Dimensions: 9"L x 8"W x 3"H / 2 lbs.
(10/cs.) 13"L x 13"W x 10"H / 12 lbs.

Ultrasonic Cleaner


Dentsonic UC 400 / 1000 / 1500



Stainless Steel Cover



Stainless Steel Basket



Drain kit only available on UC-1000/1500

Model	UC-400	UC-1000	UC-1500
Input 115V / 60Hz	✓	✓	✓
Input 230V / 50Hz	✓	✓	✓
Overall Dimensions L x W x H	10.6" x 6.3" x 9.5" 23cm x 16cm x 24cm	14.2" x 10.6" x 12.2" 36cm x 27cm x 31cm	15" x 13" x 12.5" 39cm x 33cm x 31cm
Tank Internal Dimensions L x W x H	9.5" x 5.5" x 4" 24cm x 14cm x 10cm	11.8" x 9.5" x 5.9" 30cm x 24cm x 15cm	13" x 11.8" x 5.9" 33cm x 30cm x 15cm
Ship Weight	10 lbs. 4.5Kg	19 lbs. 8.6Kg	24 lbs. 10.6K
Tank Capacity	3.8 Qts. / .95 Gallon 3.59 Liter	10.6 Qts. / 2.65 Gallon 10 Liter	16 Qts. / 4 Gallon 15 Liter
Timer	LCD Display 1-99 minutes	LCD Display 1-99 minutes	LCD Display 1-99 minutes
302 Stainless Steel	✓	✓	✓

Features

- Powerful 40 KHz of cavitation for thorough, quieter cleaning
- Stainless steel cover and tank with welded drain
- Durable vinyl clad steel housing is scratch resistant and easy to clean
- Digital timer for quick and easy process
- Sweep technology ensures exceptional cleaning efficiency
- 1 year warranty
- 3 tank sizes available (3.8 qt, 10.6qt, 15.8qt)
- Includes Stainless Steel basket and lid

Description	Model	Retail Price
UC-400 Ultrasonic Cleaner (3.8 qt.)	UC-400	\$698.00
UC-1000 Ultrasonic Cleaner (10.6 qt)	UC-1000	\$1,338.00
UC-1500 Ultrasonic Cleaner (15.8 qt)	UC-1500	\$2,100.00

Package includes the following:

Ultrasonic cleaner	✓	✓	✓
Stainless steel basket	✓	✓	✓
Stainless steel lid	✓	✓	✓
Timer	✓	✓	✓
Drain kit	✗	✓	✓

Options	Retail Price
UC400-BASKET	Replacement Basket (3.8 qt.) \$90.00
UC1000-BASKET	Replacement Basket (10.6 qt) \$140.00
UC1500-BASKET	Replacement Basket (15.8 qt) \$180.00

Ultrasonic Cleaner

Dentsonic UC 250 / 300



UC-250



UC-300



UC-250 Basket



UC-600-GBS
Glass Beaker Set
500 ML x 2 and holder
Fits: UC300, 600, 750



UC-300-S

Features

- Powerful 35 KHz of cavitation for thorough, quieter cleaning
- Stainless steel tank
- Durable plastic housing is scratch resistant and easy to clean
- Digital timer for quick and easy process
- Sweep technology ensures exceptional cleaning efficiency
- 2 tank sizes available 2.5qt., 3.0qt.
- Includes basket and lid
- 1 year warranty

Shipping Dimensions

UC-250	14"L x 11"W x 10"H / 8 lbs. (4/cs) 23L x 15W x 20H / 35lbs.
UC-300	20"L x 14"W x 12"H / 16 lbs. (2/cs) 20"L x 14"W x 24"H / 34Lbs

Model	UC-250	UC-300
Input 115Vac	✓	✓
Input 220Vac	✓	✓
Frequency	50/60 Hz	50/60 Hz
Overall dimensions	12.6" x 9.6" x 8.0"	17.5" x 10.8" x 8.1"
Tank Internal Dimensions	9.8" x 6.0" x 3.0"	9.8" x 7.0" x 3.0"
Shipping Weight	8 lbs	13 lbs
Tank Capacity	2.5 qt, 0.65 Gal, 2.46 L	3.0 qt, 0.8 Gal, 3.028 L
Digital Timer	LCD Display 1-30 min	LCD Display 1-30 min
302 Stainless Steel Tank	✓	✓

Description	Model	Suggested Retail Price
-------------	-------	------------------------

UC-250 Ultrasonic Cleaner 2.5 Qt	UC-250	\$300.00
UC-300 Ultrasonic Cleaner 3.0 Qt (W/ Stainless Steel Basket)	UC-300-S	\$460.00

Options	Model
---------	-------

UC-600-GBSGlass Beaker Set (500ML X 2 and holder) Fits UC 300, UC600, UC 750

Package includes the following:

Ultrasonic cleaner	✓	✓
Basket	✓	✓
Lid	✓	✓
Drain kit	✓	✓
Digital timer	✓	✓
Heater	✓	✓

Ultrasonic Cleaner

Dentsonic UC-600/750

Features

- Powerful 35 KHz of cavitation for thorough, quieter cleaning
- Stainless steel tank
- Durable plastic housing is scratch resistant and easy to clean
- Digital timer for quick and easy process
- Sweep technology ensures exceptional cleaning efficiency
- 2 tank sizes available: 1.6 Gallon/2.0 Gallon
- Includes basket and lid
- 1 year warranty



Shipping Dimensions

UC-600	23 "L x 16"W x 14 "H / 23 lbs.
UC-750	24 "L x 16"W x 14"H / 27lbs.

Model	UC-600	UC-750
Input 115Vac	✓	✓
Input 220Vac	✓	✓
Frequency	50/60 Hz	50/60 Hz
Overall dimensions	21.0" x 12.0" x 9.6"	21.3" x 13.3" x 10.9"
Tank Internal Dimensions	13" x 8.0" x 4.0"	13" x 9.4" x 4.1"
Shipping Weight	20 lbs	25 lbs
Tank Capacity	6.4 qt, 1.6 Gal, 6.056 L	8 qt, 2 Gal, 7.5 L
Digital Timer	LCD Display 1-30 min	LCD Display 1-30 min
302 Stainless Steel Tank	✓	✓

UC-600
W/Plastic
Basket



UC-600
W/S.S. Basket



UC-750
W/S.S. Basket



Description

UC-600 Ultrasonic Cleaner 1.6 Gallon Plastic Basket
UC-600 Ultrasonic Cleaner 1.6 Gallon Stainless Steel Basket
UC-750 Ultrasonic Cleaner 2.0 Gallon Stainless Steel Basket

UC-600
UC-600-S

Model

UC-750

Suggested Retail Price

\$798.00
\$818.00
\$1,090.00

Package includes the following:

Ultrasonic cleaner	✓	✓
Basket	✓	✓
lid	✓	✓
Drain kit	✓	✓
Digital timer	✓	✓
Heater	✓	✓

Options

UC600-BSKT-P	Replacement Plastic Basket
UC600-BSKT-SS	Replacement Stainless Steel Basket (70cm)
UC750-BSKT-SS	Replacement Stainless Steel Basket (70cm)

Suggested Retail Price

\$25.00
\$116.00
\$120.00

TPC Dentsonic Ultrasonic Cleaner

Cassette Cross reference

851 South Lawson Street City of Industry CA 91748 • tpcdental.com • 800-560-8222

	9.6"x5.9"x3.0"	9.8"x6.9"x3.1"	9.5"x5.5"x4"	12.7"x8.1"x3.9"	12.5"x9"x3.5"	11.8"x9.5"x5.9"	13"x11.8"x5.9"

Hu-Friedy							
Large Cassette 16-Instruments 8" x 11" x 1.25"				✓	✓	✓	✓
Large Cassette 20-Instruments 4.5" x 8" x 1"			✓	✓	✓	✓	✓
Small Cassette 8-Instruments 5.5" x 8" x 1.25"	✓	✓	✓	✓	✓	✓	✓
Small Cassette 10-Instruments 5.5" x 8" x 1.25"	✓	✓	✓	✓	✓	✓	✓
Exam Cassette 5-Instrument 3" x 8" x 1.25"	✓	✓	✓	✓	✓	✓	✓
Space Saver 7-Instrument 4.5" x 8" x 1"	✓	✓	✓	✓	✓	✓	✓
Double Decker Cassette 8-Inst. 3" x 8" x 1.5"	✓	✓	✓	✓	✓	✓	✓
Double Decker Cassette 14-Inst. 4.5" x 8" x 1.5"	✓	✓	✓	✓	✓	✓	✓
Double Decker Cassette 20-Inst. 6" x 9.7" x 1.5"		✓	✓	✓	✓	✓	✓
Signa Stat Cassette 12-Inst. 6.5" x 10.5" x 1.25"				✓	✓	✓	✓
Oral Surgery Cassette 12-Inst. 8" x 14.5" x 1.25"							
Radiography Cassette 4.5" x 8" x 1.5"		✓	✓	✓	✓	✓	✓
Handpiece Cassette 3" x 8" x 1.25"	✓	✓	✓	✓	✓	✓	✓
Ultrasonic Insert Cassette 1.5" x 8" x 1.25"	✓	✓	✓	✓	✓	✓	✓

Zirc							
E-Z Jett Cassette 10-Place 7.93" x 7" x 1.125"				✓	✓	✓	✓
E-Z Jett Cassette 8-Place 7.93" x 4.4375" x 1.125"	✓	✓	✓	✓	✓	✓	✓
E-Z Jett Cassette 5-Place 8" x 3.125" x 1.1875"	✓	✓	✓	✓	✓	✓	✓
Compact Cassette 7.125" x 3.875" x 0.625"	✓	✓	✓	✓	✓	✓	✓
Mini Matt 3.875" x 3.5" x .5"	✓	✓	✓	✓	✓	✓	✓
Instrument Mat 6.75" x 2" x 0.375"	✓	✓	✓	✓	✓	✓	✓
Standard Steri-Container 8.125" x 1.875" x 1.875"	✓	✓	✓	✓	✓	✓	✓
Compact Steri-Container 7.125" x 1.5" x 1.5"	✓	✓	✓	✓	✓	✓	✓

Warning: There should absolutely be NO metal-metal contact inside of the tank during ultrasonic cycle
This includes No direct: Loose Instrument-to-Instrument, Cassette-to-Tank, Cassette-to-Cassette, Cassette-to-Basket; metal contact. Make sure Cassettes have standard silicone insulations and feet; if no basket is used.
 Tank Dimensions are measured **without** baskets; basket sizes vary depending on unit.
 Please contact your cassette manufacturer or TPC if you have any questions.

D650N

Digital Amalgamator



D-650N

Features

- ISO 9002 / ISO 13488 quality management certified
- Input voltage: 110V or 230V
- Touch pad selector
- Pre-programmed time settings
- Quiet, high-speed mix capsules
- Built-in automatic timer
- Accepts all capsules
- Customized trituration times
- Vibration frequency $\geq 4,400$ RPMs
- 1 year warranty

Description	Model	Suggested Retail Price
D650N Digital Amalgamator	D-650N	\$450.00

Package includes the following:

D650N digital amalgamator
AC power cord

✓
✓

Shipping Dimensions: 13"L x 10"W x 10"H / 10 lbs.
(2/cs.) 20"L x 14"W x 10"H / 20 lbs.

N7

Micromotor and Handpiece



MM7

Features

- Powerful torque and variable speed motor
- No heat after long hours of operation
- Easy change of driving direction for efficiency
- Variable speed foot control
- Compact and ergonomic design
- Automatic stop with alarm beep and overload reset
- Includes control unit, motor handpiece and foot pedal
- 0–35,000 RPMs
- 110/220V ; 50/60Hz, 0.7A
- 1 year warranty

Description	Model	Suggested Retail Price
TPC Marathon N7 Micromotor and Handpiece	MM7	\$598.00
TPC Marathon N2 Micromotor and Handpiece	MM2	\$398.00

Package includes the following:

N7 micromotor and handpiece	✓
AC power cord	✓
Foot pedal	✓

Options:

		Retail
MM7-HP	Handpiece only with cable	\$358.00
MM7-MOTOR	Motor control box	\$358.00
MM7-VFC	Replacement foot control	\$69.00
MM7-BR	Replacement carbon brushes	\$15.00

Shipping Dimensions: 11"L x 11"W x 5"H / 7 lbs.
(8/cs.) 24"L x 12"W x 21"H / 58 lbs.

Disposable Micro Applicators



M6500-R



M6500-F



M6500-SF



M6500-XSL



• 36 Boxes/case

- Four individual dispenser tubes included with each box

Part No.	Description	Suggested Retail Price
M6500-R	Disposable Micro Applicators, 400 ct. (Green / Orange)	\$17.50
M6500-F	Disposable Micro Applicators, 400 ct. (Yellow / Blue)	\$17.50
M6500-SF	Disposable Micro Applicators, 400 ct. (White / Purple)	\$17.50
M6500-XSL	Disposable Micro Applicators, 400 ct. (Black)	\$17.50

Package includes the following:

Four dispenser tubes, 100 brushes per tube

Shipping Dimensions:
(36 boxes/ctn) 420 x 320 x 370 mm/16kg
17”L x 13”W x 15”H/35 lbs

Features

- Non-linting fibers
- Bendable stiff tip
- Dispensing tubes eliminate the potential for cross contamination.

Disposable Air/Water Syringe Tips



ITEM	Description	Retail
P7700	White Syringe Tip W/Metal Interior	\$40.00



ITEM	Description	Retail
P7705	White Syringe Tip W/Plastic Interior	\$30.00



ITEM	Description	Retail
P7725	Clear Syringe Tip W/ Colored interior	\$30.00

- 2 Syringe tip dispensers included (125 tips per dispenser)
- 250 tips/Box
- 36 boxes per case
- Shipping dimensions (36Bx/Case)
467 x 340 x 345 / 14.5kg 18"L x 13"W x 14"H / 32 Lbs



ITEM	Description	Retail
P7700-Bulk	White Syringe Tip W/ Metal Interior	\$200.00
P7705-Bulk	White Syringe Tip W/ Plastic Interior	\$148.00
P7725-Bulk	Clear Syringe Tip W/ Colored interior	\$148.00

- Does not include dispenser
- 200 tips /Bag
- 8 Bags / Box (1600 tips)
- 5 Boxes / case (8000 tips)
- Shipping dimensions (5 Bx/Case)

380 x 255 x 230 / 8.2 kg
18"L x 10"W x 9"H / 18 Lbs

Disposable Prophylaxis Angles



(Teal Color / White Cup)
Regular



(Pink Color / Grey Cup)
Soft

- Latex free
- Individually packaged
- PA500
- 500PCS/BOX, 8 Boxes/case
- PA144
- 144PCS/BOX, 20 Boxes/case



Shipping dimensions (8 Bx/Case)
PA500: 570 x 350 x 795mm / 16.5kg
22" x 14" x 31" / 36 Lbs

Shipping dimensions (20 Bx/Case)
PA144: 500 x 350 x 730mm / 11.75kg
20" x 14" x 29" / 26 Lbs

ITEM	Description	Retail	1 Case	3+ Case
PA-144R	Disposable prophylaxis angles, 144 ct., Regular	\$56.00 / Box	\$51.00 / Box	\$47.00 / Box
PA-500R	Disposable prophylaxis angles, 500 ct., Regular	\$178.00 / Box	\$170.00 / Box	\$160.00 / Box
PA-144S	Disposable prophylaxis angles, 144 ct., Soft	\$56.00 / Box	\$51.00 / Box	\$47.00 / Box
PA-500S	Disposable prophylaxis angles, 500 ct., Soft	\$178.00 / Box	\$170.00 / Box	\$160.00 / Box

Hygiene Prophy Handpiece



	PH-777 (Red) PH-777 (BLK)		
Description	Model	Suggested Retail Price	
Prophy Aide Hygiene Handpiece	PH-777(Red)	\$398.00	
Prophy Aide Hygiene Handpiece	PH-777(Black)	\$398.00	
Package includes the following:			
4 Hole fixed back end	✓		
5,000 RPM	✓		
Twist open chuck	✓		
360 Head Swivel Capabilities	✓		
Forward and Reverse operation	✓		

Disposable Prophylaxis Cups & Brushes



Description			Model	Retail
Latch type prophylaxis brush	flat / white	pkg of 144	PB300	\$43.90
Latch type prophylaxis brush	flat / black	pkg of 144	PB310	\$45.90
Screw-on type prophylaxis brush	flat / white	pkg of 144	PB320	\$41.20
Screw-on type prophylaxis brush	pointed / white	pkg of 144	PB330	\$41.20
Latch type prophylaxis brush	pointed / white	pkg of 144	PB340	\$43.90
Latch type prophylaxis brush	pointed / black	pkg of 144	PB350	\$45.90
Latch type prophylaxis brush	cup shade / white	pkg of 100	PB360	\$43.90
Latch type prophylaxis brush	cup shade / bristle	pkg of 100	PB370	\$48.00
Latch type prophylaxis cup	webbed / white	pkg of 144	PC505	\$25.20
Screw-on type prophylaxis cup	webbed / white	pkg of 144	PC515	\$28.00
Pointed-on prophylaxis cup	webbed / yellow	pkg of 144	PC525	\$28.00
Snap-on prophylaxis cup	webbed / white	pkg of 144	PC535	\$15.00

Packing : (100 boxes / ctn) 415 x 330 x 185 mm

Disposable Pre-Bent Tips



18G-Pink



22G-Grey



20G-Yellow



25G-Blue



20G-Black

- Blunt end, pre-bent tips
- Provides fast and accurate dispensing of etchant, resins, sealants, and flowable composite

- Universal Luer Lock design
- Available in color-coded 18G, 20G, 22G, and 25G size
- 100 PCS/BX, 100 BX/CASE
- Shipping dimensions (100 BX/CASE)


- 17" x 13" x 8" / 12 Lbs

ITEM	Description	Packaging	Retail
PBT-2210	18G-Pink	100/BX, 100BX/CS	\$4.00 / Box
PBT-2220	20G-Yellow	100/BX, 100BX/CS	\$4.00 / Box
PBT-2230	20G-Black	100/BX, 100BX/CS	\$4.00 / Box
PBT-2240	22G-Grey	100/BX, 100BX/CS	\$4.00 / Box
PBT-2250	25G-Blue	100/BX, 100BX/CS	\$4.00 / Box

High Volume Evacuation Tips (HVE)



- Vented on one side/non-vented on the other
- Designed to fit all high volume evacuation valves
- Made of rigid plastic
- Shipping dimensions (4 Bx/ Case)




HVE-1000	White Color High Volume Combo tip	100 /Bag, 10 Bags/Bx, 4 Bx/Cs	\$43.00/Box
----------	-----------------------------------	-------------------------------	--------------------

Saliva Ejectors (SE)



- Soft, non-removable tip for patient safety
- Wire reinforced tube maintains desired shape
- Shipping dimensions (4 Bx/Case)



Suggested Retail Price			
SE-1000 W/W	White with White tip	100 / Bag, 10 Bags/ Bx, 4 Bx / Cs	\$43.00/Box
SE-1000 B/C	Clear with Blue tip	100 / Bag, 10 Bags/ Bx, 4 Bx / Cs	\$43.00/Box

Surgical Aspirator Tips



- Designed to fit all high-volume evacuation valves
- 30° angle for maximum patient comfort
- Available in 1/8", 1/4"
- 25 Pcs/Bag, 20Bag/Bx, 8 Bx/Case
- Shipping dimensions: 8 Bx(160Bags)/Case

578 x 358 x438 mm / 20 Kg 23" x 14" x 17"

ITEM	Description	Packaging	Retail
SAT-SW	Small 1/8" x 6 1/4"(White)	25 Pcs/Bag, 20 Bag/Bx, 8 Bx/Cs	\$80.00 / Box
SAT-LG	Large 1/4" x 6 1/4"(Green)	25 Pcs/Bag, 20 Bag/Bx, 8 Bx/Cs	\$80.00 / Box

Tray Cover Papers



- 8 1/2" x 12 1/4"
- With contoured edges, perfectly fits size "B" Trays
- Shipping dimensions (4 BX / CASE)

18" x 13" x 10" / 42 Lbs



Blue TC-6231	Purple TC-6233	Green TC-6235	White TC-6237	Pink TC-6239
-----------------	-------------------	------------------	------------------	-----------------

ITEM	Description	Packaging	Retail
TC-6231	Blue color	1000 PCS / BX, 4 BX / CS	\$120.00
TC-6233	Purple color	1000 PCS / BX, 4 BX / CS	\$120.00
TC-6235	Green color	1000 PCS / BX, 4 BX / CS	\$120.00
TC-6237	White color	1000 PCS / BX, 4 BX / CS	\$120.00
TC-6239	Pink color	1000 PCS / BX, 4 BX / CS	\$120.00

Plastic Set Up Tray



ITEM	Description	Packaging	Retail
PT-2021	Blue color	20 / Bx, 4 Bx / Cs	\$6.00
PT-2022	White color	20 / Bx, 4 Bx / Cs	\$6.00
PT-2023	Pink color	20 / Bx, 4 Bx / Cs	\$6.00
PT-2024	Green color	20 / Bx, 4 Bx / Cs	\$6.00
PT-2025	Purple color	20 / Bx, 4 Bx / Cs	\$6.00
PT-2026	Yellow color	20 / Bx, 4 Bx / Cs	\$6.00



Stainless Steel Tray Rack



- Fit for size " B " Tray • Up to 8
- Trays • Shipping dimensions (4 Bx / Case) 310 x 290 x 425mm / 12.4 kg
- 12" x 11" x 17"/32 Lbs

ITEM	Description	Packaging	Retail
TR-1008	Stainless Steel Rack	1 Set / Bx, 4 Bx / Cs	\$70.00

Alligator X-Ray Film Holder



- Straight
- Cold sterile only
- 3 Pcs / Box

Description	ITEM	Retail
Film Holder-straight	M7000-S	\$12.50

Post Mount Utility Tray



- Available in white, beige, grey color
- Fits post 1 3/8" -2 1/2"
- Dimensions: 13" x 10"
- Contoured edge to contain spills and materials
- 20 trays / case
- Shipping dimensions (20 Trays / Case)

Description	ITEM	Retail
Post mounted utility tray (White, Beige, Grey) (Specify color)	#US100	\$70.00

Packing

1 Set/Box, 20 Sets/Case
645 x 540 x 455 mm / 16.8 Kg
25" x 22" x 18" /37 Lbs

Applier & Activator



AP311 Applier




AP322 Activator

ITEM	Description	Packaging	Dimension & Weight	Retail
AP311	Applier	1 Pcs/Bx, 100 Bx/Cs	480 x 340 x 310 mm/15.56 Kg 19" x 13.5" x 12.2" / 34.5Lbs	\$43.00
AP322	Activator	1 Pcs/Bx, 100 Bx/Cs	575 x 350 x 465 mm/ 24.96 Kg 23" x 14" x 18" / 55 Lbs	\$43.00

Universal Film Mounts




CP102U



- Universal select & cut
- Roll of 500 strips of 2 mounts
- 1000 windows / Box
- 12 Bx / Cs
- Dimensions & weight(12Bx/Cs)

410 x 380 x 380 mm /17.24 Kg
16" x 15" x 15" /38 Lbs

ITEM	Description	Retail
CP102U	1000 Windows/Bx, 12 Bx/Cs	\$48.00



CLEAR POCKET MOUNT



NOTCHED MOUNT



OPEN WINDOW MOUNT

FILM MOUNT CROSS REFERENCE CHART

TPC				CASTERM		RINN			FLOW		ADA
POCKET	NOTCHED	WINDOW		POCKET	POCKET	POCKET	NOTCHED	WINDOW	NOTCHED	WINDOW	WINDOW
CP102U				N.A.	20-0102						
CP8 PEDO				108 PEDO	N.A.				N.A.		N.A.
CP1C	N1H2	W1H2			20-1018	25-0178	30-6011		23200	601-200G	22004
	N2H0	W2H0			20-1024	25-0811	30-4021		23608	402-608G	22616
CP2H	N2H2	W2H2		102H	20-1020	25-0179	30-6021		23203	602-203G	22020
	N2H3	W2H3			20-1023	25-0180	30-7021		N.A.	702-206G	22061
	N4H0	W4H0			20-1044	25-0187	30-4040		23609	404-609G	22624
CP4C	N4H2	W4H2		104C	20-1048	25-0181	30-6040		23209	604-209G	22103
CP4H				104H	20-1040	25-0815	30-6045		23239	6045-239G	22152
CP4H2V	N4H2V2	W4H2V2		106HV	20-1062	25-1061	30-6062		23216	606-216G	22202
CP6H	N6H2	W6H2		106H	20-1060	25-0182	30-6061		23215	606-215G	22194
CP8C	N8H2	W8H2		108C	20-1088	25-0825	30-6081		23219	608-219G	22251
CP10C				110C	20-1108	N.A.	30-6101		23220	610-220G	22285
CP8H6V	N8H6V2	W8H6V2		114HV	20-1142	25-0184	30-6141		23225	614-225G	22335
CP12H6V	N12H6V2	W12H6V2		116-118HV	20-1182	25-0187	30-6180		23228	618-268G	22400
CP10C	N10BW2	W10BW2		110C	20-1108	25-0821	30-2345		23261	610-261G	22558
	N10BW3				N.A.	25-0820	30-3456		23262	N.A.	N.A.
		W20BW2			80-0138	N.A.	30-4567		N.A.	618-230G	22574



2H #2 6H #0



1H #2



2H #0



2H #2



2H #3



4H #0



4H #2



4H



4H 2V #2



6H #2



8H #2



10H #2



8H 6V #2



12H 6V #2



10BW #2




10BW #3



20BW #2

Film Positioning System Complete



Operation Manual

- Be sure the Bite Block is centered in aiming ring
- Place plain or white side of film toward X-Ray Tubehead
- Do not bend back the bite block when inserting the film or phosphor plates

Description

X ray film positioning system

ITEM

FPS3000

Retail

\$56.00

Packing


1 set/box, 10 sets/carton
520 x 370 x 280mm / 3.8kg
20.5" x 15" x 11" / 8.5 Lbs

Lab Products (Crown Boxes and Articulators)



ITEM	Description	Packaging	Retail
LAB-CRWBX-2	Plastic Crown Box W/2" Foam 500 / cs(Blue)	375 x 290 x 490 mm/7.88 Kg 15" x 12" x 20" / 18 Lbs	\$170.00
LAB-CRWBX-2BLK	Plastic Crown Box W/2" Foam 500 / cs(Black)	375 x 290 x 490 mm/7.88 Kg 15" x 12" x 20" / 18 Lbs	\$170.00
LAB-CRWBX-1	Plastic Crown Box W/1" Foam 1000 /cs(Blue)	375 x330 x 330 mm/6.8 Kg 15" x 13" x 13" / 15 Lbs	\$210.00
LAB-CRWBX-1BLK	Plastic Crown Box W/1" Foam 1000 /cs(Black)	375 x330 x 330 mm/6.8 Kg 15" x 13" x 13" / 15 Lbs	\$210.00
LAB-DART-WH	Disposable Articulators White 100 / Bag, 20 Bags / Case	483 x 432 x 305 mm/14.6 Kg 19" x 17" x 12" / 32 Lbs	\$300.00
LAB-DART-BLK	Disposable Articulators Black 100 / Bag, 20 Bags / Case	483 x 432 x 305 mm/14.6 Kg 19" x 17" x 12" / 32 Lbs	\$280.00

Disposable Evacuation Traps



- Directly replaces current evacuation trap
- Easily removed and replaced
- Made of premium plastic

ITEM	Description	Packaging	Retail
53016-S	Diameter 1 7/8", Clear with trap cover 100 / Bx, 8 Bx / Cs	673 x 345 x 426 mm / 9.7 Kg 27" x 14" x 17" / 22 Lbs	\$239.20 / case
53016-L	Diameter 2 1/8", Clear with trap cover 100 / Bx, 8 Bx / Cs	618 x 443 x 390 mm / 12.5 Kg 24" x 18" x 16" / 28 Lbs	\$239.20 / case
DT-6902	Diameter 1 7/8", Green color 144 / Bx, 8 Bx / Cs	550 x 550 x 420 mm / 11.1 Kg 22" x 22" x 17" / 25 Lbs	\$152.00 / case
DT-6901	Diameter 2 1/8", Blue color 144 / Bx, 8 Bx / Cs	550 x 550 x 460 mm / 14.2 Kg 22" x 22" x 18" / 32 Lbs	\$152.00 / case

Disposable Vacuum Pump Canisters



Small (2 3/4" x 3 5/8")
VC-3001

Large (3 1/2" x 4 3/8")
VC-3002

- Designed to replace the present vacuum canister and mesh filter in your central vacuum pump
- Clear plastic for easy viewing
- Save time on cleaning

ITEM	Description	Retail
VC-3001	Small canister (2 3/4" x 3 5/8") with mesh screen, 12 / Bx, 6 Bx / Cs	\$56.00
VC-3002	Large canister (3 1/2" x 4 3/8") with mesh screen, 8 /Bx, 6 Bx / Cs	\$59.00

Disposable Ear-Loop Face Masks



Blue



Pink



- Meets ASTM Level 2 Standard and EN 14683 Rating
- Superior breathability
- Splash resistant
- Ultra soft inner layer
- 50 /Bx, 20 Bx/Cs
- 3 Layer

- Shipping Dimensions (20 Bx/CS)

- 534 x 203 x 381 mm / 4.54 Kg



ITEM	Description	Retail
PPE-PEMBL-2	Ear loop Level 2 Blue Color, 50 Pcs / Bx, 20 Bx / Cs	\$100.00 / case
PPE-PEMPK-2	Ear loop Level 2 Pink Color, 50 Pcs / Bx, 20 Bx / Cs	\$100.00 / case

Disposable Isolation Gowns



Disposable – Knit Cuff – Latex Free

- Knit cuff, Latex free, fluid resistant
- One size fit for most people
- Extra coverage with back and neck tie
- Provides Level 3 barrier protection
- Weight 45 G
- Size 49.2" (long), 27.6" (wide), 23.62" (sleeve)
- Standard GB/T 38462-2020
- Color BLUE
- 50 pcs/box
- Shipping dimensions (50 Pcs/Box) Standard



ITEM	Description	Packaging	Retail
#PPE-IGKB	Isolation Gown, Disposable case of 50	432 x 368 x 406 mm / 6.6 Kg 17" x 14.5" x 16" / 14.5 Lbs	\$270.00

Disposable Face Shields



Features

- Anti-fog coating provides a clear view
- Comfortable elastic band holds face shield in place
- Soft foam allows use of mask and eyeglasses
- Covers: face, side of face, chin & below chin
- One size fits all
- Single use

ITEM	Description	Packaging	Retail
PPE-FS (24)	Disposable face shields, 24 / Bx	50 x 330 x 150 mm / 1.28 Kg 14" x 13" x 6" / 3 Lbs	\$28.00
PPE-FS(100)	Disposable face shields, 100 / Bx	430 x 360 x 360 mm / 5.16 Kg 17" x 14" x 14" / 12 Lbs	\$100.00
PPE-FS(240)	Disposable face shields, 240 / Bx(10 x 24 / Bx)	810 x 380 x 690 mm / 15 Kg 32" x 15" x 27" / 33 Lbs	\$224.00

Digital Thermometer



Non-Contact Digital-Infrared Features

- Non-contact design
- Hygienic & professional
- Comfortable to hold
- Fever warning
- Multiple measurement modes
- 1 year warranty

Part No.	Description	Retail Price
PPE-DT	Thermometer, Non-Contact Digital Infrared	\$70.00

Press Ring Syringe

3-Way Syringe Push Ring



Part No.	Description	Retail Price
TPC3WAY-PR	3-Way Syringe	\$70.00
#52046	Replacement Cartridge (2pcs per pack)	\$17.00
• 360 degree tip rotation • Accurate air andf water control • Autoclaveable quick disconnect aerator tip • Quick connect push ring		

Autoclavable 3-Way Syringe Tip



Part No.	Description	Retail Price
M7500	Autoclavable 3 - way syringe tips (Box of 100) (Compatable with DCI 3-way syringe)	\$290.00
• Syringe tip only • Box of 100		

Shipping Dimensions: 6"L x 4"W x 2"H / 2 lbs.

Deluxe Autoclaveable Valve

Autoclavable High Vacuum Valve - Short



Part No.	Description	Retail Price
#P2330	Deluxe Autoclavable high volume ejector-short	\$41.00
#P2166	Replacement Deluxe HVE lever & spool	\$12.00
• Autoclavable valve • Finger tip lever control for on/off flow • Swivel quick disconnect for 360° rotation • Adapts to 1/2" ID vacuum tubing • Disconnect at swivel for simple autoclaving • Anodized aluminum one piece valve body construction with anodized aluminum swivel		

Autoclavable High Vacuum Valve - Long



Part No.	Description	Retail Price
#P2335	Deluxe Autoclavable high volume ejector-long	\$41.00
#P2166	Replacement Deluxe HVE lever & spool	\$12.00
• Autoclavable valve • Finger tip lever control for on/off flow • Swivel quick disconnect for 360° rotation • Adapts to 1/2" ID vacuum tubing • Disconnect at swivel for simple autoclaving • Anodized aluminum one piece valve body construction with anodized aluminum swivel		

Autoclavable Saliva Ejector Valve



Part No.	Description	Retail Price
#P2340	Deluxe Autoclavable saliva ejector valve	\$35.00
#P2169	Replacement Deluxe saliva ejector lever & spool	\$12.00
#P2185	Replacement deluxe saliva ejector tip only (5 pcs/Bag)	\$13.00
• Autoclavable valve • Finger tip lever control for on/off flow • Swivel quick disconnect for 360° rotation • Disconnect at swivel for simple autoclaving • Anodized aluminum one piece valve body construction with anodized aluminum swivel		

Parts



Curing Bulb

- 12 Volt 75 Watt
- Power Light 100
- Replacement for PL100

#P201 \$ 39.00



Operatory Light Reflector

- Quality coated glass
- No core charge
- Original factory OEM specifications
- Replacement reflector for Belmont®, LSM® and Celux light

#P705 \$ 160.00



Scaler Tubing 30K without Connector

- 25KHz 3-wire handpiece assembly
- 7' length
- Only Compatible with TPC brand scalers

#S8002 \$ 120.00



Scaler Tubing 25K without Connector

- 25KHz 3-wire handpiece assembly
- 7' length
- Only Compatible with TPC brand scalers

#S8001 \$ 120.00



Mini HP Control Block

- Mini control block for 3 handpieces
- Solid brass construction
- Extra mounting holes allow a variety of placements
- 10-32 threaded barbs
- Longest lasting inter-woven diaphragm
- High flow design

#52003 \$110.00

#52004 Diaphragm \$13.00



Water Bottle System

- 750 ml Extra Heavy duty bottles
- Easy to read pressure gauge
- Pre-set integrated regulator
- Ideal for bactericidal solutions
- Post clamp
- Will mount to vertical surface

#W5030 \$150.00

#W5035 – 750 ml bottle only \$15.00

#W5040 - 1000 ml bottle only \$18.00



Flow Control Valve

- Chrome plated brass

#52007 \$22.00



Momentary Valve

- Chrome plated brass

#52008 \$14.00



On/Off 2-Way Valve

- Chrome plated brass

#52009 \$14.00

Parts



On / Off 3 Way Valve

- Chrome plated brass


#52010 \$15.00



Check Valve

- Chrome plated brass

#52019 \$10.00



Routing Valve

- Chrome plated brass

#53014 \$15.90



Non-Retraction Water Relay Valve

- Chrome plated brass
- Equivalent to DCI #7068

#52011 \$35.00



4-hole Adapter, Metal

- Quality plated brass

#52039 \$9.00



Plated Brass Nut 4-hole Midwest

- Brass extends HP thread life
- Standard connector

#52040 \$4.00



HP Tubing Adapter, 2 - hole HP to 4-hole Tubing

- Adapts 2 hole HP to 4-hole tubing

#52041 \$39.00



HP Tubing Adapter, 4-hole HP to 2-hole Tubing

- Adapts 4-hole HP to 2-hole tubing

#52042 \$39.00



Stainless Steel Tray

- Premium quality stainless steel
- Size 13F (9-3/4" x 13-1/2")

#52026 \$42.00

#52027 Skid Pad \$3.16

Parts



Female Quick Disconnect

- Panel mount with auto shut-off
- 1/4" poly-flow rear connection
- Chrome plated brass
- 1/4" Coupler

#53012 1/4" Coupler **\$15.00**

#53013 3/8" coupler **\$25.00**



Male Quick Disconnect w/ Shut-off

- With auto shut-off
- 1/4" poly-flow rear connection
- Chrome plated brass
- 1/8" Male QD

#53042 **\$ 13.00**



1/4" Coupler with shut off

- barb type
- Chrome plated brass
- Equivalent to DCI #0006

#53047 **\$ 15.90**



Foot Control, Wet/Dry

- Precision-touch disc foot control with wet/dry toggle
- Weighted for maximum stability
- Durable chrome steel dome
- Smooth, variable pressure control
- Non-skid foot pads

#53200 **\$ 116.00**

#53200-TS Toggle (only) **\$1.98**

#53121 Toggle assembly **\$24.00**



Foot Control, Shrouded

- Large stable base with shroud
- Precision-touch disc foot control
- Weighted for maximum stability
- Durable chrome steel dome
- Smooth, variable pressure control
- Non-skid foot pads

#53117 **\$ 78.00**



Air Electric Switch

- Normally open air switch
- 15 PSI activated Check this spec

#A800N-SW **\$ 28.50**



4 Hole HP Straight Asepsis Assembled

- Assembled with metal adapter / nut
- Silicone Jacket
- 5.5' length

#52036 **\$33.00**



Syndge Tubing

#52038 **\$18.00**

- 6' Silicone Jacket tubing



Cuspidor Strainer

#80332 **\$2.50**

- 2000-C 2.0 cuspidor bowl strainer
- 5 pcs / bag
- Comparable to PINN. #6100

Parts



53016-S

53016-L

Disposable Amalgam Traps

- 53016-S-100 1 7/8" Clear
- 53016-L-100 2 1/8" Clear
- Clear trap pkg of 100

#53016-S-100 **\$45.00**

#53016-L-100 **\$45.00**



Solids Collector Canister

- Includes: Lid, Base, o-ring
- Bottom ported 2015 unit
- 2 1/8" trap

#53008 **\$34.00**



Solids Collector Canister

- Includes: Lid, Base, o-ring
- Side ported for 3550 unit
- 2 1/8" trap

#53008-R **\$34.00**



Solids Collector Canister

- Includes: Lid, Base, o-ring
- Bottom ported, 2600 unit
- 1 7/8" trap

#53008-2600 **\$34.00**



Nylon Bushing

#LP-RI-S **\$25.00**

- OD 1 7/8"
- ID 1 3/8"



Nylon Bushing

#LP-PC-L **\$25.00**

- OD 1 7/8"
- ID 1 11/16"



Nylon Bushing

#L600N-ARM-RP24

- \$22.00**
- OD 1 7/8"
- ID 40mm
- L600N Arm Only



48" Light Post

#LP48 **\$210.00**

- 48" light post
- LP-PC-L included
- Optional: LP-RI-S



Swing Light Post Upper part

#LP-26 **\$258.00**

- Upper post for LP100



Chair Adaptor

#BELBRKT **\$238.00**

- BLMT Chair Bracket
- Fits BLMT 20 / 50 chairs
- Reversible left or right



BLMT Radius Bracket (Unit)

#2400-BRKT-BEL

\$550.00

- Fits: BLMT 20, 50



BLMT Rear Radius Bracket (Oblong Shape)

#2450-BRKT-BEL-N **\$252.00**

- Fits: BLMT 50 Rear Mount



BLMT Rear Radius Bracket (U shape)

#2450-BRKT-BEL **\$252.00**



Chair Adaptor

L2001 \$

- TPC chair bracket
- Fits Laguna 2.0 only
- Reversible left or



Chair Adaptor

M4001 \$


- TPC chair bracket
- Fits Mirage 2.0 only
- Reversible left or right



Chair Adaptor

\$8

- TPC chair bracket
- Fits Mirage 1.0 only
- Reversible left or right



Chair Adaptor

#TPCRDSBKT

- TPC swing bracket
- Fits Mirage 1.0 - 2.0
- Fits Laguna 1.0, LP105 Post



Chair Adaptor

#L2000-BKT

- Laguna chair bracket
- Fits Laguna 1.0 only

Parts PC2630 / PC 2930



Waste Container

- Fits: PC-2630,2635,3600
OMC2375-CVSL
- Includes Cap and jar
- Color: Blue

#PC-2630N-SWC

\$13.90



Water Bottle

- Water Bottle (Clear)
- Fits: PC2630, 2635, 3600,
OMC2375, 2930
- 700 ml

#PC-2630N-CWB

\$10.00



H.P. holder Automatic

- Panel Mount HP holder
- Fits: PC2930 only
- Plastic w/ Metal base
- Includes: Base, Holder

#PCP-N-RP1

\$



H.P. holder

- Panel Mount HP holder
- Fits: PC2630 only
- Plastic

#PC-2630N-RP11

\$31.00



On / Off 3 Way Valve

- Chrome plated brass
- Water bottle on/off

#52010 **\$15.90**



Routing Valve

- Chrome plated brass
- H.P Selector Valve

#53014 **\$15.90**



On/Off 2-Way Valve

- Chrome plated brass
- Suction on / off valve

#52009 **\$14.00**



Flow Control Valve

- Chrome plated brass
- Water Adjustment for PC
2630, 2930

\$22.00



1/4" Coupler with shut off

- barb type
- Chrome plated brass
- Equivalent to DCI #0006

#53047 **\$15.90**



4 Hole H.P. Tubing Silicone Jacket

- With metal adapter / nut
- Silicone Jacket
- 5.5' length

#52036 **\$33.00**



Water Relay Block

- 2 position water relay block

#PC--P2023 **\$39.00**



Foot Control, Shrouded

- Durable chrome steel dome
- Smooth, variable pressure control
- Non-skid foot pads

#53117 **\$78.00**



Syringe Tubing

- 6' Silicone Jacket tubing

#52038 **\$18.00**



Power Switch

- Main power switch

#PCN-P5001 **\$17.00**



110V Pressure relief valve

- 110V Pressure relief valve

#PC-P2031 **\$60.00**



Air Venturi Valve

- S.E. Air Venturi Valve

#PC-2630-P001 **\$99.00**

Parts MC3600CV-SL



Waste Container

- Fits: PC-2630, 2635, 3600 OMC2375-CVSL
 - Includes Cap and jar
 - Color: Blue
- #PC-2630N-SWC
\$13.90



Water Bottle

- Water Bottle (Clear)
 - Fits: PC2630, 2635, 3600, OMC2375, 2930
 - 700 ml
- #PC-2630N-CWB
\$10.00



H.P. holder Automatic

- Bar Mount HP holder
 - Fits: MC3600CV/SL
 - Plastic w/ Metal base
 - Includes: Base, Holder
- #MC3600-RP-32
\$40.00
#MC3600-RP-30 (Valve)
\$



Top Cover

- Plastic top cover
 - Fits: MC3600/3600CV-SL
- #MC3600-RP-1
\$40.00



On / Off 3 Way Valve

- Chrome plated brass
 - Water bottle on/off
- 52010 **\$15.90**



Routing Valve

- Chrome plated brass
 - H.P Selector Valve
- #53014 **\$15.90**



On/Off 2-Way Valve

- Chrome plated brass
 - Suction on / off valve
- #52009 **\$14.00**



Flow Control Valve

- Chrome plated brass
 - Water Adjustment for PC 2630, 2930
- #52007 **\$22.00**



1/4" Coupler with shut off

- barb type
 - Chrome plated brass
 - Equivalent to DCI #0006
- #53047 **\$15.90**



4 Hole H.P. Tubing Silicone Jacket

- With metal adapter / nut
 - Silicone Jacket
 - 5.5' length
- #52036 **\$33.00**



Water Relay Block

- 2 position water relay block
- #PC--P2023 **\$39.00**



110V Main Motor

- 110Vac motor with start cap
 - Quite operation
 - Included mounts, air filter
- #OMC2375CV-MTR-2 **\$650.00**



Syringe Tubing

- 6' Silicone Jacket tubing
- #52038 **\$18.00**



Power Switch

- Main power switch
- #PCN-P5001 **\$17.00**



110V Pressure relief valve

- 110V Pressure relief valve
- #PC-P2031 **\$60.00**



Tri Block

- 3 position Tri Block Complete
- #MC3600-RP-50 **\$100.00**

Parts OMC2375CV



Waste Container

- Fits: PC-2630,2635,3600 OMC2375-CVSL
 - Includes Cap and jar
 - Color: Blue
- #PC-2630N-SWC
\$13.90



Water Bottle

- Water Bottle (Clear)
 - Fits: PC2630, 2635, 3600, OMC2375, 2930
 - 700 ml
- #PC-2630N-CWB
\$10.00



H.P. holder Shell

- Bar Mount HP holder
- Fits: OMC2375CVSL
- Plastic Shell only

#OMC-RP-12-C
\$88.00



H.P. Holder Coplete

- Bar Mount HP holder
- Fits: OMC2375CVSL
- Includes: Shell, Block, Valve

#OMC-RP-12
\$42.00



N.C. Valve

- Normally Closed
 - HP activation valve
- #OMC-RP-13 **\$11.00**



Routing Valve

- Chrome plated brass
 - H.P Selector Valve
- #53014 **\$15.90**



On/Off 2-Way Valve

- Chrome plated brass
 - Suction on / off valve
- #52009 **\$14.00**



Flow Control Valve

- Chrome plated brass
 - Water Adjustment for PC 2630, 2930
- #52007 **\$22.00**



Pressure Gauge

- OMC Cabinet Pressure Gauge
 - Front Panel Mounted
 - Bracket not included
- #53047 **\$15.90**



4 Hole H.P. Tubing

Silicone Jacket

- With metal adapter / nut
 - Silicone Jacket
 - 5.5' length
- #52036 **\$33.00**



110V Main Suction Motor

- 110Vac motor with start cap
 - Quite operation
 - Included mounts, air filter
- #OMC-2375CV-MTR-1 **\$590.00**



Water Relay Block

- 2 position water relay block
- #PC--P2023 **\$33.00**



Syringe Tubing

- 6' Silicone Jacket tubing
- #52038 **\$33.00**



Power Switch

- Main power switch
- #PCN-P5001 **\$17.00**



110V Pressure relief valve

- 110V Pressure relief valve
- PC-P2031 **\$60.00**



Tri Block

- 3 position Tri Block Complete
- # OMC-RP-17N **\$105.00**

Custom Umbilical Lengths / Utility Center

2000 / 2015 Umbilical

- 1 (6 X 4) MM Air "Yellow"
 - 1 (6 X 4) MM Water "Blue"
 - 4 (4 X 2.7mm)
 - Yellow
 - Grey
 - Purple
 - Red
 - 1 (5/8") Drain Line
 - 1 (5/8") Cuspidor Drain 2000 unit
 - 1 (2") Umbilical
 - No Electrical / Low Voltage
- #2015-UMB-FT **\$ 50.00**
 #2015-UMB-10 **\$300.00** 10' length

3500 Umbilical

- 1 (6X4) MM Air "Black"
 - 9 (4 X 2.7mm)
 - Yellow
 - Grey
 - Purple
 - Red
 - Green
 - Blue
 - Maroon
 - Black (spare)
 - Red (spare)
 - No Low Voltage
- #3500-UMB-FT **\$50.00**
 #3500-UMB-10 **\$300.00** 10' length

3550 Umbilical

- 1 (6X4)mm Air "Yellow"
 - 1/2" Spiral Loc Tubing
 - Syringe tubing.
- #3550-UMB **\$130.00** 6' length
 #3550-UMB-10 **\$250.00** 10' length

MC500 / 501 Umbilical

- 1 6X4MM Air "Black "
 - (1 1/2") Umbilical
 - 8 4 X 2.7mm
 - Yellow
 - Grey
 - Red
 - Green
 - Blue
 - Maroon
 - Black (spare)
 - Red (spare)
 - Purple (spare)
 - 1/2" Spiral Lock Suction
- MC501 only Add
 \$2.62 per ft
- Low Voltage
- #MC500-UMB-FT **\$ 50.00**

OMC2375 Umbilical

- 1 6X4MM Air "Yellow"
 - 1 (1 1/2") Umbilical
 - 8 4 X 2.7mm
 - Yellow
 - Grey
 - Clear
 - 5/8" Spiral Lock Suction
 - No Electrical
- #OMC2375-UMB-FT **\$ 50.00**
 #OMC2375-UMB-10 **\$ 300.00**

MS3000 / MR4000

- 1 6X4mm Air "Black"
 - 8 4 X 2.7mm
 - Yellow
 - Grey
 - Red
 - Green
 - Blue
 - Maroon
 - Black (spare)
 - Red (spare)
 - Purple (Spare)
 - No Low Voltage
- #3000-UMB-FT **\$50.00**

Junction Box

- 12.5" X 12.5"
 - Metal frame
 - Plastic top cover
- #53120 **\$138.00**



Master Shut-Off Valve

- Includes air / water valved
 - 3/8" compression tube
 - Includes pressure regulator
 - Includes pressure guage
- #53119 **\$458.00**



Rebulid Kit for 53119

- Diaphragm
 - Spring
 - Plunger
- #53199 **\$30.00**



40 micron filter

- Use for air and water
 - 40 micron filter element
- #53107 **\$3.90**



Raw Suction Tubing per Ft.

Part Number	Description	Price per Ft
53096	1/2 Grey Spiral Lock Tubing	\$3.30 Per Ft
53093	5/8" Clear Serial Lock Tubing	\$3.90 Per Ft
53085	5/8" Grey Serial Lock Tubing	\$3.90 Per Ft
53021	2" Umbilical Tubing	\$7.90 Per Ft
2000-HVETubing	1/2" Smooth HVE tubing	\$3.90 Per Ft
2000-SETubing	Saliva Ejector Tubing	\$1.30 per Ft



851 South Lawson Street, City of Industry CA 91748
 800-560-8222 (PH) 626-810-4245 (FX)

PURCHASE ORDER NO. _____

PURCHASE ORDER

TO: Dealer		SHIP TO:	
		End User	

P.O. DATE	SHIP VIA		F.O.B.	TERMS		
BUYER	FREIGHT	REQ. DATE	CONFIRMING TO	REMARKS	TAX	
QTY. REQ.	UNITS	DESCRIPTION			UNIT COST	EXTENDED COST
		TOTAL				

Credit Card Number: _____ Date: _____

Expiration Date: _____ Code: _____

Authorized by: _____ Date: _____

Terms of Sale and Payment: Purchaser hereby authorizes TPC to ship the above equipment / supplies. The purchaser understands and agrees that the complete payment or this entire order is due prior to shipping order. Cancellatioins are subjected to a 25% handling fee. By siging this form, the customer agrees to authroize TPC to bill the credit card in the amount shown above plus all freight charges.

Standard Upholstery Color Selection



LK-03
Sand

LK-08
Taupe

LK-04
Stormy

LK-11
Greystone

LK-27
Orange



LK-31
Cobalt

LK-34
Desert Rose

LK-49
Grey Teal

LK-59
Atlantis

LK-72
Black

TPC recommends final upholstery color to be selected from actual color swatch chart.



TPC

TPC Advanced Technology, Inc.

851 S. Lawson St. City Of Industry, CA 91748

Tel. (800) 560-8222 Fax (626) 810-4245

tpcdental.com