



Mirage

OPERATORY SYSTEM



MODEL 2015



TPC Mirage Operatory System

Model 4000 Dental Chair

Standard Features

Three operator-programmable position, one preset exit position
 Electronic foot switch preposition / auto return switch
 Fuse / surge protection
 Toe guard protection
 Synchronized seat and toe board tilt
 All metal construction
 Thin back with lumbar support
 Double articulating headrest
 Safety shut-off plate
 Pivoting armrest
 Seamless asepsis upholstery
 Input voltage (110v or 220v)- specify
 5 year warranty

Mirage Chair Mounted Delivery System

Standard Features

Asepsis automatic control for three handpieces with syringe
 Individual handpiece pressure adjustments
 Non-retracting water coolant with on/off
 Individual water coolant flow controls
 Air coolant flow control
 Handpiece tubing flush system
 Ultra soft handpiece tubing
 4 hole straight gray asepsis tubing with Midwest connector
 Stainless steel tray with non-skid pad
 Disc foot control wet/dry toggle with 7" of smooth tubing
 Equipment utility center with asepsis housing
 Vacuum canister with solids collector screen
 Straight light post for unit mount dental light
 Pre-wiring for fiber optic light source
 Air-actuated counterbalanced arm lock

Mirage Post Mount Utility Center

Standard Features

Gravity drain swing away porcelain cuspidor
 Adjustable timed bowl rinse with push-button cup filler
 Self contained water system with city/ bottle selector
 Water outlet with flow control
 Air outlet with flow control
 Vacuum canister with solids collector screen

Assistant's Instrumentation

Standard Features

Telescoping arm with movable assistant's instrument holder
 One autoclavable HVE with asepsis tubing
 One autoclavable saliva ejector with tubing

Dental Operatory Light

L700-1 Advance Operatory Light, 12,000 lux / 24,000 lux, 4,500° Kelvin
L825 Lucent Operatory Light (sensor), 29,000 lux / 17,000 lux, 4,600° Kelvin
L800 Lucent Operatory Light (no sensor), 29,000 lux / 17,000 lux, 4,600° Kelvin

	4000 (Dental Chair)	2000 (Unit w/Cuspidor)	2015 (Unit w/o Cuspidor)	L700 (Standard Operatory Light)	L825 (Sensor Operatory Light)	L800 (No Sensor Light)
Three operator-programmable position, one preset exit position	•					
Electronic foot switch preposition / auto return switch	•					
Fuse / surge protection	•					
Toe guard protection	•					
Synchronized seat and toe board tilt	•					
All metal construction	•					
Thin back with lumbar support	•					
Double articulating headrest	•					
Safety shut-off plate	•					
Pivoting armrest	•					
Seamless asepsis upholstery	•					
Input voltage (110v or 220v)- specify	•					
5 year warranty	•					
Asepsis automatic control for three handpieces with syringe		•	•			
Individual handpiece pressure adjustments		•	•			
Non-retracting water coolant with on/off		•	•			
Individual water coolant flow controls		•	•			
Air coolant flow control		•	•			
Handpiece tubing flush system		•	•			
Ultra soft handpiece tubing		•	•			
4 hole straight gray asepsis tubing with Midwest connector		•	•			
Stainless steel tray with non-skid pad		•	•			
Disc foot control wet/dry toggle with 7" of smooth tubing		•	•			
Equipment utility center with asepsis housing		•	•			
Vacuum canister with solids collector screen		•	•			
Straight light post for unit mount dental light		•	•			
Pre-wiring for fiber optic light source		•	•			
Air-actuated counterbalanced arm lock		•	•			
Gravity drain swing away porcelain cuspidor		•	•			
Adjustable timed bowl rinse with push-button cup filler		•	•			
Self contained water system with city/ bottle selector		•	•			
Water outlet with flow control		•	•			
Air outlet with flow control		•	•			
Vacuum canister with solids collector screen		•	•			
Telescoping arm with movable assistant's instrument holder		•	•			
One autoclavable HVE with asepsis tubing		•	•			
One autoclavable saliva ejector with tubing		•	•			
L700-1 Advance Operatory Light, 12,000 lux / 24,000 lux, 4,500° Kelvin				•		
L825 Lucent Operatory Light (sensor), 29,000 lux / 17,000 lux, 4,600° Kelvin					•	
L800 Lucent Operatory Light (no sensor), 29,000 lux / 17,000 lux, 4,600° Kelvin						•

Factory Options

Additional 3 way syringe
 Additional HVE with tubing
 Power optic system, ISO B (5-hole)
 Power optic system, ISO C (6-pin)
 Built-in piezo scaler with transformer
 Super thin x-ray viewer- unit mounted
 Standard LCD monitor bracket
 Deluxe LCD monitor bracket
 Mirage LCD monitor bracket
 15" standard multi media LCD monitor
 15" Deluxe multi media LCD monitor w/ x-ray viewer
 Plastic utility shelf- (grey, white or beige)-post mounted
 Touch pad-unit head mounted


Item Number

#TPC3WAY
 #2000HVE
 #H7000
 #H7660
 #A800-1
 #XV100/150
 #WIC800
 #WIC825
 #WIC840
 #M015
 #MTFT15
 #US100
 #TP-2000

Packing Dimensions

55 x 33 x 33	350 lbs.
50 x 5 x 22	20 lbs.
51 x 20 x 17	85 lbs.
23 x 20 x 19	36 lbs.

Upholstery Selections:



WIC800

WIC840

L700

M015

L800

A800-1

XV100/150