

TPC

Mirage

Mobile Cart





MC-500 Mirage Mobile Cart
MC-501 Mirage Mobile Cart with Assistant's Package

Standard Features

	MC-501	MC-500
Automatic control for three handpieces	•	•
Master ON/OFF	•	•
Individual handpiece pressure controls	•	•
Individual air and water coolant controls	•	•
Three-way syringe	•	•
Easy-clean work surface, 13 1/2" x 15 1/2"	•	•
Adjustable height mobile frame, 30" - 38"	•	•
Handpiece tubings	•	•
Disc type wet / dry foot control	•	•
Central vacuum assistant's instrumentation Central vacuum solids separator Central vacuum valve with easy-clean tubing Saliva ejector valve with tubing, Hve valve with tubing	•	
Standard utility center, 4 1/2" x 9" x 14" Manual shutoff valves	•	•
Air water filters, regulators	•	•
Umbilical, 7'	•	•
Non-retraction water coolants	•	•
Handpiece flush system	•	•
9" x 13" stainless steel tray with anti-skid pad	•	•

Packing dimensions: 24" x 24" x 36" / 55lbs

Optional Accessories

	MC-501	MC-500
#A800 Advance 800 piezo scaler-built-in	•	•
#ALED5B LED 5B curing light-built-in	•	•
#2000-hvet Additional Hve valve and tubing	•	
#H7000 Fiber optic light source complete-5 hole	•	•
#H7040 Additional fiber optic tubing-5 hole	•	•
#H7660 Fiber optic light source complete-ISO 6 pin	•	•
#H7045 Additional fiber optic tubing-ISO 6 pin	•	•
#HP6000 High speed handpiece-push button	•	•
#H5000 Low speed handpiece complete	•	•
#XV100/150 Super thin x-ray viewer-cart mounted	•	•
#AIC300-888 AdvanceCAM intraoral camera-built-in	•	•

