

# ***TPC***

## **Software Compatibility List**

The following imaging applications have been tested and are compatible with all AIC 5XXX intraoral cameras

- Apteryx X-Ray vision and X-ray Vision lite.
- Dexis Versions 8,9,10,1
- Patterson Eagle soft
- TPC Soft Image
- Professor Suni “Apteryx imaging”
- CliniView Quality Systems INC
- Dentimax “Apteryx Imaging”
- Kodak Imaging.
- Eva
- EZ 2000
- XDR Digital imaging.
- Romexis

We have confirmed that the following applications are not compatible with the Advance Cam intraoral camera systems. The 'freeze' button does not trigger image capture within these software applications.

- Visix Must purchase additional plug-in from the software manufacture.
- Vixwin~~AA~~
- Tiger View
- Sota Imaging (Requires Twain Driver)

If you or your customer are using an application not listed above, please feel free to contact us to discuss potential integration options.

This list will be updated as new information becomes available. Please note that certain software updates may affect camera functionality or integration. If you experience any issues with your camera following an update, we encourage you to reach out so we can assist in resolving the problem.

- System requirements
- Windows XP, Vista and windows 7, 10, 11
- 32 & 64 bit compatible hardware / software
- Minimum 512 MB ram. Only common on windows XP 32 bit OS. 64 Bit Compatible
- USB 2.0

There is currently no support for the advance cam on MAC / Linux OS.

This list is to be used as a reference only. If you have any questions please contact [service@tpcdental.com](mailto:service@tpcdental.com)