

Superb Air

Oil-less Air Compressor

The TPC Superb Air Oil-less Air Compressors meet the high standards of the dental professionals, providing the reliability you need from a brand you can trust. The new motor design delivers a reduced noise with a 5-year warranty. The membrane air drying system efficiently filters out compressed air contaminants and moisture from the compressed air. Purge tank removes residual humidity in the compressor following each cycle.



DC8124

Model	DC8112	DC8122	DC8124	DC8224	DC8226	DC8228
User Capacity	2 users	2 users	3 users	4 users	5 users	7 users
CDFM @ 80psi	4.9	4.9	9.8	9.8	14.7	19.6
# of Heads	1	1	1	2	2	2
Total HP	1	1	2	2	3	4
Voltage	110V	230V	230V	230V	230V	230V
Line Frequency	60Hz	60Hz	60Hz	60Hz	60Hz	60Hz
Amps	14	14	15	15	22	30
Breaker Amps	20	20	20	20	30	40
dBa	72	72	72	73	73	73
Tank Capacity	10 Gallon	10 Gallon	10 Gallon	16 Gallon	23 Gallon	30 Gallon
Size (HxWxD)	29x28x23	29x28x23	29x28x23	9x28x23	31x34x26	31x34x26



DC8228



Features:

- New motor design to reduce noise and increase reliability
- Quiet operation
- Membrane air drying system to remove moisture from compressed air
- Dry air prevents bacteria growth and air if filtered to 0.01 micron keeping 99.9997% of air contaminants out of your air system
- Energy efficient thermally protected motors
- Superb Air Oil-less compressors are available in single or twin heads
- 5-year warranty

Description

Model

Superb Air Compressor, 1 HP, 110V	DC8112
Superb Air Compressor, 1 HP, 220V	DC8122
Superb Air Compressor, 1 x 2 HP, 220V	DC8124
Superb Air Compressor, 1 x 1 HP + 1 x 1 HP, 220V	DC8224
Superb Air Compressor, 1 x 1 HP + 1 x 2 HP, 220V	DC8226
Superb Air Compressor, 2 x 2 HP, 220V	DC8228

Options

Model

Control Panel 3 button Vac, Air, Water	SVCP-3
Control Panel 4 button Vac, Vac, Air, Water	SVCP-4
Air Compressor and Vacuum pump stacking rack	STACK

Shipping Dimensions: L X W X H

DC8112: 25" X 26" X 35"
DC8122: 25" X 26" X 35"
DC8124: 30" X 29" X 37"
DC8224: 30" X 29" X 37"
DC8226: 39" X 29" X 37"
DC8228: 39" X 29" X 37"

SuperbAIR air compressors not to be used for CAD / 3D printers