

— MONITOR / CAMERA COMBO · AIC3IN1W

17" LCD Monitor / Camera Combo, Corded (AIC3IN1W)

17" high-definition monitor with integrated corded intraoral camera. Easy integration on a dental delivery unit with no additional cabling.



01 FEATURES

- ✓ Camera with 1/4" high definition image sensor and six LEDs
- ✓ 17" monitor, 1280 x 1024 resolution
- ✓ Full USB port for easy image transfer
- ✓ Input supports HDMI and AV
- ✓ Playback of MP3 audio, MP4 video and JPG images
- ✓ Medical grade tempered glass front panel

Specifications

02 MODELS

Model	Description
AIC3IN1W	17" LCD monitor with integrated IO camera, corded. Includes MO17DS, AIC5888

03 OPTIONS

Part	Description
AIC500	Advance Cam Barriers, 500 / Box