

— MONITOR / CAMERA COMBO · AIC3IN1WL

17" LCD Monitor / Camera Combo, Wireless (AIC3IN1WL)

17" high-definition monitor with integrated wireless intraoral camera, built-in Wi-Fi transmitter and Type-C fast charging.



01 FEATURES

- ✓ Camera with 1/4" high definition image sensor and six LEDs
- ✓ 17" monitor, 1280 x 1024 resolution
- ✓ Real-time high-definition video input
- ✓ Built-in Wi-Fi transmitter
- ✓ Wireless camera: 1200mAh lithium battery, 120 minutes of use, Type-C fast charging
- ✓ Medical grade tempered glass front panel

Specifications

02 MODELS

Model	Description
AIC3IN1WL	17" LCD monitor with integrated IO camera, wireless. Includes MO17DS, AIC5855

03 OPTIONS

Part	Description
AIC500	Advance Cam Barriers, 500 / Box