

— CHAIR MOUNTED DELIVERY SYSTEM · 2000-2.0 / 2015-2.0

Mirage2.0

Asepsis delivery that mounts directly to the patient chair, with an air-actuated counterbalanced arm and touch pad controls. The 2000-2.0 adds a porcelain cuspidor with a swing-away bowl.

2000-2.0 shown



3^{HP}

Asepsis automatic handpiece control

1,500^{ml}

Self-contained clean water bottle

6-pin

Asepsis tubing, Midwest (5-hole optional)

5^{year}

Limited warranty



Control head. Simple adjustments, ultra-soft tubing.



Porcelain cuspidor. Timed bowl rinse and cup filler.



Assistant's instrumentation. HVE, saliva ejector and touch pad.

Specifications

TPC

01 STANDARD FEATURES

- ✓ Asepsis automatic control for three handpieces
- ✓ Individual handpiece pressure adjustments
- ✓ Non-retracting water coolant with on/off
- ✓ Individual water and air coolant flow controls
- ✓ Handpiece tubing flush system
- ✓ Ultra-soft 6-pin asepsis tubing, Midwest
- ✓ Disc foot control, wet/dry toggle, 7' tubing
- ✓ Utility center with asepsis housing
- ✓ Self-contained water, city/bottle selector
- ✓ Straight light post for unit-mount light
- ✓ Pre-wired for fiber-optic light source
- ✓ Air-actuated counterbalanced arm lock
- ✓ 3-way syringe on unit head
- ✓ Touch pad controls
- ✓ Vacuum canister with solids collector screen
- ✓ Water quick disconnect

02 MODELS

Included	2000-2.0	2015-2.0
Asepsis 3-handpiece control head, flex arm with dual air brake, tri-block design, utility center, light post, 7' umbilical	✓	✓
Rotating 4-position assistant's package, HVE, SE, water QD, solids collector, clean water system	✓	✓
Porcelain cuspidor, sensor bowl rinse and cup filler, swing-away bowl	✓	—

03 OPTIONS

TPC3WAY-PRT	Additional autoclavable 3-way syringe	L550-LED	Lustrous LED light, motion sensor
P2335+TB	Deluxe autoclavable HVE with tubing	BELBRKT	BLMT chair mount bracket (20 / 50)
WIC800	Standard monitor bracket	TPCBRKT	TPC chair mount bracket
WIC840	Mirage monitor bracket, 16" extension	TP-2005	Touch pad: 3 presets, auto return
M017	17" LCD multimedia monitor	TPE-2011	Chair and electric motor touch pad
A800N-1	Advance 800 piezo scaler, built in	53016-L-100	Disposable solids trap, 2½" (100)
H7000	Illumination system, ISO B (5-hole)	2000-C-2.0	Mirage cuspidor
H7660	Illumination system, ISO C (6-pin)	EMB-901	Built-in electric motor system
ALED55BN-A	LED 5W curing light, built in	XV100	Super-thin X-ray viewer
L700-LED	Vivid LED operatory light	XV150-R	X-ray viewer bracket, round arm
L600-LED	Lucent LED operatory light		

04 PACKING DIMENSIONS

Unit	39 × 29 × 23 in · 60 lb	Unit PMU	43 × 25 × 12 in · 40 lb	Cuspidor	26 × 18 × 10 in · 20 lb
-------------	-------------------------	-----------------	-------------------------	-----------------	-------------------------