

— CHAIR MOUNTED DELIVERY SYSTEM · 5000-2.0 / 5015-2.0

Bellagio

Chair-mounted asepsis delivery with a touch pad, piezo scaler, LED curing light and timer built into the control head. The 5000-2.0 adds a porcelain cuspidor with rotating assistant's instrumentation.



5000-2.0 shown

3^{HP}

Asepsis automatic handpiece control

1,500^{ml}

Self-contained clean water bottle

3-in-1

Built-in scaler, curing light and timer

5^{year}


Limited warranty



Control head. Easy-access controls and ultra-soft tubing.



Porcelain cuspidor. Deep contoured bowl, easy to clean.



Optional monitor. Chairside imaging and patient education.

Specifications

TPC

01 STANDARD FEATURES

- ✓ Asepsis automatic control for three handpieces
- ✓ Individual handpiece pressure adjustments
- ✓ Non-retracting water coolant with on/off
- ✓ Individual water and air coolant flow controls
- ✓ Handpiece tubing flush system
- ✓ Ultra-soft 5-hole asepsis tubing, Midwest
- ✓ Disc foot control, wet/dry toggle, 7' tubing
- ✓ Utility center with asepsis housing
- ✓ Vacuum canister with solids collector screen
- ✓ Self-contained water, city/bottle selector
- ✓ Straight light post for unit-mount light
- ✓ Pre-wired for fiber-optic light source
- ✓ Air-actuated counterbalanced arm lock
- ✓ 3-way syringe on unit head
- ✓ Integrated touch pad control
- ✓ Integrated piezo scaler and LED curing light
- ✓ Integrated timer
- ✓ Water quick disconnect

02 MODELS

Included	5000-2.0	5015-2.0
Asepsis 3-handpiece control head, flex arm with dual air brake, tri-block design	✓	✓
Utility center with filters, regulators and master shut-off; light post with monitor cutout; 7' umbilical	✓	✓
Rotating 4-position assistant's package with HVE, SE and touch pad	✓	✓
Water quick disconnect, solids collector and clean water system	✓	✓
Porcelain cuspidor, sensor bowl rinse and cup filler, swing-away bowl	✓	—

03 OPTIONS

WIC800	Standard flat-screen monitor bracket	L550-LED	Lustrous LED operatory light, motion sensor
WIC840	Mirage monitor bracket, 16" extension	BELBRKT	BLMT chair mount bracket (Bel 20 / 50)
M017	17" LCD multimedia monitor	53016-L-100	Disposable solids trap, 2½" (100)
H7000	Illumination system, ISO B (5-hole)	EMB-901	Built-in electric motor system
H7660	Illumination system, ISO C (6-pin)	XV100	Super-thin X-ray viewer
L700-LED	Vivid LED operatory light	XV150-R	X-ray viewer bracket, round arm

04 PACKING DIMENSIONS

Unit	39 × 29 × 23 in · 60 lb	Unit PMU	43 × 25 × 12 in · 40 lb	Cuspidor	26 × 18 × 10 in · 20 lb
-------------	-------------------------	-----------------	-------------------------	-----------------	-------------------------

Options are not factory-installed unless specified on the purchase order. Additional charges apply.