

— LED PIEZO SCALER · A750N / A750LED

Advance 750LED

A compact piezo ultrasonic scaler and endodontic unit with an intense LED light source, thin subgingival tips and a power range for every clinical application.



32 kHz

Operating frequency

5 tips

G1, G2, G4, P4, P5 included

2.1 kg

Compact unit

1 year

Limited warranty

STANDARD FEATURES

- ✓ Intense LED light source (A750LED)
- ✓ No heat build-up at the tip
- ✓ Fine spray eliminates overspray
- ✓ Gentle, little or no patient sensitivity
- ✓ Autoclavable handpiece
- ✓ Endo file holder for root canal enlargement

Specifications

01 MODELS

Included	A750N	A750LED
Piezo ultrasonic scaler and endodontic unit	✓	✓
Tips G1, G2, G4, P4, P5; foot control; tip wrench	✓	✓
Water bottle; endo file holder and files	✓	✓
LED handpiece light source	—	✓

02 SCALER TIPS

G1 Included	G2 Included	G4 Included	P3 Optional, slim	P4 Included	P5 Included
-----------------------	-----------------------	-----------------------	-----------------------------	-----------------------	-----------------------

03 SPECIFICATIONS

Input	115 / 230 V, 50/60 Hz
Frequency	32 kHz
Water pressure	1–4 kg/cm ²
Dimensions	160 × 155 × 70 mm
Weight	2.1 kg
Fuses, 110 V	1 × F1A 250V
Fuses, 230 V	2 × F0.5A 250V
Warranty	1-year (tips excluded)

05 TIPS AND ENDO

A761	Tip #G1
A762	Tip #G2
A767	Tip #G4
A763	Tip #P3 slim
A764	Tip #P4
A766	Tip #P5
A760	Endo files and holder (21 mm #25)

04 REPLACEMENT PARTS

A746	Handpiece with LED
A745	Handpiece, non-LED
A782	HP tubing with LED source
A780	HP tubing, non-LED
A746-RP101	Replacement LED light
A778	Water line 7.5' with QC
A776	Tip wrench
A749	O-rings (5 pack)

06 PACKING DIMENSIONS

Carton	L × W × H in	Weight
Single unit	13 × 10 × 5	5 lb
Case of 8	21 × 14 × 18	42 lb

Indicate EMS or Satellec style tips when ordering.