

— PIEZO ULTRASONIC BUILT-IN SCALER · A800N / A800LED

Advance 800N

A piezo scaler that builds into the delivery unit and runs from the unit foot control. No external water line, no power cord, and a clean operatory with nothing extra in sight.



A800N built in

24-28 kHz

Automatic tuning

5 tips

G1, G2, G4, P4, P5 included

2³/₈ in

Small electronic control box

1 year

Limited warranty



Complete kit. Control box, transformer, tubing, handpiece and tips.



Autoclavable handpiece. LED handpiece on the A800LED.

Specifications

TPC

01 MODELS

Included	A800N	A800LED
Piezo scaler with handpiece, transformer and control box	✓	✓
Tips G1, G2, G4, P4, P5 and tip wrench	✓	✓
Advance 800 LED handpiece	—	✓

02 SPECIFICATIONS

Transformer	24 V, 1 A
Activation	Unit foot control
Control box	2¾" × 1.5" × 1.25"
Absorbed power	15 W max
Scaler	10 W max, 24–28 kHz
Water pressure	30–145 psi
Water use	10–50 cc/min
Warranty	1-year (tips excluded)

04 TIPS AND ENDO

A761	Tip #G1
A762	Tip #G2
A767	Tip #G4
A763	Tip #P3 slim
A764	Tip #P4
A766	Tip #P5
A760	Endo files and holder (21 mm #25)

03 REPLACEMENT PARTS

A746	Handpiece with LED
A745	Handpiece, non-LED
A782	HP tubing with LED source
A780	HP tubing, non-LED
A746-RP101	Replacement LED light
A800N-PCB	Main PCB
A800N-WH-HARNESS	Wire harness
A776	Tip wrench
A749	O-rings (5 pack)

05 PACKING DIMENSIONS

Carton	L × W × H in	Weight
Single unit	9 × 7 × 4	4 lb
Case of 10	14 × 10 × 18	40 lb

Indicate EMS or Sateltec style tips when ordering.

Factory or field installed. Order the A800N-1 as the built-in scaler option on TPC Mirage and Bellagio delivery systems.