

— LED PIEZO SCALER WITH BOTTLE SYSTEM · A850-LED

Advance 850LED

A compact LED piezo scaler and endodontic unit with its own water bottle, so it runs anywhere in the operatory. Gentle on patients, with no heat build-up at the tip.

A850-LED shown



24-28 kHz

Automatic tuning

5 tips

G1, G2, G4, P4, P5 included

10-50 cc/min

Water consumption

1 year

Limited warranty



Simple controls. General, perio and endo power ranges.



LED handpiece. Autoclavable, with fine water spray.



Bottle or external water. Run self-contained or from the water line.

Specifications

01 PACKAGE INCLUDES

- ✓ A850-LED piezo scaler with LED handpiece
- ✓ Foot control
- ✓ Water bottle
- ✓ Tips G1, G2, G4, P4 and P5
- ✓ Tip wrench
- ✓ Endo file holder and files

02 SCALER TIPS

G1 Included	G2 Included	G4 Included	P3 Optional, slim	P4 Included	P5 Included
-----------------------	-----------------------	-----------------------	-----------------------------	-----------------------	-----------------------

03 TECHNICAL DATA

Absorbed power	15 W max
Scaler	10 W max, 24–28 kHz
Tuning	Automatic
Water pressure	30–145 psi
Water use	10–50 cc/min
Approvals	ETL, per UL 2601-1
MDD class	Ila (93/42/EEC)
Warranty	1-year limited

04 REPLACEMENT PARTS

A746	Handpiece with LED
A781	LED handpiece tubing
A746-RP101	Replacement LED light
A775	Foot control
A850-RP-BTL	Water bottle
A778	Water line 7.5' with QC
A776	Tip wrench
A749	O-rings (5 pack)

05 TIPS AND ENDO

A761	Tip #G1
A762	Tip #G2
A767	Tip #G4
A763	Tip #P3 slim
A764	Tip #P4
A766	Tip #P5
A760	Endo files and holder (21 mm #25)

06 PACKING DIMENSIONS

Carton	L × W × H in	Weight
Single unit	18 × 14 × 5	6 lb
Case of 4	18 × 15 × 21	26 lb

Indicate EMS or Satellec style tips when ordering.

Scaling and endo in one unit. General, perio and endo power ranges, with the endo file holder and files included for fast canal enlargement.