

— CURING LIGHT SYSTEM · ALED-39N / ALED-39BN

LED 39N

A slim, corded LED curing pen at 1,000 mW/cm² with three curing modes. Choose the countertop ALED-39N or the ALED-39BN that runs on 24 V direct from the dental unit.

ALED-39N shown



1,000 mW/cm²

Light intensity

430–490

nm

Wavelength range

24v

Direct built-in option
(39BN)

2 year

Limited warranty

STANDARD FEATURES

- ✓ Fast, pulse and step curing modes
- ✓ Countertop stand included
- ✓ Over-temperature indicator
- ✓ Slim, lightweight handpiece
- ✓ Built-in power transformer
- ✓ Class II (IEC 60601-1)

Specifications

01 CURING MODES

Fast

F



Full power for a 20-second cycle, with a beep every 10 seconds.

Pulse

P



Interrupted output to manage heat, with a beep every 10 seconds.

Step

S



Soft start from zero to full power over 10 seconds, then full power for 10 seconds.

02 MODELS

ALED 39N Curing light with stand

ALED 39BN Built-in, 24 V direct

03 SPECIFICATIONS

Wavelength	430–490 nm
Intensity	1,000 mW/cm ²
Class	Class II (IEC 60601-1)
Warranty	2-year limited

04 PACKAGE INCLUDES

- ✓ LED 39 curing light
- ✓ 8 mm light guide
- ✓ Protective eye-shield
- ✓ Stand
- ✓ Built-in power transformer

05 OPTIONS

P109-LG	Light guide for LED 39N
P208	Protective eye-shield
ALED 39N-PS	Replacement power supply

06 PACKING DIMENSIONS

Carton	L × W × H in	Weight
Single unit	13 × 7 × 3	2 lb
Case of 20	24 × 14 × 14	30 lb

Two ways to power it.

Use the ALED 39N on its countertop stand, or choose the ALED 39BN to run on 24 V direct from the dental unit.

