

— CORDLESS CURING LIGHT SYSTEM · ALED-60N

LED 60N

Cordless convenience with corded backup: a 5 W LED at 1,400 mW/cm², an LCD display and a built-in radiometer in the charging base.

ALED-60N shown



1,400 mW/cm²

Light intensity

5-10 sec

Typical cure time

3 modes

Fast, pulse and step

2 year

Limited warranty

STANDARD FEATURES

- ✓ Built-in radiometer
- ✓ Easy-to-read LCD display
- ✓ Quiet fan extends LED life
- ✓ Corded or cordless operation, no down time
- ✓ Visual and audible low-battery alarm
- ✓ Sleep mode conserves battery

Specifications

TPC

01 CURING MODES

Fast

F



Full power for a 20-second cycle, with a beep every 10 seconds.

Pulse

P



Interrupted output to manage heat, with a beep every 10 seconds.

Step

S



Soft start from zero to full power over 10 seconds, then full power for 10 seconds.

02 SPECIFICATIONS

Light source	5 W LED
Wavelength	430–490 nm
Intensity	1,400 mW/cm ²
Battery	Lithium ion
Power supply	AC 100–240 V, 50/60 Hz
Adapter size	10" L × 2" W × 8" H
Warranty	2-year limited

03 PACKAGE INCLUDES

- ✓ LED 60N cordless curing light
- ✓ Charging base
- ✓ Countertop stand with power indicator
- ✓ Built-in power radiometer
- ✓ 8 mm light guide
- ✓ Protective eye-shield
- ✓ AC power supply

04 OPTIONS

P109-LG	Replacement light guide, 8 mm
P208	Replacement protective eye-shield
ALED 60N-BATT	Replacement rechargeable battery
P117	11 mm light guide

05 PACKING DIMENSIONS

Carton	L × W × H in	Weight
Single unit	11 × 8 × 4	3 lb
Case of 10	18 × 12 × 17	24 lb

Cordless or corded.

Work cordless, or keep curing while the battery charges. The built-in radiometer in the base verifies output in seconds.

