

— LED BLEACHING SYSTEM · ALED 2505

Advance 2505

Chairside whitening with three high-intensity blue and UV LEDs, a microprocessor-controlled timer and a compact mobile stand that rolls away for easy storage.

ALED 2505 shown



1,500 mW/cm²

Light intensity

395–490

nm

Broad-spectrum output

3 LEDs

High-power blue and UV

2 year

Limited warranty

STANDARD FEATURES

- ✓ Multi-arch whitening application
- ✓ Digital display with audio feedback
- ✓ Quiet fan operation
- ✓ Three high-intensity blue and UV LEDs
- ✓ Adjustable, microprocessor-controlled timer
- ✓ Compact mobile stand for easy storage

Specifications



01 PACKAGE INCLUDES

- ✓ Advance 2505 LED bleaching light (illuminator)
- ✓ AC power supply and power cord
- ✓ 2-year warranty
- ✓ Mobile stand with adjustable arm
- ✓ User manual

02 SPECIFICATIONS

Spectrum	395–490 nm, broad
Intensity	1,500 mW/cm ²
Light source	3 high-power blue / UV LEDs
Timer	Adjustable, microprocessor
Power	Auto-select 110–240 VAC, 50/60 Hz
Dimensions	60" H × 20" W
Weight	14 lb
Warranty	2-year limited

03 PACKING DIMENSIONS

Carton	L × W × H in	Weight
ALED 2505	21 × 21 × 6	14 lb

Whitening that moves with you. The mobile stand and adjustable arm position the illuminator over the full arch, then roll the system out of the way when you're done.