

— OPERATORY STOOL

# *Mirage* Assistant's Stool

Raised seating for the dental assistant, with an adjustable anatomical body support bar and full-ring foot rest, set up for right- or left-handed use.

AS-1101



**23-30 in**

Seat height range

**18.5 in**

Seat width

**24 in**

Five-leg base diameter

**1 year**

Warranty

## STANDARD FEATURES

- ✓ Adjustable anatomical body support bar
- ✓ Adjustable full-ring foot rest
- ✓ Adjustable for right- or left-handed use
- ✓ Large, firm cushion seat
- ✓ Five-leg, ultra-stable 24" base
- ✓ Durable tile / carpet casters
- ✓ Gas cylinder lift for seat height adjustment
- ✓ 10 upholstery colors

# Specifications

## 01 DIMENSIONS

|                        |                     |
|------------------------|---------------------|
| Backrest adjustment    | 0" – 3.25" vertical |
| Seat width             | 18.5"               |
| Seat height range      | 23" – 30"           |
| Five-leg base diameter | 24"                 |
| Packing dimensions     | 24" × 23" × 21"     |
| Packed weight          | 42 lb               |
| Warranty               | 1 year              |

## 02 UPHOLSTERY



|                                                                                   |                                                                                     |                                                                                     |                                                                                     |                                                                                     |
|-----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|  |  |  |  |  |
| Sand<br>LK-03                                                                     | Stormy<br>LK-04                                                                     | Taupe<br>LK-08                                                                      | Greystone<br>LK-11                                                                  | Orange<br>LK-27                                                                     |
|  |  |  |  |  |
| Cobalt<br>LK-31                                                                   | Desert<br>Rose<br>LK-34                                                             | Grey Teal<br>LK-49                                                                  | Atlantis<br>LK-59                                                                   | Black<br>LK-72                                                                      |

Specify color and product code when ordering. Colors are approximate.

## 03 OPTIONS

|                        |                               |
|------------------------|-------------------------------|
| <b>AS1101-16</b>       | Replacement seat cushion      |
| <b>AS1101-21</b>       | Replacement back cushion      |
| <b>AS1101-18</b>       | Replacement body bar cushion  |
| <b>AS1101-Vinyl-UC</b> | Non-standard vinyl upholstery |
| <b>AS1101-ULTR-UC</b>  | Ultraleather upholstery       |

## 04 COMPARE THE MIRAGE STOOL LINE

|                      |  |  |  |  |  |  |
|----------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|                      | <b>Saddle</b><br><b>SS-5105</b>                                                     | <b>Doctor's</b><br><b>DR-1102</b>                                                   | <b>Solace</b><br><b>DR-7000</b>                                                     | <b>Assistant's</b><br><b>AS-1101</b>                                                 | <b>Hygiene</b><br><b>HY-5107</b>                                                      | <b>Mica Surgeon</b><br><b>DRSA-7100</b>                                               |
| Seat height          | 24.5–31.8"                                                                          | 19.5–24.8"                                                                          | 19.5–24.8"                                                                          | 23–30"                                                                               | 20–25"                                                                                | 20–25"                                                                                |
| Seat width           | 16"                                                                                 | 18"                                                                                 | 18.3"                                                                               | 18.5"                                                                                | 17"                                                                                   | 18" × 15.7"                                                                           |
| Seat tilt            | Yes                                                                                 | Yes                                                                                 | Yes                                                                                 | —                                                                                    | Yes                                                                                   | Yes                                                                                   |
| Armrest              | —                                                                                   | Optional                                                                            | Optional                                                                            | —                                                                                    | —                                                                                     | Standard                                                                              |
| Foot ring & body bar | —                                                                                   | —                                                                                   | —                                                                                   | Yes                                                                                  | —                                                                                     | —                                                                                     |

Options are not factory-installed unless specified on the purchase order. Replacement cushions are Duration vinyl.