

— SONIC AIR SCALER · AS5000

AS 5000

A lightweight, autoclavable sonic scaler that runs on handpiece air, with no control box or extra foot pedal. Available for 2-hole, 4-hole and KaVo-type quick-connect couplers.

AS5000 shown



3 tips

Universal, sickle and perio

2 / 4 hole


Or KaVo-type quick connect

50 ml/min


Water flow

6 month

Limited warranty



4-hole Midwest coupling. [AS5000-M4](#)



2-hole Borden coupling. [AS5000-B2](#)



Quick connect. KaVo-type SK4 sold separately. [AS5000-K](#)




Specifications

01 MODELS

Included	AS5000-M4	AS5000-B2	AS5000-K
AS5000 air sonic scaler handpiece	✓	✓	✓
AS1, AS2 and AS3 tips	✓	✓	✓
A776 tip wrench	✓	✓	✓
6-month warranty	✓	✓	✓

AS5000-M4: 4-hole · AS5000-B2: 2-hole · AS5000-K: KaVo-type quick connect (SK4 adaptor not included).

02 INCLUDED TIPS



03 SPECIFICATIONS

Air pressure	0.25–0.35 MPa
Water flow	50 ml/min
Noise	≤ 70 dB
Sterilization	Autoclavable
Warranty	6-month limited

04 OPTIONS

AS1	Universal tip
AS2	Sickle tip
AS3	Perio tip
A776	Scaler tip wrench only
SK4	KaVo-type quick connector

05 PACKING DIMENSIONS

Carton	L × W × H in	Weight
AS5000 kit	9 × 4 × 2	1 lb

No box. No extra pedal.

The AS5000 connects to standard handpiece tubing and runs on unit air, so scaling is ready at any chair with a 2-hole, 4-hole or KaVo-type coupler.

