

— ADVANCECAM INTRAORAL CAMERA · AIC5888 + AIC5925A

AdvanceCAM Corded Camera & Docking Station System (AIC5888 + AIC5925A)

1080P corded intraoral camera with docking station offering AV, USB, VGA and HDMI output and built-in Wi-Fi image transfer.



01 FEATURES

- ✓ 1080P high definition macro lens
- ✓ Built-in 6 high-brightness white LEDs
- ✓ Docking station with AV, USB, VGA and HDMI image output
- ✓ Standard 16G USB flash drive for storage
- ✓ Built-in Wi-Fi for wireless image transfer to a computer

Specifications

02 MODELS

Model	Description
AIC5888+AIC5925A	Corded Camera & Docking Station System

03 SPECIFICATIONS

- ✓ Camera: 1080P CMOS image sensor
- ✓ Pixels: 1920×1080
- ✓ Frame rate: 30 fps max
- ✓ Light source: 6 white LEDs

04 SYSTEM REQUIREMENTS

- ✓ Intel Pentium 4 1.5 GHz or higher
- ✓ 512 MB RAM or higher
- ✓ Windows 10 or higher
- ✓ USB 2.0 port