

— ADVANCECAM INTRAORAL CAMERA · AIC5899

AdvanceCAM Corded Direct USB Camera System (AIC5899)

High-definition corded USB intraoral camera with six white LEDs. Connects directly to a PC USB 2.0 port.



01 FEATURES

- ✓ High-definition CMOS 720P imager
- ✓ Capture size 1280 x 720
- ✓ 6 white LED light source
- ✓ USB 5VDC powered, 30 fps max
- ✓ Direct connection to a USB 2.0 port

Specifications

02 MODELS

Model	Description
AIC5899	Corded Direct USB Camera System

03 SYSTEM REQUIREMENTS

- ✓ Intel Pentium 4 1.5 GHz or higher
- ✓ 512 MB RAM or higher
- ✓ Windows 10 or higher
- ✓ USB 2.0 port

04 OPTIONS AND ACCESSORIES

Part	Description
AIC500	Advance Cam Barriers, 500 / Box
AIC930	Advance Cam Camera Holder
AIC928-USB	Advance Cam USB Camera Cable