

— ADVANCECAM INTRAORAL CAMERA · AIC5855 + AIC5925B

AdvanceCAM Wireless Camera & Docking Station System (AIC5855 + AIC5925B)

1080P wireless intraoral camera with 5G Wi-Fi transmission, Type-C fast charging and a docking station with AV, USB, VGA and HDMI output.



01 FEATURES

- ✓ 1080P ultra-high definition macro lens
- ✓ Built-in 6 high-brightness white LEDs
- ✓ 5G Wi-Fi transmission for stable, lag-free images
- ✓ Type-C fast charging, large capacity lithium battery
- ✓ Docking station with AV, USB, VGA and HDMI output
- ✓ Standard 16G USB flash drive included

Specifications

02 MODELS

Model	Description
AIC5855+AIC5925B	Wireless Camera & Docking Station System

03 SPECIFICATIONS

- ✓ Camera: 1080P CMOS image sensor
- ✓ Pixels: 1920×1080
- ✓ Frame rate: 30 fps max
- ✓ Battery: 1200 mAh (120 minutes of continuous use)
- ✓ Camera charging: Type-C
- ✓ Docking station power: DC 9V
- ✓ Light source: 6 white LEDs

04 SYSTEM REQUIREMENTS (USB CONNECTION)

- ✓ Windows 10 or higher
- ✓ USB 2.0 ports
- ✓ 4 GB RAM