

— BELLAGIO OPERATORY PACKAGE · BP5000-L700LED-2.0-CB / BP5015-L700LED-2.0-CB

Bellagio

Chair Mounted Operatory System, Cast Aluminum Base

Complete operatory with the 4000-2.0-CB hydraulic chair, 5000-2.0 chair mounted delivery unit and L700-LED Vivid light. Available with cuspidor or with assistant's instrumentation.



Specifications

01 MODELS

Model	Description
BP5000-L700LED-2.0-CB	Bellagio Operatory Package with Cuspidor
BP5015-L700LED-2.0-CB	Bellagio Operatory Package with Assistant's Instrumentation

02 PACKAGE INCLUDES

Part	Description
4000-2.0-CB	Hydraulic Cast Aluminum Base patient chair with wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 3 programmable positions with automatic return and foot control, cast aluminum base
5015-2.0	Asepsis 3 Handpiece automatic control with coolant controls, flex arm with air brake, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air / water filters and regulators, master shut off valves. Includes 3 way syringe and light post with cutout for internal wires and monitor bracket Asepsis telescoping 4 position assistant's package, 3-way syringe, Hve and saliva ejector, water / air quick disconnects, solids collector clean water system, Built-in Piezo Scaler, Built-in LED Curing Light, internal timer, integrated touch pad controls
L700-LED	Vivid LED post mounted light with adjustable intensity, 10,000-35,000 Lux composite cure mode, motion sensor, manual toggle on/off
2000-C-2.0	Porcelain cuspidor with timed bowl rinse and cup filler

03 PACKING DIMENSIONS

Chair: 44 x 29 x 33 340 lbs.

Cushion: 51 x 24 x 10 35 lbs.

Unit: 39 x 29 x 23 60 lbs.

Assist PKG: 43 x 25 x 12 40 lbs.

Dental Light Arm L600: 34 x 11 x 5 30 lbs.

Cuspidor 26 x 18 x 10 20 lbs.

04 OPTIONS AND PARTS

Part	Description
WIC800	Standard flat screen monitor bracket
WIC840	Mirage flat screen monitor bracket with 16" extension Arm
M017	17" LCD multi media flat screen monitor
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply
H7660	Illumination System- ISO C (6 pin) – tubing, control box, power supply
53016-L-100	Disposable Solids Collector Trap, 2 1/8' (pack of 100)
2000-C-2.0	Mirage Cuspidor
EMB-901	Built in Electric Motor System
AS-1101	Mirage Assistants stool
DR-1102	Mirage Doctor's stool

Options are not factory-installed unless specified on the purchase order. Contact TPC for pricing.