

— BELLAGIO OPERATORY PACKAGE · BP5000-L-7601LED-2.0 / BP5015-L-7601LED-2.0

Bellagio

Chair Mounted Operatory System, Steel Base, Luminous LED

Complete operatory with the 4000-2.0 hydraulic steel base chair, 5015-2.0 delivery unit and L-7601-LED Luminous light. Available with cuspidor or with assistant's instrumentation.



01 FEATURES

✓ Built-in Scaler

✓ Built-in Curing light

Specifications

02 MODELS

Model	Description
BP5000-L-7601LED-2.0	Bellagio Operatory Package with Cuspidor
BP5015-L-7601LED-2.0	Bellagio Operatory Package with Assistant's Instrumentation

03 PACKAGE INCLUDES

Part	Description
4000-2.0	Hydraulic patient chair with steel base wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 3 programmable positions with automatic return and foot control
5015-2.0	Asepsis 3 Handpiece automatic control with coolant controls, flex arm with air brake, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves. Includes 3 way syringe and light post with cutout for internal wires and monitor bracket Asepsis telescoping 4 position assistant's package, 3-way syringe, Hve and saliva ejector clean water / air quick disconnects, solids collector with touch pad control
L-7601-LED	Luminous LED post mounted light with adjustable intensity, 3,000-35,000 Lux, manual toggle on/off, optional poly-block filter
2000-C-2.0	Porcelain cuspidor with timed bowl rinse and cup filler

04 STANDARD FEATURES

✓ Built-in Scaler

✓ Built-in Curing light

05 PACKING DIMENSIONS

Chair: 44 x 29 x 33 340 lbs.

Cushion: 51 x 24 x 10 35 lbs.

Unit: 39 x 29 x 23 60 lbs.

Assist PKG: 43 x 25 x 12 40 lbs.

Dental Light Arm L-7601: 34 x 11 x 5 30 lbs.