

— OIL-LESS DESICCANT DRYER AIR COMPRESSOR

Superb Air

Oil-less compressors with a desiccant drying system that removes moisture and contaminants from the air, and a purge tank that clears residual humidity after every cycle. Seven models for 1–9 users.



DC6112 / DC6122
1 head



DC6214 / DC6224
2 heads



DC6326
3 heads



DC6228
2 heads, 30 gal

1–9 users

Capacity across seven models

4.9–19.6

CFM at 80 psi

62 dBA

Quiet operation, 1- and 2-head models

5 year

Limited warranty

STANDARD FEATURES

- ✓ Desiccant air drying system
- ✓ Purge tank removes residual humidity
- ✓ One to four heads, 10–30 gallon tanks
- ✓ New motor design for lower noise
- ✓ Energy-efficient, thermally protected motors
- ✓ 110V and 230V models

Specifications

TPC

01 CHOOSE YOUR MODEL

Model	DC6112	DC6122	DC6214	DC6224	DC6326	DC6228	DC6428
User capacity	1-3	1-3	3-4	5	6-7	8-9	8-9
CFM @ 80 psi	4.9	4.9	9.8	9.8	14.7	19.6	19.6
Heads	1	1	2	2	3	2	4
Total HP	1.15	1.15	2.30	2.30	3.45	4.80	4.60
Voltage	110V	230V	110V	230V	230V	230V	230V
Amps	8	4	16	8	12	16	16
Breaker amps	20	20	30	20	30	30	30
Noise (dBA)	62	62	62	62	72	72	72
Tank (gallons)	10	10	16	16	23	30	30
Size L×W×H (in)	16×23×27	16×23×27	26×21×29	26×21×29	34×22×29	41×22×36	41×22×36
Shipping (in)	29×23×35	29×23×35	34×26×35	34×26×35	39×25×38	46×27×41	46×27×41

Superb Air compressors are not to be used for CAD / 3D printers.

02 CONTROL PANELS & STACKING RACK



CP-2AV
2-button: Vac / Air



CP-3AV
3-button: Vac1 / Vac2 / Air



CP-3
3-button: Vac / Air / Water



CP-4
4-button: Vac1 / Vac2 / Air / Water



CV22
Stacking rack, 33½ × 22 × 26"

03 COMPLETE THE UTILITY ROOM



Superb Vac
WV3 / WV5
Wet-ring vacuum, 1-5 users



Superb Vac Twin
WV6 / WV10
Twin wet-ring vacuum, 6-10 users



Stack them
CV22
Compressor over vacuum on one rack