

— OPERATORY STOOL

Mirage Doctor's Stool

A contoured seat and adjustable backrest on a gas-lift column and five-leg, dual-caster base, with seat tilt and an optional adjustable armrest.

DR-1102**19.5–24.8ⁱⁿ**

Seat height range

18ⁱⁿ

Seat width

24ⁱⁿ

Five-leg base diameter

1^{year}

Warranty

STANDARD FEATURES

- ✓ Five-leg, ultra-stable 24" base
- ✓ Gas cylinder lift for seat height adjustment
- ✓ Large, firm cushion seat
- ✓ Adjustable for right- or left-hand use
- ✓ Dual-wheel casters
- ✓ Seat tilt
- ✓ Optional adjustable armrest
- ✓ 10 upholstery colors

Specifications

01 DIMENSIONS

Backrest adjustment	0" – 3.25" vertical
Seat width	18"
Seat height range	19.5" – 24.8"
Five-leg base diameter	24"
Shipping dimensions	23" × 23" × 17"
Shipping weight	33 lb
Warranty	1 year

02 UPHOLSTERY


				
Sand LK-03	Stormy LK-04	Taupe LK-08	Greystone LK-11	Orange LK-27
				
Cobalt LK-31	Desert Rose LK-34	Grey Teal LK-49	Atlantis LK-59	Black LK-72

Specify color and product code when ordering. Colors are approximate.

03 OPTIONS

DR1102-9	Replacement seat cushion
DR1102-11	Replacement back cushion
DR1102-Vinyl-UC	Non-standard vinyl upholstery
DR-1102-ULTR-UC	Ultraleather upholstery

04 COMPARE THE MIRAGE STOOL LINE

						
	Saddle SS-5105	Doctor's DR-1102	Solace DR-7000	Assistant's AS-1101	Hygiene HY-5107	Mica Surgeon DRSA-7100
Seat height	24.5–31.8"	19.5–24.8"	19.5–24.8"	23–30"	20–25"	20–25"
Seat width	16"	18"	18.3"	18.5"	17"	18" × 15.7"
Seat tilt	Yes	Yes	Yes	—	Yes	Yes
Armrest	—	Optional	Optional	—	—	Standard
Foot ring & body bar	—	—	—	Yes	—	—

Options are not factory-installed unless specified on the purchase order. Replacement cushions are Duration vinyl.