

— OPERATORY STOOL

Laguna Doctor's Stool

A supportive seat and adjustable backrest with seat tilt, on a gas-lift column and 24" five-leg base with tile/carpet casters.

DR-5045

**20–25 in**

Seat height range

18 × 15.7 in

Seat width

24 in

Five-leg base diameter

1 year

Warranty

STANDARD FEATURES

- ✓ Five-leg, ultra-stable 24" base
- ✓ Gas cylinder lift for seat height adjustment
- ✓ Large, firm cushion seat
- ✓ Durable tile / carpet casters
- ✓ Seat tilt
- ✓ 10 Laguna upholstery colors

Specifications

01 DIMENSIONS

Backrest adjustment	0" – 3.25" vertical
Seat width	18" × 15.7"
Seat height range	20" – 25"
Five-leg base diameter	24"
Packing dimensions	24" × 23" × 20"
Packed weight	33 lb
Warranty	1 year

02 UPHOLSTERY

				
Sand LK-03	Stormy LK-04	Taupe LK-08	Greystone LK-11	Orange LK-27
				
Cobalt LK-31	Desert Rose LK-34	Grey Teal LK-49	Atlantis LK-59	Black LK-72

Specify color and product code when ordering. Colors are approximate.

03 COMPARE THE LAGUNA STOOL LINE



Assistant's Stool
AS-5047

Seat height	23-30"
Seat width	18.5"
Backrest	—
Seat tilt	—
Foot ring & body bar	Yes



Doctor's Stool
DR-5045

Seat height	20-25"
Seat width	18" × 15.7"
Backrest	0-3.25" adjustment
Seat tilt	Yes
Foot ring & body bar	—