

— OPERATORY STOOL

Mirage Mica Surgeon Stool

A saddle seat with adjustable backrest, seat tilt and adjustable armrest, on a gas-lift column and 24" five-leg base. Upholstered in black.

DRSA-7100-DU100

**20–25 in**

Seat height range

18 × 15.7 in

Seat width

24 in

Five-leg base diameter

1 year

Warranty

STANDARD FEATURES

- ✓ Adjustable armrest
- ✓ Seat tilt
- ✓ Large, firm cushion seat
- ✓ Black DU-100 upholstery
- ✓ Five-leg, ultra-stable 24" base
- ✓ Durable tile / carpet casters
- ✓ Gas cylinder lift for seat height adjustment

Specifications



01 DIMENSIONS

Backrest adjustment	0" – 3.25" vertical
Seat width	18" × 15.7"
Seat height range	20" – 25"
Five-leg base diameter	24"
Packing dimensions	24" × 23" × 22"
Packed weight	45 lb
Warranty	1 year


02 UPHOLSTERY



Black DU-100

The Mica surgeon stool is available in black only.

03 COMPARE THE MIRAGE STOOL LINE

						
	Saddle SS-5105	Doctor's DR-1102	Solace DR-7000	Assistant's AS-1101	Hygiene HY-5107	Mica Surgeon DRSA-7100
Seat height	24.5–31.8"	19.5–24.8"	19.5–24.8"	23–30"	20–25"	20–25"
Seat width	16"	18"	18.3"	18.5"	17"	18" × 15.7"
Seat tilt	Yes	Yes	Yes	—	Yes	Yes
Armrest	—	Optional	Optional	—	—	Standard
Foot ring & body bar	—	—	—	Yes	—	—