

— LED OPERATORY LIGHT

Lucent LED

A 3,000–35,000 lux LED light with dimmer mode, offered as a post mount or ceiling mount, each with or without a motion sensor.

L600-LED shown



3K–35K lux

Adjustable intensity at 700 mm

100×180

Illumination area (mm) at 700 mm

8 VA

Maximum absorbed power

2 year

Warranty

STANDARD FEATURES

- ✓ Dimmer mode, 3,000–35,000 lux
- ✓ Motion sensor models available
- ✓ Replaceable LED bulb group, handles and shield
- ✓ Static heat sink cooling, no fan
- ✓ 24 V AC/DC operation
- ✓ Post mount or ceiling mount

Specifications

TPC

01 SPECIFICATIONS

Light intensity	3,000 – 35,000 lux @ 700 mm
Illumination area	100 × 180 mm @ 700 mm
Absorbed power	Max 8 VA
Input voltage	110V / 220V, 60 / 50 Hz
Output voltage	24 V AC / DC
Cooling	Static heat sink, no fan
Dimmer mode	Standard on all models
Motion sensor	L625-LED and L695-LED
Warranty	2 years

03 OPTIONS

RLL-H04	Replacement LED optical bulb group
RLL-H06	Replacement handles, left and right
RLL-H01	Replacement shield
RLL-H03	Replacement reflectors, set of 2
RLL-PS	Power supply, 110–240 V in / 24 V out
LP-48	48" light post
L690-08	Ceiling post, 8' ceiling
L690-09	Ceiling post, 9' ceiling
L690-010	Ceiling post, 10' ceiling
L690-011	Ceiling post, 11' (special order)

05 COMPARE THE TPC LED LIGHT LINE

					
	Vivid L700-LED	Lustrous L550-LED	Lucent L600-LED	Luminous L-7601-LED	Radiant R6101-LED
Intensity (lux)	3,000–35,000	8,000–30,000	3,000–35,000	3,000–35,000	3,000–35,000
Motion sensor	Standard	Standard	L625 / L695	L-7605	R6105
PoliBlock filter	Yes	Yes	—	Optional	—
Ceiling mount	L795-LED	L570-LED	L690 / L695	L-7901-LED	—

02 MODELS

	
Post mount L600-LED L625-LED motion sensor	Ceiling mount L690-LED L695-LED motion sensor

Specify 8', 9' or 10' ceiling for ceiling models. The LP-48 light post must be ordered if the light is ordered separately.

04 PACKING DIMENSIONS

Carton	L × W × H in	Weight
L600 post-mount arm	34 × 11 × 5	10 lb
L690 ceiling arm	32 × 18 × 6	18 lb
Light head	14 × 9 × 14	7 lb
L900-08 post, 8' ceiling	15 × 3 × 3	3 lb
L900-09 post, 9' ceiling	30 × 3 × 3	5 lb
L900-10 post, 10' ceiling	44 × 4 × 4	7 lb