

— LCD MONITOR BRACKETS · WIC800 / WIC840 / WIC855

LCD Monitor Brackets (WIC800 / WIC840 / WIC855)

Standard and Mirage monitor brackets with 4-way pivot adjustment and standard 75mm spacing.



01 FEATURES

- ✓ Easily position the monitor angle for comfort
- ✓ 4-way pivot adjustment
- ✓ Standard 75mm spacing
- ✓ Mirage bracket: 16" extension arm with horizontal travel
- ✓ 1 year warranty

Specifications

02 MODELS

Model	Description
WIC800	Standard LCD Monitor Bracket, fits 45-50mm post
WIC840	Mirage LCD Monitor Bracket with 16" Extension Arm, fits 45-50mm post
WIC855	Mirage LCD Monitor Bracket with 16" Extension Arm, fits 55-60mm post
M100MM	Expansion Plate for 100 mm spacing pattern
WIC875	AdvanceCam Intraoral Camera & Docking Station Monitor Adapter

03 CLAMP SIZES

- ✓ WIC800 / WIC840: 1.77" - 1.96" (45mm - 50mm)
- ✓ WIC855: 2.16" - 2.36" (55mm - 60mm)

04 PACKING DIMENSIONS

Standard bracket: 11" x 7" x 6" / 6 lbs.

Mirage bracket: 15" x 10" x 5" / 9 lbs.