

— CURING LIGHT METER · LM-300

LM300

A pocket-size digital meter that confirms your curing light is delivering the output your restorations need. Hold the light guide to the sensor and read the result.

LM-300 shown



LCD

Digital intensity readout

Battery

Powered, battery included

1lb

Shipping weight

1 year

Limited warranty



Quick check. Place the light guide over the sensor and activate the light.



Works with TPC lights. Verify output on every LED model we offer.

Specifications

01 PACKAGE INCLUDES

- ✓ LED light meter
- ✓ Battery
- ✓ 1-year warranty

02 CHECK AGAINST RATED OUTPUT

Curing light	Model	Rated intensity
LED 70N	ALED-70	2,000 mW/cm²
LED 60N	ALED-60N	1,400 mW/cm²
LED 55N	ALED55N	1,250 mW/cm²
LED 55BN	ALED55BN	1,250 mW/cm²
LED 39N	ALED 39N	1,000 mW/cm²
LED 50N	ALED-50N	800–1,000 mW/cm²

03 WHY MEASURE

- ✓ Light guides chip, scratch and collect composite over time
- ✓ LED output can drift as lights age
- ✓ Undercured restorations risk sensitivity and failure
- ✓ A quick reading documents your equipment check

04 PACKING DIMENSIONS

Carton	L × W × H in	Weight
LM-300	8 × 4 × 2	1 lb

Make it part of the morning routine. A ten-second reading before the first patient tells you every curing light in the office is ready to work.