

— LAGUNA 2.0 OPERATORY PACKAGE · LP2100-550LED-2.0 / LP2115-550LED-2.0

Laguna 2.0

Chair Mounted Operatory System

Complete operatory with the L2000-2.0 electromechanical chair, 2015-2.0 delivery unit and L550-LED Lustrous light. Available with cuspidor or with assistant's instrumentation.



Specifications

01 MODELS

Model	Description
LP2100-550LED-2.0	Laguna 2.0 Operatory Package with Cuspidor
LP2115-550LED-2.0	Laguna 2.0 Operatory Package with Assistant's Instrumentation

02 PACKAGE INCLUDES

Part	Description
L2000-2.0	Electromechanical patient chair with wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 3 programmable positions with automatic return and foot control
2015-2.0	Asepsis 3 Handpiece automatic control with coolant controls, flex arm with air brake, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves. Includes 3 way syringe and light post with cutout for internal wires and monitor bracket Asepsis rotating 4 position assistant's package, Hve and saliva ejector, water disconnect solids collector and clean water system. With Touch pad control. on both Dr. & Assistant instrument
L550-LED	Lustrous LED post mounted light with adjustable intensity, 8,000-30,000 Lux
2000-C-2.0	Porcelain cuspidor with timed bowl rinse and cup filler

03 PACKING DIMENSIONS

Chair: 44 x 29 x 33 340 lbs.

Cushion: 51 x 24 x 10 35 lbs.

Unit: 39 x 29 x 23 60 lbs.

Unit: 43 x 25 x 12 40 lbs.

Dental Light: 34 x 11 x 5 Light 30 lbs.

Cuspidor 20 x 13 x 14 20 lbs.

04 OPTIONS AND PARTS

Part	Description
TPC3WAY-PRT	Additional autoclavable 3 way syringe with tubing
P2335+TB	Deluxe autoclavable Hve Long valve with tubing
WIC800	Standard flat screen monitor bracket
WIC840	Mirage flat screen monitor bracket with 16" extension Arm
M017	17" LCD multi media flat screen monitor
A800N-1	Advance 800 piezo scaler – built in
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply
H7660	Illumination System- ISO C (6 pin) – tubing, control box, power supply
ALED55BN-A	LED 5W curing light – built in to dental unit
53016-L-100	Disposable Solids Collector Trap, 2 1/8' (pack of 100)
TPE-2011	Touch pad Control - Chair & Electric Motor System touch pad
EMB-901	Built in Electric Motor System
XV100	Super thin X-ray film viewer
XV150-R	X-ray film viewer bracket (curved) for TPC Mirage unit only
AS-5047	Laguna Assistant's stool
DR-5045	laguna Doctor's stool

Options are not factory-installed unless specified on the purchase order. Contact TPC for pricing.