

Laguna 2.0

Electromechanical Patient Chair, L3000-2.0

The Laguna 2.0 electromechanical patient chair with two programmable positions, a preset exit position and an electronic foot control with auto return. All-metal construction and seamless upholstery.



01 FEATURES

- ✓ Two programmable positions plus one preset exit position
- ✓ Electronic foot control with auto return function
- ✓ Built-in fuse / surge protection
- ✓ Toe guard protection to prevent scuffing
- ✓ All-metal construction
- ✓ Double articulating headrest

Specifications

02 SPECIFICATIONS

Model	L3000-2.0
Chair type	Electromechanical patient chair
Positions	Two programmable positions, one preset exit position
Foot control	Electronic, with auto return function
Headrest	Double articulating, adjustable for height and wheelchair positioning

Model L3000-2.0

Input voltage	110V or 220V
Packing: chair (L x W x H, in)	44 x 29 x 33, 340 lbs
Packing: cushion (L x W x H, in)	51 x 24 x 10, 35 lbs

03 PRODUCT INCLUDES

- ✓ Two programmable positions
- ✓ One preset exit position
- ✓ Electronic foot control with auto return function
- ✓ Built-in fuse / surge protection
- ✓ Toe guard protection to prevent scuffing
- ✓ All-metal construction
- ✓ Double articulating headrest, adjustable for height and wheelchair positioning
- ✓ Seamless upholstery
- ✓ Input voltage 110V or 220V

04 OPTIONS

TP-2005	Touch pad control, chair mounted, 3 preset, auto return, last position
PC-5000-DU	Upcharge: dental chair upholstery, wide back, LK vinyl (specify color and product code)
PC-5000-Vinyl-UC	Upcharge: non-standard dental chair upholstery, wide back, vinyl
PC-5000-ULTR-UC	Upcharge: Ultra leather dental chair upholstery, wide back
LAG.2.0-BRKT	Laguna L2000-2.0 post mount chair bracket

Options are not factory-installed unless specified on the purchase order. Additional charges apply. Contact TPC for pricing.