

Laguna 2.0

Electromechanical Patient Chair

An electromechanical patient chair on a solid steel base. Hands-free foot control with two programmable positions plus an entry/exit position, a wide backrest and pivoting armrests.



01 FEATURES

- ✓ Independent backrest
- ✓ Articulating dual axis headrest
- ✓ Solid steel base
- ✓ Programmable foot control

Specifications

02 SPECIFICATIONS

Model	L2000-2.0
Lifting / lowering	Electromechanical, cantilever type
Initial height / stroke	15.5" / 15"
Backrest reclining angle	0° (horizontal) to 70°
Seat tilting	10° to 20°, synchronized with backrest
Preset positions	Microcomputer control; three preset positions, auto return, last position

Model L2000-2.0

Chair rotation	Fixed
Headrest	Dual axis type
Safety shut-off plate	Under the main link
Emergency stop	Any switch can cancel chair movement
Input voltage	110V or 220V
Chair size (L x W x H, in)	44 x 29 x 33
Warranty	5 Year Limited

03 PRODUCT INCLUDES

- ✓ Two programmable positions
- ✓ One preset exit position
- ✓ Electronic foot control with auto return function
- ✓ Built-in fuse / surge protection
- ✓ Toe guard protection to prevent scuffing
- ✓ All-metal construction
- ✓ Double articulating headrest, adjustable for height and wheelchair positioning
- ✓ Seamless upholstery

04 OPTIONS

TP-2005	Touch pad control, chair mounted, 3 preset, auto return, last position
PC-5000-DU	Upcharge: dental chair upholstery, wide back, LK vinyl (specify color and product code)
PC-5000-Vinyl-UC	Upcharge: non-standard dental chair upholstery, wide back, vinyl
PC-5000-ULTR-UC	Upcharge: Ultra leather dental chair upholstery, wide back
LAG.2.0-BRKT	Laguna L2000-2.0 post mount chair bracket

Options are not factory-installed unless specified on the purchase order. Additional charges apply. Contact TPC for pricing.