

— LAGUNA 2.0 HYGIENE PACKAGE · LHP1001-2.0

Laguna 2.0

Hygiene Operator System

Hygiene operator with the L2000-2.0 electromechanical chair and H1001 rear radius delivery unit.



Specifications

01 MODELS

Model	Description
LHP1001-2.0	Laguna 2.0 Hygiene Operator System (without light)

02 PACKAGE INCLUDES

Part	Description
L2000-2.0	Electromechanical patient chair with wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 4 programmable positions with automatic return and foot control Rear radius asepsis 3 Handpiece automatic control with coolant
H1001	controls, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves, Unit Head Height Adjustment from 27" to 39" Includes 3 way syringe and clean water system, with touch pad control. Hve and saliva ejector, and solids collector, piezo scaler with light

03 OPTIONS AND PARTS

Part	Description
LP-100	Mirage swing mount light post (ordered w/unit)
LP-105	Mirage swing mount monitor post with monitor mounting bracket (ordered w/unit)
L570-LED	Lustrous LED Ceiling Mounted Operatory Light
L550-LED	Lustrous LED Operatory Light
LP48	Lustrous LED Operatory Light
WIC800	Standard flat screen monitor bracket
WIC840	Mirage flat screen monitor bracket with 16" extension
M017	17" LCD multi media flat screen monitor
A800N-1	Advance 800 piezo scaler – built in
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply
ALED55BN-A	LED 5W curing light – built in to dental unit TPCRDSBRKT TPC Radius Bracket
TPE-2011	Touch pad Control - Chair & Electric Motor System Touch pad
EMB-901	Built in Electric Motor System
AS-5047	Laguna Assistants stool
DR5045	Laguna Doctor's Stool
M4001	Mirage chair Bracket
2000-PMU	PMU with cuspidor & cup filler

Options are not factory-installed unless specified on the purchase order. Contact TPC for pricing.