

— LAGUNA 2.0 ORTHODONTIC PACKAGE · LOP2000-550LED-2375 / CV / CV-SL

Laguna 2.0

Orthodontic Package

Orthodontic package with the L2000-2.0 electromechanical chair, OMC2375 mobile delivery cabinet and L550-LED Lustrous light, available centrally plumbed or self-contained.



01 FEATURES

- ✓ L2000-2.0 electromechanical chair, solid steel base
- ✓ OMC2375 mobile delivery cabinet
- ✓ L550-LED Lustrous operatory light
- ✓ Central plumbing or self contained (compressor + vacuum)
- ✓ 1 year warranty

Specifications

02 MODELS

Model	Description
LOP2000-550LED-2375	Laguna 2.0 Orthodontic Package, OMC2375 cabinet
LOP2000-550LED-2375CV	Laguna 2.0 Orthodontic Package, OMC2375CV: cabinet + compressor + vacuum (self contained)
LOP2000-550LED-2375CV-SL	Laguna 2.0 Orthodontic Package, OMC2375CV-SL: cabinet + compressor + vacuum + LED curing light + scaler (self contained)

03 PACKAGE INCLUDES

Part	Description
L2000-2.0	Electromechanical patient chair with wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 4 programmable positions with automatic return and foot control
OMC-2375	Orthodontic mobile delivery cabinet with a pullout writing surface, two instrument drawers, one medium and two large storage drawers that fully extend and quickly remove for easy cleaning, and a large open space for accessories. Oversized locking dual-casters, central plumbing with junction box.
L550-LED	Lustrous LED post mounted light with dual intensity, 8,000 – 30,000 Lux, dual axis light head rotation
L2001-LP48	The Laguna chair bracket and LP48" light post are included at no charge when purchasing the Orthodontic Package

04 CHAIR FEATURES

- ✓ Independent backrest
- ✓ Articulating dual axis headrest
- ✓ Solid steel base
- ✓ Programmable foot control
- ✓ Built in fuse / surge protection
- ✓ Electric motor operated on chair base and backrest

05 CABINET FEATURES

- ✓ Five easy to remove and clean full extension drawers and one pull-out writing surface
- ✓ High quality hardware throughout
- ✓ Oversized caster wheels with locks
- ✓ A removable counter-top plastic tray
- ✓ Standard 12" x 12" junction box
- ✓ 7' length of umbilical

06 LIGHT FEATURES

- ✓ Composite cure mode
- ✓ 24V DC power supply
- ✓ Bulb type: 4 white LEDs, 2 yellow
- ✓ Light intensity (Lux): high 30,000, low 8,000
- ✓ Focal length: 27" (70 cm)
- ✓ Illumination area: 100mm x 180mm
- ✓ Sensor mode, manual toggle on/off

07 CONTROL FEATURES

- ✓ 3 HP automatic selection
- ✓ 3-way syringe
- ✓ Vacuum package with saliva ejector and HVE
- ✓ Self contained water bottle system
- ✓ Variable speed disc-type foot control
- ✓ Water coolant on/off switch with flow control
- ✓ 1 year warranty
- ✓ Requires central air / vacuum systems. Includes utility center and junction box.
- ✓ Optional: self contained compressor and vacuum (no plumbing required)

08 PACKING DIMENSIONS

Chair: 44 x 29 x 33 / 340 lbs.

Cushion: 51 x 24 x 10 / 35 lbs.

Light arm L600: 26 x 14 x 10 / 30 lbs.; light head L550: 17 x 15 x 7

Cabinet: 32 x 28 x 41 / 140 lbs.

09 OPTIONS AND PARTS

Part	Description
WIC840	Mirage flat screen monitor bracket with 16" extension
M017	17" LCD multi media flat screen monitor
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply
OMC2200	LCD monitor bracket
OMC-ST81	Optional Stone Surface on OMC 2375, Color Artic
OMC-ST82	Optional Stone Surface on OMC 2375, Color White
OMC-ST85	Optional Stone Surface on OMC 2375, Color Cream
OMC-ST86	Optional Stone Surface on OMC 2375, Color Bamboo
L2001	Laguna L2000-2.0 Chair Bracket
AS-5047	Laguna Assistants stool
DR-5045	Laguna Doctor's stool
PTH-101	Pony Tail Headrest

Options are not factory-installed unless specified on the purchase order. Contact TPC for pricing.