

— MOBILE DOCTOR'S CART · MC-500-2.0 / MC-501-2.0

Mirage 2.0

Asepsis delivery on a painted H-frame cart that adjusts from 27" to 39" high. Self-contained water, 10 feet of umbilical, and an optional on-board vacuum package on the MC-501-2.0.



MC-501-2.0

With vacuum package: HVE and saliva ejector



MC-500-2.0

Without vacuum package

27–39 in

Height-adjustable H-frame

3 HP

Asepsis automatic handpiece control

10 ft


Umbilical, more available by the foot

5 year

Limited warranty



Doctor's cart only. Without vacuum package, [MC-500-2.0](#)



Compact footprint. 19½" top, height adjusts 27"–39".

Specifications

TPC

01 PACKAGE INCLUDES

Included	MC-500-2.0	MC-501-2.0
Asepsis 3-handpiece automatic selection with bottle system	✓	✓
3-way syringe with quick-change autoclavable tip	✓	✓
Master on/off toggle	✓	✓
Air/water pressure adjustment for each handpiece	✓	✓
Pressure gauge	✓	✓
Non-retracting water/air coolant adjustment and purge	✓	✓
Asepsis handpiece tubing	✓	✓
Painted H-frame, height adjusts 27" to 39"	✓	✓
Vacuum canister with removable screen, 2 $\frac{1}{8}$ "	–	✓
1 autoclavable HVE and saliva ejector with tubing	–	✓
Variable-speed foot control	✓	✓
10 feet of umbilical	✓	✓

02 OPTIONS

H7000	Illumination system, ISO B (5-hole)
H7660	Illumination system, ISO C (6-pin)
A800N-1	Advance 800 N piezo scaler, built in
53016-L-100	Disposable solids trap, 2 $\frac{1}{8}$ " (100)
MC500-UMB-FT	Additional umbilical, per foot
MC501-UMB-FT	Additional umbilical, per foot

03 PACKING DIMENSIONS

Carton	L × W × H in	Weight
MC-500-2.0	29 × 23 × 21	50 lb
MC-501-2.0	29 × 23 × 21	50 lb

Roll it where you need it. The MC-501-2.0 brings its own HVE and saliva ejector, so one cart can serve any operator, from hygiene rooms to surgical suites.