

— MIRAGE MOBILE DELIVERY CART

MC3600CV

A self-contained mobile delivery cart for treatment wherever central plumbing isn't practical. Dual oil-less compressors, built-in suction and clean and waste water bottles, all in one rolling metal cabinet.



MC3600CV



MC3600CV-SL

2

Handpiece positions, plus syringe, SE and HVE

115 psi

Max air pressure, 40 L/min

68 dB


Sound level

1 year

Limited warranty



Clean water bottle. Self-contained supply, no plumbing.



Waste bottle. 750 ml canister, simple to empty.



Optional monitor. 17" LCD with AIC5888 camera.

Specifications

TPC

01 PACKAGE INCLUDES

Included	MC3600CV	MC3600CV-SL
3-way syringe, saliva ejector and HVE	✓	✓
2 handpiece positions	✓	✓
Variable-speed disc-type foot control	✓	✓
Self-contained dual compressor	✓	✓
Clean water bottle	✓	✓
Waste water bottle	✓	✓
Advance 800 built-in piezo scaler	–	✓
LED 39 built-in curing light	–	✓
Dual 110V outlets, rear of cart	✓	✓
1-year warranty	✓	✓

02 CART SPECIFICATIONS

Electrical	(2) 550W heads, 115 VAC, 50/60 Hz, 6.0 A
Compressor	Oil-less, self-contained
Max pressure	115 psi / 40 L per min
Vacuum	80" of water (6" Hg)
Waste canister	750 ml
Sound	68 dB
Construction	Metal case
Dimensions	31" H × 20" W × 16.5" D
Weight	96 lb
Shipping	Freight only

03 OPTIONS

OMC-2260	LCD monitor bracket with post
AIC3N1N1W	17" LCD monitor with AIC5888 intraoral camera
US100	Utility tray

04 PACKING DIMENSIONS

Carton	L × W × H in	Weight
Cart	39 × 27 × 23	160 lb

Moves with ease. No plumbing required. Efficient, compact and simple to maintain: ideal for satellite rooms, temporary locations or any space where connecting to central utilities isn't feasible.