

## — WALL / CABINET MOUNTED SYSTEM

# Mirage 2.0

Mirage 2.0 delivery that mounts to the wall or cabinetry: a side delivery unit on a flex arm, a compact rear delivery unit, and a matching vacuum package for the assistant.




**MS3000-2.0**

Side delivery unit



**MR4000-2.0**

Rear delivery unit



**MVP5000-2.0**

Vacuum package

**3 models**

Side, rear and vacuum package

**3 HP**


Asepsis automatic handpiece control

**9 in**


Rigid arm with pneumatic brake (MS3000)

**5 year**


Limited warranty



**MS3000-2.0.** 16 1/4" wide, 18" deep head.



**MR4000-2.0.** About 20" wide under the counter.



**MVP5000-2.0.** 26 1/2" pivoting arm reach.

# Specifications

**TPC**

## 01 PACKAGE INCLUDES

| Included                                                   | MS3000 | MR4000 | MVP5000 |
|------------------------------------------------------------|--------|--------|---------|
| Asepsis 3-handpiece automatic selection with bottle system | ✓      | ✓      | —       |
| 3-way syringe with quick-change autoclavable tip           | ✓      | ✓      | ✓       |
| Master on/off toggle                                       | ✓      | ✓      | —       |
| Air/water pressure adjustment for each handpiece           | ✓      | ✓      | —       |
| Pressure gauge                                             | ✓      | ✓      | —       |
| Non-retracting water/air coolant adjustment and purge      | ✓      | ✓      | —       |
| Asepsis handpiece tubing                                   | ✓      | ✓      | —       |
| Asepsis flex arm, 9" rigid arm, pneumatic brake            | ✓      | —      | —       |
| 2" diameter wall bracket                                   | ✓      | —      | —       |
| Pivoting arm-mounted vacuum canister with screen           | —      | —      | ✓       |
| Water quick disconnect                                     | —      | —      | ✓       |
| Vacuum canister with removable screen                      | —      | —      | ✓       |
| 4-position holder, autoclavable HVE and SE with tubing     | —      | —      | ✓       |
| Variable-speed foot control                                | ✓      | ✓      | ✓       |

## 02 OPTIONS

|                    |                                     |
|--------------------|-------------------------------------|
| <b>H7660</b>       | Illumination system, ISO C (6-pin)  |
| <b>H7000</b>       | Illumination system, ISO B (5-hole) |
| <b>A800N</b>       | Advance 800 piezo scaler, built in  |
| <b>P2335+TB</b>    | Deluxe autoclavable HVE with tubing |
| <b>53016-L-100</b> | Disposable solids trap, 2½" (100)   |

## 03 PACKING DIMENSIONS

| Carton      | L × W × H in | Weight |
|-------------|--------------|--------|
| MS3000-2.0  | 17 × 21 × 51 | 60 lb  |
| MR4000-2.0  | 19 × 15 × 20 | 25 lb  |
| MVP5000-2.0 | 13 × 24 × 11 | 15 lb  |

**Mix and match.** Pair a side or rear delivery unit with the MVP5000-2.0 vacuum package for a complete wall- or cabinet-mounted operator.