

— SWING MOUNT DELIVERY SYSTEM · MSP2400-2.0

Mirage2.0 – BLMT

A Mirage 2.0 swing-mount delivery unit and rear-radius assistant's arm built to fit BLMT 7, 10, 20, 20N and 50 patient chairs. Upgrade the operatory without replacing the chair.



5 chairs

Fits BLMT 7, 10, 20, 20N and 50

3 HP

Asepsis automatic handpiece control

2 arms

Swing-mount unit and rear-radius assistant's arm

5 year

Limited warranty



Optional light post. LP-100 swing-mount post with an LED operatory light.



Mirage 2.0 control head. Asepsis 3-handpiece control, clean water system.

Specifications

01 PACKAGE INCLUDES · MSP2400-2.0

2400-2.0	Rotating asepsis 3-handpiece automatic control with coolant controls, flex arm with air brake, diaphragm block design, wet/dry variable-speed foot control, deluxe utility center with air/water filters and regulators, master shut-off valves, 3-way syringe and clean water system. Fits 7, 10, 20, 20N, 50.
2450-2.0	Rear-radius asepsis assistant's arm system with HVE and saliva ejector, water quick disconnect and solids collector. Specify chair (7, 10, 20, 20N or 50) when ordering.

02 OPTIONS

2400-2.0	Swing mount delivery system	L550-LED	Lustrous LED operatory light
2450-2.0-0	Rear radius vacuum pkg (20N, 50)	P2335+TB	Deluxe autoclavable HVE, long
2450-2.0-U	Rear radius vacuum pkg (7, 10, 20)	WIC800	Standard monitor bracket
2450-BRKT-BEL-N	Rear radius bracket only (20N, 50)	WIC840	Monitor bracket, 16" extension
2450-BRKT-BEL	Rear radius bracket only (7, 10, 20)	M017	17" LCD multimedia monitor
BELBRKT	BLMT chair mount bracket	A800N-1	Advance 800 piezo scaler, built in
2400-BRKT-BEL	Radius chair bracket only (all)	H7000	Illumination system, ISO B
BEL20N-FC	BLMT foot control (20N, 50)	H7660	Illumination system, ISO C
BEL20-FC	BLMT foot control (7, 10, 20)	ALED55BN-A	LED 5W curing light, built in
LP-100	Swing mount light post*	XV100	Super-thin X-ray film viewer
LP-105	Swing mount monitor post*	XV150-R	X-ray viewer bracket, Mirage
L600-LED	Lucent LED operatory light	AS-1101	Mirage assistant's stool
L690-LED	Lucent LED ceiling light	DR-1102	Mirage doctor's stool

03 PACKING DIMENSIONS

Carton	L × W × H in	Weight
2400-2.0 unit	36 × 28 × 39	100 lb
2450-2.0 arm	29 × 12 × 24	25 lb
LP-100 light post	45 × 26 × 44	30 lb

*LP-100 and LP-105 must be ordered with the unit. Options are not factory-installed unless specified on the purchase order. Additional charges apply.

Keep the chair, upgrade the delivery. MSP2400-2.0 brings Mirage 2.0 asepsis delivery and a rear-radius assistant's arm to existing BLMT 7, 10, 20, 20N and 50 chairs. Specify the chair model when ordering.