

— PATIENT CHAIRS · MODEL 3000-2.0-CB

Mirage 2.0

Hydraulic Patient Chair with Cast Aluminum Base, 3000-2.0-CB

Hands-free foot control positioning with two programmable positions plus an entry/exit position and last position. Cast aluminum base and a wide, reinforced seat, without armrests.



01 FEATURES

- ✓ Synchronous backrest
- ✓ Articulating dual axis headrest
- ✓ Cast aluminum base
- ✓ Programmable foot control
- ✓ Long pump cover

Specifications

02 SPECIFICATIONS

Model	3000-2.0-CB
Lifting / lowering	Hydraulic, cantilever type
Initial height / stroke	15.5" / 15"
Backrest reclining angle	0° (horizontal) to 70°
Seat tilting	10° to 20°, synchronized with backrest
Preset positions	Microcomputer control; three preset positions, auto return, last position

Model 3000-2.0-CB

Chair rotation	30° left / 30° right of centerline
Headrest	Dual axis type
Safety shut-off plate	Under the main link
Emergency stop	Any switch can cancel chair movement
Input voltage	110V or 220V
Chair size (L x W x H, in)	44 x 29 x 34
Warranty	5 Year Limited

03 PRODUCT INCLUDES

- ✓ Two programmable positions
- ✓ One preset exit position
- ✓ Electronic foot control with auto return function
- ✓ Built-in fuse / surge protection
- ✓ Toe guard protection to prevent scuffing
- ✓ All-metal construction
- ✓ Safety shut-off base plate
- ✓ Double articulating headrest, adjustable for height and wheelchair positioning
- ✓ Chair top swivel with lock, 60 degree
- ✓ Seamless upholstery
- ✓ Wide chair back
- ✓ Coordinated leg rest tilt

04 OPTIONS

TP-2005	Touch pad control, chair mounted, 3 preset, auto return, last position
TPC-BRKT-2.0	Mirage M4000-2.0 post mount chair bracket
PC-1000-DU	Replacement dental chair upholstery, LK vinyl, wide back
PC-1000-Vinyl-UC	Upcharge: non-standard dental chair upholstery, vinyl, wide back
PC-1000-ULTR-UC	Upcharge: Ultra leather dental chair upholstery, wide back
PC-1000-PR-UC	Upcharge: NaguaSoft dental chair upholstery, vinyl, wide back (PR50, 49, 47, 40, 37)

Options are not factory-installed unless specified on the purchase order. Additional charges apply. Contact TPC for pricing.