

— HYDRAULIC PATIENT CHAIR

Mirage2.0

Smooth hydraulic cantilever lift, a seat that swivels 30° each way for open access to the oral cavity, and a choice of cast aluminum or steel base.



4000-2.0-CB shown

15in

Hydraulic lift stroke from a 15.5" initial height

0-70°

Backrest recline, synchronized seat tilt

30° L/R

Seat swivel with lock

5year

Limited warranty



Cast aluminum base. 4000-2.0-CB



Steel base, durable cover. 4000-2.0



Dual-axis headrest. Rotates for oral-cavity access without repositioning the chair.

Specifications

TPC

01 STANDARD FEATURES

- ✓ Hydraulic cantilever lift for smooth, quiet travel
- ✓ Two programmable positions plus preset exit position
- ✓ Electronic foot control with auto return
- ✓ Built-in fuse and surge protection
- ✓ Safety shut-off plate stops seat-down and recline
- ✓ All-metal construction, toe guard against scuffing
- ✓ Cast aluminum (4000-2.0-CB) or steel (4000-2.0) base
- ✓ Seat swivels 30° left or right, with lock
- ✓ Double-articulating dual-axis headrest, height and wheelchair adjustable
- ✓ Pivoting armrests for easy entry and exit
- ✓ Synchronous backrest with coordinated leg rest tilt
- ✓ Seamless, reinforced wide-back upholstery
- ✓ Long pump cover
- ✓ Input voltage 110V or 220V

02 CHAIR SPECIFICATIONS

Lift mechanism	Hydraulic, cantilever
Height / stroke	15.5" initial / 15" stroke
Backrest recline	0° (horizontal) to 70°
Seat tilt	10°–20°, synced with backrest
Positions	2 programmable + exit
Chair rotation	30° left / 30° right of center
Headrest	Dual-axis type
Safety plate	Under the main link
Emergency stop	Any switch stops motion
Warranty	5-year limited

03 OPTIONS

TP-2005	Touch pad: 3 presets, auto return, last position
TPC-BRKT-2.0	M4000-2.0 post-mount chair bracket
PC-1000-DU	Wide-back upholstery, LK vinyl
PC-1000-Vinyl-UC	Non-standard vinyl upholstery
PC-1000-ULTR-UC	Ultraleather upholstery
PC-1000-PR-UC	Nauga Soft vinyl (PR50, 49, 47, 40, 37)

04 PACKING DIMENSIONS

Carton	L × W × H in	Weight
Chair	44 × 29 × 34	340 lb
Cushion	51 × 24 × 10	35 lb
Chair bracket	11 × 5 × 25	35 lb

05 STANDARD UPHOLSTERY

									
Sand LK-03	Taupe LK-08	Stormy LK-04	Greystone LK-11	Orange LK-27	Cobalt LK-31	Desert Rose LK-34	Grey Teal LK-49	Atlantis LK-59	Black LK-72

Options are not factory-installed unless specified on the purchase order. Colors are approximate.