

— MIRAGE 2.0 ORTHODONTIC PACKAGE · MOP3000-550LED-2375 / CV / CV-LS

Mirage 2.0

Orthodontic Package, Steel Base

Orthodontic package with the 3000-2.0 hydraulic chair, OMC2375 mobile delivery cabinet and L550-LED Lustrous light, available centrally plumbed or self-contained.



01 FEATURES

- ✓ 3000-2.0 hydraulic chair with steel base
- ✓ OMC2375 mobile delivery cabinet with pull-out writing surface
- ✓ L550-LED Lustrous operatory light
- ✓ Central plumbing or self contained (compressor + vacuum)
- ✓ 1 year warranty

Specifications

02 MODELS

Model	Description
MOP3000-550LED-2375	Mirage 2.0 Orthodontic Package, OMC2375 cabinet, central plumbing with junction box
MOP3000-550LED-2375CV	Mirage 2.0 Orthodontic Package, OMC2375CV: cabinet + compressor + vacuum (self contained)
MOP3000-550LED-2375CV-LS	Mirage 2.0 Orthodontic Package, OMC2375CV-LS: cabinet + compressor + vacuum + LED curing light + scaler (self contained)

03 PACKAGE INCLUDES

Part	Description
3000-2.0	Hydraulic patient chair with steel base, wide back and full lumbar support, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 4 programmable positions with automatic return and foot control
OMC2375	Orthodontic mobile delivery cabinet with a pullout writing surface, two instrument drawers, one medium and two large storage drawers for materials and supplies that fully extend and quickly remove for easy cleaning, and convenient access from the front. Oversized locking dual-caster wheels, central plumbing with junction box.
L550-LED	Lustrous LED operatory light with Poliblock mode, 8,000 – 30,000 Lux, dual axis light head rotation

04 CHAIR FEATURES

- ✓ Synchronous backrest
- ✓ Articulating dual axis headrest
- ✓ Laser cut steel base
- ✓ Long pump cover
- ✓ Hydraulic pump for quiet operation
- ✓ Programmable foot control

05 LIGHT FEATURES

- ✓ Composite cure mode
- ✓ 24V DC power supply
- ✓ Bulb type: 4 white LEDs, 2 yellow
- ✓ Light intensity (Lux): high 30,000, low 8,000
- ✓ Focal length: 27" (70 cm)
- ✓ Illumination area: 100mm x 180mm
- ✓ Sensor mode, manual toggle on/off

06 CONTROL FEATURES

- ✓ 3 HP automatic selection
- ✓ 3-way syringe
- ✓ Vacuum package with saliva ejector and HVE
- ✓ Self contained water bottle system
- ✓ Variable speed disc-type foot control
- ✓ Water coolant on/off switch with flow control
- ✓ 1 year warranty
- ✓ Requires central air / vacuum systems. Includes utility center and junction box.
- ✓ Optional: self contained compressor and vacuum (no plumbing required)

07 PACKING DIMENSIONS

Chair: 44 x 29 x 33 / 340 lbs.

Cushion: 51 x 24 x 10 / 35 lbs.

Light arm L600: 26 x 14 x 8 / 30 lbs.; light head L570: 17 x 15 x 7

Cabinet: 32 x 28 x 41 / 140 lbs.

08 OPTIONS AND PARTS

Part	Description
WIC840	Mirage flat screen monitor bracket with 16" extension
M017	17" LCD multi media flat screen monitor
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply
ALED55BN-A	LED 5W curing light – built in to dental unit

Options are not factory-installed unless specified on the purchase order. Contact TPC for pricing.

Part	Description
TP-2005	Touch pad control – dental chair mounted 3 preset, auto return, last position
OMC2200	LCD monitor bracket
OMC-ST81	Optional Stone Surface on OMC 2375, Color Artic
OMC-ST82	Optional Stone Surface on OMC 2375, Color White
OMC-ST85	Optional Stone Surface on OMC 2375, Color Cream
OMC-ST86	Optional Stone Surface on OMC 2375, Color Bamboo
M4000-09	Optional Magnetic Headrest
AS-1101	Mirage Assistants stool
DR-1102	Mirage Doctor's stool
PTH-101	Pony Tail Headrest