

PATIENT CHAIRS · MODEL 3000-2.0

Mirage 2.0

Orthodontic Patient Chair, 3000-2.0

A hydraulic patient chair with steel base, wide back and full lumbar support, double articulating headrest with sliding height adjustment, and 4 programmable positions with automatic return. Shown here as supplied with the Mirage 2.0 Orthodontic Package.



01 FEATURES

- ✓ Synchronous backrest
- ✓ Articulating dual axis headrest
- ✓ Laser cut steel base
- ✓ Long pump cover
- ✓ Hydraulic pump for quiet operation
- ✓ Programmable foot control

Specifications

02 SPECIFICATIONS

Model	3000-2.0
Base	Laser cut steel base, hydraulic pump
Backrest	Synchronous backrest, wide back with full lumbar support
Headrest	Double articulating headrest with sliding height adjustment; optional pony tail headrest
Chair swivel	30° left and right of center

Model 3000-2.0

Positions	4 programmable positions with automatic return and foot control
-----------	---

Chair size (L x W x H, in)	44 x 29 x 33
----------------------------	--------------

03 PRODUCT INCLUDES

- | | |
|---|--|
| ✓ Hydraulic patient chair with steel base | ✓ 30 degree swivel from left and right of center |
| ✓ Wide back with full lumbar support | ✓ 4 programmable positions with automatic return |
| ✓ Double articulating headrest with sliding height adjustment | ✓ Programmable foot control |

04 OPTIONS

TP-2005	Touch pad control, chair mounted, 3 preset, auto return, last position
---------	--

M4000-09	Optional magnetic headrest
----------	----------------------------

PTH 101	Pony tail headrest
---------	--------------------

Options are not factory-installed unless specified on the purchase order. Additional charges apply. Contact TPC for pricing.